



PCR clean-up Gel extraction

User Manual

NucleoTrap[®]

NucleoTrap[®]CR

July 2010/Rev.04

PCR clean-up, Gel extraction

Protocol-at-a-glance (Rev. 04)















	NucleoTraP®CR PCR clean-up		NucleoTrap® Gel extraction											
1 NucleoTrap®: Excise DNA fragment / Solubilize gel slice NucleoTraP®CR: Adjust binding conditions	4 vol NT2 / 1 vol sample 	 300 µl NT1 / 100 mg gel												
2 Bind DNA	10 µl silica matrix / 100 µl sample RT 10 min 10,000 x g 30 s  	4 µl silica matrix / µg DNA 50°C 5 – 10 min 10,000 x g 30 s  												
3 Wash silica matrix	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">1st</td> <td>400 µl NT2</td> </tr> <tr> <td>2nd</td> <td>400 µl NT3</td> </tr> <tr> <td>3rd</td> <td>400 µl NT3</td> </tr> </table> 10,000 x g 30 s 10,000 x g 30 s 10,000 x g 30 s  	1st	400 µl NT2	2nd	400 µl NT3	3rd	400 µl NT3	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">1st</td> <td>500 µl NT2</td> </tr> <tr> <td>2nd</td> <td>500 µl NT3</td> </tr> <tr> <td>3rd</td> <td>500 µl NT3</td> </tr> </table> 10,000 x g 30 s 10,000 x g 30 s 10,000 x g 30 s  	1st	500 µl NT2	2nd	500 µl NT3	3rd	500 µl NT3
1st	400 µl NT2													
2nd	400 µl NT3													
3rd	400 µl NT3													
1st	500 µl NT2													
2nd	500 µl NT3													
3rd	500 µl NT3													
4 Dry silica matrix	RT or 37°C 10 – 15 min	RT or 37°C 10 – 15 min												
5 Elute DNA	25 – 50 µl NE RT 10 – 15 min 10,000 x g 30 s  	25 – 50 µl NE RT 10 – 15 min 10,000 x g 30 s  												

Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Consumables and equipment to be supplied by the user	5
1.3	About this User Manual	5
2	Product description	6
2.1	The basic principle	6
2.2	Kit specifications	6
2.3	Elution procedures	7
3	Storage conditions and preparation of working solutions	8
4	Safety instructions – risk and safety phrases	9
5	Protocol for direct purification of PCR products using NucleoTraP®CR	10
6	Protocol for DNA extraction from agarose gels using NucleoTrap®	12
7	Support protocol for concentration, desalination, removal of enzymes, etc.	14
8	Appendix	15
8.1	Troubleshooting	15
8.2	Ordering information	17
8.3	References	17
8.4	Product use restriction/warranty	18

1 Components

1.1 Kit contents

NucleoTraP®CR		
Cat. No.	10 preps 740587.10	100 preps 740587
NucleoTraP®CR Suspension	100 µl	1000 µl
Buffer NT2	10 ml	2 x 50 ml
Wash Buffer NT3 (Concentrate)*	4 ml	20 ml
Elution Buffer NE**	5 ml	15 ml
User Manual	1	1

NucleoTrap®		
Cat. No.	10 preps 740584.10	100 preps 740584
NucleoTrap® Suspension	100 µl	1000 µl
Buffer NT1	6 ml	2 x 30 ml
Buffer NT2	10 ml	2 x 50 ml
Wash Buffer NT3 (Concentrate)*	4 ml	20 ml
Elution Buffer NE**	5 ml	15 ml
User Manual	1	1

* For preparation of working solutions and storage conditions see section 3.

**Composition of Elution Buffer NE: 5 mM Tris/HCl, pH 8.5

1.2 Consumables and equipment to be supplied by the user

Consumables:

- 96 – 100% ethanol
- 1.5 ml microcentrifuge tubes

Equipment:

- Centrifuge for microcentrifuge tubes
- Manual pipettors and disposable tips
- Vortex mixer
- Heating-block
- Personal protection equipment (lab coat, gloves, goggles)

1.3 About this User Manual

It is strongly recommended that first-time users of the **NucleoTrap®CR/NucleoTrap®** kit read the detailed protocol sections of this User Manual. Experienced users, however, may refer to the Protocol-at-a-glance instead. The Protocol-at-a-glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure.

All technical literature is available on the internet at www.mn-net.com.

2 Product description

2.1 The basic principle

With the **NucleoTraP[®]CR/Trap[®]** method, DNA binds in the presence of chaotropic salts (Buffer NT1 and Buffer NT2) to specially activated silica particles (matrix). Buffer NT1 contains additional components in order to dissolve agarose gel slices. Afterwards, the **NucleoTraP[®]CR/Trap[®]** matrix is added to the binding mixtures. Contaminations like salts and soluble macromolecular components are removed by a simple washing step with ethanolic Wash Buffer NT3. Pure DNA is finally eluted under low ionic strength conditions with slightly alkaline Elution Buffer NE (5 mM Tris-Cl, pH 8.5).

2.2 Kit specifications

- The **NucleoTraP[®]CR** kit is designed for direct purification of PCR products.
- The **NucleoTrap[®]** kit is designed for the purification of DNA from TAE/TBE agarose gels.
- In contrast to the **NucleoTrap[®]** matrix, the **NucleoTraP[®]CR** matrix will not bind DNA fragments < 100 bp due to a larger pore size of the silica matrix.
- Standard as well as low melting agarose gels can be used.
- The prepared DNA fragments can be used directly in applications like automated fluorescent DNA sequencing, PCR, or any kind of enzymatic manipulation.

Table 1: Kit specifications at a glance

Parameter	NucleoTraP [®] CR	NucleoTrap [®]
DNA fragments from agarose gels	-	++
Desalination, removal of enzymes, nucleotides and/or labeling reagents like biotin or radioactive ATP etc.	++	+
Direct purification of amplified DNA	++	-
Elution volume	20 – 50 µl	20 – 50 µl
Binding capacity	6 µg/10 µl matrix	6 µg/10 µl matrix
Time/prep	45 min for 6 preps	60 min for 6 preps
- not recommended	+ possible	++ optimal

2.3 Elution procedures

- For the elution of DNA one of the following solutions can be used: Buffer NE (supplied) / TE buffer, pH 8.5 / distilled water, pH 8.5.
- If water is used, the pH should be checked and adjusted to pH 8-8.5 since deionized water usually exhibits a pH below 7. Furthermore, absorption of CO₂ leads to a decrease in pH of unbuffered solutions.
- Note: EDTA in TE buffer may cause problems in subsequent reactions. See Table 2 for correlation between fragment size and typical recoveries for purification of 1 - 5 µg of PCR fragments (for gel extraction, recoveries are approximately 10% lower).

Table 2: DNA recovery with NucleoTraP®CR/NucleoTrap®

Fragment length	NucleoTraP®CR	NucleoTrap®
20 bp	0%	50%
40 bp	0%	68%
120 bp	68%	78%
200 bp	76%	85%
520 bp	80%	87%
2.5 kbp	81%	88%
5.3 kbp	80%	86%
8.7 kbp	76%	80%
19.4 kbp	74%	74%

3 Storage conditions and preparation of working solutions

Attention:

Buffers NT1 and NT2 contain chaotropic salts! Wear gloves and goggles!

- The **NucleoTraP®CR/NucleoTrap®** kits should be stored at room temperature and are stable for up to one year.

Before starting any **NucleoTraP®CR/NucleoTrap®** protocol prepare the following:

- **Wash Buffer NT3:** Add the indicated volume of 96 – 100% ethanol to **Wash Buffer NT3 Concentrate**.





NucleoTraP®CR		
Cat. No.	10 preps 740587.10	100 preps 740587
Wash Buffer NT3 (Concentrate)	4 ml Add 16 ml ethanol	20 ml Add 80 ml ethanol

NucleoTrap®		
Cat. No.	10 preps 740584.10	100 preps 740584
Wash Buffer NT3 (Concentrate)	4 ml Add 16 ml ethanol	20 ml Add 80 ml ethanol

4 Safety instructions – risk and safety phrases

The following components of the **NucleoTraP®CR/NucleoTrap®** kits contain hazardous contents.

Wear gloves and goggles and follow the safety instructions given in this section.

Component	Hazard contents	Hazard symbol		Risk phrases	Safety phrases
NT1	Sodium perchlorate	 O**  Xn*	Explosive when mixed with combustible material - Harmful if swallowed	R 9-22	S 13-27
NT2	Sodium perchlorate	 O**  Xn*	Explosive when mixed with combustible material - Harmful if swallowed	R 9-22	S 13-27

Risk phrases

- R 9 Explosive when mixed with combustible material
R 22 Harmful if swallowed

Safety phrases

- S 13 Keep away from food, drink, and animal feedstuffs
S 27 Take off immediately all contaminated clothing

* Hazard labeling not necessary if quantity per bottle below 125 g or ml (certificate of exemption according to 67/548/EEC Art. 25, 1999/45/EC Art. 12 and German GefStoffV § 20 (3) and TRGS 200 7.1). For further information see Material Safety Data Sheet.

**Hazard labeling not necessary if quantity per bottle below 50 g or ml (certificate of exemption according to 67/548/EEC Art. 25, 1999/45/EC Art. 12 and German GefStoffV § 20 (3) and TRGS 200 7.1). For further information see Material Safety Data Sheet.

5 Protocol for direct purification of PCR products using NucleoTraP[®]CR

Before starting the preparation:

- Check if Wash Buffer NT3 was prepared according to section 3.

1 Adjust DNA binding conditions

Add **4 volumes** of **Buffer NT2** to **1 volume** of **sample** (e.g., 400 μ l Buffer NT2 and 100 μ l PCR reaction mixture).

For sample volumes < 100 μ l adjust the volume of the reaction mix to 100 μ l using TE buffer (pH 7.5).

Note: If the volume of the PCR reaction mixture is > 100 μ l, the volumes of Buffer NT2 and NucleoTraP[®]CR Suspension must be increased proportionally. Example: a volume of 150 μ l reaction mixture needs 600 μ l of Buffer NT2, and 15 μ l NucleoTraP[®]CR Suspension to adjust proper binding conditions.



**+ 4 vol NT2
per
1 vol sample**

2 Bind DNA

Vortex the NucleoTraP[®]CR Suspension thoroughly until a homogeneous mixture results. Add **10 μ l** of **NucleoTraP[®]CR Suspension** to each **100 μ l** of **reaction mixture**.

Incubate the mixture for **10 min** at **room temperature** and vortex briefly every 2 – 3 min.

Centrifuge the sample at **10,000 x g** for **30 s** and discard the supernatant.



**10 μ l silica
matrix**

**RT
10 min**



**10,000 x g
30 s**

3 Wash silica matrix

1st wash

Add **400 μ l Buffer NT2** to the pelleted silica matrix and vortex briefly for resuspension of the pellet. Centrifuge for **30 s** at **10,000 x g** and remove the supernatant completely.



+ 400 μ l NT2



**10,000 x g
30 s**

2nd wash

Add **400 µl Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 s** at **10,000 x g** and remove the supernatant completely.



+ 400 µl NT3

10,000 x g
30 s
3rd wash

Add **400 µl Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 s** at **10,000 x g**. Remove the supernatant and centrifuge the pellet again briefly. Remove residual Buffer NT3 completely.



+ 400 µl NT3

10,000 x g
30 s
4 Dry silica matrix

Dry the pelleted silica matrix at **room temperature** or at **37°C** for **10 – 15 min**.

It is not recommended to dry the sample by vacuum since over-dried pellets lead to lower recoveries.

Residual ethanol from Buffer NT3 would inhibit subsequent reactions and has to be removed in this step.

RT or 37°C
10 – 15 min
5 Elute DNA

Add **25-50 µl Buffer NE** to the silica matrix. Resuspend the pellet by vortexing. Incubate the mixture at **room temperature** for **10-15 min**. Vortexing the mixture 2 – 3 times during incubation is recommended. Centrifuge the sample at **10,000 x g** for **30 s** and transfer the DNA containing supernatant to a clean tube (not provided). Repeating this step will increase the yield by approximately 10%.

+ 25 – 50 µl
NE
10,000 x g
1 min

Yield of larger fragments (>5 – 20 kbp) can be increased by performing the incubation at 55°C.


6 Protocol for DNA extraction from agarose gels using NucleoTrap®

Before starting the preparation:

- Check if Wash Buffer NT3 was prepared according to section 3.
- Set heating block to 50°C.

1 Excise DNA fragment / Solubilize gel slice

Take a clean scalpel to excise the DNA fragment from agarose gel. Excise gel slice containing the fragment carefully to minimize the gel volume. Determine the weight of the gel slice and transfer it to a clean tube (not provided).



For each **100 mg** agarose gel add **300 µl NT1**.

For gels containing >2% agarose, double the volume of Buffer NT1.

Note: If the weight of the gel slice is >100 mg, the volume of Buffer NT1 must be increased proportionally. Example: a 150 mg gel slice (<2% agarose) needs 450 µl Buffer NT1.



+ 300 µl NT1

2 Bind DNA

Vortex the NucleoTrap® Suspension thoroughly until a homogeneous mixture results. **For each µg of DNA** add **4 µl** of the **NucleoTrap® Suspension**, but at least 10 µl.



4 µl silica matrix / µg DNA

Incubate sample at 50°C until the gel slice is dissolved (5 – 10 min). Vortex the sample briefly every 2 – 3 min until the gel slice is dissolved completely.

**50°C
5 – 10 min**

Centrifuge for **30 s** at **10,000 x g** and discard supernatant.



**10,000 x g
30 s**

3 Wash silica matrix

1st wash

Add **500 µl Buffer NT2** to the pelleted silica matrix and vortex briefly for resuspension of the pellet. Centrifuge for **30 s** at **10,000 x g** and remove the supernatant completely.



+ 500 µl NT2



**10,000 x g
30 s**

2nd wash

Add **500 µl Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 s** at **10,000 x g** and remove the supernatant completely.



+ 500 µl NT3

10,000 x g
30 s**3rd wash**

Add **500 µl Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 s** at **10,000 x g**. Remove the supernatant and centrifuge the pellet again briefly. Remove residual Buffer NT3 completely.



+ 500 µl NT3

10,000 x g
30 s**4 Dry silica matrix**

Dry the pelleted silica matrix at **room temperature** or at **37°C** for **10 – 15 min**.

It is not recommended to dry the sample by vacuum since over-dried pellets lead to lower recoveries.

Residual ethanol from Buffer NT3 would inhibit subsequent reactions and has to be removed in this step.

RT or 37°C
10 – 15 min**5 Elute DNA**

Add **25 – 50 µl Buffer NE** to the silica matrix. Resuspend the pellet by vortexing. Incubate the mixture at **room temperature** for **10 – 15 min**. Vortexing the mixture 2 – 3 times during incubation is recommended. Centrifuge the sample at **10,000 x g** for **30 s** and transfer the DNA containing supernatant to a clean tube (not provided). Repeating this step will increase the yield by approximately 10%.



+ 20 – 50 µl NE

10,000 x g
30 s

Yield of larger fragments (>5 – 20 kbp) can be increased by performing the incubation at 55°C.

7 Support protocol for concentration, desalination, removal of enzymes, etc.

Before starting the preparation:

- Check if Wash Buffer NT3 was prepared according to section 3.

1 Adjust DNA binding conditions

Add **4 volumes Buffer NT2** to **1 volume** of DNA containing **sample** (e.g., 400 µl Buffer NT2 and 100 µl reaction mixture).



**+ 4 vol NT2
per
1 vol sample**

2 Bind DNA

Vortex the NucleoTraP®CR/NucleoTrap® Suspension thoroughly until a homogeneous mixture results. For **each µg of DNA** add **4 µl of silica matrix** but at least 10 µl.



**4 µl silica
matrix /
µg DNA**

Incubate the mixture for **10 min** at **room temperature** and vortex briefly every 2 – 3 min.

**RT
10 min**

Centrifuge for **30 s** at **10,000 x g** and discard supernatant.



**10,000 x g
30 s**

Important note: Be aware of the NucleoTrap® Suspension binding fragments down to 20 bp (see Table 2, section 2.3).

Continue with **section 5, step 3**.

8 Appendix

8.1 Troubleshooting

Problem	Possible cause and suggestions
Incomplete lysis of agarose slices	<p data-bbox="309 347 613 367"><i>High concentration of agarose</i></p> <ul data-bbox="309 379 981 432" style="list-style-type: none"> <li data-bbox="309 379 981 432">• Use doubled volumes of Buffer NT1 for highly concentrated agarose gels.
	<p data-bbox="309 472 443 491"><i>Wrong buffer</i></p> <ul data-bbox="309 504 807 528" style="list-style-type: none"> <li data-bbox="309 504 807 528">• Buffer NT2 cannot be used for gel dissolution.
	<p data-bbox="309 568 533 587"><i>Time and temperature</i></p> <ul data-bbox="309 600 981 735" style="list-style-type: none"> <li data-bbox="309 600 981 735">• Check incubation temperature. Depending on the weight of gel slice, incubation (section 6, step 2) can be prolonged up to 20 min. Vortex every 2 min and check integrity of the gel slice. Heavy weight gel slices may be quenched or crushed before addition of Buffer NT1.
No DNA yield	<p data-bbox="309 775 613 794"><i>Reagents not applied properly</i></p> <ul data-bbox="309 807 981 860" style="list-style-type: none"> <li data-bbox="309 807 981 860">• Add indicated volume of 96 – 100% ethanol to Wash Buffer NT3 Concentrate and mix well before use.
	<p data-bbox="309 900 981 919"><i>Insufficient drying of the NucleoTraP®CR/NucleoTrap® silica matrix</i></p> <ul data-bbox="309 932 981 1038" style="list-style-type: none"> <li data-bbox="309 932 981 1038">• Ethanolic Wash Buffer NT3 has to be removed quantitatively before elution. Prolong the drying time up to 30 min. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots).
	<p data-bbox="309 1078 636 1098"><i>Isolation of large DNA fragments</i></p> <ul data-bbox="309 1110 981 1163" style="list-style-type: none"> <li data-bbox="309 1110 981 1163">• Add room-temperature Elution Buffer NE and incubate at 55°C for 10 – 15 min.

Problem	Possible cause and suggestions
Suboptimal performance of DNA in sequencing reactions	<i>Carry-over of ethanol/ethanolic Buffer NT3</i>
	<ul style="list-style-type: none">• Make sure to dry the silica matrix in order to achieve complete removal of ethanolic Buffer NT3 after the washing step. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots).• Buffers other than Buffer NE, for example TE buffer (Tris/EDTA), were used for elution of DNA. Note: EDTA may inhibit sequencing reactions. In this case it is recommended to re-purify DNA and elute in Buffer NE or water.
	<i>Not enough DNA used for sequencing reaction</i>
	<ul style="list-style-type: none">• Quantitate DNA by agarose gel electrophoresis before setting up sequencing reactions.
	<i>NucleoTraP®CR or NucleoTrap® particles were not removed quantitatively</i>
	<ul style="list-style-type: none">• Centrifuge the eluate again and transfer the supernatant to a new tube.

8.2 Ordering information

Product	Cat. No.	Pack of
NucleoTraP®CR	740587.10	10 preps
	740587	100 preps
NucleoTrap®	740584.10	10 preps
	740584	100 preps
NucleoTraP®CR Suspension	740564	100 preps
NucleoTrap® Suspension	740569	100 preps
Buffer NT1	740596.100	2 x 50 ml
Buffer NT2	740597	2 x 50 ml
Buffer NT3 Concentrate (for 100 ml Buffer NT3)	740598	20 ml
Collection Tubes (2 ml)	740600	1000

Visit www.mn-net.com for more detailed product information.

8.3 References

Vogelstein B., and D. Gillespie. 1979. Preparative and analytical purification of DNA from agarose. Proc. Natl. Acad. Sci. USA 76: 615-619.

8.4 Product use restriction/warranty

NucleoTraP®CR/NucleoTrap® kit components were developed, designed, distributed, and sold **FOR RESEARCH PURPOSES ONLY**. They are suitable **FOR IN-VITRO USES ONLY**. No claim or representation is intended for its use to identify any specific organism or for clinical use (diagnostic, prognostic, therapeutic, or blood banking).

It is rather the responsibility of the user to verify the use of the **NucleoTraP®CR/NucleoTrap®** kit for a specific application range as the performance characteristic of this kit has not been verified to a specific organism.

This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish an extra copy.

MACHEREY-NAGEL does not warrant against damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; against defects in products or components not manufactured by MACHEREY-NAGEL, or against damages resulting from such non-MACHEREY-NAGEL components or products.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, REPRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications. This warranty is exclusive and MACHEREY-NAGEL makes no other warranty expressed or implied.

The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications.

Please contact:
MACHEREY-NAGEL Germany
Tel.: +49 (0) 24 21 969 270
e-mail: TECH-BIO@mn-net.com

Last updated: 12/2006, Rev.02

Patents/Trademarks:

PCR is patented by Roche Diagnostics
NucleoTrap is a trademark of MACHEREY-NAGEL GmbH & Co KG

All used names and denotations can be brands, trademarks or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information, i.e. it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment. Regarding these products or services we can not grant any guarantees regarding selection, efficiency or operation.