

PCR clean-up Gel extraction

User Manual

NucleoTrap[®]

NucleoTraP[®]CR

January 2008/Rev. 03

Protocol-at-a-glance (Rev. 03)

NucleoTrap[®]/NucleoTraP[®]CR





	Gel extraction NucleoTrap [®]	PCR clean-up NucleoTraP [®] CR
1 NucleoTrap[®] - Excise DNA fragment / Solubilize gel slice NucleoTraP[®]CR - Adjust binding conditions	 300 µl NT1 / 100 mg	 4 vol. NT2 / 1 vol. sample
2 Bind DNA	4 µl silica matrix / µg DNA 50°C 5-10 min 30 sec 10,000 x g	10 µl silica matrix / 100 µl sample RT 10 min 30 sec 10,000 x g
3 Wash silica matrix	1st 500 µl NT2 2nd 500 µl NT3 3rd 500 µl NT3 30 sec, 10,000 x g 30 sec, 10,000 x g 30 sec, 10,000 x g	1st 400 µl NT2 2nd 400 µl NT3 3rd 400 µl NT3 30 sec, 10,000 x g 30 sec, 10,000 x g 30 sec, 10,000 x g
4 Dry silica matrix	RT or 37°C 10-15 min	RT or 37°C 10-15 min
5 Elute DNA	25-50 µl NE RT 10-15 min 30 sec 10,000 x g	25-50 µl NE RT 10-15 min 30 sec 10,000 x g

Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Consumables and equipment to be supplied by user	5
2	Product description	6
2.1	The basic principle	6
2.2	About this User Manual	6
2.3	Kit specifications	6
2.4	Elution procedures	7
3	Storage conditions and preparation of working solutions	9
4	Safety instructions – risk and safety phrases	10
5	Protocol for DNA extraction from agarose gels using the NucleoTrap [®] kit	11
6	Protocol for direct purification of PCR products using the NucleoTraP [®] CR kit	13
7	Support protocol for concentration, desalination, removal of enzymes etc.	15
8	Appendix	16
8.1	Troubleshooting	16
8.2	Ordering information	18
8.3	References	18
8.4	Product use restriction / warranty	18

1 Components

1.1 Kit contents

NucleoTrap®		
Cat. No.	10 preps 740584.10	100 preps 740584
NucleoTrap® Suspension	100 µl	1000 µl
Binding Buffer NT1	6 ml	2 x 30 ml
Binding Buffer NT2	10 ml	2 x 50 ml
Wash Buffer NT3 (Concentrate)*	4 ml	20 ml
Elution Buffer NE**	5 ml	15 ml
User Manual	1	1

NucleoTraP® CR		
Cat. No.	10 preps 740587.10	100 preps 740587
NucleoTraP® CR Suspension	100 µl	1000 µl
Binding Buffer NT2	10 ml	2 x 50 ml
Wash Buffer NT3 (Concentrate)*	4 ml	20 ml
Elution Buffer NE**	5 ml	15 ml
User Manual	1	1

* For preparation of working solutions and storage conditions see section 3.

** Composition of Elution Buffer NE: 5 mM Tris/HCl, pH 8.5

1.2 Consumables and equipment to be supplied by user

Consumables:

- 96-100% ethanol
- 1.5 ml microcentrifuge tubes

Equipment:

- Centrifuge for microcentrifuge tubes
- Manual pipettors and disposable tips
- Vortex mixer
- Heating-block

Personal protection equipment (lab coat, gloves, goggles)

2 Product description

2.1 The basic principle

With the **NucleoTrap®/TraP®CR** method, DNA binds in the presence of chaotropic salts (Buffer NT1 and Buffer NT2) to specially activated silica particles (matrix). Buffer NT1 contains additional components in order to dissolve agarose gel slices. Afterwards, the **NucleoTrap®/TraP®CR** matrix is added to the binding mixtures. Contaminations like salts and soluble macromolecular components are removed by a simple washing step with ethanolic Wash Buffer NT3. Pure DNA is finally eluted under low ionic strength conditions with slightly alkaline Elution Buffer NE (5 mM Tris-Cl, pH 8.5).

2.2 About this User Manual

Experienced users who are performing the purification of PCR* products or DNA fragments from agarose gels using **NucleoTrap®** or **NucleoTraP®CR** isolation kits may refer to the Protocol-at-a-glance instead of this User Manual. The Protocol-at-a-glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure. First-time users are strongly advised to read this User Manual.

2.3 Kit specifications

- The **NucleoTrap®** kit is designed for the purification of DNA from TAE/TBE agarose gels.
- The **NucleoTraP®CR** kit is designed for direct purification of PCR* products.
- In contrast to the **NucleoTrap®** matrix, the **NucleoTraP®CR** matrix will not bind DNA fragments < 100 bp due to a larger pore size of the silica matrix.
- Standard as well as low melting agarose gels can be used.
- The prepared DNA fragments can be used directly in applications like automated fluorescent DNA sequencing, PCR*, or any kind of enzymatic manipulation.

* PCR is patented by Roche Diagnostics

Kit specification at a glance		
Parameters	NucleoTrap [®] CR	NucleoTrap [®]
DNA fragments from agarose gels	-	++
Desalination, removal of enzymes, nucleotides and /or labeling reagents like biotin or radioactive ATP etc.	++	+
Direct purification of amplified DNA	++	-
Elution volume	20-50 µl	20-50 µl
Binding capacity	6 µg/10 µl matrix	6 µg/10 µl matrix
Time/prep	45 min for 6 preps	60 min for 6 preps

- not recommended

+ possible

++ optimal

2.4 Elution procedures

- For the elution of DNA one of the following solutions can be used: Buffer NE (supplied) / TE buffer, pH 8.5 / distilled water, pH 8.5.
- If water is used, the pH should be checked and adjusted to pH 8-8.5 since deionized water usually exhibits a pH below 7. Furthermore, absorption of CO₂ leads to a decrease in pH of unbuffered solutions.
- Note: EDTA in TE buffer may cause problems in subsequent reactions. See Table 1 for correlation between dispensed elution buffer volumes and typical recoveries for purification of 1-5 µg of PCR* fragments (for gel extraction, recoveries are approximately 10% lower).

Table 1: DNA recovery with NucleoTrap[®]/NucleoTraP[®]CR

Fragment length	NucleoTrap[®]	NucleoTraP[®]CR
20 bp	50 %	0 %
40 bp	68 %	0 %
120 bp	78 %	68 %
200 bp	85 %	76 %
520 bp	87 %	80 %
2.5 kbp	88 %	81 %
5.3 kbp	86 %	80 %
8.7 kbp	80 %	76 %
19.4 kbp	74 %	74 %

3 Storage conditions and preparation of working solutions

Attention:

Buffers NT1 and NT2 contain chaotropic salts. Wear gloves and goggles!

- The **NucleoTrap®/NucleoTraP®CR** kits should be stored at room temperature and are stable for up to one year.

Before starting any **NucleoTrap®/NucleoTraP®CR** protocol prepare the following:

- Add the indicated volume of 96-100% ethanol to Wash Buffer NT3 Concentrate.





NucleoTrap®		
	10 preps	100 preps
Cat. No.	740584.10	740584
Wash Buffer NT3 (Concentrate)	4 ml add 16 ml ethanol	20 ml add 80 ml ethanol

NucleoTraP® CR		
	10 preps	100 preps
Cat. No.	740587.10	740587
Wash Buffer NT3 (Concentrate)	4 ml add 16 ml ethanol	20 ml add 80 ml ethanol

4 Safety instructions – risk and safety phrases

The following components of the NucleoTrap[®]/NucleoTraP[®]CR kits contain hazardous contents.

Wear gloves and goggles and follow the safety instructions given in this section

Buffer	Hazard Contents	Hazard Symbol	Risk Phrases	Safety Phrases
NT1	Sodium perchlorate	 O*  Xn	Explosive when mixed with combustible material. Harmful if swallowed	R 9-22 S 13-27
NT2	Sodium perchlorate	 O*  Xn	Explosive when mixed with combustible material. Harmful if swallowed	R 9-22 S 13-27

Risk Phrases

R 9 Explosive when mixed with combustible material

R 22 Harmful if swallowed

Safety Phrases

S 13 Keep away from food, drink and animal feedstuffs


S 27 Take off immediately all contaminated clothing

* Label not necessary, if quantity below 50 g or ml (according to 67/548/EEC Art. 25, 1999/45/EC Art. 12 and German GefStoffV § 42 and TRGS 200 7.1)

5 Protocol for DNA extraction from agarose gels using the NucleoTrap® kit

1 Excise DNA fragment / Solubilize gel slice


Take a clean scalpel to excise the DNA fragment from agarose gel. Excise gel slice containing the fragment carefully to minimize the gel volume. Determine the weight of the gel slice and transfer it to a clean tube (not provided).



For each **100 mg** agarose gel add **300 µl NT1**.

For gels containing > 2% agarose, double the volume of Buffer NT1.


Remark: If the weight of the gel slice is > 100 mg, the volume of Buffer NT1 must be increased proportionally. Example: a 150 mg gel slice (< 2% agarose) needs 450 µl Buffer NT1.



+ 300 µl NT1

2 Bind DNA

Vortex the NucleoTrap® Suspension thoroughly, resulting in a homogeneous mixture. For each **µg of DNA** add **4 µl** of the **NucleoTrap® Suspension**, but at least **10 µl**.



4 µl silica matrix / µg DNA

Incubate sample at **50°C** until the gel slices are dissolved (**5-10 min**). Vortex the sample briefly every 2-3 min until the gel slices are dissolved completely.

**50°C
5-10 min**

Centrifuge for **30 sec** at **10,000 x g** and discard supernatant.




**30 sec
10,000 x g**

3 Wash silica matrix

1st wash

Add **500 µl Buffer NT2** to the pelleted silica matrix and vortex briefly for resuspension of the pellet. Centrifuge for **30 sec** at **10,000 x g** and remove the supernatant completely.




+ 500 µl NT2




**30 sec
10,000 x g**

2nd wash

Add **500 µl Buffer NT3** and vortex briefly. Centrifuge for **30 sec** at **10,000 x g** and remove the supernatant completely.

**+ 500 µl NT3****30 sec**
10,000 x g**3rd wash**

Add **500 µl buffer NT3** and vortex briefly. Centrifuge for **30 sec** at **10,000 x g**. Remove the supernatant and centrifuge the pellet again briefly. Remove residual Buffer NT3 completely.

**+ 500 µl NT3****30 sec**
10,000 x g**4 Dry silica matrix**

Dry the pelleted silica matrix at **room temperature** or at **37°C** for **10-15 min**.

It is not recommended to dry the sample by vacuum since over-dried pellets lead to lower recoveries.

Residual ethanol from Buffer NT3 would inhibit subsequent reactions and has to be removed in this step.

37°C
10-15 min**5 Elute DNA**

Add **25-50 µl Elution Buffer NE** to the silica matrix. Resuspend the pellet by vortexing. Incubate the mixture at **room temperature** for **10-15 min**. Vortexing the mixture 2-3 times during incubation is recommended. Centrifuge the sample for **30 sec** at **10,000 x g** and transfer the DNA containing supernatant to a clean tube (not provided). Repeating this step will increase the yield by approximately 10%.

Yield of larger fragments (> 5-10 kbp) can be increased by performing incubation at 55°C.

**+ 25-50 µl NE****RT**
10 min**30 sec**
10,000 x g


6 Protocol for direct purification of PCR products using the NucleoTraP®CR kit

1 Adjust DNA binding conditions

Mix **4 volumes** of **Buffer NT2** with **1 volume** of **sample** (e.g. 400 µl Buffer NT2 and 100 µl PCR reaction mixture).

For sample volumes < 100 µl adjust the volume of the reaction mix to 100 µl using TE buffer (pH 7.5).


Remark: If the volume of the PCR reaction mixture is > 100 µl, the volumes of Buffer NT2 and NucleoTraP®CR Suspension must be increased proportionally. Example: a volume of 150 µl reaction mixture needs 600 µl of Buffer NT2, and 15 µl NucleoTraP®CR Suspension to adjust proper binding conditions.



**mix 4 vol. NT2
with
1 vol. sample**

2 Bind DNA

Vortex the NucleoTraP®CR Suspension thoroughly, resulting in a homogeneous mixture. Add **10 µl of NucleoTraP®CR Suspension** to each **100 µl of reaction mixture**. Incubate the mixture for **10 min** at **room temperature** and vortex briefly every 2-3 min.



**10 µl silica
matrix

RT
10 min**

Centrifuge the sample at **10,000 x g** for **30 sec** and discard the supernatant.




**30 sec
10,000 x g**

3 Wash silica matrix

1st wash

Add **400 µl Buffer NT2** to the pelleted silica matrix and vortex briefly for resuspension of the pellet. Centrifuge the sample for **30 sec** at **10,000 x g** and remove the supernatant completely.




400 µl NT2



**30 sec
10,000 x g**

2nd wash

Add **400 µl Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 sec** at **10,000 x g** and remove the supernatant completely.




400 µl NT3



30 sec
10,000 x g

3rd wash

Add **400 µl Buffer NT3** to the pelleted silica matrix, vortex briefly, and centrifuge for **30 sec** at **10,000 x g**. Remove the supernatant and centrifuge the pellet again briefly. Remove residual Buffer NT3 completely.



400 µl NT3



30 sec
10,000 x g

4 Dry silica matrix

Dry the pelleted silica matrix at **room temperature** or at **37°C** for **10-15 min**.

It is not recommended to dry the sample by vacuum since over-dried pellets lead to lower recoveries.

Residual ethanol from Buffer NT3 would inhibit subsequent reactions and has to be removed in this step.

37°C
10-15 min

5 Elute DNA

Add **25-50 µl Elution Buffer NE** to the silica matrix. Resuspend the pellet by vortexing. Incubate the mixture at **room temperature** for **10-15 min**. Vortexing the mixture 2-3 times during incubation is recommended. Centrifuge the sample at **10,000 x g** for **30 sec** and transfer the DNA containing supernatant to a clean tube (not provided). Repeating this step will increase the yield by approximately 10%.

Yield of larger fragments (> 5-10 kbp) can be increased by performing the incubation at 55°C.



+ 25-50 µl NE

RT
10 min



30 sec
10,000 x g

7 Support protocol for concentration, desalination, removal of enzymes etc.

DNA from samples like e.g. reaction mixtures can be purified with the **NucleoTraP[®]CR** kits as well as with the **NucleoTrap[®]** kits.

1 Adjust DNA binding conditions

Add **4 volumes Buffer NT2** to **1 volume** of DNA containing **sample** (e.g. 400 µl Buffer NT2 and 100 µl reaction mixture).



mix 4 vol. NT2

with

1 vol. sample

2 Bind DNA

Vortex the NucleoTrap[®] / NucleoTraP[®]CR Suspension thoroughly until a homogeneous mixture results. For every µg of DNA add 4 µl of silica matrix (**at least 10 µl**). Incubate the mixture for **10 min** at **room temperature** and vortex briefly every 2-3 min.

Continue with **section 6, step 3**.

Important note: Be aware of the NucleoTrap[®] Suspension binding fragments down to 20 bp (also see Table 1, section 2.4).

8 Appendix

8.1 Troubleshooting

Problem	Possible cause and suggestions
Incomplete lysis of agarose slices	<p><i>High concentration of agarose</i></p> <ul style="list-style-type: none"> Use doubled volumes of Buffer NT1 for highly concentrated agarose gels.
	<p><i>Wrong buffer</i></p> <ul style="list-style-type: none"> Buffer NT2 can not be used for gel dissolution.
	<p><i>Time and temperature</i></p> <ul style="list-style-type: none"> Check incubation temperature. Depending on the weight of gel slice, incubation (section 5, step 2) can be prolonged up to 20 min. Vortex every 2 min and check integrity of the gel slice. Heavy weight gel slices may be quenched or crushed before addition of Buffer NT1.
No DNA yield	<p><i>Reagents not applied properly</i></p> <ul style="list-style-type: none"> Add indicated volume of 96-100% ethanol to Wash Buffer NT3 concentrate and mix well before use.
	<p><i>Insufficient drying of the NucleoTrap®/NucleoTraP®CR silica matrix</i></p> <ul style="list-style-type: none"> Ethanol Wash Buffer NT3 has to be removed quantitatively before elution. Prolong the drying time up to 30 min. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots).
	<p><i>Isolation of large DNA fragments</i></p> <ul style="list-style-type: none"> Add room-temperature Elution Buffer NE and incubate at 55°C for 10-15 min.

Problem	Possible cause and suggestions
Suboptimal performance of DNA in sequencing reactions	<i>Carry-over of ethanol/ethanolic Buffer NT3</i>
	<ul style="list-style-type: none">• Make sure to dry the silica matrix in order to achieve complete removal of ethanolic Buffer NT3 after the washing step. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots)• Buffers other than Buffer NE e.g. TE buffer (Tris/EDTA) were used for elution of DNA. Note: EDTA may inhibit sequencing reactions. In this case it is recommended to re-purify DNA and elute in Buffer NE or water.
	<i>Not enough DNA used for sequencing reaction</i>
	<ul style="list-style-type: none">• Quantitate DNA by agarose gel electrophoresis before setting up sequencing reactions.
	<i>NucleoTrap[®] or NucleoTraP[®]CR particles were not removed quantitatively</i>
	<ul style="list-style-type: none">• Centrifuge the eluate again and transfer the supernatant to a new tube.

8.2 Ordering information

Product	Cat. No.	Pack of
NucleoTrap [®]	740584.10	10 preps
NucleoTrap [®]	740584	100 preps
NucleoTraP [®] CR	740587.10	10 preps
NucleoTraP [®] CR	740587	100 preps
NucleoTrap [®] Suspension	740589	100 preps
NucleoTraP [®] CR Suspension	740564	100 preps
Buffer NT1	740596.100	2 x 50 ml
Buffer NT2	740597	2 x 50 ml
Buffer NT3 Concentrate (for 100 ml Buffer NT3)	740598	20 ml
Collection Tubes (2 ml)	740600	1000

8.3 References

Vogelstein B., and D. Gillespie. 1979. Preparative and analytical purification of DNA from agarose. Proc. Natl. Acad. Sci. USA **76**: 615-619.

8.4 Product use restriction / warranty

NucleoTrap[®]/NucleoTraP[®]CR kit components were developed, designed and sold **for research purposes only**. They are suitable **for in vitro uses only**. No claim or representation is intended for its use to identify any specific organism or for clinical use (diagnostic, prognostic, therapeutic, or blood banking).

It is rather the responsibility of the user to verify the use of the **NucleoTrap[®]/NucleoTraP[®]CR** kit for a specific application range as the performance characteristic of this kit has not been verified to a specific organism.

This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish an extra copy.

MACHEREY-NAGEL does not warrant against damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; against defects in products or components not manufactured by MACHEREY-NAGEL, or against damages resulting from such non-MACHEREY-NAGEL components or products.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, REPRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications. This warranty is exclusive and MACHEREY-NAGEL makes no other warranty expressed or implied.

The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications.

Please contact:

MACHEREY-NAGEL Germany
Tel.: +49 (0) 24 21 969 270
e-mail: TECH-BIO@mn-net.com

Last updated 12/2006, Rev.02