# DiaQuant®

## Peroxide S

de

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



en

### Description:

DiaQuant® Peroxide S are test strips for the semi-quantitative detection of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) and peroxides in rinse solutions.

### Measuring range:

0.5 – 25 mg/L H<sub>2</sub>O<sub>2</sub>

### Color gradation:

0 · 0.5 · 2 · 5 · 10 · 25 mg/L H<sub>2</sub>O<sub>2</sub>

### Pack content:

100 test strips

### Reaction principle:

Hydrogen peroxide reacts with the peroxidase (POD) and the organic redox indicator in the test field to form a blue-colored oxidation product.

### General information:

Only remove the required number of test strips from the container, and then immediately close the container securely afterwards. Ideally, the color evaluation should be carried out by diffuse daylight. Artificial light can lead to difficulties or an incorrect color reading.

### Hazard warnings:

This test does not contain hazardous substances that must be labelled.

### Instructions for use:

1. Insert the test strip into the test solution for 1 sec.
2. Shake off any excess liquid.
3. Wait 15 s.
4. Then immediately compare the test field with the color scale. Take the value which matches closest with the colored test field (reading accuracy: ± ½ colored field of the scale). The reaction color of the test field may change after the value has been taken. It is therefore crucial to evaluate the coloration within the prescribed timescale in order to achieve correct results.

### Quality control:

To check the correct functioning of the test strips, use a hydrogen peroxide solution with a concentration of 5 mg/L. For this purpose, first prepare a stock solution of 500 mg/L H<sub>2</sub>O<sub>2</sub>, by diluting 1.5 mL of a 30% hydrogen peroxide solution in 1000 mL distilled water. Then, dilute 5 mL of this stock solution in 500 mL distilled water (= 5 mg/L H<sub>2</sub>O<sub>2</sub>). Immediately perform the measurement with the test strip. If the control solution produces a negative result even after repeating the process, then the remaining unused test strips must be discarded. Even during a negative control (inserting a test strip into distilled water) no blue coloration may occur. Possible reasons for both faults may be that the use-by-date has been exceeded, the test strip container was left open for too long or the test strips were stored incorrectly.

### Interferences:

If the solution to be tested has a pH value of between 2 and 12, the reaction will take place without interferences. Strong acid solutions have to be buffered using sodium acetate; alkaline solutions have to be set to a pH value of 5–7 using citric acid. The presence of other strong oxidants such as chlorine will also lead to positive results.

### Storage:

Avoid exposing the strips to sunlight and moisture. Keep container cool and dry (storage temperature 4 – 30 °C).

If correctly stored, the test strips may be used until the use-by-date printed on the packaging.

### Additional information:

The test strip container stopper contains a non-toxic drying agent. If swallowed, drink plenty of water.

### Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from [www.mn-net.com/SDS](http://www.mn-net.com/SDS).

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