NANOCONTROL BSB₅-Standard

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NANOCOLOR BOD₅-Standard

for analytical quality assurance of the Tests 8-22 / 8-25

Method:

The **NANOCOLOR BOD₅-Standard** is used as a control system for the **NANOCOLOR[®]** tube tests BOD₅ - Test 8-22 (REF 985822) and BOD₅-TT - Test 8-25 (REF 985825) for the purposes of analytical quality assurance. In freeze-dried form, the standard contains D(+)-glucose and L-glutamic acid in accordance with the requirements stipulated in the German Standard Method **DIN ISO 5815-1** and is dissolved in **inoculated diluting water**. The ready BOD₅ standard solution can directly be used as test sample in the form of a 1:25 or 1:50 dilution.

Contents of reagent set:

10 test tubes with standard reagent, freeze-dried

Hazard warning:

This reagent set does not contain any harmful substances which must be specially labelled as hazardous.

Storage:

The reagent set is to be stored in cool and dry conditions at +2 °C to +8 °C. The expiry date printed on the kit is to be observed.

Procedure:

Always use a **freshly prepared standard!** At the end of the day discard the BOD₅ standard solution that is no longer needed! Carry out the analysis as a **double determination**, wherever possible.

Preparation of the BOD₅ standard solution for tube test BOD₅-TT - Test 8-25 (REF 985825):

Requisite accessories:

Reagent set **NANOCOLOR[®]** BOD₅-TT - Test 8-25 (REF 985825), BOD₅-TT Accessories Set (REF 916925), measuring cylinder (nominal volume 500 mL), glass beaker or Erlenmeyer flask (minimum volume 500 mL)

| | |
|----|--|
| 1. | Dissolve the contents of a BOD ₅ standard test tube in 500 mL inoculated diluting water . (For production and handling of inoculated diluting water, see BOD ₅ -TT Accessories Set, REF 916925). Then dilute the resulting sample 1+1 with diluting water. The BOD ₅ -Standard solution is now ready for direct use as a test preparation in 1:50 dilution . No additional, separate preliminary dilution of the BOD ₅ -Standard solution is necessary. |
| 2. | Pour 20 mL of the ready BOD ₅ standard solution into reaction vessel (component of BOD ₅ - TT Accessories Set) and stir to enrich oxygen content by shaking vigorously for 30 seconds. |
| 3. | Open test tube with reagent BOD ₅ -TT R0 and fill to the brim with BOD ₅ standard solution <i>without letting any air bubbles</i> in. |
| 4. | Close test tube <i>without letting any air bubbles</i> in, label and incubate in water bath or incubator for 5 days at (20 ± 1) °C in the dark . |

Preparation of BOD₅ standard solution for tube test BOD₅ - Test 8-22 (REF 985822):

Requisite accessories:

NANOCOLOR[®] BOD₅ - Test 8-22 (REF 918822), BOD₅-accessories Set (REF 916918), BOD₅ Nutrient Mixture (REF 918994) or BOD₅ Nutrient Mixture PLUS (REF 918995), measuring cylinder (nominal volume 500 mL)

| | |
|----|---|
| 1. | Dissolve the contents of one BOD ₅ standard test tube in 500 mL diluting water . (For production and handling of diluting water for dilution see BOD ₅ -Accessories Set, REF 916918) and pour this solution into the 1 L laboratory flask of the BOD ₅ -Accessories Set. The BOD ₅ -Standard solution is now ready for direct use as a test preparation in 1:25 dilution . No additional, separate preliminary dilution of the BOD ₅ -Standard solution is necessary. |
| 2. | Add 1.25 mL inoculated nutrient solution R1 and R2 of the test kit BOD ₅ Nutrient Mixture (REF 918994) or BOD ₅ Nutrient Mixture PLUS (REF 918995), and stir to enrich oxygen content by shaking vigorously for a few seconds. |
| 3. | Open 1 oxygen Winkler flask and 1 test tube, and fill to the brim with ready BOD ₅ standard solution <i>without letting any air bubbles</i> in. |
| 4. | Close the oxygen Winkler flask <i>without letting any air bubbles</i> in, label and incubate in water bath or incubator for 5 days at (20 ± 1) °C in the dark . Close also the test tube <i>without letting any air bubbles</i> in, label and immediately perform oxygen measurement. |

Measurement and reporting of results:

The measurements and the evaluation of the measurement results are carried out in accordance with the instructions for use of reagent set **NANOCOLOR[®]** BOD₅ - Test 8-22 (REF 985822) or **NANOCOLOR[®]** BOD₅-TT - Test 8-25 (REF 985825). Give due consideration to the **dilution factor** when calculating the final results!

Required value for the BOD₅ standard: 198 ± 40 mg/l O₂

The required value of 198 + 40 mg/l O₂ corresponds to the generally recognised variation range of the biochemical oxygen demand BOD₅ in accordance with the German Standard Method **DIN ISO 5815-1**.

Disposal:

Place used test tubes back in the original package. **MACHEREY- NAGEL** takes back all **NANOCOLOR BOD₅-Standard** reagent sets free of charge for proper disposal. Please observe local regulations concerning transportation of waste.

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