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User's manual

Mode d'emploi

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NANOCOLOR[®] T-Set

MACHEREY-NAGEL

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1 Technical Data

Type:	electronic temperature sensor for temperature control and calibration and creation of an inspection certificate for the monitoring of test resources
Detector:	PT 1000 (95 x ø 4 mm)
Resolution:	± 0.1 °C
Display:	via LED displays of the heating block and T-Set software
Operation:	via foil-covered symbol keys on the heating block and T-Set software
Temperature range:	0 – 200 °C
Accuracy:	± 0.2 °C
Long-term stability:	± 0.1 °C
Interface:	bi-directional serial RS 232
Operating range:	0 – 50 °C, up to 90 % rel. humidity
Power supply:	via RS 232
Power consumption:	max. 20 mW
Dimensions:	112 cm (length)
Weight:	approx. 60 g
Approval symbol:	CE
Certificate:	calibrated with standard calibration thermometer
Guarantee:	2 years

CE Declaration of conformity:

This instrument conforms to the following directives:
89/336 EC dated 03.05.1989 (including amendment directive 92/31/EC) – EMC directive

2 Operation

2.1 Operating Instructions

Insert the RS 232 plug of the *NANOCOLOR*[®] T-Set in the RS 232 interface connector of the switched off heating block. Fully insert the temperature sensor of the *NANOCOLOR*[®] T-Set in the small hole provided on the safety cover of the heating block.

Whilst holding the -key pressed (for the *NANOCOLOR*[®] *VARIO 3*, the key on the side required), switch the heating block on with the mains switch on the rear of the unit. Keep the -key pressed until a long signal tone is heard. In the upper display, „tEST“ is shown. Now, by pressing the -key, a change-over is possible between the test program („tEST“) and calibration („CAL“).

2.2 Program Selection

The program required is started by pressing the -key.

Before running a test program or a calibration, the temperature of the heating block should be less than the lowest programmed temperature in the heating block. Otherwise, the selected program will not start until the temperature has fallen below the lowest programmed temperature.

2.2.1 Test Program (tEST)

All temperatures programmed in the heating block are measured by T-Set, from the lowest to the highest, registered and then saved in the heating block (this takes about 40 min). A signal tone is heard when the heating block has reached a programmed setpoint temperature (shown in the upper display). The actual temperature measured by T-Set is shown in the lower display and saved in the heating block.

While the test program is running, the PROG, WAIT and RUN LED's flash and symbolise the on-going temperature response. The temperatures measured by T-Set can be logged by hand while the test program is running, or output and documented with the T-Set software after the test

program has finished.

After testing all the temperatures programmed in the heating block, „tEST“ is shown in the upper display and „End“ in the lower display. The heating block must then be switched off.

Safety Notice:

The temperature sensor is hot!

When removing the sensor from the heating block, grasp the black insulation.

2.2.2 Calibration (CAL)

All temperatures programmed in the heating block are measured by T-Set, from the lowest to the highest, registered and then saved in the heating block (this takes about 40 min). When the heating block has reached a programmed setpoint temperature (shown in the upper display) and calibrated the value, a signal tone is initiated. The actual temperature measured and calibrated by T-Set is shown in the lower display and saved in the heating block.

While the test program is running, the PROG, WAIT and RUN LED's flash and symbolise the on-going temperature response. The temperatures measured and calibrated by T-Set can be logged by hand during calibration, or output and documented with the T-Set software after the calibration has finished.

After calibration of all temperatures programmed in the heating block, „CAL“ is shown in the upper display and „End“ in the lower display. The heating block must then be switched off.

Safety Notice:

The temperature sensor is hot!

When removing the sensor from the heating block, grasp the black insulation.

3 T-Set Software

3.1 Installation

Insert the *NANOCOLOR*[®] software CD in the CD-drive of your computer. The Start menu is displayed on the screen. Alternatively, the installation can be started via the Windows Explorer.

Select the language preferred, then the function „Install *NANOCOLOR*[®] T-Set“. Follow the installation instructions shown on the screen. By default, the *NANOCOLOR*[®] T-Set software is copied to the folder C:\Program\MACHE-REY_NAGEL\T-Set Software. It is recommended to retain this as the target folder. An icon is automatically created on the screen. An installation on another drive does not usually result in any malfunction.

3.2 Program Start

Connect the heating block to the computer using the serial null-modem cable (Cat. No. 919 680). Switch on the heating block. Start the *NANOCOLOR*[®] T-Set software by clicking the icon shown on the screen. The heating block connected, is automatically recognised by the software. In the top area of the start window, a menu bar is displayed (File, Setup, Help) together with 4 buttons for the various functions in the *NANOCOLOR*[®] T-Set software: „Graph“, „Test“, „Calibration“ and „Settings“.

3.3 Menu

File

The only option available in the File menu, is „Exit“. If this item is selected, the program is ended and the window is closed. The closing of the window has no effect on the processes that are on-going in a heating block that may be connected.

Setup

This option allows the user access to various program parameters and to make adjustments.

Connection

COM (serial interface) Setup, settings ex-works:

Interface:	COM 1
Baud rate:	2400
Data bits:	8
Stop bits:	1
Parity:	None
Flow control:	None

Language

Select the preferred language. Changes are immediately accepted.

Configuration

User

The name, address and user name entered here, appear later on reports and certificates.

Connection

The automatic search for an equipment can be switched ON or OFF.

Help

Help for various functions of the T-Set and the T-Set software

3.4 Graph

This function is used to record, display and save heating curves. The operating panel on the heating block is also displayed on the screen (temperature, time and keyboard). The heating program required can be selected via the *NANOCOLOR*[®] T-Set software with the symbol key icons or directly at the heating block via the symbol keys. Start the heating program with the corresponding -key. The on-going temperature curve is displayed on the screen and can be saved at any required instant of time by pressing the „Save“ button. The heating curve is automatically

formatted. The temperature is displayed on the vertical axis and the time on the horizontal axis. The two horizontal lines that can be seen after starting the recording of a heating curve, show the permissible range of temperature, $\pm 2\text{ }^{\circ}\text{C}$. Yellow vertical lines show the start and end of the heating process. The timer starts as soon as the temperature has reached the programmed value $\pm 2\text{ }^{\circ}\text{C}$.

The „Open“ button is used to call up temperature curves that have previously been saved. When a saved file is opened, the „Off-line mode“ is initiated. In order to continue the measurements and to close the opened heating curve, this option must be deactivated.

After selecting the „Report“ button, a recorded heating curve can be edited (input of details, e.g. title, commentaries, etc.). When the „Report“ button is pressed, the *NANOCOLOR*[®] T-Set software creates a protocol of the recorded heating curve together with all information that has been input, which can then be saved and/or printed.

3.5 Test

After completing the run of a test program, the test data measured by the *NANOCOLOR*[®] T-Set and saved in the heating block, can be read out, displayed and edited.

At the completion of a test, information can be entered in the entry boxes (date of the test, number of the *NANOCOLOR*[®] T-Set, the certificate number applicable and the date of issue). In the „Notes box“, any notes that may be relevant can be entered.

Note:

If the heating block does not contain any data or the information has already been entered, it is not possible to input any data or amend any existing data (the fields are greyed). The input of data or any modification is possible, only after a new test run has been completed.

The temperature table is divided into two columns. In the „nominal“ column, all temperatures are listed that have been programmed in the heating block. The „actual“ column shows the temperatures last measured with the

NANOCOLOR® T-Set.

When the „Report“ button is activated, the **NANOCOLOR®** T-Set software creates a certificate of the test run, that is then automatically printed.

3.6 Calibration

After completing a calibration run, the calibration data measured by the **NANOCOLOR®** T-Set and saved in the heating block, can be read out, displayed and edited.

After calibration, information can be entered in the entry boxes (date of the calibration, number of the **NANOCOLOR®** T-Set, the certificate number applicable and the date of issue). In the „Notes box“, any notes that may be relevant can be entered.

Note:

If the heating block does not contain any calibration data or the data has already been entered, it is not possible to input any data or amend any existing data (the fields are greyed). The input of data or any modification is possible, only after a new calibration run has been completed.

The temperature table is divided into three columns. In the „nominal“ column, all temperatures are listed that have been programmed in the heating block. The „actual“ column shows the temperatures measured with the **NANOCOLOR®** T-Set. In the „Calibrated to“ column, all calibrated temperatures are listed.

When the „Report“ button is activated, the **NANOCOLOR®** T-Set software creates a certificate of the calibration run, that is then automatically printed.

3.7 Settings

With this program item, existing settings on the equipment can be checked and edited.

The temperature and time tables are divided into two sections. In the grey-shaded cells the temperatures and heating times are set by the manufacturer and cannot be changed. In the white cells, the user can input his own settings, or the existing values can be amended. The temperature values must be entered as whole numbers; the time values in the form of hh:mm.

The „Factory Data“ box contains the serial number and the software version number.

In the „Sound“ box, the periodic signal tone that is sounded at the completion of a heating program, can be switched On or Off. If the signal tone is switched off, at the end of a heating program only 2 long signal tones are heard.

In the box „Error Log“, faults or errors are recorded. This is a read-only field. The information given can be of help in locating faults when problems arise. Using the „Write“ button, the settings on the screen can be written into the memory of the unit.

Disposal in accordance with EU directive 2002/96/EC

In compliance with EU directive 2002/96/EC, MACHEREY-NAGEL will accept back and dispose of the old instrument at no cost.

ATTENTION!

It is not allowed to dispose of the instrument using municipal waste disposal services. Please talk to your local MACHEREY-NAGEL contact.

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MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren