MACHEREY-NAGEL

Photometer PF-3

- Handbuch
- Manual
- Manuel

Deutsch	4
English	26
Français	47
Español	70

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Tel.: +49 24 21 969-0

sales@mn-net.com · www.mn-net.com

1. General information and introduction	28
2. CE mark	28
3. Identification plate.....	28
4. Explanation of symbols used	29
5. Safety information.....	30
6. Chemical and biological safety.....	30
7. Technical description	31
7.1 Technical data.....	31
8. Setting up and connecting the device.....	32
8.1 Scope of delivery*	32
8.2 Connecting the device	32
8.3 Battery operation	33
9. Device views.....	33
9.1 Interface/power supply.....	34
10. Initial operation	34
10.1 Switching on.....	34
10.2 Operation and user guidance	35
11. Measure	35
11.1 Test selection.....	35
11.2 Additional tests	36
11.3 <i>VISOCOLOR</i> [®] and <i>NANOCOLOR</i> [®] tests.....	36
11.4 Performing <i>VISOCOLOR</i> [®] <i>ECO</i> , Powder Pillow and <i>NANOCOLOR</i> round cuvette tests	38
12. Settings.....	38
12.1 Enabling/disabling tests	39
12.2 Date.....	39
12.3 Time.....	39
12.4 Reaction time.....	39
12.5 Define switch-off time.....	40
12.6 Update.....	40
12.7 Language.....	40
12.8 Calibration.....	40
12.9 20–80% bar	41
12.10 0 set function	41
12.11 Calibration adjustment	41
13. Quality control	42
14. Memory	42
14.1 Display memory	42

14.2 Transferring data	43
15. User-defined special method creation.....	43
15.1 Maintenance and cleaning of the device.....	43
15.2 Problems, causes, and correction	44
16. Service	45
16.1 Replacement parts, accessories, and consumables.....	45
17. Disposal	45
17.1 Warranty, liability, and claims.....	46
17.2 Contact information	46

1. General information and introduction

The compact photometer PF-3 is a filter photometer from MACHEREY-NAGEL for mobile use and use in the laboratory. It is used for the evaluation of colorimetric test kits in 16 mm round cuvettes from MACHEREY-NAGEL. The illuminated display and intuitive menu navigation make working with the device fast, convenient and pleasant.

2. CE mark

The CE mark indicates that the device meets the harmonization legislation of the European Community listed below:

European directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)

European directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

European directive 2014/30/EU on electromagnetic compatibility (EMC)

3. Identification plate

Fig 1: Identification plate

4. Explanation of symbols used

The device identification plate and these operating instructions contain the symbols and terms which are listed below and have the following meaning:

Term/symbol	Meaning
Photometer PF-3	Device name
SN	Serial number of the device
3.6-4.5 V, 500 mA	Power supply
	DC voltage
IP68	Dust- and waterproof (test conditions: 1 m water depth, 1 hour)
	In accordance with 2012/19/EU, disposal of the device via public waste disposal systems is not permissible. See also the section on disposal.
	CE denotes compliance with the applicable harmonization legislation of the European Community.
MACHEREY-NAGEL	Manufacturer identification
	This symbol indicates that there is a hazard due to electrical shock.
	This symbol indicates that the chemicals used have a corrosive effect. Follow the safety measures in the laboratory and wear the required protective equipment. Observe the information in the current safety data sheets (SDS) of the products used.
	This symbol indicates that there is a hazard due to the use of flammable materials.
	Explanations about the text. Tips and tricks to work better.

5. Safety information

BE SURE TO READ THE FOLLOWING SAFETY INFORMATION CAREFULLY BEFORE YOU USE THE DEVICE.

Failure to comply with the information can lead to severe injuries to the operator, to malfunctions, or to damage to the device.

Keep these operating instructions so you can also consult them later.

Follow the safety information and directions in the operating instructions and take note of the stickers and information on the device.

Work may not be performed on inner parts of the device. Noncompliance invalidates any warranty claims.

Meaning of the safety information:

HAZARD

Identifies an impending or possible hazardous situation which, if not avoided, results in death or serious injuries.

WARNING

Identifies an impending or possible hazardous situation which, if not avoided, may result in death or serious injuries.

CAUTION

Identifies a possible hazardous situation which, if not avoided, may result in minor or moderate injuries.

NOTE

Identifies a situation which, if not avoided, can lead to malfunctions or damage to the device. Information which explains the circumstances of the text and which must be observed in particular.

6. Chemical and biological safety

WARNING

Possible hazards upon contact with chemical and biological substances.

Working with samples, reagents, and the corresponding accessories is associated with hazards.

When working with the cuvettes, it must be ensured that suitable protective equipment is worn. Observe the safety data sheet (SDS) of the test kit used.

Normal operation of this device may require the use of chemicals which are hazardous to health or biologically harmful samples.

Before handling these substances, observe the safety instructions and safety information printed on all containers of the original solution and in the safety data sheet. All solutions used must be disposed

of in accordance with the national regulations and laws. The type of protective equipment must be selected in accordance with the concentration and quantity of the hazardous substance at the respective workstation.

7. Technical description

The compact photometer PF-3 enables measurements at 3 predefined wavelengths. The light generated by the respective LED is directed through an interference filter into the cuvette shaft of the device. After passing through the sample, the proportion of absorbed light is calculated by detecting the attenuated light. The proportion of absorbed light is converted with the aid of a stored calibration into the concentration of the analyte examined.

7.1 Technical data

Type:	LED photometer with microprocessor control, self-test, and autocalibration
Optics:	LED + interference filter Not sensitive to extraneous light for fast measurements without covering the cuvette shaft
Wavelengths:	Depends on version, 3 wavelengths possible
Wavelength accuracy:	± 2 nm, full width at half maximum 10–12 nm
Number of filters:	3
Light source:	LED
Detector:	Silicon photovoltaic cell
Measurement modes:	Preprogrammed tests (<i>NANOCOLOR</i> [®] round cuvette tests, <i>VISOCOLOR</i> [®] Powder Pillow tests and <i>VISOCOLOR</i> [®] <i>ECO</i> tests)
Cuvette slot:	Round cuvettes, 16 mm OD
Memory:	50 measured values
Display:	Illuminated graphic display, 64 × 128 pixels
Operation:	Intuitive, icon-based menu navigation, membrane keypad
Interface:	Mini USB
Update:	Free of charge, via Internet/PC
Operating range:	5–50°C at 90% rel. humidity Battery operation: 10–40°C To 3000 m
Power supply:	3 AA disposable batteries, rechargeable batteries, USB interface; optional internal battery pack
Protection class according to IEC 60529:	IP68 (device is dust- and waterproof.) IP6X ≙ Dustproof IPX8 ≙ Protection against prolonged submersion.
Overvoltage category:	Device: I Mains adapter: II
Degree of contamination:	2
Dimensions:	170 × 95 × 68 mm
Weight:	0.5 kg
Warranty:	2 years

8. Setting up and connecting the device

8.1 Scope of delivery*

Check the delivery for completeness. If parts are missing from the delivery, contact MACHEREY-NAGEL or your distribution partner directly.

- Compact photometer PF-3
- 4 × AA batteries
- Manual
- Screwdriver
- Certificate

* The scope of delivery refers to the product itself and may differ if the photometer is contained in a different product.

8.2 Connecting the device

CAUTION

Multiple hazards: Only qualified staff may carry out the instructions described in this section.

WARNING

Exposure of the eye to light. Do not look directly into the light rays of the lamps.

WARNING

Risk of electrical shock: When used for operation, ensure that the USB cable does not show any signs of damage. Check the suitability of the voltage source for the device.

NOTE

Damage due to condensation and violation of the temperature limits: The device is designed for indoor and outdoor use. Do not expose the device to direct sunlight. During use outdoors, protect the interface from exposure to humidity and moisture.

WARNING

Risk of fire: When using an external power supply, ensure that the power outlet is not overloaded. There is a risk of overload and fire. Ensure that the power cable does not show any signs of damage. Check the suitability of the voltage source used for the device.

WARNING

Risk of injury: Familiarize yourself with the device before working with it and read this document carefully. Use the device only if you have been trained in how to handle it.

NOTE

Malfunction due to defects at the power supply unit and housing: If the device has an obvious crack in the housing or a damaged power supply unit, the device should be taken out of service.

8.3 Battery operation

WARNING

Risk of corrosion due to battery leaks: Take the batteries out of the device if the PF-3 is not used for longer periods.

NOTE

Damage due to leaking batteries: MACHEREY-NAGEL recommends the use of high-quality batteries.

9. Device views

- ① Graphic display 128 × 64 pixels
- ② Cuvette shaft for round cuvettes with a 16 mm outer diameter
- ③ 4-part keypad
- ④ Battery compartment

9.1 Interface/power supply

The device has a USB interface for power supply and data export. The interface is located underneath the cover of the battery compartment. The battery compartment is at the lower rear of the device and can be opened with a Phillips screwdriver.

5 USB interface

6 Battery compartment

The device is normally operated with three AA disposable/rechargeable batteries. Optionally, the device can also be equipped with a battery pack (REF 919391) which is charged via the USB interface. Operation via the USB connection is also possible.

10. Initial operation

10.1 Switching on

Switch on the device using the I/O switch on the front. The device name and firmware version appear on the display. After a few seconds, the start screen with the most recently selected method appears.

10.2 Operation and user guidance

The device is designed for the simplest possible operation. The user guidance is provided entirely by means of menu symbols which can be selected through the corresponding action symbols with the four action buttons below the screen.

The screen is divided into three sections. The upper status bar shows the date, battery status, and time.

The middle of the screen contains the actual display, for example, of measurement results or menu entries and the corresponding menu symbols.

The action symbols are arranged in the bottom line of the display. These symbols always stand for a specific action which can be selected directly with the action button beneath it.

After exiting a menu using , your changes are made and stored. When exiting a menu using , no changes are made.

Meaning of the symbols

 / 	Scroll left/right		Cancel
	Measure the zero value		Go back
	Measure		Change number values
	Confirmation/selection		Enable/disable

11. Measure

11.1 Test selection

To start a measurement, you must first select the corresponding test. This is done using  / . By default, all tests available for the corresponding PF-3 version are stored as favorites on the top menu level. You may disable or enable these as desired (see 12.1, page 39).

Some tests are stored multiple times in the device in various units, known as submethods. This means that the unit relevant for the measurement can always be selected.

NOTE

Upon initial use, we recommend adding desired favorites to the top menu level so that they can be easily and conveniently used for daily work.

11.2 Additional tests

MACHEREY-NAGEL regularly updates the tests available for the PF-3. Newly available parameters can also be subsequently loaded onto the device free of charge. The corresponding program can be found at www.mn-net.com/nanocolor-geraete-updates.

11.3 VISOCOLOR® and NANOCOLOR® tests

NOTE

Only the associated mains adapter GTM46101-1005-USB (REF 919220) may be used.

CAUTION

Multiple hazards: Only qualified staff may carry out the instructions described in this section.

WARNING

Risk of electrical shock: When used for operation, ensure that the USB cable and mains adapter used do not show any signs of damage. Check the suitability of the voltage source for the device.

NOTE

Damage due to condensation and violation of the temperature limits: The device is designed for indoor and outdoor use. Do not expose the device to direct sunlight. During use outdoors, protect the interfaces from exposure to humidity and moisture.

WARNING

Risk of fire: When using an external power supply, ensure that the power outlet is not overloaded. There is a risk of overload and fire. Ensure that the mains adapter including cable does not show any signs of damage. Check the suitability of the voltage source used for the device. Use only the mains adapter recommended by MACHEREY-NAGEL.

WARNING

Risk of injury: Familiarize yourself with the device before working with it and read this document carefully. Use the device only if you have been trained in how to handle it.

NOTE

Malfunction due to defects at the power supply unit and housing: If the device has an obvious crack in the housing or a damaged power supply unit, the device should be taken out of service.

WARNING

Possible hazards upon contact with chemical and biological substances.

Working with samples, reagents, and the corresponding accessories is associated with hazards.

When working with the cuvettes, it must be ensured that suitable protective equipment is worn. Observe the safety data sheet (SDS) of the test kit used.

NOTE

Soiling of the cuvette shaft and incorrect measurements: Clean the cuvettes before each measurement using a lint-free cotton cloth in order to avoid soiling of the cuvette shaft and incorrect measurements. Keep the cuvette shaft free of soiling.

The PF-3 can evaluate tests from our *VISOCOLOR*[®] and *NANOCOLOR*[®] programs. There is a manual for the *VISOCOLOR*[®] tests which describes the implementation of the tests in pictogram form. This manual is included in the device case or it can be downloaded from www.mn-net.com. The implementation of the *NANOCOLOR*[®] tests is described in the corresponding package leaflets.

In addition to the parameter designation, *VISOCOLOR*[®] and *NANOCOLOR*[®] tests are always identified in the device by the addition of "NANO" or "VISO" and a 4-digit test number which is unique for this test.

11.4 Performing VISOCOLOR® ECO, Powder Pillow and NANOCOLOR round cuvette tests


Instead of the VISOCOLOR® ECO measuring glass, a clean, dry round cuvette (16 mm OD) is used. A second round cuvette is used for the zero setting (an untreated sample or a reagent blank value, depending on the test).

The blank values of many tests are already stored in the PF-3. A zero measurement is therefore necessary only in the case of individual tests.

To measure a test, first prepare the sample (and the zero value, if applicable) according to the test instructions. If no zero measurement is necessary (NANOCOLOR® tests / 0 set function), step 3 can be skipped.

1. Switch on photometer
2. Select desired parameter
3. Perform the zero measurement (not necessary for all tests) using “NULL/ZERO”
4. Perform the sample measurement using “MESSEN” [measure]
5. Read measured value

For detailed implementation of the individual tests, follow the corresponding test instructions or the separately available pictogram instructions for VISOCOLOR® ECO on the PF-3. These manual pages in the pictogram depiction are available online at www.mn-net.com and on the enclosed software DVD.

If the measurement result is outside of the measuring range of the test, the device displays “less than” (<) or “greater than” (>) and the measurement range limit of the test. If a measured value was measured with a test for which a calibration adjustment is stored, an exclamation point is shown to the right of the measured value. To delete the calibration adjustment, you must reset the parameter properties to factory settings (see 12.11, page 41).

12. Settings

The settings menu can be reached at the top menu level alongside the tests using  / . To get to the settings menu, scroll over the wrench symbol and press the action button .

12.1 Enabling/disabling tests

In this menu, you can decide which tests should be directly selectable (enabled) on the top menu level. Using / you can select between the various tests. You can enable or disable the test using the action button . The symbol in the upper right-hand part of the screen indicates whether a test is enabled or disabled (/). All enabled tests will be displayed directly in the top level after the device is started. Nonetheless, all disabled tests are still stored in the device.

12.2 Date

You can adjust the date in the date menu. Once you select the menu, you will see the date currently set. The first number is highlighted and can be changed directly using the button. Using the right action button you can select the date digit to be edited.

After you have set the desired date, you can change the date format using / . The following options are available:

dd.mm.yyyy

mm/dd/yyyy

yyyy-mm-dd

Press the action button to save your changes. If you exit the menu using , your changes will be discarded.

12.3 Time

You can adjust the time in the time menu. Once you select the menu, you will see the time currently set. The first number is highlighted and can be changed directly. You can set a different number value using the left action button . Using the right action button you can select the time digit to be edited.

After you have set the desired time, you can change the time format using the action buttons / . You can switch between a 24h and a 12h display.

Press the action button to save your changes. If you exit the menu using , your changes will be discarded.

12.4 Reaction time

In the reaction time menu, you can set whether the device should display the reaction time for the test to be measured and should count down prior to the measurement. If the reaction time is disabled, you can start the measurement directly. You can enable or disable the reaction time using the action button .


NOTE


When the reaction time is enabled, you can start a measurement early by pressing the action button twice.

12.5 Define switch-off time

In the “Define switch-off time” menu, you can define whether and when the device automatically switches off after inactivity. Using  /  you can select from four different switch-off times (5, 10, 15, 20 min) and continuous operation (∞).


Press the action button  to save your changes. If you exit the menu using , your changes will be discarded.


12.6 Update

If a new update or new parameters are to be installed on the device, it may be necessary to switch the device into update mode. Available updates and new parameters for the device are found on our website at www.mn-net.com/de/nanocolor-device-updates. Detailed information on installing an update can be found in the corresponding update instructions.

12.7 Language

You can select various languages for the PF-3. However, language changes always relate only to the names of the individual parameters and tests; all other functions and settings are shown using language-neutral icons. You can select the various languages using  / .

Press the action button  to save your changes. If you exit the menu using , your changes will be discarded.

12.8 Calibration

MACHEREY-NAGEL recommends device calibration at regular intervals, however at least once annually. To perform a calibration, press the  button in this menu to start the calibration process. After a few seconds, the prompt to insert the calibration cuvette appears. Insert the cuvette and confirm once again with . After the calibration, you can exit the menu using .

NOTE

Ensure that the calibration cuvette used is free of scratches and clouding. Clean the outside of the cuvette with a soft cotton cloth.

NOTE

Risk of an incorrect calibration: Do not insert the calibration cuvette until you have been prompted to do so.

12.9 20–80% bar

The device can optionally display a 20–80% measurement range bar to check whether the result is in the middle measurement range. In this menu, you can enable and disable the 20–80% bar using the action button . The measurement range bar is disabled as a factory preset.

12.10 0 set function

A zero measurement must be performed for all *VISOCOLOR*[®] *ECO* tests. Therefore, when switching between tests, it is necessary to measure a new zero. If you conduct a series of tests for various parameters from the same sample, you can use the 0 set function to store a zero for all *VISOCOLOR*[®] tests so that the zero measurement must be performed only once.

When the 0 set function is switched on, you can call up a new test after measuring a zero, insert the measurement cuvette directly, and perform the measurement using . If you measure a new zero, this is stored in the device as a new, overriding zero.

Once you switch off the device, it is necessary to measure a new zero. This function is not enabled by default. You can switch the 0 set function on or off in the settings.

NOTE


Retaining a single zero measurement as the basis for further tests can be recommended only if measurements for various parameters are performed from one sample or the samples are very similar.

If you are not sure whether there is sufficient similarity, we always recommend measuring a zero for our *VISOCOLOR*[®] tests.

12.11 Calibration adjustment

The PF-3 has already been factory-calibrated for direct use without user adjustment. The calibration adjustment can be used to meet legal requirements. To do this, you can adjust the MACHEREY-NAGEL factory programming for the *VISOCOLOR*[®] and *NANOCOLOR*[®] tests in the “Adjust calibration” menu with the aid of a standard. You can select the test to be adjusted using . Confirm with , insert the clean cuvette with the zero solution in the photometer, and press the **NULL ZERO** button. Insert the cuvette with the standard with a known concentration and press the action button . To change the value displayed, you can set a different number value using the left action button . Using the right action button , you can select the digit or decimal place to be edited. The standard adjustment can be performed only within the measurement range of the respective test. Entering the number value zero is not permitted. Confirm your entry with the action button

. The new action button  appears in the menu window. To reset the calibration adjustment to factory settings, confirm the action button .

The action button disappears from the menu window and the calibration is reset to factory settings. After performing or removing the calibration adjustment, you can exit the menu using . After performing the calibration adjustment, an exclamation point is added to the values measured for the adjusted test in the display and also in the memory.

NOTE

Incorrect results due to calibration adjustment possible: We recommend using the calibration specified by MACHEREY-NAGEL stored in the device and not adjusting it. The calibration adjustment can be used to meet legal requirements. To perform the adjustment, we recommend selecting a standard in the upper third of the measurement range of the test to be adjusted in order to minimize calibration errors.

13. Quality control

Each institution must implement its own quality control guideline.

To check the performance of the test kit and device used, we recommend measuring several standard solutions with a known concentration within the measurement range.

The photometric accuracy of the measurement of the PF-3 can be checked with the aid of the *NANOCONTROL NANOCHECK 2.0* color solutions (REF 925703) provided by MACHEREY-NAGEL.

14. Memory

The PF-3 has an internal memory for up to 50 measurement results. The device also stores the date and time of the measurement, in addition to the measured value. If all of the free memory space is full, the device overwrites the oldest entry after a new measurement.

NOTE

No warning is given when the oldest result is overwritten.

14.1 Display memory

In the “Display memory” menu, you will find all stored measured values. You can select the individual entries using  / . In the upper right-hand part of the screen, you will see the number of the entry in the memory as well as the total number of available entries. If a value in the memory is shown with an exclamation point, it was determined using a calibration adjustment of the corresponding test (see 12.11, page 41).

14.2 Transferring data

You can transfer the data from the memory of the PF-3 to a PC. This option is shown in the memory menu only if the device is actually connected to a computer. To do this, connect the device to your computer via the USB interface (see 9.1, page 34).

To transfer the data, you need the free *NANOCOLOR*[®] data export software. This can be downloaded from the MACHEREY-NAGEL website. All other information can be found in the corresponding software instructions.

15. User-defined special method creation

In addition to the preprogrammed *VISOCOLOR*[®] and *NANOCOLOR*[®] tests, up to 10 user-defined special methods can be stored in the device. To create and manage special methods, you need the *NANOCOLOR*[®] data export software from MACHEREY-NAGEL which can be downloaded free of charge from our website. All other information can be found in the corresponding software instructions.

15.1 Maintenance and cleaning of the device

NOTE

The device is to be switched off and disconnected from the power supply during all cleaning work.

WARNING

Risk of fire: Do not use any flammable cleaning agents or organic solvents to clean the device.

WARNING

Incorrect measurements: A dirty cuvette shaft can lead to soiling of the cuvettes used and to incorrect results. Check the cleanliness of the cuvette shaft at regular intervals and during maintenance. Clean the cuvette shaft if it is soiled (see cleaning).

WARNING

Discard chemicals, waste, and used cuvettes in accordance with the national, regional, and local regulations.

NOTE

Do not clean the device with acetone or similar products.

NOTE

Ensure that there is no cuvette in the cuvette shaft.

The housing of the photometer can be cleaned with a damp cloth. Immediately wipe off liquid contaminants from the device and clean the housing as previously described. The cuvette shaft should be cleaned using a soft, lint-free cotton cloth and then dried. The photometer and cuvettes must always be kept clean. The cuvettes are cleaned with a lint-free cloth.

NOTE

Do not use any solvents to clean the housing and display.

15.2 Problems, causes, and correction

NOTE

Malfunction of or damage to the device due to incorrect maintenance: The maintenance of the device may be performed only by qualified staff.

Various messages may be displayed, depending on the operating status. The source of the error may either be incorrect operation or a malfunction of the device.

If errors recur, contact MACHERY-NAGEL.

Error message	Reason	Solution
Display remains dark	No power supply	Check whether disposable/rechargeable batteries are inserted/USB cable is connected.
Device does not work when switched on	Device has crashed	Switch off device, wait 10 s, switch device back on
One or more buttons cannot be operated.	The keyboard is defective.	Contact MACHERY-NAGEL or your responsible distributor.
The device displays incorrect measured values.	The device calibration must be renewed. The device filters are soiled or defective.	Perform a device calibration (see 12.8, page 40) and check the device (see 13, page 42). Contact MACHERY-NAGEL or your responsible distributor.

Error message	Reason	Solution
<p>After starting, the device shows the following image on the display:</p> 	<p>The device calibration is no longer working properly.</p>	<p>Clean the cuvette shaft and press "CAL" the next time the symbol appears. When doing so, no cuvette should be in the cuvette shaft.</p>

16. Service

16.1 Replacement parts, accessories, and consumables

WARNING

Risk of injury: Injury to persons and damage to or malfunction of the instrument or equipment if parts not authorized by the manufacturer are used. The replacement parts in this section are authorized by the manufacturer.

Description	REF
Sturdy case with foam inlay, for PF-3, accessories, and VISOCOLOR® tests	934102
Calibration cuvette	916908
USB cable for data export and power supply	919390
Optional rechargeable battery pack	919391
Reaction tubes (20)	91680
Manual	919392
VISOCOLOR® test instructions for PF-3	934001

17. Disposal

NOTE

Disposal of the device via public waste disposal systems is not permissible. Contact your local MACHEREY-NAGEL contact person.

 Disposal according to EU directive 2012/19/EU. In accordance with EU directive 2012/19/EU, MACHEREY-NAGEL takes back old equipment and disposes of it free of charge

17.1 Warranty, liability, and claims

The warranty for this device is for 24 months from the date of purchase. The original invoice serves as proof and must be submitted when asserting a claim. Improper handling and/or maintenance of the device invalidates the warranty. It does not cover any defects which can be attributed to any external power supply unit other than the one supplied.

The warranty is limited to the repair of defective parts or – at MACHEREY-NAGEL's discretion – to the delivery of a defect-free replacement device. The 24-month warranty period is not affected by a warranty claim. There is no right of withdrawal. Other claims are excluded. These include, in particular, all claims for damages resulting from consequential damages or indirect damages. Our general conditions of sale and delivery as printed on all price lists additionally apply, as amended.

17.2 Contact information

If you still have questions after reading the operating instructions or need technical assistance, contact

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MACHEREY-NAGEL GmbH & Co. KG