

Special Instruction

Version 3

Simplification of NANOCOLOR® standard tests

General information

The measuring sample is not prepared in a 25-mL-volumetric flask, but in a beaker or in an Erlenmeyer flask. There is no need to fill up to the final volume of 25 mL.

Overview	
Parameters	Chlorine, Iron, Hydrazine, Silica, Copper, Manganese, Nickel, Nitrite, <i>ortho</i> -Phosphate
Required tests	NANOCOLOR® Chlorine (REF 91816) NANOCOLOR® Iron (REF 91836) NANOCOLOR® Iron LR (REF 918128) NANOCOLOR® Hydrazine (REF 91844) NANOCOLOR® Silica (REF 91848) NANOCOLOR® Copper (REF 91853) NANOCOLOR® Manganese (REF 91860) NANOCOLOR® Manganese LR (REF 918126) NANOCOLOR® Nickel (REF 91862) NANOCOLOR® Nitrite (REF 91867) NANOCOLOR® <i>ortho</i> -Phosphate (REF 91877 or REF 91878)
Accessories	Piston pipette with tips, beaker or Erlenmeyer flask 50 mL, 50-mm-rectangular cuvette of glas
Measurement and calculation	Call up absorbance measurement program and measure. Consider the factors in the below given table for calculation of measuring results. The result can be calculated with this formula: $F4 \times A^4 + F3 \times A^3 + F2 \times A^2 + F1 \times A + F0 = \dots \text{ mg/L}$ F = Factor; A = Absorbance

Remarks

- Procedure for blank value and sample is described in the table below
- Pour blank value and sample in 50-mm-rectangular cuvette separately
- Clean cuvette from outside
- Consider reaction time
- Consider device-depending measuring wavelength λ
- The given factors are for measurements in 50-mm-rectangular cuvettes
- Consider device-depending factors while calculating the results
- Consider the respective photometer manual for programming as a special method

Procedure

Parameter REF, measuring range	Wavelength	Reaction time	Sample Procedure	Blank value Procedure	Factor 400 D, 500 D	Factor VIS, ^{UV} VIS	Factor VIS II, ^{UV} VIS II
Chlorine, free REF 91816 0.02 – 2.00 mg/L Cl ₂	$\lambda = 540 \text{ nm}$	30 s	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix • 1 spoon R2, mix 	20 mL sample solution	F1: 0.92	F1: 0.92	F1: 0.92
Chlorine, total REF 91816 0.02 – 2.00 mg/L Cl ₂	$\lambda = 540 \text{ nm}$	3 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix • 1 spoon R2, mix • 5 drops R3, mix 	• 20 mL sample solution	F1: 0.93	F1: 0.93	F1: 0.93
Iron REF 91836 0.01 – 2.00 mg/L Fe	$\lambda = 470 \text{ nm}$	5 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix • 1 spoon R2, mix • 1 mL R3, mix • 1 mL R4, mix 	<ul style="list-style-type: none"> • 20 mL dest. water • 1 mL R1, mix • 1 spoon R2, mix • 1 mL R3, mix • 1 mL R4, mix 	F1: 1.30	F1: 1.30	F1: 1.30
Iron LR REF 918128 0.005 – 0.500 mg/L Fe	$\lambda = 540 \text{ nm} /$ 563 nm	3 min	<ul style="list-style-type: none"> • 20 mL sample solution • 3 mL R1, mix • 1 spoon R2, mix 	<ul style="list-style-type: none"> • 20 mL dest. water • 3 mL R1, mix • 1 spoon R2, mix 	F1: 0.546 ($\lambda = 540 \text{ nm}$)	F1: 0.455 ($\lambda = 563 \text{ nm}$)	F1: 0.455 ($\lambda = 563 \text{ nm}$)

Parameter REF, measuring range	Wavelength	Reaction time	Sample Procedure	Blank value Procedure	Factor 400 D, 500 D	Factor VIS, ^{UV} /VIS	Factor VIS II, ^{UV} /VIS II
Hydrazine REF 91844 0.002 – 0.250 mg/L N ₂ H ₄	λ = 436 nm	10 min	<ul style="list-style-type: none"> • 20 mL sample solution • 2 mL R1, mix • 2 mL R2, mix 	<ul style="list-style-type: none"> • 20 mL dest. water • 2 mL R1, mix • 2 mL R2, mix 	F1: 0.17	F1: 0.17	F1: 0.17
Silica REF 91848 0.01 – 1.40 mg/L Si 0.02 – 3.00 mg/L SiO ₂	λ = 690 nm	15 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix, wait 3 min • 1 mL R2, mix, wait 1 min • 1 mL R3, mix 	<ul style="list-style-type: none"> • 20 mL silica-free water • 1 mL R1, mix, wait 3 min • 1 mL R2, mix, wait 1 min • 1 mL R3, mix 	F1: 0.74 for mg/L Si F1: 1.57 for mg/L SiO ₂	F1: 0.74 for mg/L Si F1: 1.57 for mg/L SiO ₂	F1: 0.74 for mg/L Si F1: 1.57 for mg/L SiO ₂
Silica high sensitive REF 91848 0.002 – 0.100 mg/L Si 0.005 – 0.200 mg/L SiO ₂	λ = 800 nm	15 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix, wait 3 min • 1 mL R2, mix, wait 1 min • 1 mL R3, mix 	<ul style="list-style-type: none"> • 20 mL silica-free water • 1 mL R1, mix, wait 3 min • 1 mL R2, mix, wait 1 min • 1 mL R3, mix 	F1: 0.294 for mg/L Si F1: 0.629 for mg/L SiO ₂	F1: 0.303 for mg/L Si F1: 0.648 for mg/L SiO ₂	F1: 0.303 for mg/L Si F1: 0.648 for mg/L SiO ₂
Copper REF 91853 0.01 – 2.00 mg/L Cu ²⁺	λ = 585 nm	15 min	<ul style="list-style-type: none"> • 20 mL sample solution • 2 mL R1, mix • The pH-value must be 8.5- 9.5, otherwise adjust with R1. • 2 mL R2, mix 	<ul style="list-style-type: none"> • 20 mL sample solution 	F1: 0.95	F1: 0.95	F1: 0.95
Manganese REF 91860 0.01 – 2.00 mg/L Mn	λ = 470 nm	5 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix • 1 mL R2, mix, wait 1 min • 1 mL R3, mix 	<ul style="list-style-type: none"> • 20 mL sample solution 	F1: 1.24	F1: 1.24	F1: 1.24

Parameter	Wavelength	Reaction time	Sample Procedure	Blank value Procedure	Factor 400 D, 500 D	Factor VIS, UV_{VIS}	Factor VIS II, UV_{VIS} II
Manganese LR REF 918126 0.005 – 0.700 mg/L Mn	$\lambda = 436 \text{ nm} / 450 \text{ nm}$	6 min	<ul style="list-style-type: none"> • 20 mL sample solution • 0.5 mL R1, mix • 1 mL R2, mix, wait 5 min • 2 mL R3, mix, wait 1 min 	<ul style="list-style-type: none"> • 20 mL sample solution 	F1: 0.406 F2: 0.046 ($\lambda = 436 \text{ nm}$)	F1: 0.333 F2: 0.030 ($\lambda = 450 \text{ nm}$)	F1: 0.333 F2: 0.030 ($\lambda = 450 \text{ nm}$)
Nickel REF 91862 0.01 – 2.00 mg/L Ni ²⁺	$\lambda = 436 \text{ nm}$	5 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1 • 1 mL R2, mix • 1 mL R3, mix • 1 mL R4, mix 	<ul style="list-style-type: none"> • 20 mL sample solution 	F1: 1.04	F1: 1.04	F1: 1.04
Nitrite REF 91867 0.002 – 0.100 mg/L NO ₂ -N 0.005 – 0.250 mg/L NO ₂ ⁻	$\lambda = 520 \text{ nm}$	10 min	<ul style="list-style-type: none"> • 20 mL sample solution • 2 mL R1, mix • 2 mL R2, mix 	<ul style="list-style-type: none"> • 20 mL dest. water • 2 mL R1, mix • 2 mL R2, mix 	F1: 0.08 for mg/L NO ₂ -N F1: 0.263 for mg/L NO ₂ ⁻	F1: 0.08 for mg/L NO ₂ -N F1: 0.263 for mg/L NO ₂ ⁻	F1: 0.08 for mg/L NO ₂ -N F1: 0.263 for mg/L NO ₂ ⁻
ortho-Phosphate REF 91877 0.04 – 1.70 mg/L PO ₄ -P 0.1 – 5.0 mg/L PO ₄ ³⁻	$\lambda = 690 \text{ nm}$	10 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix • 1 mL R2, mix 	<ul style="list-style-type: none"> • 20 mL dest. water • 1 mL R1, mix • 1 mL R2, mix 	F1: 1.88 for mg/L PO ₄ -P F1: 5.76 for mg/L PO ₄ ³⁻	F1: 1.88 for mg/L PO ₄ -P F1: 5.76 for mg/L PO ₄ ³⁻	F1: 1.88 for mg/L PO ₄ -P F1: 5.76 for mg/L PO ₄ ³⁻
ortho-Phosphate REF 91878 0.2 – 6.6 mg/L PO ₄ -P 0.5 – 20.0 mg/L PO ₄ ³⁻	$\lambda = 436 \text{ nm}$	10 min	<ul style="list-style-type: none"> • 20 mL sample solution • 1 mL R1, mix • 1 mL R2, mix 	<ul style="list-style-type: none"> • 20 mL dest. water • 1 mL R1, mix • 1 mL R2, mix 	F1: 6.07 for mg/L PO ₄ -P F1: 18.6 for mg/L PO ₄ ³⁻	F1: 6.16 for mg/L PO ₄ -P F1: 18.9 for mg/L PO ₄ ³⁻	F1: 5.98 for mg/L PO ₄ -P F1: 18.3 for mg/L PO ₄ ³⁻

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service

Tel.: +49 24 21 969-333
E-Mail: support@mn-net.com

For orders (international)

E-Mail: sales@mn-net.com