


# Cross-contamination tests – MagnetaPure 32+





The MagnetaPure 32+ is a bench-top robot device for automated nucleic acid extraction from various sample materials such as tissue, blood, food, plants, cell-free body fluids or environmental samples. It is ideally suited for medium sample throughput processing using the MACHEREY-NAGEL NucleoMag® kits.

The NucleoMag® procedure is based on reversible adsorption of nucleic acids to paramagnetic beads under appropriate buffer conditions. The magnetic rods of the MagnetaPure 32+ device transfer the paramagnetic particles through the various purification steps of binding, mixing, washing and elution, thus guaranteeing minimal hands-on time. The MagnetaPure 32+ allows parallel processing of up to 32 samples per run. The purified nucleic acids are of high quality and purity.

In the application data below, cross-contaminations and reagent crossover were evaluated for the MagnetaPure 32+ instrument using the NucleoMag® Pathogen kit as an example.

MagnetaPure 32+	
Description	Automated nucleic acid extraction instrument
Technology	Magnetic rods
Capacity	Up to 32 samples/run
Processing volume	Up to 1000 µL
Features	<ul style="list-style-type: none"> <li>▪ Compact Bench-top robot</li> <li>▪ Easy-to-use touch screen</li> <li>▪ Ready-to-use NucleoMag® scripts</li> <li>▪ Built-in UV lamp for decontamination</li> <li>▪ Built-in heating block for lysis and elution</li> <li>▪ Acoustic signal at program completion</li> <li>▪ Pause function for emergency stop</li> <li>▪ Open and flexible programming</li> </ul>

The data below indicate a cross-contamination free workflow for the processing of MACHEREY-NAGEL NucleoMag® kits on the MagnetaPure 32+ automated extraction robot. No qPCR signals were detected in the wells containing negative controls (RNase-free H<sub>2</sub>O).



### NucleoMag® Pathogen cross-contamination test

Positive (human saliva plus MS2 phage RNA spikes) and negative control (RNase-free H<sub>2</sub>O) samples (200 µL each) were arranged in an alternating pattern in the appropriate position on a 96-well deepwell plate and subjected to the NucleoMag® Pathogen procedure on the MagnetaPure 32+ automated extraction robot. Presence of RNA in the eluates (row 6 and 12) was examined by a qPCR Taqman assay for MS2-phage RNA. All positive samples were successfully amplified and detected (average mean C<sub>t</sub> = 27.06 ± 0.15). Absence of qPCR signal (ND, C<sub>t</sub> undetermined) in the negative control samples indicates a cross contamination-free workflow. Cross-contamination tests were performed in two independent experiments (left side: run 1; right side: run 2) on two different plates each. The data indicate a cross-contamination free workflow for processing NucleoMag® kits on the MagnetaPure 32+ automated extraction robot.