


MACHEREY-NAGEL

User manual



■ NucleoSpin® DNA Trace

November 2024 / Rev. 10

Contact MN

Germany and international

MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0
Toll-free: 0800 26 16 000 (Germany only)
E-mail: info@mn-net.com

Technical Support Bioanalysis

Tel.: +49 24 21 969-333
E-mail: support@mn-net.com

USA

MACHEREY-NAGEL Inc.
924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA
Toll-free: 888 321 6224 (MACH)
E-mail: sales-us@mn-net.com

France

MACHEREY-NAGEL SAS
1, rue Gutenberg – BP135 · 67720 Hoerdt Cedex · France
Tel.: +33 388 68 22 68
E-mail: sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €
Siret 379 859 531 00020 · RCS Strasbourg B379859531 · N° intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG
Hirsackerstr. 7 · 4702 Oensingen · Switzerland
Tel.: +41 62 388 55 00
E-mail: sales-ch@mn-net.com

Genomic DNA from forensic samples

Protocol at a glance (Rev.10)











		Funnel	NucleoSpin® DNA Trace	
1 Lyse sample			4–8 mL FLB 50 µL Proteinase K 56 °C, 1 h	
2 Clarify sample			≥ 5,000 x g, 10 min	
3 Adjust DNA binding conditions			3.5 mL ethanol Vortex	
4 Bind DNA			Load sample 3,000 x g, 3 min	
5 Wash silica membrane			1 st wash 2.5 mL BW 2 nd wash 5 mL B5 3 rd wash 5 mL B5 1 st , 2 nd , 3 rd	3,000 x g, 3 min
6 Dry silica membrane			3,000 x g, 3 min	
7 Elute DNA			100 µL BE (70 °C) RT, 2 min 3,000 x g, 3 min	

Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Reagents, consumables, and equipment to be supplied by user	4
1.3	About this user manual	5
2	Product description	6
2.1	The basic principle	6
2.2	Kit specifications	6
3	Storage conditions and preparation of working solutions	8
4	Safety instructions	9
4.1	Disposal	9
5	Protocols	10
5.1	Isolation of genomic DNA from solid samples, for example small amounts of cells or tissue (forensic samples)	10
5.2	Isolation of genomic DNA from human bones	12
6	Appendix	14
6.1	Troubleshooting	14
6.2	Ordering information	15
6.3	Product use restriction/warranty	16

1 Components

1.1 Kit contents

NucleoSpin® DNA Trace		
REF	4 preps 740942.4	25 preps 740942.25
Lysis Buffer FLB	50 mL	250 mL
Wash Buffer BW	13 mL	75 mL
Wash Buffer B5 (Concentrate)*	12 mL	100 mL
Elution Buffer BE**	13 mL	13 mL
NucleoSpin® DNA Trace F Columns (plus Collection Tubes)	4	25
Proteinase K (lyophilized)*	6 mg	30 mg
Proteinase Buffer PB	1.8 mL	1.8 mL
Collection Tubes (50 mL)	4	25
Elution Tubes (0.5 mL)	4	25
User manual	1	1

1.2 Reagents, consumables, and equipment to be supplied by user

Reagents

- Ethanol (96 – 100 %) to prepare Buffer B5 and to adjust DNA binding conditions)

Consumables

- Disposable pipette tips
- 15 mL and 50 mL centrifugation tubes

Equipment

- Manual pipettors
- Centrifuge with swing-out rotor, suitable for 15 mL and 50 mL tubes
- Suitable homogenization device (e.g., mortar and pestle, rotor-stator)
- Personal protection equipment (e.g., lab coat, gloves, goggles)

* For preparation of working solutions and storage conditions see section 3.

** Composition of Elution Buffer BE: 5 mM Tris/HCl, pH 8.5

1.3 About this user manual

It is strongly recommended that first-time users of the **NucleoSpin® DNA Trace** kit read the detailed protocol sections of this user manual. Experienced users, however, may refer to the Protocol at a glance instead. The Protocol at a glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure.

All technical literature is available on the internet at **www.mn-net.com**.

Please contact Technical Service regarding information about changes of the current user manual compared to previous revisions.

2 Product description

2.1 The basic principle

NucleoSpin® DNA Trace allows DNA isolation from cells, tissue, and many other sources. Lysis is achieved by incubation of homogenized samples in a solution containing chaotropic ions and Proteinase K. Appropriate conditions for binding of DNA to the silica membrane in the **NucleoSpin® DNA Trace F Columns** are created by chaotropic salt and ethanol. The binding process is reversible and specific to nucleic acids. Contaminations are removed by repeated washing with 2 different ethanolic buffers. Pure genomic DNA is finally eluted under low ionic strength conditions in a slightly alkaline elution buffer.

2.2 Kit specifications

- NucleoSpin® DNA Trace kit is designed for the preparation of highly pure genomic DNA from small amounts of any tissue, cells and, forensic samples, for example dried blood spots. The NucleoSpin® DNA Trace F Columns included in the kit are ideally suited for collecting small amounts of nucleic acids from large volumes because these columns are shaped like a funnel combining a large volume capacity with a small diameter of the binding membrane (F means funnel). The DNA isolated by NucleoSpin® DNA Trace F Columns can be used directly for PCR or other enzymatic reactions.
- Age, storage conditions, quantity, and consistency of samples can affect DNA quality, and therefore the protocol may be adapted accordingly (e.g., increasing incubation time). For successful DNA preparation, it is essential that the sample is lysed well and separated afterwards – only clear lysates should be loaded onto **NucleoSpin® DNA Trace F Columns** in order to avoid clogging of the silica membrane.
- The **NucleoSpin® DNA Trace** kit allows purification of up to 20 µg of pure genomic DNA with an A_{260}/A_{280} ratio of between 1.70 and 1.90. Some samples (especially forensic samples) may contain only traces of DNA. However, the amount will be sufficient for amplification and detection reactions.
- Additional enzymes, which are not included in the kit, may be necessary for lysis of certain bacteria (e.g., lysozyme, lysostaphine).
- Support protocol for the isolation of genomic DNA from human bones. For this application additional Buffer T1, Buffer B3, and Proteinase K are necessary. Therefore MACHEREY-NAGEL offers the **NucleoSpin® DNA Trace Bone Buffer Set** (see ordering information). This buffer set is especially designed for completion of the **NucleoSpin® DNA Trace** kit. It is suited for 25 preparations of genomic DNA from human bones in conjunction with the NucleoSpin® DNA Trace kit (REF 740942.25).

Table 1: Kit specifications at a glance

Parameter	NucleoSpin® DNA Trace
Technology	Silica membrane technology
Format	Funnel columns
Sample material	Forensic samples, buccal swabs, blood spots

Table 1: Kit specifications at a glance

Sample size	Forensic samples which can be extracted with up to 8 mL Lysis Buffer FLB (in general 10 mg tissue, < 10 ⁵ cells)
Fragment size	200 bp–approx. 50 kbp
Typical recovery	> 70 % for amounts > 10 ng
A ₂₆₀ /A ₂₈₀	1.7 – 1.9
Elution volume	100 µL
Preparation time	60 min/prep (without Proteinase K incubation time which needs > 1 h)
Binding capacity	20 µg
Use	For research use only

• **Forensic quality product:**

NucleoSpin® DNA Trace is certified as forensic quality product. Consumables used in forensics need to be treated carefully to prevent DNA contamination. MACHEREY-NAGEL therefore has a stringently controlled production process to avoid DNA contamination of consumables. Further, MACHEREY-NAGEL uses ethylene oxide (EO) treatment to remove amplifiable DNA, which might still be introduced during the manufacturing process. MACHEREY-NAGEL products carrying the forensic quality seal, contain plastic materials that are EO treated. This means, DNA of any kind, which might still be introduced into plastic consumables during the production process, is inactivated by means of the treatment with ethylene oxide, in order to prevent the generation of accidental human profile by PCR amplification. Ethylene oxide treatment has been shown to be the method of choice to prevent DNA profiles due to DNA contamination. (Shaw et al., 2008).




Figure 1 According to Shaw et al., 2008, Comparison of the effects of sterilization techniques on subsequent DNA profiling. *Int J Legal Med* 122: 29–33.

3 Storage conditions and preparation of working solutions

Attention: Buffers FLB and BW contain chaotropic salts! Wear gloves and goggles!

CAUTION: Buffers FLB and BW contain guanidine hydrochloride which can form highly reactive compounds when combined with bleach (sodium hypochlorite). DO NOT add bleach or acidic solutions directly to the sample-preparation waste.

- All kit components can be stored at room temperature (18–25 °C) and are stable for at least one year.
- Upon storage, especially at low temperatures, a white precipitate may form in Lysis Buffer FLB. Such precipitates have to be dissolved by incubating at 45–50 °C for 10 min before use.

Before starting any **NucleoSpin® DNA Trace** protocol, prepare the following:

- Wash Buffer B5: Add the indicated volume of ethanol (96–100 %) to **Buffer B5 Concentrate**. Mark the label of the bottle to indicate that ethanol was added. Store Wash Buffer B5 at room temperature (18–25 °C) for at least one year.
- Before first use of the kit, add the indicated volume (see table below or on the bottle) of Proteinase Buffer PB to dissolve lyophilized **Proteinase K**. Proteinase K solution is stable at -20 C for at least 6 months.

NucleoSpin® DNA Trace		
REF	4 preps 740942.4	25 preps 740942.25
Wash Buffer B5 (Concentrate)	12 mL Add 48 mL ethanol	100 mL Add 400 mL ethanol
Proteinase K	6 mg Add 300 µL Proteinase Buffer	30 mg Add 1.5 mL Proteinase Buffer

4 Safety instructions

When working with the **NucleoSpin® DNA Trace** kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles). For more information consult the appropriate Material Safety Data Sheets (MSDS available online at www.mn-net.com/msds).



Caution: Guanidine hydrochloride in buffer BW, can form highly reactive compounds when combined with bleach! Thus, do not add bleach or acidic solutions directly to the sample preparation waste.

The waste generated with the **NucleoSpin® DNA Trace** kit has not been tested for residual infectious material. A contamination of the liquid waste with residual infectious material is highly unlikely due to strong denaturing lysis buffer and Proteinase K treatment but it cannot be excluded completely. Therefore, liquid waste must be considered infectious and should be handled and discarded according local safety regulations.

4.1 Disposal

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

5 Protocols

5.1 Isolation of genomic DNA from solid samples, for example small amounts of cells or tissue (forensic samples)

Before starting the preparation:

- Check that Wash Buffer B5 was prepared according to section 3.
- Preheat Elution Buffer BE to 70 °C.



1 Lyse sample

Place the sample in a 15 mL centrifuge tube (not provided) and add **4–8 mL Buffer FLB**. The sample should be covered completely with Buffer FLB.

Solid samples should be homogenized by commercial tools (pestle and mortar, rotor-stator homogeniser). In general, 10 mg tissue, < 10⁵ cells or any DNA-containing solid sample can be used. Forensic samples (dried blood spots, chewing gum, swabs, etc.) should be covered completely with lysis buffer.


Add **50 µL Proteinase K** stock solution, mix by vortexing, and incubate at **56 °C** in a (shaking) water bath until complete lysis is obtained (**1–3 h or overnight**).

Vortexing every 15 min (3–4 times) leads to shorter lysis times if no shaking water bath/incubator is available. Final incubation at 70–100 °C for 5 min may be recommended for optimal denaturation and lysis of difficult samples (e.g., dried, old or clotted blood samples).




2 Clarify sample

Afterwards, any insoluble particles remaining in the sample have to be removed by centrifugation for **10 min** at **≥ 5.000 x g** in order to avoid clogging of the NucleoSpin® DNA Trace membrane.



3 Adjust DNA binding conditions

Add **3.5 mL ethanol** (96–100 %) to **4 mL cleared FLB-lysate** and vortex the mixture. Use proportionally up scaled volumes of ethanol, if more FLB-lysate has been prepared in step 1.




4 Bind DNA

Pipette mixture onto the **NucleoSpin® DNA Trace F Column**.

Centrifuge for **3 min** at **3,000 x g**. Discard flowthrough with Collection Tube. Put the NucleoSpin® DNA Trace F Column into a fresh Collection Tube (provided).

**Load sample****3,000 x g,
3 min****5 Wash silica membrane****1st wash**

Add **2.5 mL Buffer BW** to the NucleoSpin® DNA Trace F Column. Centrifuge for **3 min** at **3,000 x g**.

**+ 2.5 mL BW****3,000 x g,
3 min****2nd wash**

Add **5 mL Buffer B5** to the NucleoSpin® DNA Trace F Column. Centrifuge for **3 min** at **3,000 x g**, discard flowthrough and reuse Collection Tube.

**+ 5 mL B5****3,000 x g,
3 min****3rd wash**

Add **5 mL Buffer B5** to the NucleoSpin® DNA Trace F Column. Centrifuge for **3 min** at **3,000 x g**, discard flowthrough and reuse Collection Tube.

+ 5 mL B5**3,000 x g,
3 min****6 Dry silica membrane**

Centrifuge additional **10 min** at **3,000 x g** in order to remove **Buffer B5** completely.

**3,000 x g,
10 min****7 Elute DNA**

Attach the supplied Elution Tube (0.5 mL) with adaptor to the NucleoSpin® DNA Trace F Column and insert assembly into a new 50 mL centrifuge tube (not provided). Pipette **100 µL Buffer BE** (preheated to 70 °C) onto the NucleoSpin® membrane and incubate for **2 min** at **room temperature**.

**+ 100 µL BE
(70 °C)****RT,
2 min**

Centrifuge for **3 min** at **3,000 x g** to collect the nucleic acid-containing fraction.

**3,000 x g,
3 min**

Remove the elution tube containing the nucleic acids and keep it for further use

5.2 Isolation of genomic DNA from human bones

Before starting with the preparation, please read remarks below.

- Before starting with the preparation, set incubators or water baths to 56 °C and 70 °C, respectively. Before elution, equilibrate Elution Buffer BE to 70 °C.

Attention:

- The list numbers in this support protocol do not correspond with the list numbers in section 5.1 and protocol at a glance.
- Additional Buffer T1, Buffer B3 and Proteinase K is necessary. The **NucleoSpin® DNA Trace Bone Buffer Set** (REF 740943.25) is especially designed for completion of the NucleoSpin® DNA Trace kit. It is suited for 25 preparations of genomic DNA from human bones in conjunction with the NucleoSpin® DNA Trace kit (REF 740942.25).
- For each prep, 2 mL additional buffer is necessary (0.5 M EDTA / 0.25 M NaPO₄, pH 8, not included in the NucleoSpin® DNA Trace Bone Buffer Set).

1 Prepare sample

Mill 1 g bone to a fine powder.

2 Pre-lyse sample

Add **2 mL buffer** (0.5 M EDTA/0.25 M NaPO₄, pH 8) and **7 mL Buffer T1** and **100 µL Proteinase K solution**. Vortex to mix. Be sure that the samples are completely covered with lysis solution.

If processing several samples, Proteinase K and Buffer T1 may be premixed directly before use. Do never mix Buffer T1 and Proteinase K more than 10–15 min before addition to the sample: Proteinase K tends to self-digestion in Buffer T1 without substrate.

Incubate at **56 °C overnight**.

Afterwards incubate sample for **48 h at 4 °C** on a shaking incubator.

3 Lyse sample

Vortex the samples. Add **8 mL Buffer B3**, vortex vigorously and incubate at **70 °C for 10 min**. Vortex briefly.

Centrifuge for **10 min at 5,000 x g** and transfer the supernatant to a new microcentrifuge tube.

4 Adjust DNA binding conditions

Add **8.4 mL ethanol** (96–100 %) to the sample and vortex vigorously.

5 Bind DNA

For each sample, take one **NucleoSpin® DNA Trace F Column** placed in a Collection Tube (50 mL). Apply the sample successively to the column. Centrifuge for **3 min at 3,000 x g**. Discard the flowthrough and place the column back into the Collection Tube.

6 Wash silica membrane

1st wash

Add **3 mL Buffer BW**. Centrifuge for **3 min** at **3,000 x g**. Discard the flowthrough and place the column back into the Collection Tube.

2nd wash

Add **3 mL Buffer B5** to the column and centrifuge for **3 min** at **3,000 x g**. Discard the flowthrough and place the column back into the Collection Tube.

3rd wash

Add **3 mL Buffer B5** to the column and centrifuge for **3 min** at **3,000 x g**. Discard the flowthrough and place the column back into the Collection Tube.

7 Dry silica membrane

Centrifuge the column for **10 min** at **3,000 x g**.

Residual ethanol is removed during this step.

8 Elute highly pure DNA

Attach the supplied Elution Tube with adaptor to the NucleoSpin® DNA Trace F Column and insert assembly into a new 50 mL centrifuge tube (not provided). Add **60 µL Buffer BE** (preheated to 70 °C). Incubate at **room temperature** for **2 min**.

Centrifuge for **3 min** at **3,000 x g** to collect the nucleic acid-containing fraction.

Remove the elution tube containing the nucleic acids and keep it for further use.

6 Appendix

6.1 Troubleshooting

Problem	Possible cause and suggestions
No or poor DNA yield, poor DNA quality	<i>Incomplete sample lysis</i>
	<ul style="list-style-type: none"> Sample was not thoroughly homogenized and mixed with Buffer FLB / Proteinase K. The mixture has to be shaken continuously. Alternatively, prolong incubation time with Proteinase K.
No or poor DNA yield, poor DNA quality	<i>Reagents not applied properly</i>
	<ul style="list-style-type: none"> Prepare Buffer B5 and Proteinase K solutions according to instructions (section 3). Add ethanol to lysates before loading them on NucleoSpin® DNA Trace F Columns.
No or poor DNA yield, poor DNA quality	<i>Suboptimal elution of DNA from the column</i>
	<ul style="list-style-type: none"> Apply Elution Buffer BE (70 °C) directly onto the center of the silica membrane and incubate for 2 min. Elution efficiencies decrease dramatically, if elution is done with other buffers at $\text{pH} \leq 7.0$.
Poor DNA quality and/or suboptimal performance of genomic DNA in enzymatic reactions	<i>RNA in sample</i>
	<ul style="list-style-type: none"> If RNA-free DNA is desired, add 20 μL of RNase A solution (20 mg/mL) to Lysis Buffer FLB.
Poor DNA quality and/or suboptimal performance of genomic DNA in enzymatic reactions	<i>Carry-over of ethanol</i>
	<ul style="list-style-type: none"> Be certain to centrifuge ≥ 5 min at $3,000 \times g$ in order to remove all of ethanolic Buffer B5 before eluting the DNA. If for any reason, the level of Buffer B5 has reached the column outlet after the second wash, discard flowthrough. Place the NucleoSpin® DNA Trace F Column back into the Collection Tube, and centrifuge again.

6.2 Ordering information

Product	REF	Pack of
NucleoSpin® DNA Trace	740942.4/.25	4/25 preps
NucleoSpin® DNA Forensic	740840.10/.50/.250	10/50/250 preps
NucleoSpin® Funnel Column	740959	30 columns
Buffer FLB	740322.500	500 mL
Buffer BW	740922	100 mL
Buffer B5 Concentrate (for 100 mL Buffer B5)	740921	20
NucleoSpin® DNA Trace Bone Buffer Set	740943.25	1 set
Proteinase K	740506	100 mg
RNase A	740505.50 740505	50 mg 100 mg
NucleoSpin® Forensic Filters	740988.10/.50/.250	10/50/250 pieces
NucleoSpin® Forensic Filters (Bulk)	740988.50B/.250B/1000B	50/250/1000 pieces

Visit www.mn-net.com for more detailed product information.

6.3 Product use restriction/warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety instructions from the respective Material Safety Data Sheet of the product.


This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY-NAGEL literature. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY-NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.

Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:
MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-333
support@mn-net.com



MACHEREY-NAGEL

www.mn-net.com

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

