

MACHEREY-NAGEL

OPTIMA[®] 1701 MS

Chromatography



Lowest bleed and inert mid-polarity column

- Optimiert für die Umweltanalytik (PAHs, PCBs, Pestizide, ...)
- 100 % Ion Trap- und Quadrupol-MS Kompatibilität
- Referenzsäule für Strukturaufklärung, z. B. in Kombination mit OPTIMA[®] 5 MS


MACHEREY-NAGEL

www.mn-net.com




Outstanding features and excellent performance

OPTIMA® 1701



OPTIMA® 1701 MS



The combination of a novel silarylene synthesis strategy with symmetrically substituted cyanopropylsilanes led to improved stability and substantially reduced column bleed.

Silarylene phase with selectivity analog to 14 % cyanopropyl-phenyl – 86 % dimethylpolysiloxane · USP G46

- Reference column for structure identification, e.g., in combination with OPTIMA® 5 MS

Excellent deactivation


- Reliable quantification even for critical samples at ultra trace levels

Mid-polar ultra low bleed silarylene phase

- 100 % Ion-Trap and Quadrupol-MS-compatibility
- It is possible to inject aqueous samples, water stable phase


Special selectivity due to high cyanopropyl content

- Ideal for mid-polar analytes



Inert column · bleed comparison test


MN Appl.-Nr. 215050



| | |
|-----------------------|---|
| Column | OPTIMA® 1701 MS, 0.25 µm film, 30 m x 0.25 mm ID |
| Injection temperature | 280 °C |
| Carrier gas | Helium, 0.8 bar |
| Temperature | 200 °C (20 min) – 225 °C (20 min), 8 °C/min – 250 °C (20 min), 8 °C/min – 275 °C (20 min), 8 °C/min – 300 °C (20 min), 8 °C/min |
| Detector | FID, 280 °C |

Separation of 16 EPA PAH

MN Appl.-Nr. 215070




Peaks:

1. naphthalene, 2. acenaphthylene, 3. acenaphthene, 4. fluorene, 5. phenanthrene, 6. anthracene, 7. fluoranthene, 8. pyrene, 9. benz[a]anthracene, 10. chrysene, 11. benzo[b]fluoranthene, 12. benzo[k]fluoranthene, 13. benzo[a]pyrene, 14. indeno[1,2,3-cd]pyrene, 15. dibenz[a,h]anthracene, 16. benzo[ghi]perylene

| | |
|-----------------------|---|
| Column | OPTIMA® 1701 MS, 0.25 µm film, 30 m x 0.25 mm ID |
| Sample | PAH test mixture acc. to EPA (REF 722314) (20 µg/mL each in toluene) |
| Injection volume | 1 µL, splitless (for 1 min) |
| Injection temperature | 300 °C |
| Carrier gas | Helium, lin. velocity 34 cm/sec |
| Temperature | 90 °C for 1 min – 220 °C, 10 °C/min – 300 °C, 4 °C/min |
| Detector | MSD |

Herbicide mix

MN Appl.-Nr. 215100



Peaks:


1. 3,5-Dichlorobenzoic acid methyl ester (SS), 2. 4-Nitroanisoole, 3. Dicamba methyl ester, 4. Dichlorprop, methyl ester, 5. 2,4-D methyl ester, 6. Pentachloroanisoole, 7. 2,4,5-TP, methyl ester, 8. 2,4,5-T, methyl ester, 9. Chloramben, methyl ester, 10. 2,4-DB methyl ester, 11. Dinoseb methyl ester, 12. Bentazon methyl ester, 13. DCPA methyl ester, 14. Picloram methyl ester, 15. Acifluorfen methyl ester

| | |
|-----------------------|---|
| Column | OPTIMA® 1701 MS, 0.25 µm film, 30 m x 0.25 mm ID |
| Sample | Herbicide mix (400 ng/mL in hexane) |
| Injection volume | 1 µL, splitless (for 1 min) |
| Injection temperature | 300 °C |
| Carrier gas | Helium, lin. velocity 42 cm/sec |
| Temperature | 80 °C – 200 °C, 20 °C/min – 245 °C, 8 °C/min – 260 °C, 20 °C/min |
| Detector | MSD |

OPTIMA® 1701 MS

Triazine pesticide mix

MN Appl.-Nr. 215080



| | |
|-----------------------|---|
| Column | OPTIMA® 1701 MS, 0.25 µm film, 30 m x 0.25 mm ID |
| Sample | Triazine pesticide mix |
| Injection volume | 1 µL, split 1:100 |
| Injection temperature | 250 °C |
| Carrier gas | Helium, lin. velocity 42 cm/sec |
| Temperature | 160 °C for 1 min – 180 °C, 15 °C/min – 220 °C, 2 °C/min |
| Detector | MSD |

Peaks:

1. prometon, 2. atraton, 3. propazine, 4. atrazine, 5+6. simazine + tertbutylazine, 7. secbumetone, 8. prometryn, 9. ametryn, 10. simetryn, 11. terbutryn

Optima® 1701 MS

Max. temperature for isothermal operation 280 °C, max. temperature for short isotherms in a temperature program 300 °C; 0.53 mm ID columns: max. temperatures 260 and 280 °C, resp.

Similar phases:

VF-1701ms, TG-1701MS, OV-1701, DB-1701, HP-1701, Rtx-1701, SPB-1701, CP Sil 19 CB, 007-1701, BP10, ZB-1701

Ordering information


| ID | Length | Film thickness | Agilent J&W VF-1701ms | Restek Rtx-1701 | Phenomenex ZB-1701 | MACHEREY-NAGEL OPTIMA® 1701 MS |
|---------|--------|----------------|-----------------------|-----------------|--------------------|--------------------------------|
| 0.25 mm | 30 m | 0.25 µm | CP9151 | 12023 | 7HG-G006-11 | 726630.30 |
| 0.25 mm | 60 m | 0.25 µm | CP9154 | 12026 | 7KG-G006-11 | 726630.60 |
| 0.25 mm | 30 m | 0.50 µm | – | 12038 | – | 726631.30 |
| 0.25 mm | 60 m | 0.50 µm | – | 12041 | – | 726631.60 |
| 0.25 mm | 30 m | 1.00 µm | CP9152 | 12053 | 7HG-G006-22 | 7266631.30 |
| 0.25 mm | 60 m | 1.00 µm | CP9156 | 12056 | – | 726632.60 |
| 0.32 mm | 30 m | 0.25 µm | CP9162 | 12024 | 7HM-G006-11 | 726633.30 |
| 0.32 mm | 60 m | 0.25 µm | CP9165 | 120240 | 7KM-G006-11 | 726633.60 |
| 0.32 mm | 30 m | 0.50 µm | – | 12039 | – | 726634.40 |
| 0.32 mm | 60 m | 0.50 µm | – | 12042 | – | 726634.60 |
| 0.32 mm | 30 m | 1.00 µm | CP9163 | 12054 | 7HM-G006-22 | 726635.30 |
| 0.32 mm | 60 m | 1.00 µm | CP9166 | 12057 | – | 726635.60 |

All used names and denotations can be brands, trademarks or registered labels of their respective owner – also if they do not have a special denotation.

www.mn-net.com

MACHEREY-NAGEL






TÜVRheinland
CERTIFIED

Management System
EN ISO 13485:2016
ISO 9001:2015

www.tuv.com
ID 0000056401



MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com