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Description:

DiaQuant[®] Chlorine are test strips for the semi-quantitative detection of total chlorine in rinse solutions.

Measuring range:

0.1–10 mg/L Cl₂

Color gradation:

0 · 0.1 · 0.5 · 1 · 3 · 10 mg/L Cl₂

Pack content:

100 test strips

Reaction principle:

Total chlorine reacts with the potassium iodide and the organic redox indicator in the test field, forming a green-colored oxidation product.

General information:

Only remove the required number of test strips from the container, and then securely close the container immediately afterwards. Do not touch the test field. Ideally, the color evaluation should be carried out by diffuse daylight. Artificial light can lead to difficulties or an incorrect color reading.

Instruction for use:

1. Insert the test strip into the sample solution, moving it back and forth for 15 s.
2. Shake off any excess liquid.
3. Compare the test field with the color scale. Take the value which matches closest with the colored test field (reading accuracy: ± ½ colored field on the scale).

The reaction color of the test field may change after the value has been taken. It is therefore crucial to evaluate the coloration within the prescribed time scale in order to achieve correct results.

Quality control:

To check the correct functioning of the test strips, use a total chlorine solution with a concentration of 1 mg/L. For this purpose, first prepare a stock solution of 100 mg/L chlorine by adding 0.1 g calcium hypochlorite to 1000 mL distilled water while continuously stirring. Then, filter the solution and subject it to photometric analysis (i.e. with *NANOCOLOR[®]* Chlorine/Ozone 2, REF 985017). Dilute 1 mL from the stock solution (100 mg/L) in 100 mL distilled water and then add one spatula tip of ammonium chloride (= 1 mg/L Cl₂). Immediately perform the measurement with the test strip. If the control solution produces a negative result even after repeating the process, then the remaining unused test strips must be discarded. Even during a negative control (inserting a test strip into distilled water), no positive coloration may occur. Possible reasons for incorrect functioning of the test strips may be that the use-by-date has been exceeded, the container was left open for too long or was stored incorrectly.

Interferences:

The presence of other strong oxidants such as bromine, iodine, and hydrogen peroxide will also lead to false positive results.

Storage:




Avoid exposing the test strips to sunlight or moisture. Keep container cool and dry (storage temperature between 4 °C and 30 °C). If stored correctly, the test strips may be used until the use-by-date printed on the packaging.

Additional information:

The test strip container stopper contains a non-toxic drying agent. If swallowed, drink plenty of water.

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

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