REF 985674

Robot

de

Test 6-74

03.23

**NANOCOLOR® Phenol-Index 5**

#### Methode:

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**NANOCOLOR® Phenolic index 5**

#### Method:

Photometric determination of phenols and other compounds capable of oxidative coupling, which form antipyrene dyes with 4-aminoantipyrene

Range:	0.2 – 5.0 mg/L phenolic index
Wavelength:	520 nm

#### Contents of reagent set:

20 test tubes Phenolic index 5

1 bottle with 40 mL Phenolic index 5 Reagent A

1 test tube with 200 mg Phenolic index 5 Reagent B1

1 test tube with 5 mL Phenolic index 5 Reagent B2

*Note: The content of the 20 test tubes Phenolic index 5 may be colored lightly yellow, which has no influence to the color reaction.*

#### Hazard warning:

Test tubes contain 4-aminoantipyrene 25 – 100 %, reagent A contains ammonia 1 – 5 %, reagent B1 contains potassium peroxodisulfate 20 – 100 %.

H317, H334 May cause an allergic skin reaction. May cause allergy or asthma symptoms or breathing difficulties if inhaled.

P261, P272, P280, P302+352, P304+340, P333+313, P342+311, P363 Avoid breathing dust. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection.

IF ON SKIN: Wash with plenty of water / ... IF INHALED: Remove person to fresh air and keep comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER/doctor / ... Wash contaminated clothing before reuse. For further information ask for a safety data sheet.

#### Interferences:

- oxidizing substances
- reducing substances
- cyanides

The method can not be applied for the analysis of sea water.

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**MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Allemagne**
**France : MACHEREY-NAGEL SAS** · 1, rue Gutenberg – BP135 · 67720 Hoerdt · France
MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €