MACHEREY-NAGEL

URYXXON® Relax VET

- Bedienungsanleitung
- User manual
- Mode d'emploi
- Manual de instrucciones

Inhalt / Contents / Contenu / Contenido

Bedienungsanleitung.....	3
User manual	55
Mode d'emploi	103
Manual de instrucciones	155

MACHEREY-NAGEL

URYXXON® Relax VET

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URYXXON® Relax VET

User manual


1. Quickstart	60
2. Introduction	62
2.1 Intended Use	62
2.2 Intended User.....	62
2.3 Indications for use	62
2.4 Patient population	62
2.5 Type of specimen	62
2.6 System description	63
2.6.1 Measuring principle	63
2.6.2 Functional principle URYXXON® Relax VET	63
2.7 Safety warnings	64
2.8 Labels on the device and the box	65
2.9 Symbols explanation	66
3. Unpacking and set up	67
3.1 List of delivered parts	67
3.2 Consumables	67
3.3 Note on ambient surrounding	68
3.4 Setting up the instrument	68
3.5 Description of instrument parts	68
3.6 How to plug the instrument in	69
3.7 How to start the instrument	69
3.8 How to switch off the instrument	69
3.9 How to load the printer paper	70
3.10 How to install batteries (optional)	70
3.11 Instrument self test	71
3.12 Calibration	71
3.13 How to prevent biological contamination	71
3.14 Use of the instrument	71
3.14.1 Buttons	71
3.14.2 Scroll Menu	72
3.14.3 Round buttons	72
4. User menu	73
4.1 Flow-chart of the menu structure	73
4.2 Description of the menu items	74
5. Analysis of test strips	75
5.1 How to perform a measurement	75
5.2 Display of results	76
5.3 Measurement errors	76
5.4 Entering the patient identification.....	76
5.5 Changing the sequence number ("SN")	77
5.6 Transferring data to a PC	77

6. Enter the main menu	78
7. Recall results	79
7.1 How to scroll through memory.....	79
7.2 How to find specific results (filtering)	79
7.2.1 Selecting the date	80
7.2.2 Selecting search criteria.....	80
7.2.3 Display suitable matches	81
7.3 How to delete results from memory	81
8. Quality control testing	82
8.1 How to review old QC measurements	82
8.2 How to perform a QC measurement	82
9. Equipment Settings	83
9.1 How to modify strip settings.....	83
9.1.1 Type	83
9.1.2 Units.....	83
9.1.3 Order of Parameters	83
9.1.4 Test Strip LOT	83
9.2 How to protect settings from unauthorized access	83
9.3 How to turn the printer on and off	84
9.4 How to enable and disable acoustic signals	84
9.4.1 Acoustic confirmation of user inputs	84
9.4.2 Acoustic warning on positive results	84
9.5 How to deactivate and activate the autostart.....	84
9.6 How to set energy saving options for the battery mode	84
9.7 How to change the language	84
9.8 How to set time and date	84
9.9 How to activate data transfer	85
9.10 How to change the text of the printout header	85
9.11 How to print the settings	85
10. Cleaning and maintenance	86
10.1 How to clean the housing	86
10.2 How to clean the strip holder.....	86
11. Service menu	87
11.1 How to reset the system (Load Default)	87
11.2 How to control the LOT-control (LOT activate).....	87
11.2.1 Entry of test strips' LOT	87
11.2.2 Warning on expired test strips	88
11.3 How to change the sensitivity	88
11.3.1 Settings Reset	88
11.3.2 Appointing new Sensitivity Settings	89
11.3.3 Reported value – (original threshold) – modification.....	89

11.4 How to update the instrument	89
11.5 Control strip	90
12. Interface description	91
12.1 Serial interface	91
12.2 USB 1.1-interface	91
12.3 Transmission protocol	91
12.4 Barcode scanner, PC-keyboard	91
12.5 Ethernet-Interface	91
13. Error Messages and Fault Clearance	92
14. Warranty	93
14.1 General usage comments	93
15. Technical information	94
15.1 Technical data	94
15.2 CE Certification	94
15.3 Analytical Performance Characteristics	95
15.3.1 Table of results	95
15.3.2 Precision	96
15.3.3 Linearity	97
15.3.4 Cut-off concentration	97
15.3.5 Comparative study	98
15.3.6 Interference study	100
15.5 Electromagnetic Compatibility	101
15.6 Waste disposal	101
15.7 CAN ICES	101
16. Additional Information	102
16.1 Manufacturer Information	102
16.2 Technical support	102

1. Quickstart

Unpack the instrument and place it on an even, hard surface. Connect the power supply and turn the equipment on with the On/Off-switch (Pic 10-⑧). After the self test the start screen will appear on the display.


Display 1: Start menu

- Dip a test strip into the urine sample for approx. one second (Pic 1).
- Blot by touching the edge of the strip to a paper towel to remove excess urine (Pic 2).

Pic 1: Dip test strip

Pic 2: Remove excess urine

CAUTION

Biological Hazard: Urine residues and used test strips bare the danger of infection. Always use protective gloves during handling and disposal. The disposal of used test strips should be performed according to the regulations of the handling of potentially infectious material.

CAUTION

Please make sure to remove access of urine from the test strip to avoid a carryover into the instrument.

- Place the strip on the strip holder (Pic 3)

Pic 3: Correct

Pic 4: Wrong

- Slide or push the strip to the end of the strip holder. Do not touch the reagent pads on the test strip.

The instrument will automatically detect an applied strip. The measurement cycle will be started. A progress bar on the display shows the remaining measurement time.

NOTICE

If "Autostart" (Chapter 9.5 "How to deactivate and activate the autostart") is deactivated, the measurement must be started using the start control panel .

NOTICE

The strip will be drawn into the instrument after 30 sec.

At the end of the measurement the result will be displayed on the screen and transferred to the printer and interfaces.

RESULT		Combi 10 VET	
Seq. No:	0001		
Pat. ID:	102358		
2022-06-26	17:19		
BLD	NEG		
UBG	NORM		
BIL	NEG		
PRO	NEG		
NIT	NEG		
KET	NEG		
GLU	NEG		
pH	5		
SG	1.030		
LEU	NEG		

Display 2: Result

By pressing the printer symbol  the result can be printed again. Choosing the return panel  will lead back to the start screen.

Another analysis may be started by applying the next test strip.

NOTICE

To start a new measurement it is not necessary to go back to the start screen. A new strip is detected at any time and the measurement is then started automatically.

2. Introduction

The URYXXON® Relax VET is a portable and easy to use urine analyzer for professional use. It is a *veterinary* diagnostic device designed to read test strips for urine analysis produced by MACHEREY-NAGEL only. The URYXXON® Relax VET provides semi-quantitative values, which may be displayed, printed and transferred. The tests performed with the URYXXON® Relax VET are intended for *veterinary* diagnostic use only.

2.1 Intended Use

The URYXXON® Relax VET is a portable reflection photometer for the evaluation of urine test strips. The URYXXON® Relax VET is intended for use with compatible test strips from MACHEREY-NAGEL. The *veterinary* analysis system is an aid to diagnoses in the areas of kidney function, urinary tract infections, metabolic diseases (e.g. diabetes mellitus) and liver function. The URYXXON® Relax VET is intended for use by professionals. The URYXXON® Relax VET enables integration into laboratory information and management systems (LIMS).

2.2 Intended User

The URYXXON® Relax VET is intended for use by trained professionals with experience in human urine sampling and is not suitable for self-testing.

2.3 Indications for use

For indication regarding the indications for use refer to the instructions for use of the used strips.

2.4 Patient population

There are no restrictions regarding the patient population with MACHEREY-NAGEL urine test strips.


2.5 Type of specimen

Testing of fresh, native, well-mixed and non-centrifuged urine is recommended. Samples should be protected from light. Only urine samples which have stood for a maximum of 4 hours should be used.

2.6 System description

2.6.1 Measuring principle

The test strip moves below a fixed measuring head on a sled with an embedded reference pad. The reflectometric analysis of the test strip and the reference field take place during withdrawal and release of the sled.

Pic 5: Measuring Principle

The strip is illuminated with an LED and a detector registers the intensity of light reflected by the test strip at three different wavelengths. Using an internal calibration, the intensity of the reflected light is correlated to the concentration of the analyte. Whenever samples are strongly alkaline, a density correction is automatically conducted.

2.6.2 Functional principle URYXXON® Relax VET

A measurement is started by placing a strip on the holder. If the Autostart-feature is turned off, the measurement is started by pressing the start panel  in the display. The result is shown on the display, printed out and released via the interfaces after the measurement has been completed. After three minutes the instrument will go to stand-by. Touching the screen will reactivate the instrument.

All user inputs are performed via the touch-screen (Chapter 3.14 "Use of the instrument").

2.7 Safety warnings

BE SURE TO READ THE FOLLOWING SAFETY INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE.

Failure to follow these instructions may result in serious injury to the operator, malfunction or damage of the device.

Keep this manual in a safe place for future reference.

Follow the safety notes and instructions in the operating manual and observe the stickers and notices on the device.

In case of a test strip sled blockage, switch off the device and turn it back on after 10 seconds.

If the blockage persists, do not continue using the device and contact technical support.

If liquids have entered the interior of the device, do not use it any further and contact technical support.

If there are any cracks or holes in the housing, stop using the device and contact technical support.

If the software freezes, switch off the device and turn it back on after 10 seconds. If the software

remains unresponsive, do not continue using the device and contact technical support.

Do not work on internal parts of the unit. Non-compliance will invalidate any warranty claim.

The following safety warnings are used throughout the manual:

DANGER

Indicates an imminent or potentially hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Indicates an imminent or potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may result in malfunction or damage of the instrument. Information that clarifies facts in the text and that requires special attention.

2.8 Labels on the device and the box

Pic 6: Exemplary carton label

Pic 7: Exemplary type plate label

2.9 Symbols explanation

	Batch code
	Biological risk
	Catalogue number
	Caution
	CE Mark
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for < n > tests
	Date of Manufacture
	Device for near-patient testing
	Device not for self-testing
	Direct current
	Do not dispose of the device in common household waste
	Do not reuse
	Do not use if package is damaged and consult instructions for use
	Keep away from sunlight
	Keep dry
	Manufacturer
	Non-sterile
	Serial number
	Temperature limits
	Use-by date

3. Unpacking and set up

3.1 List of delivered parts

- ① URYXXON® Relax VET reflectometer
- ② Power pack 100 – 240 V, 47/63 Hz, 9 V, incl. adapter
- ③ Printer paper
 - + User manual (this booklet)
 - + Quick guide

Pic 8: Content

Check all delivered parts for visible damages. In case a part is damaged, please contact your local distributor or MACHEREY-NAGEL.

Read the operating manual for URYXXON® Relax VET carefully before the first startup in order to ensure an error free operation.

NOTICE

Position the device with at least 5 cm of free space at the back to ensure safe access to the power switch.

3.2 Consumables

- Printer paper, pack of five rolls, REF 93065
- Check solution for quality controls, Medi-Test Control, 1 test tube with 15 mL reagent solution Medi-Test Control N and 1 test tube with 15 mL reagent solution Medi-Test Control P, REF 93038
- 6 AA type batteries (optional)

The following test strips can be evaluated with the URYXXON® Relax VET:

- Medi-Test Combi 10 VET REF 930870

NOTICE

Only use consumables within their expiration date.

3.3 Note on ambient surrounding

If the device is exposed to higher temperature fluctuations (e.g. after transport or distribution), it must be switched on not before sufficient acclimatization is given. The device should not be used close to electrical fields (e.g. by microwaves, radio units et cetera). In worst case the measurement results can be affected.

3.4 Setting up the instrument

Place the instrument on a hard, even surface where humidity and temperature are fairly constant. Make sure that the instrument is allowed to acclimate to room temperature prior to use.

Make sure that you

- Do not place the instrument near strong electromagnetic fields
- Do not place the instrument near heating plates, ovens or radiators
- Do not expose the instrument to strong light sources (i. e. direct sunlight)

3.5 Description of instrument parts

Pic 9: Front view

Pic 10: Backside view

Actuator	Function
1. Touch-Screen	Control of equipment functions
2. Test Strip Slide	Test strip retainer and autonomous start of analysis
3. Printer Flap	Opening the printer flap for paper replacement
4. Serial Interface	Connection of a computer (cable length ¹ max. 3 m)
5. USB - Interface	Connection of a computer (cable length ¹ max. 3 m)
6. USB - Interface	Connection of a keyboard or a bar code scanner ²
7. Mains Connection	Contact for the provided power pack
8. On/Off Switch (I/O)	Turning the equipment on and off
9. LAN-Interface	Connection to network (serial over LAN)

¹ We recommend to use shielded cables.

² We recommend to make use of the hand scanner from MACHEREY-NAGEL

CAUTION

Do only connect the intended devices to the corresponding interface. If you connect another device to the corresponding interface, the analyzer or the connected device may be damaged. Check all cables prior to use and verify the proper connection.

3.6 How to plug the instrument in

WARNING

Risk of electric shock. Only use the included power supply. Use of other power supplies can damage the device and harm the user.

WARNING

Risk of electric shock. Cable break can lead to exposure to conduction cables. Stop using the device immediately and remove the power supply from the outlet if possible.

Pic 11: Power pack

Pic 12: DC in

Four adapters are provided for adapting the power pack to the available mains connection. The sealing protecting the connector pins is removed and the adapter matching the mains connection is plugged on to the power pack (Pic 11). After plugging the power pack cable into the jack "DC IN" (Pic 10-⑦) and connecting the power pack to the power socket the URYXXON® Relax VET is ready for operation. To ensure safe functioning of the instrument, only use the included power supply!

3.7 How to start the instrument

You can find the On/Off switch on the back of the instrument. Turn on the instrument by pressing the switch. Right after switching on, the instrument performs a self-test. The touch screen lights up and the start screen appears. If a serious error is detected during the self-test, no measurements can be performed. In this case, contact the MACHEREY-NAGEL technical service.

3.8 How to switch off the instrument

Always switch off the device after use by pressing the on/off switch. If you do not use the device for a longer period of time, disconnect the device from the power outlet.

3.9 How to load the printer paper

Pic 13: Printer A

Pic 14: Printer B

Open the printer flap by pressing the rectangular key next to the printer flap (Pic 13).

Pic 15: Printer C

Pic 16: Printer D

Unroll the paper roll by 5 cm and place the roll in the paper compartment with the end on the lower side. Fix the end of the paper to the housing with your finger while closing the flap (Pic 15 + Pic 16).

3.10 How to install batteries (optional)

The URYXXON® Relax VET can be operated with 6 type AA batteries independent of the mains supply. The battery compartment is on the underside of the equipment. Notice the designated polarity (+/-) marked on the battery compartment while inserting the batteries.

Pic 17: Battery compartment

3.11 Instrument self test

The instrument will perform an automatic self test each time it is turned on. If an error message appears, the instrument will not start measurements. In this case, please contact your local service provider.

3.12 Calibration

The instrument will perform an automatic calibration each time a test is run.

3.13 How to prevent biological contamination

Using the instrument makes contact to urine inevitable. Therefore, use protective gloves whenever measuring with the instrument and be cautious when handling urine samples.

CAUTION

Biological Hazard: Urine residues and used test strips bare the danger of infection.

Always use protective gloves during handling and disposal. The disposal of used test strips should be performed according to the regulations of the handling of potentially infectious material.

3.14 Use of the instrument

All user inputs are done via a touch-screen (touch-display). All functions are activated directly by slight pressure with the finger on explicit pictograms or text representing the menu items.

3.14.1 Buttons


Framed areas react to pressure and trigger the action linked to it. The caption of an area describes its function.

Examples:

	Confirm / perform action
	Cancel action
	Equipment settings
	Check mode
	Memory
	Enter main menu

3.14.2 Scroll Menu

Press the up-and-down arrows on the right side of the screen to scroll through a list of information on the left side of the screen. The desired information on the left side is highlighted.

Display 3: Selective lists

Pressing will  select the highlighted line. You can leave the menu by pressing .

3.14.3 Round buttons


These buttons typically appear on screens that require a selection among serial items. The button with a filled circle is the current selection.

Pressing the circle will activate a selection. Save your selection by pressing . Pressing  will quit the menu without performing any changes.

4. User menu

4.1 Flow-chart of the menu structure

4.2 Description of the menu items

- SN:
Chapter 5.5 "Changing the sequence number ("SN")"
- ID:
Chapter 5.4 "Entering the patient identification"
- (P):
Standby
- ▼:
Main menu
- Memory:
Chapter 7 "Recall results"
- Check mode:
Chapter 8 "Quality control testing"
- Settings:
Chapter 9.1 "How to modify strip settings"
Chapter 9.2 "How to protect settings from unauthorized access"
Chapter 9.3 "How to turn the printer on and off"
Chapter 9.4 "How to enable and disable acoustic signals"
Chapter 9.5 "How to deactivate and activate the autostart"
Chapter 9.6 "How to set energy saving options for the battery mode"
Chapter 9.7 "How to change the language"
Chapter 9.8 "How to set time and date"
Chapter 9.9 "How to activate data transfer"
Chapter 9.10 "How to change the text of the printout header"
Chapter 9.11 "How to print the settings"

5. Analysis of test strips

5.1 How to perform a measurement

The URYXXON® Relax VET is very easy to use. In order to start the measurement, the test strip is placed on the strip holder. The instrument automatically detects a new strip and starts the measurement. A progress bar appears, that indicates the remaining analysis time. After 30 seconds the test strip is drawn into the instrument.

NOTICE

Make sure to remove excess urine by blotting the test strip carefully on a lint-free cloth (Chapter 9.5 "How to deactivate and activate the autostart").
For details on sample handling please refer to Chapter 1 "Quickstart".

NOTICE

If auto mode (Chapter 9.5 "How to deactivate and activate the autostart") is deactivated, the analysis needs to be started by pressing  on the touch-screen.

After the measurement, the instrument will release the analyzed test strip which can now be discarded. The result is displayed on the screen and is transferred via the interfaces and/or printed according to equipment settings.

For additional information on the test strip, please read the instruction leaflet that comes with the strips.

CAUTION

Biological Hazard: Urine residues and used test strips bare the danger of infection. Always use protective gloves during handling and disposal. The disposal of used test strips should be performed according to the regulations of the handling of potentially infectious material.

5.2 Display of results

The sequence number (Seq.No.) as well as the patient identification (ID) will be displayed with the results.

RESULT		Combi 10 VET	
Seq. No:	0001		
Pat. ID:	102358		
2022-06-25	17:19		
BLD	NEG		
UBG	NORM		
BIL	NEG		
PRO	NEG		
NIT	NEG		
KET	NEG		
GLU	NEG		
pH	5		
SG	1.030		
LEU	NEG		

Display 4: Result

Positive findings are clearly marked by an asterisk (*) on the printout and on the display. Additionally, positive findings are marked red on the display. It is also possible to enable an acoustic signal on positive findings. Additionally, it is possible to enable an acoustic signal for the findings "Repeat with new sample", "abnormal" and "high abnormal".

The printout is light-sensitive and may turn yellow when exposed to light during storage. For archiving purposes the printouts should be kept in a dark place (patient file) or as a photocopy.

The result displayed may be printed again by pressing . The return panel  will lead back to the start screen.

5.3 Measurement errors

If the display shows a measurement error (e.g. wrong strip) instead of a result please read the instructions in Chapter 13 "Error Messages and Fault Clearance".

Repeat the measurement. In case of permanent errors please contact MACHEREY-NAGEL service.

Measurement errors will be printed and send according to the instrument settings.

5.4 Entering the patient identification

The patient identification needs to be entered before starting the analysis. This can be done as follows:

- Directly on the equipment: Pressing  in the start menu brings up an alphanumeric keypad. Enter the ID using the keys. To enter characters (i. e. "Miller") press  to change the character entry. Repeated pressing on the same field within 0.5 seconds switches through the characters displayed on the key. Wrong entries may be erased by pressing 

Pat. ID:			
123578901			
1	2	3	ABC
4	5	6	←
7	8	9	0
✗		✓	

Display 5: Entering ID (numeric)

Pat. ID:			
ABCDEF			
—	ABC	DEF	abc
GHI	JKL	MNO	←
PQRS	TUV	WXYZ	-(#)
✗		✓	

Display 6: Entering ID (letter)

- Using a standard PC-keyboard: Connect the keyboard to the instrument. User inputs on the keyboard will automatically be interpreted as patient identifications.
- Using a bar code reader: Connect the barcode reader to the instrument. Barcode readings will automatically be interpreted as patient identifications.

After entering the patient identification start the measurement. The patient identification is saved together with the diagnostic findings.

NOTICE

A new ID cannot be entered before the present analysis has been completed.

5.5 Changing the sequence number ("SN")

Pressing **SN:** in the start menu brings up a numerical pad. Enter a new sequence number using the keys on the pad. All following measurements will now be counted from this number on.

Seq. No:			
1234			
1	2	3	
4	5	6	←
7	8	9	0
✗		✓	

Display 7: Seq.-Input


5.6 Transferring data to a PC

The results may be transferred to a PC via the USB- or RS232-interface.


A detailed description of the interface can be found in Chapter 12 "Interface description".

6. Enter the main menu

Pressing  on the start screen will bring up the main menu.

Display 8: Start menu


Display 9: Main menu

From here the other functions e.g. memory (Chapter 7 "Recall results"), the test mode (Chapter 8 "Quality control testing") as well as the settings (Chapter 9 "Equipment Settings") can be reached.


7. Recall results

The URYXXON® Relax VET has memory for 200 measurements. Every result is automatically saved after the analysis. After 200 measurements, new data will overwrite the eldest saved dataset.

Access the memory by pressing  in the main menu .


Display 10: Main menu

Display 11: Memory mode

7.1 How to scroll through memory

Pressing  will bring up Display 12. Scrolling through the memory is possible by pressing the arrows on the right side. The next  or previous  result will be displayed.


Display 12: Memory contents

It is possible to print  and send  the dataset displayed. The memory menu will reappear upon pressing return .

7.2 How to find specific results (filtering)


To find a result, you may select the date of the measurement and a specific parameter (Display 13).

Display 13: Filtering

7.2.1 Selecting the date

By pressing  you will reach the menu displayed below (Display 14).


Display 14: Select day

Set the day with the buttons. Selecting "Date" will bring up a screen with the list of available dates (only days with measurements are shown on the screen). Select the desired date with the up-and-down arrows and confirm your selection by pressing . After confirmation your selection will be displayed on the screen "SET FILTER PAR."

7.2.2 Selecting search criteria

Pressing  in Display 13 brings up Display 15.


Display 15: Select parameter

Use the arrow keys to select the desired criteria and confirm with . The filter criteria will be displayed on the filter settings screen (Display 16).

7.2.3 Display suitable matches


After setting "Day" and "Parameter" the search may be started by pressing  (Display 16).

Display 16 is a filtering interface titled "SET FILTER PAR.". It contains two input fields: "DAY" with the value "Today" and "PAR" with the value "Positive". At the bottom, there are two buttons: a red "X" button on the left and a green checkmark button on the right.

Display 16: Filtering

When suitable matches are found, an option will appear, allowing to print the datasets , send them to a PC  or display them on the screen  (Display 17).

Display 17 is a search result interface titled "1 ENTRIES FOUND". It contains three options, each with an icon and a label: "Print" with a printer icon, "Send" with a PC icon, and "Display" with a monitor icon. At the bottom, there is a red "X" button.

Display 17: Search result

If no matching results are found, the equipment returns to the memory menu.

7.3 How to delete results from memory


Pressing  will delete all data in the memory. You need to confirm this again on a further screen. The quality control measurements are not affected by this action.

8. Quality control testing


Quality control measurements should be performed regularly with check solutions in order to ensure the correct functioning of the combination of equipment and test strips.

Every facility has to implement its own QC policy.

Upon pressing  on the start menu the "Check mode" screen will appear (Display 18 + Display 19).

Display 18: Main menu

Display 19: Check mode


8.1 How to review old QC measurements

The equipment saves the results of the last 20 QC measurements in a separate memory. They may be displayed by pressing  (Display 19) and may be printed for documentation purposes.

8.2 How to perform a QC measurement

Prepare the urine controls as described in the instruction leaflet and test them in test mode. Do not use water as negative control.

Handle the check solutions exactly as the patient samples. By pressing  (Display 19) the sample selection will appear. You may select what type of control sample you want to analyze either positive or negative (Display 20).


Display 20: Check mode


In case several results do not match the expected results (indications in the checking solutions' package inserts), please contact the MACHEREY-NAGEL service.

9. Equipment Settings

Enter the main menu by pressing the menu key  to reach the "SETTINGS" display press  (Display 21 + Display 22).

Display 21: Main menu

Display 22: Settings

9.1 How to modify strip settings

In the "Strip" menu, settings for result displaying and the test strip type can be chosen.

9.1.1 Type

The URYXXON® Relax VET is in principle able to measure several different strip types. When the equipment has data for various approved strip types, the strip type can be chosen in a selective list.

9.1.2 Units

The URYXXON® Relax VET can report the results in different units:

- Conventional (e.g. 10 mg/dL)
- SI (e.g. 56 mmol/L)
- ARB, Plus-System (e.g. +++)
- Conventional + ARB
- SI + ARB

Choose the desired unit from the selective list.

9.1.3 Order of Parameters

The output order of the parameters may be customized via a selective list. The parameters need to be selected in the desired order and confirmed by pressing . After the last parameter the equipment will ask whether the setting is to be saved. Save by pressing  or go back to the preprogrammed order by pressing .

9.1.4 Test Strip LOT

The LOT-administration of URYXXON® Relax VET is deactivated in the default setting. See Chapter 11.2 "How to control the LOT-control (LOT activate)".

9.2 How to protect settings from unauthorized access

Select "Password" in the settings menu if you want to protect the instrument settings with a PIN. An option panel with the possibilities "ON" and "OFF" will be displayed. Selecting "ON" will enable the PIN-protection.

After enabling the PIN-protection a numerical pad will appear. Enter a 4-digit PIN and confirm by pressing . The PIN will be asked for upon the next attempt to change the settings.

CAUTION

A forgotten PIN cannot be reconstructed. Only a complete reset of the instrument will delete the PIN-protection. This will result in loss of all results and settings!

9.3 How to turn the printer on and off

Selecting "Printer" in the settings menu will bring up an option panel. Choose the desired option and confirm.

9.4 How to enable and disable acoustic signals

Select "Sound" in the settings menu to enter the settings for acoustic signals.

9.4.1 Acoustic confirmation of user inputs

In the preprogrammed settings, all user inputs are confirmed with an acoustic signal. Disable or enable these signals by choosing "ON" or "OFF" in the box "Touch".

9.4.2 Acoustic warning on positive results

In the preprogrammed settings, an acoustic signal will be given on positive findings. Disable or enable this signal by choosing "ON" or "OFF" in the box "Positive".

9.5 How to deactivate and activate the autostart

Select "Autostart" in the settings menu. In basic mode, URYXXON® Relax VET automatically detects an applied test strip and starts the measurement. This function may be deactivated via an option panel.

If autostart is deactivated, the analysis must be triggered by pressing a panel in the Start menu.

9.6 How to set energy saving options for the battery mode

Select "Battery" in the settings menu. Settings in this menu will only apply when the instrument is operated with batteries.

To increase the lifetime of the batteries, the LCD backlight and the printer can be turned off using the option panels.

9.7 How to change the language

Select "Language" in the settings menu. The language of the URYXXON® Relax VET menu can be switched to the following languages using the respective selective list:

English, Deutsch, Espanol, Francais, Italiano, Portugues, Polski, Türkce, Nederlands, Norsk, Svenska, Suomi, Dansk, Indonesia

9.8 How to set time and date

Select "Date/Time" in the setting menu. To change, press on the respective number. A numerical pad appears. Enter the correct number and confirm by pressing .


The date may be formatted in three ways. The active format is shown on a button  on the right hand side of the date (Display 23). Select a date format by pressing this button.

Displayed Format Abbreviation	Meaning	Example
YMD	Year - Month - Day	2007-12-17
DMY	Day . Month . Year	17.12.2007

MDY	Month / Day / Year	12/17/2007
-----	--------------------	------------

The time format may be changed to 12 or 24 hours with the button displayed next to the time

24

Display 23: Date / time

9.9 How to activate data transfer

Select "Interface" on the settings menu. The data transfer via the interfaces can be activated or deactivated via an option panel.

9.10 How to change the text of the printout header

Select "Customization" in the settings menu. The first two lines of the printout may be filled with a user-specific identifier. Each line contains 23 characters.

To enter the text an external keyboard or the alphanumerical pad on the display can be used. The keys on the touch-screen are linked to several letters. Repeated pressing within half a second switches through the letters displayed on the key.

9.11 How to print the settings

Select "Print Settings" in the settings menu to print the equipment settings for documentation purposes. Thermo printings fade with time. Therefore, please store it in a dark place or copy the printout.

10. Cleaning and maintenance

CAUTION

Biological Hazard: Urine residues and used test strips bare the danger of infection. Always use protective gloves during handling and disposal. The disposal of used test strips should be performed according to the regulations of the handling of potentially infectious material.

10.1 How to clean the housing

Use commercially available disinfection wipes to clean and disinfect the instrument from the outside.

CAUTION

Do not use liquids that can enter the instrument as these may cause permanent damage of the instrument.

Clean the surface of the instrument with a disinfection wipe once daily.

10.2 How to clean the strip holder

Wipe off urine residues from the strip holder with a lint-free cloth after each measurement. This prevents carry-over and drying of urine residues.

To remove the strip holder for more comprehensive cleaning, make sure to turn the instrument off. Pull the strip holder from the metal sled below (Pic 18). Wipe down the strip holder with disinfection wipes. When removed from the instrument, you may also use water with a soft brush and/or liquid disinfectants.

After cleaning, insert the holder on the metal sled (Pic 19). The rectangular notches of the transport mechanism and the strip holder must be placed on top of each other (Pic 20). Use reasonable force to push the strip holder completely onto the metal sled. Doing so the metal slit and the strip holder will be pushed back completely into the housing. If you get an error message after turning the instrument on: Turn the instrument off and push the strip holder again on the metal slit with more force.

Pic 18: test sled
(bottom view)

Pic 19: Notch A

Pic 20: Notch B


CAUTION

Ensure that the strip holder is completely clean and dry before inserting.

11. Service menu

The URYXXON® Relax VET has a password protected service menu. To enter the service menu press the touch-screen three times during the self test after turning the equipment on (Display 24). Upon request input the PIN "1111".

A selective list with different menu items appears.

Display 24: Service menu

11.1 How to reset the system (Load Default)

Select "Load default" from the service menu. The instrument will be reset to delivery status.

NOTICE

All settings modified by the user will be cleared!

11.2 How to control the LOT-control (LOT activate)


Select "LOT activate" from the service menu to activate the LOT-control. This will minimize the risk of using expired test strips. You will be asked for the LOT of strips currently used. The system will give a warning when the test strips are expired. After 100 measurements the system will ask for the LOT-number of the next tube.

11.2.1 Entry of test strips' LOT


If LOT-control is activated, the actual LOT-number of the test strips can be entered using "Strip" in the menu "Settings" (Display 25).

When selecting "LOT number" a numerical pad will appear, which can be used for entering the LOT printed on the strip packaging. If the entry does not have the format expected for a LOT number, an error message will appear.

After entering the LOT, information on the number of strips from that LOT is requested (Display 26).

Display 25: LOT code


Display 26: LOT size

NOTICE


The instrument counts the number of strips. If the number of measurements reaches the previously entered number of strips of that LOT the instrument requests to enter a new LOT number. Then, please enter the LOT number of the strips you want to use.

11.2.2 Warning on expired test strips

The date of expiry of the test strips is calculated from the LOT-number. If the expiry date has passed, a warning will appear (Display 27). If you choose to continue without entering a new LOT, Display 28 appears. Please choose the number of measurements you would like to perform without additional warnings.

Display 27: Expire date A


Display 28: Expire date B

NOTICE

MACHEREY-NAGEL is not responsible for wrong measurement results due to expired test strips.

11.3 How to change the sensitivity

Select "Sensitivity" to adjust the sensitivity settings. The sensitivity of the URYXXON® Relax VET may be adjusted within specified borders for all parameters except the pH value (Display 29).

Display 29: Sensitivity A

NOTICE

Faulty measurement results due to manipulated sensitivity are sole responsibility of the operator of the equipment.

11.3.1 Settings Reset


By selecting the option panel "Default" and confirming with  all sensitivity values are reset to delivery status of the URYXXON® Relax VET.

11.3.2 Appointing new Sensitivity Settings

NOTICE


Technical and medical knowledge is required to appoint new sensitivity settings. If you do not feel confident, leave this menu without changes.

By selecting the option panel "New" and confirming with  the parameter selection will appear (Display 30).

Display 30: Sensitivity B

Use the arrows to select the parameter that should be modified. Confirm your choice by pressing . Display 31 will appear.

Display 31: Sensitivity C

11.3.3 Reported value – (original threshold) – modification

Press  to select the value that you would like to modify. Three informations are in each line:

The information is always defined "X(Y)Z", e.g. "NEG (650) 0"

- X: Measurement value "NEG"
- Y: Specific border "650"
- Z: Change "0"

Press the buttons  and  to increase or decrease the threshold value. Press  to return to the original threshold value.

Confirm your modifications by pressing .

11.4 How to update the instrument


Select "Program update" from the service menu. The URYXXON® Relax VET now expects the upload of a new firmware. Follow the instructions that come with the update-file to finalize the update process.

11.5 Control strip

The check mode for control strips allows you to check the correct functioning of the instrument using special color control strips. You can order these control strips from MACHEREY-NAGEL (REF 930988).

Choose the round button Control strip in the check mode measurement menu (Display 24). The instrument asks for insertion of the colored control strip (Display 32). The measurement of the control strip starts automatically as soon as the strip is placed on the sled, if the autostart is activated. Once the measurement has been completed, the measured values are immediately displayed and/or printed. The result of the measurement is not stored. The values are so called remission values for the pre-colored test pads on the control strip. These values have to be compared to a set of standard values from the control strip package instructions. For additional information on the control strips or in case the values do not fit the should values in the control sheet, please consult the package insert, or contact MACHEREY-NAGEL directly.

ATTENTION: Do not touch the test fields of the control strip with your fingers.

Display 32: Insert control strip

12. Interface description

The URYXXON® Relax VET may be connected to a computer via the RS232- or the USB-interface (work station or laboratory information system).

NOTICE

A comprehensive description of the interface data protocols can be obtained from your local product representative or MACHEREY-NAGEL directly.

12.1 Serial interface

Protocol RS232, 19200 Baud, 8 bit, no parity (Pic 10-④)

Connection plug arrangement:

PIN	Signal	Description	Direction
1	Nc	Not wired	
2	RxD	Data reception	Input
3	TxD	Send	Output
4	Nc	Not wired	
5	GND	Signal ground	--
6	Nc	Not wired	
7	Nc	Not wired	
8	Nc	Not wired	
9	Nc	Not wired	

12.2 USB 1.1-interface

USB-jack Type B (Pic 10-⑤). The instrument will be identified as a serial interface. The driver for the interface module may be downloaded from the MACHEREY-NAGEL homepage (www.mn-net.com).

12.3 Transmission protocol

The data is released via the interfaces as plain text. The received dataset corresponds to the format of the printout.

12.4 Barcode scanner, PC-keyboard

 **WARNING**

Be sure to read the safety instruction of any device before using additional devices (barcode reader/keyboard).

NOTICE

Make sure that the barcode scanner is set to the same language as the device to avoid character misinterpretation.

A USB-jack Type A (Pic 10-⑥) is provided for connection of a keyboard or bar code scanner.

12.5 Ethernet-Interface

A RJ45 jack (Pic 10-⑨) is provided for a wired network connection. The interface is not supported.

13. Error Messages and Fault Clearance

The instrument differentiates between different types of error messages. Advisory error messages are of less importance and are displayed on the screen during handling the instrument (e.g. Please remove strip). When the corrective action has been performed, the analyzer removes the message from the screen.

If an error message occurs during testing it will be displayed on the screen or in the result view instead of the result. Messages are displayed in plaintext and are self-explanatory.

Error Message / Error	Cause	Solution
"Dry Strip"	The test strip wasn't dipped completely	Repeat measurement with a new strip
"Wrong Strip"	A wrong test strip has been detected (wrong type)	Use correct test strips
"Wrong Position"	The strip hasn't been pushed into the strip retainer far enough	New measurement, place strip in right position
"No Paper"	Paper roll empty or printer flap open	Replace paper and close printer flap
"Battery Low"	Batteries are low	Exchange batteries or use power pack
"Instrument doesn't start"	Power supply not installed or defect	Check whether all connections are plugged in and whether the power socket is functioning

In case a fault cannot be cleared by the aid of the instructions above, please contact your local distributor or the MACHERY-NAGEL Service.

14. Warranty

The warranty for this equipment has a duration of 24 months from the date of purchase. The original copy of the bill serves as a certificate and must be submitted in case of assertion of a warranty claim. The warranty expires in case of improper handling and/or maintenance of the equipment; it does not comprise defects due to the external power supply.

The warranty is limited to the repair of faulty parts or – at our sole discretion – to the delivery of a faultless substitute. The warranty period of 24 months is not affected by claiming on the warranty during this period. There is no right of withdrawal.

Further claims are excluded. Hereunto we count in particular all claims for damages evolving from consequential damages or indirect damages.

Additionally the relevant version of our general sales terms and delivery conditions apply as printed on all price lists.

14.1 General usage comments

NOTICE

Due to the reasons stated above, MACHEREY-NAGEL declines any responsibility for the correct functioning of the instrument if it has been opened or used improperly and objects to any warranty claims made in such a case. MN's warranty and liability obligations as to product(s) or to any party on account thereof shall be void and of no effect if the product(s) or any part thereof have been subject to unauthorized change, abuse, misuse, accident or alteration, unauthorized application or installation not in compliance with the product(s) specifications as set forth in the product leaflets/instructions hereof ("unauthorized action") and such unauthorized action resulted in the non-conformance or defect of the product(s).

15. Technical information

15.1 Technical data

General info	Size	7.5 cm × 16 cm × 20 cm (H × W × D)
	Weight	710 g (without batteries and power pack)
	Power supply	External power pack 100–240 V, Output 9 V, 1.5 A Alternative: battery operation with 6 mignon batteries 1.5 V (AA)
Set-up location	Temperature	10 °C – 40 °C
	Humidity	20 % – 80 % (non-condensing)
	Altitude	Up to 3000 m
	Pollution degree	2
Protection Class	Equipment	Protection class III
	Power supply	Protection class I
	similar to UL 50E*	Type 1*
Overvoltage category	Equipment	Category II
	Power supply	Category I
Transport condition	-10 °C – 45 °C, max. 80 % humidity (non-condensing)	
Storage condition	10 °C – 40 °C, max. 80 % humidity (non-condensing)	
Application area	Commercial indoor laboratory	

* This classification for the enclosure corresponds to typical Type 1 enclosures in terms of design and protective effect, but it has not been tested in accordance with UL 50 / UL 50E and is not UL-certified.

The following test-strips from MACHEREY-NAGEL are programmed for evaluation:

Medi-Test Combi 10 VET REF 930870

The maximum sound level and maximum sound pressure are subject to the limits specified in ISO 3746 or ISO 9614-1.

15.2 CE Certification

The CE mark declares that the product complies with the harmonization legislation of the European Community listed below:

European Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS 2)

This device complies with the following regulations:

IEC 61010-1

IEC 61010-2-101

IEC 61326-1

IEC 61326-2-6

FCC47 CFR Part 15B, ICES-003 Issue 6

ISED ICES-001 Issue 5 - CAN/CSA CISPR 11:19

15.3 Analytical Performance Characteristics

15.3.1 Table of results

Param.	Conventional	International	ARB
BLD	NEG	NEG	NEG
	10 Ery/ μ L	10 Ery/ μ L	+
	50 Ery/ μ L	50 Ery/ μ L	++
	250 Ery/ μ L	250 Ery/ μ L	+++
UBG	NORM	NORM	NORM
	2 mg/dL	35 μ mol/L	+
	4 mg/dL	70 μ mol/L	++
	8 mg/dL	140 μ mol/L	+++
	12 mg/dL	200 μ mol/L	++++
BIL	NEG	NEG	NEG
	1 mg/dL	17 μ mol/L	+
	2 mg/dL	35 μ mol/L	++
	4 mg/dL	70 μ mol/L	+++
PRO	NEG	NEG	NEG
	trace (optional)	trace (optional)	trace (optional)
	30 mg/dL	0.3 g/L	+
	100 mg/dL	1 g/L	++
	500 mg/dL	5 g/L	+++
NIT	NEG	NEG	NEG
	POS	POS	+
KET	NEG	NEG	NEG
	25 mg/dL	2.5 mmol/L	+
	100 mg/dL	10 mmol/L	++
	300 mg/dL	30 mmol/L	+++
ASC	NEG	NEG	NEG
	10 mg/dL	10 mg/dL	+
	20 mg/dL	20 mg/dL	++
GLU	NEG	NEG	NEG
	NORM	NORM	NORM
	50 mg/dL	2.8 mmol/L	+
	150 mg/dL	8.3 mmol/L	++
	≥ 500 mg/dL	≥ 27.8 mmol/L	+++
pH	5	5	5
	6	6	6
	6.5	6.5	6.5
	7	7	7
	8	8	8
	9	9	9
SG	1.000	1.000	1.000
	1.005	1.005	1.005
	1.010	1.010	1.010
	1.015	1.015	1.015
	1.020	1.020	1.020
	1.025	1.025	1.025
	1.030	1.030	1.030

Param.	Conventional	International	ARB
LEU	NEG	NEG	NEG
	25 Leu/μL	25 Leu/μL	+
	75 Leu/μL	75 Leu/μL	++
	500 Leu/μL	500 Leu/μL	+++

Table 1: Results for Medi-Test parameters

Meaning of the used abbreviations:

BLD - Blood, UBG - Urobilinogen, BIL - Bilirubin, PRO - Protein, NIT - Nitrite, KET - Ketone, GLU - Glucose, SG - specific gravity of the urine (density), LEU - Leukocytes, ASC - Ascorbic acid.

15.3.2 Precision

The precision was evaluated by using the commercially available control solution Medi-Test Control, manufactured by MACHEREY-NAGEL, Düren, Germany. This Control is used for daily routine system check of the URYXXON® urinalysis system. It is available as a 2 level set (Medi-Test Control N = level 1 (N): negative-normal, Medi-Test Control P = level 2 (P): positive values). Repeatability and reproducibility amount to 100 % in all cases. The results are summarized in the following table.

Precision / Repeatability within series

	BIL	BLD	GLU	KET	LEU	NIT	PRO	UBG	pH	SG	ASC
Total measurements	40	40	40	40	40	40	40	40	40	40	40
True positive	20	20	20	20	20	20	20	20	20*	20*	20
False positive	0	0	0	0	0	0	0	0	n/a	n/a	0
False negative	0	0	0	0	0	0	0	0	n/a	n/a	0
True negative	20	20	20	20	20	20	20	20	20*	20*	20

Precision / Repeatability batch-to-batch

	BIL	BLD	GLU	KET	LEU	NIT	PRO	UBG	pH	SG	ASC
Total measurements	60	60	60	60	60	60	60	60	40	40	60
True positive	30	30	30	30	30	30	30	30	20*	20*	30
False positive	0	0	0	0	0	0	0	0	n/a	n/a	0
False negative	0	0	0	0	0	0	0	0	n/a	n/a	0
True negative	30	30	30	30	30	30	30	30	20*	20*	30

Precision / Reproducibility day-to-day

	BIL	BLD	GLU	KET	LEU	NIT	PRO	UBG	pH	SG	ASC
Total measurements	100	100	100	100	100	100	100	100	100	100	200
True positive	50	50	50	50	50	50	50	50	50*	50*	100
False positive	0	0	0	0	0	0	0	0	n/a	n/a	0
False negative	0	0	0	0	0	0	0	0	n/a	n/a	0
True negative	50	50	50	50	50	50	50	50	50*	50*	100

* pH and density do not provide negative / positive interpretation. Obtained values were in accordance to target values of Control N and P.

15.3.3 Linearity

Blood: neg. · 10 · 50 · 250 erythrocytes/ μ L

Urobilinogen: norm. · 2 · 4 · 8 · 12 mg/dL

Bilirubin: neg. · 1 · 2 · 4 mg/dL

Protein: neg. · Trace (optional) · 30 · 100 · 500 mg/dL

Nitrite: neg. · pos.

Ketones: neg. · 25 · 100 · 300 mg/dL

Ascorbic acid: neg. · 10 · 20 mg/dL

Glucose: neg. · norm. · 50 · 150 · $\geq$ 500 mg/dL

pH: 5 · 6 · 6.5 · 7 · 8 · 9

Density: 1.000 · 1.005 · 1.010 · 1.015 · 1.020 · 1.025 · 1.030 g/mL

Leukocytes: neg. · 25 · 75 · 500 leukocytes/ μ L

15.3.4 Cut-off concentration

The Cut-Off-Concentration (Medical decision point) has been defined as the concentration, where minimum 50 % of the results are positive.

Parameter	Cut-off concentration
Blood	approx. 4 Ery/ μ L
Urobilinogen	approx. 1.5 mg/dL
Bilirubin	approx. 1.0 mg/dL
Protein	approx. 15 mg/dL
Nitrite	approx. 0.1 mg/dL
Ketones	approx. 15 mg/dL
Ascorbic acid	approx. 5 mg/dL
Glucose	approx. 35 mg/dL
Leucocytes	approx. 10 Leu/ μ L

15.3.5 Comparative study

In a comparative study 240 urine samples were evaluated with urine test strips on the URYXXON® Relax VET instrument and similar competitor instruments.

The results of the study are summarized in the following table:

Parameter	Comparison URYXXON® Relax VET – Competitor 1***		Comparison URYXXON® Relax VET - Competitor 2***	
	exact agreement	agreement ± 1 bin	exact agreement	agreement ± 1 bin
Blood	80 %	99 %	87 %	100 %
Urobilinogen	88 %	100 %	91 %	100 %
Bilirubin	95 %	100 %	95 %	100 %
Protein	95 %	100 %	88 %	100 %
Nitrite	96 %	N/A	96 %	N/A
Ketones	98 %	98 %	91 %	98 %
Glucose	92 %	100 %	92 %	100 %
pH	84 %	100 %	88 %	98 %
Density	96 %*	100 %**	65 %*	90 %**
Leukocytes	88 %	99 %	85 %	100 %
overall average	91 %	100 %	88 %	98 %

* this value represents agreement ± 1 bin

** this value represents agreement ± 2 bins

*** similar competitor devices with similar intended purpose (reflectometric device for automatic read-out of urine test strips)

The exact agreement for parameters blood, urobilinogen, bilirubin, protein, nitrite, ketones, glucose, pH and leukocytes amounts to ≥ 80 % in every case. For parameter density agreement ± 1 bin amounts to 96 % and 65 % compared to two different competitor products. Agreement ± 2 bins amounts to 100 % and 90 %, respectively.

Overall average exact agreement is 91 % and 88 % when compared to competitor products. Overall agreement ± 1 bin amounts to 100 % and 98 %.

A summary of the comparative measurements for ascorbic acid with comparative product QUANTOFIX® Ascorbic acid on QUANTOFIX® Relax is shown in the following table. Agreement is 100 % in all cases.

Parameter	Comparison URYXXON® Relax VET – QUANTOFIX® Relax	
	exact agreement	agreement ± 1 bin
Ascorbic acid	100 %	100 %

Diagnostic specificity, diagnostic sensitivity, negative predictive value and positive predictive value for parameters blood, urobilinogen, bilirubin, protein, nitrite, ketones, ascorbic acid, glucose and leukocytes are summarized in the following table. Overall, diagnostic specificity is 97 % and diagnostic sensitivity is 98 %. The average negative predictive value is 98 %; the average positive predictive value is 91 %. Values for the parameters pH and density cannot be calculated because these parameters do not allow negative/positive interpretation.

Parameter	Diagnostic specificity	Diagnostic sensitivity	NPV	PPV
Blood	100 %	89 %	84 %	100 %
Urobilinogen	95 %	100 %	100 %	71 %
Bilirubin	95 %	100 %	100 %	86 %
Protein	97 %	97 %	98 %	96 %
Nitrite	96 %	98 %	99 %	86 %
Ketones	100 %	98 %	100 %	100 %
Ascorbic acid	100 %	100 %	100 %	100 %
Glucose	100 %	97 %	99 %	100 %
pH	N/A	N/A	N/A	N/A
Density	N/A	N/A	N/A	N/A
Leukocytes	88 %	100 %	100 %	76 %
overall average	97 %	98 %	98 %	91 %

Likelihood ratios for parameters URO, BIL, PRO, NIT, KET, ASC and GLU were calculated as LR+ > 10 and LR- < 0.1, which demonstrates convincing diagnostic evidence. For the parameter LEU, likelihood ratio was calculated as LR+ 5–10, which demonstrates high diagnostic evidence. LR- was < 0.1 for this parameter. For the parameter BLD, likelihood ratios were calculates as LR- 0.1–0.2, which demonstrates high diagnostic evidence. LR+ was > 10 for this parameter. Calculation of LR+ and LR- is not possible for pH and SG as diagnostic sensitivity and diagnostic specificity are not applicable for these parameters.

15.3.6 Interference study

The influence of possible interfering substances is described for each test parameter:

Blood: Normal concentrations of ascorbic acid (≤ 40 mg/dL) do not influence the test result. In the case of Combi 11, ascorbic acid concentrations > 2.5 mg/dL lead to false-negative results. False positive reactions can be caused by residues of peroxide-containing or other cleaning agents, as well as menstrual blood.

Urobilinogen: The detection is inhibited by higher concentrations of formaldehyde (> 60 mg/dL). Nitrite concentrations > 2.5 mg/dL and prolonged exposure of the urine to light can lead to low or false-negative values. Overly high or false-positive results can be caused by dyes (e.g. betanin) or medications excreted in the urine.

Bilirubin: The detection is inhibited by higher concentrations of ascorbic acid (> 10 mg/dL) and nitrite (> 2.5 mg/dL). Prolonged exposure of the urine to light can lead to low or false-negative values. Urine indican at a concentration of > 40 mg/dL may cause false-positive results.

Protein: False-positive findings can occur in the case of extremely alkaline urine ($\text{pH} > 9$), disinfectant residues (e.g. benzalkonium chloride > 12.5 mg/dL) in the urine container or in the presence of quinine (> 50 mg/dL).

Nitrite: False-negative results can occur in the case of antibiotic therapy, and in the case of an overly low nitrate level in the urine as a result of low-nitrate food or severe dilution (diuresis). Microbes without the ability to form nitrite can also be present. A false-positive reaction color can be caused by dyes (e.g. betanin) excreted in the urine.

Ketone: Phthalein compounds up to a concentration of 125 mg/dL (highest tested concentration) do not interfere with the result.

Ascorbic acid: Lower results can occur due to oxidising cleaning agents in sample containers.

Glucose: Normal concentrations of ascorbic acid (≤ 40 mg/dL) do not influence the test result. False-positive reactions can be caused by oxidising cleaning agents in the sample container.

pH: The dye betanin has no influence on the evaluation of the pH test field up to the highest tested concentration of 500 mg/dL.

Density: In the case of elevated protein excretion (> 500 mg/dL), the density values determined are too low.

Leukocytes: An attenuated reaction can be expected if preparations with nitrofurantoin (> 2 mg/dL) or phenazopyridine (> 0.2 mg/dL) are taken. Formaldehyde (as a preservative, > 20 mg/dL) and dyes (e.g. betanin) can lead to a false-positive reaction. In the case of specimens from female patients, a false-positive reaction can be simulated by vaginal secretion.

15.4 Literature

M. J. C. Schot, S. van Delft, A. M. J. Kooijman-Buiting, N. J. de Wit, R. M. Hopstaken, *BMJ Open*, 2015, 5, e006857.

K. Ekblom, A. Petersson, *Scand. J. Clin. Lab. Invest.*, 2020, 1, 39–45.

K. Durkalec-Michalski, P. M. Nowaczyk, K. Siedzik, J., *J. Int. Soc. Sport Nutr.*, 2019, 16,16.

15.5 Electromagnetic Compatibility

- Target environment of device use: Application in the commercial indoor laboratory environment.
- This device has been developed and tested in accordance with CISPR 11 Class A. It can cause radio interference in a domestic environment. In this case, you may need to take action to correct the interference.
- Use of this equipment in a dry environment, especially if synthetic materials (synthetic fiber clothing, carpets, etc.) are present, can cause destructive static discharges that can cause erroneous results.
- Do not use this device near sources of strong electromagnetic radiation, as they may interfere with proper operation.
- This *veterinary* diagnostic device is in accordance with the applicable parts of the requirements for emission and immunity described in the IEC 61326 series of standards.
- Before use, the electromagnetic environment must be assessed.

EMC Class A according CISPR 11:

The device can be used in all institutions excluding residential areas and areas which are directly connected to the public power supply, if the public power supply supplies buildings for residential purpose.

15.6 Waste disposal

Disposal according to EU Directive 2012/19/EU. In accordance with the EU Directive 2012/19/EU MACHEREY-NAGEL takes back the old device and disposes it free of charge.

DANGER

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

15.7 CAN ICES

This equipment is classified as industrial, scientific, and medical (ISM) equipment and complies with the requirements of the Canadian Interference-Causing Equipment Standard ICES-001, issue 5, Class A.

16. Additional Information

16.1 Manufacturer Information

MACHEREY-NAGEL GmbH & Co. KG

Valencienner Str. 11

52355 Düren · Germany

Phone: +49 2421 969-0

E-mail: info@mn-net.com

Internet: www.mn-net.com

16.2 Technical support

If you have any questions about the function or use of the device, or if you suspect a defect, contact the technical support of MACHEREY-NAGEL:

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