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Rev. 10/2025

## Summary

The test is suitable for determining the acute toxicity to liquid-dried luminous bacteria in accordance with DIN EN ISO 11348 2:2023-12.

The test is suitable for surface water, groundwater, seepage water, pore water, sea water, brackish water, wastewater, and aqueous extracts and eluates from sediments.

- Measuring range: 0–100 % inhibition
- Reporting of results:
  - % inhibition of light intensity in the sample compared to an uninhibited control  $G_0$  (according to DIN EN ISO 11348-2: LID); Reciprocal of the first dilution level of a sample at which the inhibition of light intensity is less than 20 %.
  - $EC_{20}$  /  $EC_{50}$ : Concentration of a sample that causes an inhibition of light intensity of exactly 20 %/50 %
- Number of tests : 200 (REF 945023) / 400 (REF 945024)
- Shelf life: 24 months
- Storage temperature: -21 – -15 °C

## Method

Determination of acute toxicity to liquid-dried luminous bacteria. The measured variable is the natural light output (luminescence) of the microorganism *Aliivibrio fischeri* NRRL B-11177 used. The inhibition of light intensity by the sample is determined in comparison to an uninhibited control sample

## Content

### REF 945023

- 10 vials BioFix® Lumi Luminous bacteria
- 10 vials BioFix® Lumi Reactivation solution
- 50 mL BioFix® Lumi Standard solution

### REF 945024

- 20 vials BioFix® Lumi Luminous bacteria
- 20 vials BioFix® Lumi Reactivation solution
- 50 mL BioFix® Lumi Standard solution

## Accessories

- Luminometer (e.g., BioFix® Lumi-10, REF 940008), with suitable cuvettes (e.g., REF 916912) and cuvette holder (e.g., REF 945013)
- Heating block/ water bath for temperature control at 15 °C
- Piston pipette 1 – 5 mL (REF 916909) with pipette tips (REF 916916)
- Piston pipette 100 – 1000 µL (REF 91677) with pipette tips (REF 91676)
- BioFix® Lumi osmotic pressure adjustment solution (REF 945602)
- BioFix® Lumi dilution solution (REF 945601)

## Quality control

The enclosed test certificate guarantees the validity criteria required by DIN EN ISO 11348 2. The sensitivity of the luminous bacteria must also be checked in your own laboratory in accordance with the standard. The necessary information, such as standard substances and test concentrations, can be found in the test certificate.

## Use of the BioFix® Lumi Standard Solution

The enclosed NaCl-based BioFix® Lumi Standard Solution is used for non-standard testing. 0.5 mL of resuspended luminous bacteria (see Procedure: Dosage – Step 1) achieve a 40–60 % inhibition of luminescence intensity with 0.5 mL of BioFix® Lumi Standard after 30 minutes at 15 °C.

## Sampling and preparation

See DIN EN ISO 5667-3-A21.

Sediment cloudy samples, centrifuge them, or filter them with CHROMAFIL® PES (REF 729242).

Check the oxygen content ( $\geq 3$  mg/L).

Adjust to pH 6.0–8.5 with NaOH or HCl solution.

Check salt content and adjust to 2 % NaCl. To do this, use one part BioFix® Lumi solution for adjusting the osmotic pressure (REF 945602) to ten parts sample.

If necessary, prepare dilution steps (for details, see DIN EN ISO 11348-2, Annex B).

## Procedure

### Reactivation

- Bring 1 vial reactivation solution to 15 °C
- Add 1 mL of reactivation solution to 1 vial of frozen luminous bacteria
- Mix by shaking
- Defrost for 2 min at room temperature
- Incubate for 15 min at 15 °C

### Dosage

- Mix with the remaining 11 mL of reactivation solution (15 °C)
- Add 0.5 mL of luminous bacteria solution to each cuvette
- Incubate for 15 min at 15 °C

### Alternative dosage with G1 stage

- Mix with 4 mL reactivation solution (15 °C)
- Add 0.2 mL luminous bacteria solution to each cuvette for G1 determination
- Mix remaining luminous bacteria solution with the remaining 7 mL reactivation solution (15 °C)
- Add 0.5 mL luminous bacteria solution to each other cuvette
- Incubate for 15 min at 15 °C

### Measurement (for details, see the manual for the luminometer used)

- Measure initial luminescence  $I_0$
- Add 0.5 mL sample/control
- Mix by swirling
- Incubate for  $t$  min at 15 °C
- Measure final luminescence  $I_t$

## Interferences

According to DIN EN ISO 11348-2, no pH correction is necessary if the pH value of the sample is between 6.0 and 8.5. pH values  $< 6.0$  or  $> 8.5$  can lead to pH-related luminescence inhibition.

Undissolved or volatile substances can falsify the test result. Turbidity in the sample can cause light loss due to light scattering.

Samples with high oxygen consumption (e.g., samples containing chlorine) or low oxygen content ( $< 3$  mg/L) may show toxicity-independent inhibition of light emission.

A salt concentration of  $> 35$  g/L NaCl or other substances of equal osmolarity should not be exceeded, as hyperosmotic effects may occur (see DIN EN ISO 11348-2, Annex D).

## Storage

The test can be stored at -21 – -15 °C until the specified expiration date. Only thaw the luminous bacteria immediately before reactivation.

Use reactivated luminous bacteria within 4 hours.

Reactivated luminous bacteria should only be stored diluted (see Procedure: Dosage - Step 1) at 2–8 °C, whereby a decrease in luminescence intensity and sensitivity may occur. Refreezing reactivated luminous bacteria is not recommended.

## Disposal

No hazardous substances requiring labeling are contained.

Luminous bacteria and test preparations can be disposed of down the drain. Care must be taken to ensure that toxic samples are disposed of properly.

Information on hazards can be found on the outer label and in the safety data sheet. The safety data sheet can be downloaded at [www.mn-net.com/SDS](http://www.mn-net.com/SDS).

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MACHERY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany