NUCLEOSIL® CHIRAL-2 / CHIRAL-3

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NUCLEOSIL® CHIRAL-2 / CHIRAL-3

**Note:** All HPLC columns from MACHEREY-NAGEL are supplied with a certificate, which contains specifications and test results of the column. NUCLEOSIL® CHIRAL-2 and CHIRAL-3 columns are quality products based on the robust silica NUCLEOSIL® with covalently bonded D or L configuration of N-(3,5-Dinitrobenzoyl)-phenylglycine as chiral selector for enantiomeric separations. They are specifically developed for HPLC analysis. If carefully and properly used excellent chromatographic results and long column lifetime can be achieved. HPLC columns are designed for qualitative and quantitative analysis of mixtures of substances and single components. They must exclusively be used in accordance with universally accepted laboratory regulations and HPLC working methods. Before running the column the entire analytical system (column and equipment) has to be carefully checked by the operator. Chromatographic conditions (mobile phase, flow, temperature etc.) must be adapted to the analytical task. MACHEREY-NAGEL does not give any warranty and is not liable for the success of a separation or application. If you have any questions after reading this manual, please call our service / technical support.

Table of contents

- Safety indication
- Description of the column
- Installation
- Precolumn filter and guard columns
- Sample
- Eluent
- Flow rate and pressure
- Temperature
- Detection
- Equilibration
- Column storage
- Application note
- Troubleshooting
- Column regeneration
- Abstract

Safety indication

Follow the general safety instructions for handling of HPLC solvents used as mobile phases (e.g., acetonitrile, methanol) and take precautions against any kind of injuries or damage to health (e.g., skin and eye protection in case of broken capillaries). Disposal of used HPLC columns must follow international, national and local environmental protection regulations. The use of HPLC columns is only permitted to staff members, who are qualified in their field. Keep HPLC columns away from children. MACHEREY-NAGEL disclaims and excludes all warranties of any kind or nature whatsoever and MN shall not be liable for any damages (whether direct, indirect, foreseeable, incidental, compensatory, consequential or special), whether based upon warranty, contract, tort or strict liability, if damages and /or losses occur caused by improper use, maintenance, neglect or improper treatment (especially opening of the column and exposure of the column bed).

Description of the column

The columns NUCLEOSIL® CHIRAL-2 and CHIRAL-3 have been developed for control of optical purity of compounds under normal phase conditions. They are recommended for the analysis of stereoisomers (enantiomers and diastereomers), control of optical purity (e.g. plant protectives, pharmaceuticals) and for product control in chiral organic synthesis. As chiral selector *N*-(3,5-dinitrobenzoyl)-phenylglycine is covalently bonded to NUCLEOSIL® silica. Both adsorbents differ only in the configuration of the chiral selector. CHIRAL-2 and CHIRAL-3 thus feature, for any pair of enantiomers, equal selectivity and chromatographic resolution, but a reversed elution sequence of the optical antipodes. For control of optical purity the user now has the option to select conditions such that the minor component, which is present as an impurity, is eluted before the main peak. Thus, overlapping peaks are avoided.

Installation

The column should be installed in the flow direction indicated on the column label. It is connected with 1/16" capillaries and fittings, typical for HPLC instruments.

Precolumn filter and guard columns

A precolumn filter containing 0.5–2.0 µm porosity stainless steel frits is recommended between sample injector and column to remove particulates from the eluent stream. For protection and an extension of column lifetime the columns should always be used with guard columns. The filter elements and the adsorbent in the guard column retain contaminants from the sample or the eluent. Connection of the guard column with the separation column is made using a suitable guard column holder (see [www.mn-net.com](http://www.mn-net.com) or MN chromatography catalog). Replacement of the guard column is required when increased column pressure and /or loss of performance is observed.

Sample

Samples, generally dissolved in the starting eluent, should be passed through a syringe filter (e.g., CHROMAFIL® Xtra PET, 0.45 µm, 25 mm, REF 729220) before entering the column. If injected sample solutions are still turbid even after filtration, the lifetime of the column may be significantly reduced. The sample volume should be as small as possible and not exceed 50 µL to achieve an optimal resolution.

Eluent

Columns are used in Normal Phase mode. Therefore they are supplied with an eluent of n-heptane – 2-propanol – trifluoroacetic acid (100:0.05:0.05, v/v/v). For separation of enantiomers nonpolar organic eluents (n-heptane, *i*-octane) are recommended. Retention time can be influenced by adding polar organic compounds (e.g., tetrahydrofuran, alcohols, chlorinated hydrocarbons). Small amounts of a strong acid like trifluoroacetic acid in the eluent (up to 0.1 vol %) can considerably improve separation of enantiomers. Pay attention to the pH stability of 2–8 and avoid any traces of water. Eluents should be filtered through a 0.2–0.45 µm membrane and degassed.

Flow rate and pressure

Flow rate (recommended: 0.5–1.0 mL/min) influences the time needed, the resolution and the column lifetime. It is limited by the back pressure, which should not exceed 400 bar. If a high pressure results from the use of the column at regular flow rates, this usually indicates contamination of the packing material, which must be removed (see troubleshooting).

Temperature

For operation of the columns temperatures up to 60 °C are possible. However, they should be at least 30 °C below the boiling temperature of the eluent, in order to ensure proper detection. Variation of the temperature influences retention times and especially the peak shape. Optimum temperatures for successful separations should be determined empirically.

Detection

Spectrophotometers, refractometers and electrochemical detectors can be used with the NUCLEOSIL® CHIRAL-2 and CHIRAL-3 columns. If electrochemical detectors are used, please note that high temperatures may be incompatible with some working electrodes. If a higher sensitivity is required, post-column derivatizations with an appropriate detector for the reaction product can be used.


Equilibration

Prior to measurement of samples the column must be rinsed with the eluent at the same flow rate and temperature as the method to be applied. Column equilibration is finished, when the baseline of the detector no longer shows a drift.

Column storage

The original eluent (see eluent) is recommended for storage. For column storage be sure the end fittings are tightly sealed using column end plugs, because storage without these seals can result in drying of the packing material.

Application note

Troubleshooting

The following outline describes the symptoms of performance loss and their cause. All columns are subject to the strict regulation and control of our quality assurance system. Columns based on silica are robust and hold their separation efficiency for long periods by correct maintenance and treatment. According to experience, column failures are mostly a result of injection of contaminants to the sorbent bed. Use of a guard column, as well as an appropriate sample pretreatment will help to minimize these risks. All NUCLEOSIL® CHIRAL-2 and NUCLEOSIL® CHIRAL-3 columns are thoroughly tested prior to shipment and are supplied with a sample chromatogram illustrating performance of that particular column.

Use the outline below to help determine the cause of a possible performance loss:

| Symptom / Error / Cause                                                                                                                                                                                                                                                                                         | Prevention / Repair                                                                                                                                                                               |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Baseline drift</b> <ul style="list-style-type: none"><li>• insufficient period for equilibration with the eluent</li><li>• contaminated eluent</li><li>• temperature</li></ul>                                                                                                                               | longer or better equilibration<br>usage of freshly prepared solvents and buffers<br>column temperature control                                                                                    |
| <b>Broad peaks</b> <ul style="list-style-type: none"><li>• mixing and /or diffusion before /behind the column</li><li>• too large sample volume</li></ul>                                                                                                                                                       | keep length and ID of capillaries at a minimum<br>smaller injection volume                                                                                                                        |
| <b>Peak interference; too fast elution</b><br>too fast elution and /or insufficient separation by: <ul style="list-style-type: none"><li>• improper column temperature or flow rate</li><li>• elution power of eluent is too high</li></ul>                                                                     | optimize concerned parameter<br>optimize eluent system                                                                                                                                            |
| <b>Increasing back pressure; degradation of the separation performance</b><br>contamination of sorbent by: <ul style="list-style-type: none"><li>• particulate accumulation on frit or sorbent bed from sample, eluent or system</li></ul>                                                                      | prefilter of samples and eluent, use an in-line filter /<br>rinse LC system, clean the sorbent<br>(see column regeneration)                                                                       |
| <b>Insufficient separation; degradation of the separation with regular column pressure</b><br>contamination with: <ul style="list-style-type: none"><li>• fats, oils, lipids from sample (coating of sorbent surface) and other organic substances from im-properly prepared eluent or matrices</li></ul>       | remove organic substances by sample preparation /<br>clean the sorbent (see column regeneration)                                                                                                  |
| <b>Double peaks (dead volume)</b> <ul style="list-style-type: none"><li>• faulty fittings (capillaries, ferrules, nuts)</li></ul><br><ul style="list-style-type: none"><li>• compression of column bed by too high flow rate of eluent</li><li>• dissolution of silica by too high pH value of eluent</li></ul> | use "PEEK Fingertight Fittings", REF 718770 or REF 718778 / replace fittings<br>consider maximum flow rate and allowed eluent /<br>replace column<br>consider pH range of column / replace column |

Column regeneration

In some cases the performance of the column can be restored by removing contaminants from the sorbent bed. It is important, however, to locate the source of contamination before using the column for the analysis of samples again.

1. **Prepare fresh eluent:** In some cases the performance loss is traced to eluent contamination. Therefore, prepare fresh eluent and flush all liquid lines before using the column again. The eluent should be filtered through a 0.2–0.45 µm membrane and degassed prior to use.
2. **Cleaning of sorbent:** To remove contamination 50 mL 2-propanol are pumped with a flow of 0.1 mL/min through the inverted column. Then return to the original eluent in the original flow direction.
3. **Regeneration:** By rinsing with n-heptane – 2-propanol – trifluoroacetic acid (100:0.05:0.05, v/v/v) at a flow rate of 0.2 mL/min the phase can be regenerated.
4. **Column replacement:** Above procedures will restore performance only in certain cases. Some organic contaminants are particularly refractory and may not respond to treatment. Under these circumstances, column replacement is necessary. It is highly advisable to locate the cause of the problem before installing a new column.

Abstract

To extend column lifetime, please keep in mind the following:

1. Recommended eluents are nonpolar organic eluents (n-heptane, *i*-octane). Polar solvents (e.g., tetrahydrofuran, alcohols, chlorinated hydrocarbons) and small amounts of acids (e.g. trifluoroacetic acid) can be added. They should be filtered through a 0.2–0.45 µm membrane and degassed.
2. Filter samples through a 0.2–0.45 µm CHROMAFIL® Xtra PET syringe filter before injection.
3. Use an in-line filter and /or a guard column for protection against impurities.
4. The recommended flow rate is 0.5–1.0 mL/min.
5. Adjust flow rate to keep column pressure below 400 bar.
6. Store the column in n-heptane – 2-propanol – trifluoroacetic acid (100:0.05:0.05, v/v/v).
7. Use analytical grade reagents and HPLC grade solvents for all work. Discard any solutions that show evidence of bacterial growth.

Please check the full range of MACHEREY-NAGEL chromatography products: [www.mn-net.com/chromatography](http://www.mn-net.com/chromatography)

... for applicative support please visit our application database  
with more than 3000 chromatography applications: [ChromaAppDB.mn-net.com](http://ChromaAppDB.mn-net.com)