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Methode:

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BOD₅-Nutrient Mixture *PLUS* (with *N*-allylthiourea)

Method:

The special chemical set *NANOCOLOR*[®] BOD₅-Nutrient Mixture *PLUS* is required for the photometric determination of BOD₅ with tube test *NANOCOLOR*[®] BOD₅ - Test 8-22 (REF 985822) and **contains the nitrification inhibiting compound *N*-allylthiourea**. The specified concentrations of *N*-allylthiourea and nutrients are in accordance with the recommendations of the German and European Standard Method **DIN EN 1899 - H51**. Additional accessory for the determination of BOD₅ is the BOD₅ Accessories Set (REF 916918).

MACHEREY-NAGEL additionally offers a nutrient mixture without the nitrification inhibiting compound *N*-allylthiourea, called *BOD₅-Nutrient Mixture* (REF 918994).

Content of reagent set:

10 test tubes, clear glass with reagent R1, lyophilized
10 test tubes, brown glass with reagent R2, lyophilized

Hazard warning:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Procedure:

Always use **freshly prepared** BOD₅-Nutrient Solutions *PLUS*!

Requisite accessories: piston pipette with tips

Preparation of BOD₅-Nutrient Solutions *PLUS* R1 and R2:

1.	Open a clear glass vial with reagent R1 , add 10 mL aerated diluting water (see BOD ₅ -Accessories Set, REF 916918), close and shake to dissolve. Now you have BOD₅-Nutrient Solution <i>PLUS</i> R1 .
2.	Open a brown glass vial with reagent R2 , add 10 mL aerated diluting water (see BOD ₅ -Accessories Set, REF 916918), close and shake to dissolve. Now you have BOD₅-Nutrient Solution <i>PLUS</i> R2 .
3.	In connection with the BOD ₅ determination, add 0.25 mL BOD₅-Nutrient Solution <i>PLUS</i> R1 and 0.25 mL BOD₅-Nutrient Solution <i>PLUS</i> R2 for every 100 mL sample dilution. <i>For further information see package inserts of tube test NANOCOLOR[®] BOD₅ - Test 8-22 (REF 985822).</i>

*At the end of the day discard all BOD₅-Nutrient Solutions *PLUS* that are no longer needed!*

Preparation of inoculated BOD₅-Nutrient Solutions *PLUS* R1 and R2:

If the sample itself does not contain a sufficient amount of appropriate micro-organisms, so-called **inoculating water** is added to the BOD₅ test preparation via the BOD₅-Nutrient Solutions *PLUS*. This inoculating water can be obtained in the following ways:

- precipitated urban waste water from presettling tank outflow or mechanical inflow
- urban waste water from sewage disposal with max. 300 mg/L COD or 100 mg/L TOC
- surface water containing urban waste water

1.	Open a clear glass vial with reagent R1 , add 10 mL inoculating water , close and shake to dissolve. Now you have inoculated BOD₅-Nutrient Solution <i>PLUS</i> R1 .
2.	Open a brown glass vial with reagent R2 , add 10 mL inoculating water , close and shake to dissolve. Now you have inoculated BOD₅-Nutrient Solution <i>PLUS</i> R2 .
3.	In connection with the BOD ₅ determination, add 0.25 mL inoculated BOD₅-Nutrient Solution <i>PLUS</i> R1 and 0.25 mL inoculated BOD₅-Nutrient Solution <i>PLUS</i> R2 for every 100 mL sample dilution. <i>For further information see package inserts of tube test NANOCOLOR[®] BOD₅ - Test 8-22 (REF 985822).</i>

*At the end of the day discard all inoculated BOD₅-Nutrient Solutions *PLUS* that are no longer needed!*

*Note: The oxygen consumption in control solutions with inoculated BOD₅-Nutrient Solutions *PLUS* shall not exceed a value of 1.5 mg/L O₂ after 5 days of incubation at (20 + 1) °C.*

Examples of samples for which inoculated BOD₅-Nutrient Solutions need to be used:

- industrial, biological inactive waste water
- surface water containing no urban waste water
- deep frozen samples!
- control tests using a glucose/glutamine acid standard solution according to DIN EN 1899-1 - H51 (*NANOCONTROL* BOD₅-Standard, REF 92582)!

Storage:

The reagent set is to be stored in dry conditions at room temperature. The expiry date printed on the kit is to be observed.

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

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