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FREESTYLE™

LC Tech FreeStyle - Report on Methods: SPE Date: 26.06.2017 Time: 13:11:02

Name: CAP-M02-Honey_Agilent.spe	SPE Column: LCTech_3ml.col
Extension cannula: yes	Standard (organic solvents)
Processing speed selection: Rinsing intensity: Use pressure limitation function during loading and washing:	Extended rinsing cycle no
Step: Conditioning_MeOH	Basic type: Conditioning Step: - ID: 795
Volume: 3 ml Suction Speed: 20 ml/min Repetitions: 0 Waiting Time after Dosage: 10 sec.	Dispensing Speed: 1 ml/min Port: 1 methanol Dispense: into Waste
Step: Conditioning_Water	Basic type: Conditioning Step: - ID: 796
Volume: 5 ml Suction Speed: 20 ml/min Repetitions: 0 Waiting Time after Dosage: 10 sec.	Dispensing Speed: 1 ml/min Port: 8 water Dispense: into Waste
Step: Load	Basic type: Load - Transfer Sample-Aliquot over sample loop Step: - ID: 797
Volume: 9 ml Vial Type: Type1 @ 16 without rinsing of vial Suction Speed: 10 ml/min Waiting Time after Dosage: 0 sec.	Dispensing Speed: 3 ml/min Waiting Time after Step: 150 sec. Dispense: into Waste
Step: Washing_water	Basic type: Washing Step: - ID: 798
Volume: 5 ml Suction Speed: 10 ml/min Repetitions: 1 Waiting Time after Dosage: 0 sec.	Dispensing Speed: 3 ml/min Port: 8 water Dispense: into Waste
Step: Drying_Nitrogen	Basic type: Drying - Drying by defined air volume Step: - ID: 799
Air volume: 100 ml Suction Speed: 100 ml/min	Dispensing Speed: 100 ml/min Dispense: into Waste
Step: Eluting_EE-MeOH_8:2	Basic type: Eluting Step: - ID: 800
Volume: 5 ml Suction Speed: 10 ml/min Repetitions: 0 Waiting Time after Dosage: 0 sec.	Dispensing Speed: 3 ml/min Port: 11 MeOH:EE 8:2 Dispense: into vials Number of vials: 1 Vial Type: Type1 @ 16
Step: Drying_after_Elution	Basic type: Drying - Drying by defined air volume Step: - ID: 801
Air volume: 20 ml Suction Speed: 100 ml/min	Dispensing Speed: 10 ml/min Dispense: stay on actual position

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