
MACHEREY-NAGEL

NANOCOLOR[®]

USB T-Set

- Handbuch
- Manual

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Handbuch NANOCOLOR® USB T-Set DE, V1.00, 06.2014, Juni 2014

Handbuch NANOCOLOR® USB T-Set DE, V2.00, 01.2023, Januar 2023

1. General information and introduction	29
2. CE labeling.....	29
3. Type plate	29
4. Explanation of symbols used	30
5. Safety	31
6. Technical description	32
6.1 Specification	33
7. Installing and connecting the device	34
7.1 Scope of delivery	34
7.2 Connecting the device	34
7.3 Transport.....	35
8. Device.....	35
9. Commissioning	36
10. Display	36
11. Programs	37
11.1 Test program	37
11.2 Calibration	38
11.3 External temperature measurements	39
12. T-Set software.....	40
12.1 Program start	40
12.2 Menu.....	40
12.3 Record heating curve	42
12.4 Read test data	43
12.5 Read calibration data	44
12.6 Settings.....	45
13. Maintenance and cleaning of the device	46
13.1 Cleaning the device.....	46

14. Malfunctions, causes and correction	47
15. Service	48
15.1 Disposal	48
15.2 Warranty, liability and complaints	49
15.3 Contact.....	49
15.4 Version history	49

1. General information and introduction

The *NANOCOLOR*[®] USB T-Set is an electronic thermosensor that is suitable for temperature control and automatic calibration of all current *NANOCOLOR*[®] heating block models. The heating blocks can be checked independently with the *NANOCOLOR*[®] USB T-Set for internal quality control.

2. CE labeling

The CE marking declares that the product complies with the following harmonization legislation of the European Community:

European Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)

European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

European Directive 2014/30/EU on electromagnetic compatibility (EMC)

3. Type plate

Figure 1: Type plate

4. Explanation of symbols used

The following symbols or terms appear on the device type plate and in this operating manual, with the meanings described below:

Term / Symbol	Meaning
NANOCOLOR® USB T-Set	Device
SN	Serial number of the USB T-Set
	In accordance with 2012/19/EU, disposal of the appliance via public disposal systems is not permitted. Please also refer to the chapter on Disposal.
	CE indicates compliance with the applicable harmonization legislation of the European Community.
MACHEREY-NAGEL	Manufacturer identification

5. Safety

PLEASE READ THE FOLLOWING SAFETY INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE.

Failure to observe these instructions may result in serious injury to the operator, malfunction, or damage to the device.

Keep this operating manual for future reference.

Follow the safety instructions and guidelines in the operating manual and pay attention to the labels and notices on the device.

Do not attempt to work on internal parts of the device. Any violation will void all warranty claims.

Meaning of the safety instructions:

DANGER

Indicates an imminent or potentially hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Indicates an imminent or potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may result in malfunction or damage of the instrument. Information that clarifies facts in the text and that requires special attention.

6. Technical description

For all temperatures programmed in the checked heating block, a target/actual comparison is measured, recorded and stored in the heating block using the *NANOCOLOR*[®] USB T-Set. Once the temperature check or calibration has been carried out, the data obtained can be transferred to the computer via the RS232 or USB interface. The free *NANOCOLOR*[®] T-Set PC software then allows GLP-compliant documentation and the creation of the corresponding test certificates.

The current temperature of the *NANOCOLOR*[®] USB T-Set is shown on the LED display. This means that temperature measurements, e.g. of samples, can also be carried out independently of the heating block.

6.1 Specification

Type	Electronic thermal sensor for temperature control and calibration as well as issuing a test certificate
Detector	PT 1000 (95 x Ø 4 mm)
Display	LED display of the USB-T-Set and the T-Set software
Resolution	± 0.01 °C (0.001 °C PC-Software)
Operation	Via the heating block's touchscreen and the T-Set software
Temperature range	0–200 °C
Accuracy	± 1 °C
Long-term stability	± 0.1 °C
Interface	USB
Operating environment	0–50 °C, up to 90 % relative humidity
Power input	max. 20 mW
Dimensions	73 cm (length)
Weight	~ 60 g
Certificate	Calibrated against standard thermometer
Warranty	2 years
Certification mark	CE, UKCA

7. Installing and connecting the device

7.1 Scope of delivery

Check that the delivery is complete. If parts of the scope of delivery are missing, please contact MACHEREY NAGEL or your local distributor immediately.

- NANOCOLOR® USB T-Set
- Software CD
- Certificate of conformity
- Manual
- Transport box

7.2 Connecting the device

CAUTION

Hazards: Only qualified personnel are permitted to carry out the instructions described in this step.

WARNING

Risk of electric shock: Ensure that the USB cable is not damaged when used for operation. Check that the voltage source is suitable for the device.

WARNING

Risk of injury: Familiarize yourself with the device before working with it and read this document carefully. Only use the device if you have received instruction on its operation.

NOTICE

Malfunction due to defects in the power supply and housing: If the housing is damaged or the power supply is defective, the device must be taken out of operation.

7.3 Transport

The product packaging, including its protective interior, is best suited for transporting the device. Therefore, please keep the packaging. If you no longer have the original packaging, use a suitable box as outer packaging and soft padding material to ensure that the device cannot move during transport.

NOTICE

Select suitable packaging for shipping.

8. Device


Figure 2: NANOCOLOR® USB T-Set

- ① Graphic display
- ② USB connection
- ③ Temperature sensor
- ④ Touch sensor for display rotation

9. Commissioning

NOTICE

The use of the USB-T set with the *NANOCOLOR*[®] heating blocks is supported from firmware version 4.0.0.10 and above.

Insert the *NANOCOLOR*[®] USB T-Set into the USB interface of the switched-off *NANOCOLOR*[®] *VARIO C2*, *VARIO 4*, or *VARIO HC* heating block (for use of the USB T-Set with the *VARIO Mini*, see the device manual).

Fully insert the temperature sensor of the *NANOCOLOR*[®] USB T-Set into the designated small borehole on the heating element of the heating block and switch on the heating block.

10. Display

The USB T-Set features an integrated display that shows the current measured temperature. By gently tapping the side of the display of the *NANOCOLOR*[®] USB T-Set, the display can be rotated by 180°.

11. Programs

NOTICE

Before running a test program or a calibration, the temperature of the heating block should be lower than the lowest temperature programmed in the heating block. Otherwise, the selected program will only start once the lowest programmed temperature has been undershot.

Pressing the “Menu” icon opens the menu of the heating block. The “T-Set” icon leads to the options available for the USB T-Set. This menu allows you to check the calibration of the heating block by pressing the “Test” icon or to recalibrate the heating block using the “CAL” icon.

11.1 Test program

NOTICE

When removing the temperature sensor from the heating block, please hold it by the black insulation. The temperature sensor may be hot!

The test program is started by pressing the “Start” icon. All temperatures programmed in the heating block are then measured by the USB T-Set, from the lowest to the highest, recorded, and stored in the heating block (duration approx. 40 minutes). When the heating block reaches a programmed target temperature, an acoustic signal sounds. The actual temperature measured by the USB T-Set is displayed, and the difference from the heating block temperature is calculated. The temperatures measured by the USB T-Set can be manually recorded during the test or exported and documented with the T-Set software after completion. Once the test is finished, the data can be sent to a connected USB

stick using the “Paper” icon. The “Up Arrow” icon returns you to the previous menu. The “Factory” icon can be used to reset the calibration of the heating block to the factory settings.

11.2 Calibration

NOTICE

When removing the temperature sensor from the heating block, please hold it by the black insulation. The temperature sensor may be hot!

The calibration program is started by pressing the “Start” icon. All temperatures programmed in the heating block are measured by the USB T-Set from the lowest to the highest, recorded, and stored in the heating block (duration approx. 40 minutes). When the heating block reaches a programmed target temperature, an acoustic signal sounds. The actual temperature measured by the USB T-Set is displayed. The temperatures measured by the USB T-Set can be manually recorded during calibration or exported and documented with the T-Set software after completion. Once all the temperatures programmed in the heating block have been calibrated, the “Checkmark” icon appears on the display. Confirming the “Checkmark” icon brings up a prompt to save the new calibration. Confirming with the “Checkmark” saves the new calibration, while pressing “X” discards it. Once calibration is complete, the data can be sent to a connected USB stick using the “Paper” icon. The “Up Arrow” icon returns you to the previous menu. The “Factory” icon can be used to reset the calibration of the heating block to the factory settings.


11.3 External temperature measurements

The *NANOCOLOR*[®] USB T-Set can be used not only for temperature control of heating blocks, but also for external temperature measurements. To do this, simply connect the USB T-Set into any USB interface. The current ambient temperature of the temperature sensor is shown on the display of the USB T-Set. With the help of the PC software, a temperature curve can be recorded (see Chapter 12. T-Set software on page 40).

Quality control

Each facility must implement its own QC policy.

To verify the accuracy of the temperature sensor, we recommend measuring it against a certified thermometer at several values within the measuring range. Alternatively, the USB T-Set can be sent to MACHEREY-NAGEL or a certified laboratory for calibration. If you wish to check the temperature of different positions in the heating block, a suitable adapter for the USB T-Set is required

(see: [qr.mn-net.com/qr/\(241\)919924](https://qr.mn-net.com/qr/(241)919924)).

12. T-Set software

NOTICE

System requirements for using the software: Microsoft Windows 10 oder 11.

Insert the *NANOCOLOR*[®] software DVD into the DVD drive of your computer. Open the “T-Set” folder and start the software by running the file “kalib.exe”. Alternatively, you can download the latest version of the software directly from our website:

<https://www.mn-net.com/media/archive/06/7a/8f/Kalib.zip>

12.1 Program start

Connect the heating block to the computer using the serial cable (REF 919680) or the included USB cable (REF 919687). Switch on the heating block. Start the *NANOCOLOR*[®] T-Set software as previously described. The connected heating block will be automatically recognized by the software. The start window displays a menu bar at the top (File, Setup, Help) and four buttons for the various functions of the *NANOCOLOR*[®] T-Set software: “Record heating curve,” “Read test data,” “Read calibration data,” and “Settings.”

12.2 Menu

12.2.1 File

The File menu contains only the “Exit” option. Selecting this option will close the program and the window. Closing the window has no effect on the processes running in a connected heating block.

12.2.2 Setup

This option allows the user to make individual settings for various program parameters.

12.2.3 Connection

COM (serial interface) setup, factory settings:

Port: COM 1

Baud rate: 2400

Data bits: 8

Stop bits: 1

Parity: none

Flow control: none

12.2.4 Language

Select the desired language. Changes are applied immediately.

12.2.5 Configuration

User

The information entered under name, address and user name will later appear on reports and certificates.

Connection

The automatic device search can be switched ON or OFF.

12.2.6 Help

Help on various functions of the T-Set and the T-Set software

12.3 Record heating curve

With this function, heating curves can be recorded, displayed, and saved. The control panel of the heating block (temperature and time displays as well as the keypad) is also shown on the screen. The desired heating program can be selected either via the symbol button icons in the *NANOCOLOR*[®] T-Set software or directly on the heating block using the symbol buttons. Start the heating program by pressing the “Start” icon. The current temperature curve is displayed on the screen and can be saved at any time using the “Save” button. The heating curve is automatically formatted. The vertical axis displays the temperature, and the horizontal axis displays the time. The two horizontal lines that become visible after starting the recording of a heating curve indicate the permissible temperature limit range of ± 2 °C. Yellow vertical lines indicate the start and end of the heating process. The timer starts as soon as the temperature reaches the programmed value ± 2 °C. The

“Open” button allows you to access previously saved temperature curves. When a saved file is opened, the software automatically switches to “offline mode.”

To continue measurements and close the opened heating curve, this option must be deactivated again.

After selecting the “Report” button, you can enter information for a recorded heating curve (title, comment, etc.). When the “Report” button is selected, the *NANOCOLOR*[®] T-Set software generates a protocol of the heating curve recording with all entered comments, which can then be saved and/or printed.

12.4 Read test data

NOTICE

If the heating block does not contain any test data or the information has already been entered, there is no longer any possibility to enter or change anything (the fields are greyed out). This is only possible if a new test procedure is carried out beforehand.

With this function, test data that was stored in the heating block after running a test program with the *NANOCOLOR*[®] USB T-Set can be read out, displayed, and edited. In the input fields, after the test has been performed, you can enter the information

“Date of test execution,” “*NANOCOLOR*[®] USB T-Set serial number,” the corresponding “Certificate number,” and the “Date of issue.” Comments can be entered in the “Comments” field.

The temperature table is divided into two columns. The “Target” column shows all temperatures programmed in the heating block. The “Actual” column displays the most recently measured temperatures with the *NANOCOLOR*[®] USB T-Set.

When the “Report” button is selected, the *NANOCOLOR*[®] USB T-Set software generates a certificate of the test execution, which is automatically sent to the PC’s default printer.

12.5 Read calibration data

NOTICE

If the heating block does not contain any calibration data or the information has already been entered, there is no longer any possibility to enter or change anything (the fields are greyed out). This is only possible if a new test procedure is carried out beforehand.

With this function, calibration data that was stored in the heating block after performing a calibration with the *NANOCOLOR*[®] USB T-Set can be read out, displayed, and edited. In the input fields, after the calibration has been performed, you can enter the

information “Date of calibration,” “*NANOCOLOR*[®] USB T-Set serial number,” the corresponding “Certificate number,” and the “Date of issue.” Comments can be entered in the “Comments” field.

The temperature table is divided into two columns. The “Target” column shows all temperatures programmed in the heating block. The “Actual” column displays the temperatures measured with the *NANOCOLOR*[®] USB T-Set. The “Calibrated to” column shows the calibrated temperature values.

When the “Report” button is selected, the *NANOCOLOR*[®] USB T-Set software generates a certificate of the calibration, which is automatically sent to the PC’s default printer.

12.6 Settings

In this menu, the current device settings can be reviewed and edited. The time table is divided into two sections. The gray cells contain the heating times preset by MACHEREY-NAGEL, which cannot be changed. In the white cells, user settings can be entered or modified. Temperature values must be entered as whole numbers, and time values in the hh:mm format. The “Device Data” field contains the serial number and software version. In the “Sound” field, the periodic signal tone after a heating program can be switched on or off. The “Error Log” field records errors and is read-only; its content can be used for troubleshooting if problems occur. The “Save” button allows the settings on the screen to be saved to the device memory.

13. Maintenance and cleaning of the device

NOTICE

For all cleaning tasks, the device must be switched off and disconnected from the power supply.

WARNING

Risk of fire: Do not use flammable cleaning agents or organic solvents to clean the device.

NOTICE

Do not clean the device with acetone or similar products.

NOTICE

Use only accessories and spare parts approved by the manufacturer. If non-approved parts are used, the operator is responsible for the conformity of the device and the warranty will be void.

13.1 Cleaning the device

The casing of the USB T-Set can be cleaned with a damp cloth. Wipe off splashes immediately. The temperature sensor can be cleaned with a cloth and cleaning agent, but should always be kept clean.

14. Malfunctions, causes and correction

NOTICE

Malfunction or damage to the device due to improper maintenance: Maintenance of the device may only be carried out by qualified personnel.

Depending on the operating status, different messages may be displayed. The source of the error can be either incorrect operation or a malfunction of the device.

In the event of recurring errors, please contact MACHEREY-NAGEL.

Error	Cause	Solution
Display remains dark.	No power supply	Check whether the USB T-Set is connected correctly and the interface is supplying sufficient voltage.
USB T-Set does not fit into the bore of the heating block.	Contamination of the T-Set bore	Check the bore for the T-Set for contamination and clean it if necessary.
	Incorrectly installed safety cover of the heating block	Check whether the safety cover of the heating block is correctly installed. If necessary, loosen the screw in the middle of the safety cover, insert the USB T-Set and lock the screw again.
	Defect USB T-Set	The thermal sensor of the USB T-Set is damaged or bent. The USB T-Set must be repaired.

Error	Cause	Solution
The heating block or the USB T-Set are not recognized in the PC software.	Drivers for connecting the heating blocks are not installed	Check in the computer's device manager whether the heating block or the USB T-Set is recognized. If necessary, install the required driver via Windows or download it from the MN website.
	Wrong cable used for the connection or cable defective	Check whether the cable used is damaged or whether the correct cable has been used (REF 919687).

15. Service

15.1 Disposal

NOTICE

Disposal via public waste disposal systems is not permitted. Please contact your local MACHEREY-NAGEL representative. Dispose of the device in accordance with EU Directive 2012/19/EU. In compliance with this directive, MACHEREY-NAGEL will take back the old device and dispose of it free of charge.

15.2 Warranty, liability and complaints

The warranty for this device is valid for 24 months from the date of purchase. The original invoice serves as proof and must be presented when making a claim. The warranty becomes void in the event of improper handling and/or maintenance of the device. It does not cover defects resulting from the use of any external power supply other than the one provided.

The warranty is limited to the repair of defective parts or - at the discretion of MACHEREY-NAGEL - the delivery of a fault-free replacement device. Making a warranty claim does not affect the 24-month warranty period. There is no right of withdrawal. Further claims are excluded. This especially includes all claims for damages resulting from consequential or indirect damages. In addition, our general terms and conditions of sale and delivery in their current version, as printed on all price lists, apply.

15.3 Contact

MACHEREY-NAGEL GmbH & Co. KG

If you still have questions or require technical assistance after reading the operating instructions, please contact:

MACHEREY-NAGEL GmbH & Co. KG

Valenciennes Str. 11 · 52355 Düren · Germany

Phone: +49 2421 969-333

E-Mail: support@mn-net.com

www.mn-net.com

15.4 Version history

Manual **NANOCOLOR**[®] USB T-Set EN, V1.00, 06.2014, June 2014

Manual **NANOCOLOR**[®] USB T-Set EN, V2.00, 04.2025, April 2025

Water Analysis

MACHEREY-NAGEL GmbH & Co. KG
Valenciennes Str. 11
52355 Düren · Germany

Tel.: +49 24 21 969-0
E-mail: info@mn-net.com

A052302/0550.05

MACHEREY-NAGEL

www.mn-net.com