

Analysis of 25-Hydroxyvitamin D in serum by automated sample cleanup and LC-MS/MS

MACHEREY-NAGEL application department

Dr. H. R. Wollseifen, Dr. A. von Bohl, T. Kretschmer, L. Schaefer

Application benefits

- Automated sample cleanup of 25-Hydroxyvitamin D in serum
- Analysis without time-consuming, costly evaporation and reconstitution steps
- Facilitates high sample throughput utilizing automation in a multi-well plate format

MN products

REF 763632.20

EC HPLC column (analytical), NUCLEOSHELL® Biphenyl, 2.7 µm, 50x2 mm

REF 738921.010M

SPE MULTI 96-well plate CHROMABOND® HLB, 96x 10 mg, 30 µm

REF 702107

Screw closure, N 9, PP, yellow, center hole, Silicone white/PTFE red, 1.0 mm

REF 702079

Screw neck vial, N 9, 11.6x32.0 mm, 1.5 mL, label, flat bottom, amber, silanized

Eppendorf product

epMotion® 5075vt –
Catalog no. 5075000044
(www.eppendorf.com/epMotion5075)

MN application numbers

SPE: 306990
HPLC: 129580

Keywords

25-Hydroxyvitamin D2,
25-Hydroxyvitamin D3, serum,
automated analysis, SPE MULTI
96-well plate CHROMABOND® HLB,
NUCLEOSHELL® Biphenyl

Introduction

Vitamin D is the generic term for a group of fat-soluble vitamins, the so-called calciferols, which belong to the secosteroids. The most important forms include vitamin D2 (ergocalciferol) and vitamin D3 (cholecalciferol). They are mainly known for their function in calcium metabolism and therefore play a key role in bone mineralization. Additionally evidence of a link between vitamin D supply and a variety of diseases (hypertension, type 2 diabetes mellitus, and cardiovascular and cancer diseases) was found [1].

The determination of vitamin D level is therefore increasingly necessary and is carried out by measuring 25-Hydroxyvitamin D, 25-OH Vit Dx (X = 2 and 3), in blood serum. In recent years, liquid chromatography combined with mass spectrometry (LC-MS/MS) has become established for the measurement of 25-OH Vit D in serum [2]. In order to ensure the highest possible sample throughput, an automated SPE methodology using a 96-well plate, packed with CHROMABOND® HLB sorbent, was developed in this application note. Automated solid phase extraction procedure was performed using the epMotion® 5075vt robot.

High recovery rates with very good reproducibility are achieved for serum samples. Finally, the extracts are analyzed using HPLC-MS/MS on a NUCLEOSHELL® Biphenyl column.




Figure 1: Molecular formula of 25-Hydroxyvitamin D (A: 25-hydroxyvitamin D2, B: 25-hydroxyvitamin D3, C: d6-25-hydroxyvitamin D3)

Analysis of 25-Hydroxyvitamin D in serum by automated sample cleanup and LC-MS/MS

Sample pretreatment

MN Appl. No. 306990

Sample preparation:


1. Put 150 μL of serum sample in 96-deep-well plate.
2. Add 20 μL internal standard-solution ($\beta = 250 \text{ ng/mL}$ in methanol).
3. Add 150 μL aqueous zinc sulphate solution (0.2 M).
4. Add standard-solution ($x \text{ } \mu\text{L}$ in methanol*).
5. Add 600 μL ($- x \text{ } \mu\text{L}$ standard-solution**) of methanol.
6. Centrifugate for 5 min, 750 x g.

Procedure of the addition of standard solution before precipitation of serum

| Concentration sample solution (ng/mL) | * Volume of standard solution (μL) | Concentration standard solution (ng/mL) | ** Volume of methanol (μL) |
|---------------------------------------|---|---|---|
| 0 | 0 | 0 | 580 |
| 2 | 28.1 | 8 | 552 |
| 6 | 28.1 | 32 | 552 |
| 9 | 28.1 | 48 | 552 |
| 21 | 28.1 | 112 | 552 |
| 51 | 28.1 | 272 | 552 |
| 99 | 28.1 | 528 | 552 |
| 201 | 28.1 | 1072 | 552 |
| 300 | 28.1 | 1600 | 552 |

Automated solid phase extraction procedure on the epMotion[®] 5075vt robot

| Parameter | Description |
|--------------|--|
| Column | SPE MULTI 96-well plate CHROMABOND [®] HLB, 96x 10 mg, 30 μm |
| Conditioning | 500 μL methanol, vacuum pressure |
| Loading | 600 μL pretreated sample (low vacuum pressure, e. g. -100–200mbar) |
| Wash 1 | 500 μL 5 % methanol (aq); vacuum pressure |
| Wash 2 | 500 μL 60 % methanol (aq); vacuum pressure |
| Drying | / |
| Elution 1 | 240 μL 95:5 (v:v) methanol/ isopropanol |
| Elution 2 | 150 μL LC-MS-water |



Analysis of 25-Hydroxyvitamin D in serum by automated sample cleanup and LC-MS/MS

Analysis by HPLC-MS / MS

MN Appl. No. 129580

Chromatographic conditions

| | | |
|-------------|--|----|
| Column | NUCLEOSHELL® Biphenyl, 2.7 µm, 50x2 mm | |
| Eluent A | 0.1 % formic acid in water | |
| Eluent B | 0.1 % formic acid in acetonitrile | |
| Gradient | Time % B | |
| | 0 min | 65 |
| | 2 min | 65 |
| | 3 min | 95 |
| | 6 min | 95 |
| | 7 min | 65 |
| | 10 min | 65 |
| Flow rate | 0.4 mL/min | |
| Temperature | 40 °C | |

| | |
|--|---------------|
| Injection volume | 5 µL |
| MS conditions for Shimadzu 8050 – Triple Quadrupole MS | |
| Acquisition mode | MRM |
| Interface | ESI |
| Polarity | positive |
| Nebulizing gas flow | 3 L/min |
| Heating gas flow | 10 L/min |
| Interface temperature | 300 °C |
| DL Temperature | 250 °C |
| Heat Block Temperature | 400 °C |
| Drying gas Flow | 10 L/min |
| CID gas | 270 kPa Argon |


MRM transitions

| Analyte | Abbreviation | CAS number | Q1 mass [Da] | Q3 mass [Da] | Q3 mass [Da] | Q3 mass [Da] | Retention time [min] |
|-------------------------|-----------------|-------------|--------------|--------------|--------------|--------------|----------------------|
| 25-hydroxyvitamin D2 | 25-OH Vit D2 | 21343-40-8 | 395.20 | 105.10 | 135.30 | 137.10 | 4.26 |
| 25-hydroxyvitamin D3 | 25-OH Vit D3 | 140710-94-7 | 383.15 | 365.30 | 120.90 | 271.30 | 3.96 |
| d6-25-hydroxyvitamin D3 | d6-25-OH Vit D3 | 78782-98-6 | 389.20 | 371.25 | 133.30 | 45.10 | 3.92 |

Table 2: MRM transitions and retention times of 25-hydroxyvitamin D.

Chromatograms

Figure 2: a



Analysis of 25-Hydroxyvitamin D in serum by automated sample cleanup and LC-MS/MS

Figure 2: b

Q 383.15>365.30 (+)

7.43e4




Figure 2: c

ISTD 389.20>371.25 (+)

2.19e4




Figure 2: Chromatograms of a matrix matched standard solutions (concentration $\beta = 21.0$ ng/mL, A: 25-OH Vit D2, B: 25-OH Vit D3, C: d6-25-OH Vit D3)

Analysis of 25-Hydroxyvitamin D in serum by automated sample cleanup and LC-MS/MS

Calibration curves




Figure 3: The matrix matched calibration range for 25-OH Vit D2 and 25-OH Vit D3 between 1.5 – 301.0 ng/mL

Recovery rates




Figure 4: Recovery rates for the presented automated SPE method using SPE MULTI 96-well plate CHROMABOND® HLB, 96x 10 mg, 30 µm.

Analysis of 25-Hydroxyvitamin D in serum by automated sample cleanup and LC-MS/MS

| Standard addition concentration (ng/mL) | Recovery rate [%] 25-OH Vit D2 | Recovery rate [%] 25-OH Vit D3 |
|---|--------------------------------|--------------------------------|
| 1.5 | 108 % ± 13 % | 51 % ± 8 % |
| 6 | 110 % ± 9 % | 69 % ± 3 % |
| 9 | 117 % ± 4 % | 107 % ± 3 % |
| 21 | 103 % ± 9 % | 93 % ± 4 % |
| 51 | 108 % ± 3 % | 96 % ± 2 % |
| 99 | 110 % ± 10 % | 99 % ± 7 % |
| 201 | 112 % ± 8 % | 100 % ± 2 % |
| 300 | 108 % ± 2 % | 102 % ± 2 % |

Recovery rates for the presented automated SPE method using SPE MULTI 96-well plate CHROMABOND® HLB, 96x 10 mg, 30 µm.

Conclusion

This application note presents a reliable and successful determination of 25-Hydroxyvitamin D2 and 25-Hydroxyvitamin D3 in serum. By using a SPE MULTI 96-well plate CHROMABOND® HLB (96x 10 mg, 30 µm) on the epMotion® 5075vt robot, it was possible to achieve high recovery rates and good reproducibility. As a result for the tested standard addition values, most recovery rates are between 90 % and 110 %. Furthermore, the presented SPE methodology eliminates the need for time-consuming solvent evaporation and reconstitution steps.

The identification and the quantification of 25-Hydroxyvitamin D in serum was finally carried out by ESI mass spectrometry on a NUCLEOSHELL® column. Using core-shell particle technology allows, highest column efficiency and resolution at a short run time with much lower back pressure compared to fully porous particles could be achieved with common HPLC systems. The Biphenyl modification successfully helps to separate analytes from matrix with different retention mechanisms: π - π interactions and hydrophobic interactions.

The chromatography application database

- Free access to more than 3,000 application examples from HPLC, GC, TLC and SPE
- Free access to the MN application database:
<https://chromaapdb.mn-net.com>



www.mn-net.com


MACHEREY-NAGEL

References

- [1] Teresa Kulie et al, Vitamin D: An Evidence-Based Review, JABFM, 2009, vol. 22 no 6.
- [2] Kornelia Galior, Hemamalini Ketha, Stefan Grebe, Ravinder J Singh: 10 years of 25-hydroxyvitamin-D testing by LC-MS/MS-trends in vitamin-D deficiency and sufficiency Bone Reports (2018), 268–273

MACHEREY NAGEL shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. Information, descriptions, and specifications in this publication are subject to change without notice.

Registered trademarks
Eppendorf, the Eppendorf Brand Design, epMotion® and epT.I.P.S.® are registered trademarks of Eppendorf SE, Germany.
CHROMABOND® MACHEREY-NAGEL GmbH & Co. KG (Germany)
NUCLEOSHELL® MACHEREY-NAGEL GmbH & Co. KG (Germany)



© Vacilando (p. 1) – stock.adobe.com



Management System
EN ISO 13485:2016
ISO 9001:2015

www.tuv.com
ID 0000056401

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

