MACHEREY-NAGEL

Photometer PF-12^{Plus}

- Handbuch
- Manual
- Manuel

Deutsch

3

English

37

Français

69

Español

105

MACHEREY-NAGEL

Photometer PF-12^{Plus}

Handbuch

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Photometer PF-12^{Plus}

Manual

1. General information and introduction	40
2. CE Marking	40
3. Type label	40
4. Explanation of symbols	41
5. Safety instructions.....	42
6. Chemical and biological safety	42
7. Technical description	43
7.1 Technical data	44
8. Setting up and connecting the device	45
8.1 Scope of delivery	45
8.2 Connecting the device	45
8.3 Battery operation	46
8.4 Battery pack operation	47
8.5 Transport.....	47
9. Device views	48
10. Initial operation	48
10.1 Turning on the device	48
11. Operation and user guidance	49
11.1 Keypad	49
11.2 Configuration	50
12. Measurements	53
12.1 Performance of <i>VISOCOLOR® ECO</i> tests.....	54
12.2 Performance of <i>NANOCOLOR®</i> tube tests	55
12.3 Changing the submethod	55
12.4 Photometric basic functions	55
12.5 Correction value determination	56
12.6 Determination of correction values for <i>NANOCOLOR®</i> tube tests.....	57
13. NTU-Check	60
14. Quality control	60
15. Memory	61
15.1 Display memory	61
15.2 Search in memory	61
15.3 Delete memory.....	61
15.4 Output from memory	61
16. Special methods.....	62
16.1 User defined methods.....	62
16.2 Sample numbering	63
16.3 Identification of sample location.....	63

16.4 Sample dilution 64

16.5 Preprogrammed special methods 64

17. Maintenance and cleaning 64

 17.1 Cleaning of the device..... 65

18. Errors, sources and solutions 65

19. Service 66

 19.1 Spare parts, accessories and consumables 66

 19.2 Disposal 67

 19.3 Warranty, liability and complaints 67

 19.4 Contact..... 67

 19.5 UKCA Authorized Representative 67

 19.6 Version history 68

1. General information and introduction

The compact photometer PF-12^{Plus} is a filter photometer from MACHEREY-NAGEL for mobile use and for use in laboratories. It is used for the evaluation of photometric test kits in 16 mm tubes from MACHEREY-NAGEL. Besides the measurement of absorbed light it offers the possibility to determine the nephelometric turbidity units by measuring the light scattering at a 90° angle. The illuminated display and the intuitive menu navigation allow a fast, comfortable and pleasant working with the instrument and make it the ideal compact photometer in the whole spectrum of water analysis.

2. CE Marking

The CE marking declares that the product complies with the following European Community harmonisation legislation:

European Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS 2)

IEC 61010-1

IEC 61010-2-081

IEC 61326-1

FCC47 CFR Part 15B, ICES-003 Issue 6

IP68

MIL 810H, 514.8

3. Type label


Figure 1: Type label

4. Explanation of symbols

The device type plate and these operating instructions contain the symbols or terms listed below, which have the following meaning

Term / Symbol	Meaning
Photometer PF-12 ^{Plus}	Device type designation
SN	Serial number of the device
4.8–6 V, 500 mA	Power supply
	Direct voltage
IP68	Dust and water proof (Test conditions: 1 m water depth, 1 hour)
	Identification USB interface
	Attention (observe documentation)! Observe the safety instructions in the operating instruction of the device
	According to 2012/19/EU, it is prohibited to dispose of the device through public waste disposal systems. Note also the information in the "Disposal" section.
MACHEREY-NAGEL	Identification of the manufacturer
	This symbol indicates that there is a risk of electric shock.
	This symbol indicates that the chemicals used have a corrosive effect. Observe the safety measures in the laboratory and wear the prescribed protective equipment. Observe the instructions in the current safety data sheets (SDS) of the products used.
	This symbol indicates that there may be a danger by using flammable substances.
CE	CE Mark
UKCA	UKCA Mark

5. Safety instructions

BE SURE TO READ THE FOLLOWING SAFETY INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE.

Failure to follow these instructions may result in serious injury to the operator, malfunction or damage to the equipment.

Keep this operating instruction in a safe place for future reference.

Follow the safety notes and instructions in the operating instruction and observe the stickers and notices on the device.

Do not work on internal parts of the unit. Non-compliance will invalidate any warranty claim.

Use of the hazard warnings:

DANGER

Indicates an imminent or potentially hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Indicates an imminent or potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may result in malfunction or damage of the instrument. Information that clarifies facts in the text and that requires special attention.

6. Chemical and biological safety

WARNING

Possible dangers of contact with chemical and biological substances.

Working with samples, reagents and corresponding accessories is associated with dangers.

Wear suitable protective equipment when working with the cuvettes. Observe the safety data sheet (SDS) of the test kit used.

During normal operation of this device, it may be necessary to use chemicals that are hazardous to your health or biologically harmful samples.

Before handling these substances, read all hazard statements and safety information printed on the containers with the original solutions and in the safety data sheets. All used solutions must be disposed of in accordance with the national laws and regulations. The type of protective equipment must be selected according to the concentration and amount of the hazardous substance present at the corresponding workplace.

7. Technical description

The PF-12^{Plus} compact photometer enables measurements in the wavelength range from 345–690 nm. The light generated by the xenon lamp is passed through an interference filter into the cuvette slot of the instrument. After passing through the sample, the proportion of absorbed light is calculated by detecting the attenuated light. The proportion of the absorbed light is converted into the concentration of the investigated analyte using a stored calibration. The PF-12^{Plus} allows the measurement of nephelometric turbidity (NTU) in the range 1–1000 NTU due to the specially placed 860 nm LED. The NTU-Check automatically controls the turbidity with each measurement.

7.1 Technical data

Type:	Single-beam filter photometer with microprocessor control, sefittest and autocalibration
Light source:	Xenon-High pressure lamp 860 nm LED
Optical system:	Automatic filter wheel with 7 interference filters
Wavelength range:	345 – 690 nm
Wavelengths:	345 / 436 / 470 / 540 / 585 / 620 / 690 nm plus 1 empty position for an additional filter 860 nm LED for NTU-measurement
Wavelength accuracy:	± 2 nm
Bandwidth at half transmission:	10 – 12 nm
Wavelength selection:	manual
Photometric range:	± 3,0 E
Photometric accuracy:	± 1 %
Height of the light beam	12 mm over bottom of cuvette slot
Long time stability:	< 0,002 E/h
Measurement modes:	More than 100 preprogrammed tests (<i>NANOCOLOR</i> ® Tube tests and <i>VISOCOLOR</i> ® ECO Test kits); Absorbance, Transmission, Factor, Standard, Turbidity; 50 free programmable methods (Starting with firmware version > 3.0.1: 20 freely programmable methods)
Turbidity measurement:	Nephelometric turbidity measurement 1 – 1000 NTU
Cuvette slot:	Test tubes 16 mm OD
Data memory:	1000 measured data sets; GLP-conform
Display:	Backlit graphic display, 64 × 128 Pixel
Operation:	Intuitive menu navigation, foil keypad, Call up tests via test number or from parameter list.
Languages:	DE / MY / TR / SI / IND / CZ / PT-BR / PL / HU / ES / IT / FR / NL / EN
External light:	Insensitive; open cuvette slot
Interfaces:	USB B 2.0
Protection class:	IP 68 (The device is dust and waterproof) IP6X ▲ Dust-tight IPX8 ▲ Protection against permanent immersion
Overvoltage category:	Device: I Power supply: II
Degree of pollution:	2
Update:	via PC
Operating conditions:	10 – 40 °C, up to 80 % rel humidity (non-condensing), up to max. 3000 m
Power supply:	via USB cable, USB power adapter, batteries / rechargeable batteries, battery pack
Dimensions:	215 × 100 × 65 mm
Weight:	0.7 kg
Warranty	2 years

The maximum sound level and maximum sound pressure are subject to the limits specified in ISO 3746 or ISO 9614-1.

8. Setting up and connecting the device

8.1 Scope of delivery*

Check the delivery for completeness. If parts of the delivery are missing, contact MACHEREY-NAGEL or your local distributor directly.

- Compact photometer PF-12^{Plus}
- USB cable AB
- 4 × AA Mignon batteries
- 4 × empty tubes 16 mm OD
- Calibration cuvette
- Funnel
- Beaker
- Syringe
- Software-DVD
- Certificate
- Screwdriver
- Operating instruction
- Packed in a rugged case

* The scope of delivery refers to the product itself and may differ if the photometer is included in another product.

8.2 Connecting the device

WARNING

Exposure of the eye to light. Do not look directly into the light beam of the lamp and the turbidity LED.

CAUTION

Multiple hazards: Only qualified personnel must perform the instructions described in this step.

WARNING

Risk of electric shock: Make sure that the USB cable is not damaged when used for operation. Check the suitability of the power source for the equipment.

NOTICE

Damage due to condensation and exceeding the temperature limits:
The unit is designed for indoor and outdoor use. Do not expose the device to direct sunlight. Protect the interface from the effects of moisture and dampness when used outdoors.

WARNING

Fire hazard: When using the external power supply, take care not to overload the power outlet. There is a risk of overload and fire.

Make sure that the power supply including cable is not damaged. Check the suitability of the voltage source used for the device. Use only the power supply recommended by MACHEREY-NAGEL.

WARNING

Risk of injury: Familiarize yourself with the device before working with it and read this document carefully. Do not use the device unless you have received instruction in its use.

NOTICE

Defects of the power supply and housing can lead to malfunction of the device. If the instrument shows an apparent breakage of the housing or a damaged power supply, it must be taken out of operation.

8.3 Battery operation

WARNING

Risk of burns: The batteries may leak. Remove the batteries from the device if the PF-12^{Plus} is not used for a longer period of time.

NOTICE

Damage due to leaking batteries: MACHEREY-NAGEL recommends the use of high quality batteries.

NOTICE

During use or charging the battery pack, the PF-12^{Plus} must be positioned in such a way that it can be safely switched off and the USB connection safely disconnected at any time.

The photometer PF-12^{Plus} works with batteries as voltage source. Four AA Mignon batteries (included in the scope of delivery) are inserted into the instrument. Rechargeable batteries, e.g. 4 × NiMH AA 2500 mAh, are also available as optional accessories. The charger and rechargeable batteries must be ordered separately. The charge level of the batteries is shown in the upper part of the display in the form of a battery symbol after switching on. When connected with a USB cable, the info "U" also appears next to the battery symbol.

When changing the batteries, always replace the complete set. The battery compartment is located on the bottom of the instrument. The cuvette slot must be empty. The instrument is placed upside down on a soft surface and the batteries are inserted as follows

The battery compartment cover may only be removed when inserting and replacing the batteries!

1. Loosen the screws of the battery compartment and remove the cover.
2. Insert 4 AA batteries (see scope of delivery) into the holder of the instrument. Make sure that the polarities on the batteries match the marking in the battery compartment.
3. Replace the battery compartment cover, tighten the screws again and turn around the device again.

NOTICE

Do not bend or crush the connecting cable of the battery compartment and battery pack.

Figure 2: Opened battery compartment of the PF-12^{Plus}

8.4 Battery pack operation

The PF-12^{Plus} can also be equipped with an internal battery pack, which is charged via the USB interface. After switching on, the charge status of the battery pack is shown in the upper part of the display in the form of a battery symbol with an "A". When connected with a USB cable, the info "L" appears next to the battery symbol to indicate that the battery is being charged. You can order this battery pack separately (REF 919201) and insert it into the photometer. The battery pack is inserted into the instrument as follows:

1. Loosen the screws of the battery compartment and remove the cover.
2. Remove the holder of the four AA batteries and carefully disconnect the plug connection to the instrument. Connect the plug of the battery pack to the connector of the PF-12^{Plus} and insert the battery pack into the battery compartment. The connecting cable must not be kinked or crushed during this process.

Replace the battery compartment cover, tighten the screws and turn the device around.

8.5 Transport

To transport the device, the delivery box with its protective interior is best suited. Therefore, keep the packaging. If you no longer have the packaging in hand, use a suitable cardboard box as outer packaging and soft padding material so that the device cannot slip during transport.

NOTICE

Choose a suitable packaging for shipment.

9. Device views

Figure 3: Top view PF-12^{Plus}


Figure 4: Rear view PF-12^{Plus}

- ❶ Grafic display 128 × 64 Pixel
- ❷ Cuvette slot for tubes with 16 mm outer diameter (OD)
- ❸ 20-part keypad
- ❹ USB B-interface (2.0) for connection of power supply, charging the optional battery pack, data transfer and update


10. Initial operation

10.1 Turning on the device

Switch the unit on using the I/O switch on the keypad. The display shows the device name and firmware version. After a few seconds the start screen for entering the method appears.

11. Operation and user guidance

The PF-12^{Plus} is designed for easy handling and operation. The photometer PF-12^{Plus} has a 20-part keypad: 10 alphanumeric keys and 10 keys with special functions.

The screen is divided into three areas. The top status bar shows the date, battery status and time. The actual display, e.g. of measurement results or menu items and the corresponding menu symbols, is in the middle of the screen, depending on the display. The bottom line, which is separated from the screen, is the taskbar. It shows which options are available in the currently selected menu. By pressing the displayed option with the corresponding keyboard command, the action is executed (e.g. **M** = Measure).

By pressing the **ESC** key you leave the respective menu and go back one level. Pressing the **ESC** key several times returns to the method selection menu.

11.1 Keypad


Key	On/Off switch
Keys ...	For entering digits, letters, special characters
Keys / /	Double-function keys, for sample number, location, dilution
Key	For entering decimal points / select menu for special methods
Key	Abort, return to method selection
Key	Scroll up
Key	Scroll down
Key	Enter, call up last result
Key	Call up memory functions / Memory on / off
Key	Absorbance program, recall the absorbance after measurement, estimated value and NTU value after measurement / selection of prefix
Key	Manual zero adjustment
Key	Measurement, consecutive measurements are numbered continuously, call up last test

11.2 Configuration

To get to the configuration mode enter „900“ in the method selection mode

or use / to select the entry „Settings“ in the list. In addition to calling up the settings via method “900”, it is also possible to call up individual settings quickly.

Each setting is assigned a number. Entering this three-digit number as the method in the start menu leads to the corresponding settings menu.

The configuration menu is based on icons. You can navigate through the menu via /, confirm by pressing and follow the instructions on the display. The icons have following meanings and functions:

11.2.1 Language (Method 954)

Select the desired language via / and press to confirm.

11.2.2 Time (Method 951)

Select the digits via / and adjust with the keypad. Confirming with opens the selection of the time format. Select the desired time format via / and press to confirm.

11.2.3 Date (Method 950)

Select the digits via / and adjust with the keypad. Confirming with opens a dialog to select the date format. The following format types are available:

23.01.2014

01/23/2014

2014-01-23

Select the desired date format via / and press to confirm.

11.2.4 Contrast (Method 959)

The contrast of the display can be changed via / in several steps. Press to store the setting.

11.2.5 Brightness (Method 960)

The brightness of the display can be changed via / in several steps. Press to store the setting.

11.2.6 Calibration (Method 956)

Press to start the calibration. The instrument asks to make sure there is no cuvette in the slot. Press to start the calibration with reference to air. Afterwards the PF-12^{Plus} asks you to insert the calibration tube (clean test tube filled with distilled water). Confirming with starts the water calibration.

11.2.7 Reaction time (Method 952)

Use / to turn the reaction time on or off and press to confirm.

11.2.8 Signal tone (Method 961)

Use / to turn the signal tone on or off and press to confirm.

11.2.9 20–80 % range (Method 957)

Use / to turn the 20–80 % range bar on or off and press to confirm. When enabled, the 20–80 % range will be displayed at every measurement and provides information on where the result lies in relation to the measurement range. If the measured value is located in the thickly marked area of the measuring range bar, measured result is outside the 20–80 % range. In this case the sample either has to be diluted or a test kit with a different measurement range has to be used.

11.2.10 Turbidity control (NTU-Check) (Method 958)

Use / to turn the turbidity control (NTU-check) on or off and press to confirm. An input mask for entering the turbidity warning limit appears. Enter the desired numerical value using the keyboard and confirm with . When this function is enabled, the nephelometric turbidity of the sample is determined parallel to each measurement. If the limit value entered is exceeded, an exclamation mark is displayed to the left of the measured value as a warning. The measured turbidity can be displayed after a measurement by pressing the key.

11.2.11 Filter 8

Enter the wavelength of the special filter using the keypad and press to confirm. This configuration can only be performed in the service mode. To calibrate the special filter insert the calibration tube (see scope of delivery) and press . Confirm with to store the entry of the special filter. Press to reset the settings.

11.2.12 AUTO-OFF (Method 953)

Use / to adjust the AUTO-OFF time option, which affects the automatic switch off of the instrument after the last action. The following options are possible: Off (the instrument is not turned off automatically), 5 min (default setting), 10 min, 15 min, 20 min, 60 min. Press to confirm.

11.2.13 SYSTEM-RESET

This configuration can only be made in the service mode. After an additional security confirmation, the photometer performs a SYSTEM-RESET, thus restoring the state of delivery.

NOTICE

When the SYSTEM RESET is executed, all measured data is deleted. Save your measurement data and your special methods if necessary.

11.2.14 PROG. UPDATE (Method 955)

In order to perform a program update, you must first download the update software and the new program file (www.mn-net.com) on your computer. Before you start the update procedure the photometer has to be connected to your computer via the USB-interface. The photometer is recognized automatically. If not, activate the update-mode manually. Therefore, confirm the program update with . After a successful update the photometer restarts automatically.

NOTICE

WARNING: Risk of incorrect measurement results and damage to the device.

For devices with serial number SN: NPF12P1XXXX, installing firmware earlier than version 4.0.0 will result in incorrect measurements and the “AD overflow” error.

Before performing the update, check your device’s serial number and the firmware version loaded in the software. Before starting a measurement, make sure your device is running the latest firmware version.

NOTICE

Your data may be lost during a program update. Save your measurement data and your special methods if necessary.

11.2.15 Calibration adjustment (Method 962)

The PF-12^{Plus} device is provided with the MACHEREY-NAGEL factory calibration and can be readily used without any need for user calibration. The calibration adjustment tool can be used to meet regulatory requirements. Therefore, the factory-provided calibration of the VISOCOLOR® ECO and NANOCOLOR® methods can be adjusted with a standard in the menu “calibration adjustment”. The parameter including submethod to be adjusted is chosen via the four-digit entry mask (e.g. 0291 for COD 1500 100 – 1500 mg/L). Place the tube with the zero value or the reagent blank into the cuvette slot of the PF-12^{Plus} and press  to start the zero measurement. Subsequently place the tube with the standard into the device and press  to start the measurement. Enter the desired value for the standard with the keypad and confirm your entry with  (Note: The entered value must be within the measurement range of the respective submethod, otherwise the entry will not be adopted). The changed calibration will be stored automatically and will be executed the next time the concerned method will be chosen. The changed calibration will be indicated in form of an asterisk behind the respective test name. When stored, the results measured with an adjusted calibration will be indicated with an asterisk, too. The calibration adjustment can be deleted and set back to factory calibration by pressing  in the same menu. The measurement values in the memory will not be affected by the reset.

NOTICE

Wrong results possible by calibration adjustment: We recommend to use the factory-provided calibration and not to change this one. The calibration adjustment can be used to meet legal requirements. When performing the adjustment, we recommend selecting a standard that is in the upper third of the measurement range of the test to be adjusted to minimize calibration errors.

NOTICE

Only the associated power supply unit GTM46101-1005-USB (REF 919220) may be used.

CAUTION

Multiple hazards: Only qualified personnel may perform the instructions described in this step.

WARNING

Risk of electric shock: Make sure that the USB cable is not damaged when used for operation. Check the suitability of the power source for the equipment.

NOTICE

Damage due to condensation and exceeding the temperature limits:
The unit is designed for indoor and outdoor use. Do not expose the device to direct sunlight. Protect the interface from the effects of moisture and dampness when used outdoors.

WARNING

Fire hazard: When using the external power supply, take care not to overload the power outlet. There is a risk of overload and fire.

Make sure that the power supply including cable is not damaged. Check the suitability of the voltage source used for the device. Use only the power supply recommended by MACHEREY-NAGEL.

WARNING

Risk of injury: Familiarize yourself with the device before working with it and read this document carefully. Do not use the device unless you have received instruction in its use.

NOTICE

Defects of the power supply and housing can lead to malfunction of the device. If the instrument shows an apparent breakage of the housing or a damaged power supply, it must be taken out of operation.

WARNING

Possible dangers of contact with chemical and biological substances.

Working with samples, reagents and corresponding accessories is associated with dangers. Wear suitable protective equipment when working with the cuvettes. Observe the safety data sheet (SDS) of the test kit used.

NOTICE

Contamination of the cuvette slot and incorrect measurement results: Clean the cuvettes with a lint-free cotton cloth before each measurement to avoid contamination of the cuvette slot and incorrect measurement results.
Also keep the cuvette slot free from dirt.

NOTICE

Note on intended use: The PF-12^{Plus} is designed exclusively for use in water analysis with the photometric test kits provided by MACHEREY-NAGEL. Any other use, in particular outside the intended areas of application or in combination with non-approved materials, constitutes misuse of the device and may lead to incorrect results and compromise the safety of the device. Any use other than that intended is expressly prohibited.

12.1 Performance of VISOCOLOR® ECO tests

A clean, dry tube (16 mm OD) is used instead of the VISOCOLOR® ECO measuring glass. A second tube (an untreated sample or reagent blank, depending on the test) is used for zero adjustment.


Switch on the photometer

03.06.2014 09:30
Method
-- --
↑↓ M: 029


The display „Method: _ _ “ appears.

A VISOCOLOR® ECO test can be performed using the test number (5-xx)

Method:
Basic functions
Tube tests
- VISOCOLOR-ECO
Bio tests
↑↓


as well as with the keys  /  via the list of VISOCOLOR® ECO tests.

03.06.2014 09:30
5081 AMMONIUM 3
↑↓ 1 2 3 6 2

Performing the zero measurement

Insert the clean tube with the zero value into the photometer and press the  key. Remove the vial with the zero value.

03.06.2014 09:30
5081 AMMONIUM 3
↑↓ 1 2 3 6 2 M

Performing the measurement

Place a clean vial with sample solution in the photometer, press the  key and read off the measured value.

For detailed information on the individual tests, refer to the corresponding test instruction or the pictogram instructions

VISOCOLOR® ECO on the PF-12^{Plus}, which are available separately. These manual pages in pictogram form can be found on the Internet at www.mn-net.com and on the enclosed software DVD.

12.2 Performance of NANOCOLOR® tube tests

The blank values of the NANOCOLOR® tube tests are prestored in the PF-12^{Plus}. A zero measurement is therefore only necessary for isolated tests.


Switch on the photometer

The display „Method: _ _ _“ appears.

A NANOCOLOR® tube test can be performed using the test number (0-xx)

As well as with the keys  /  via the list „tube tests“.

Performing the zero measurement (for isolated tests)

Place a clean vial with the zero value in the photometer and press the  key. Remove the vial with the zero value.

Performing the measurement

Place a clean vial with prepared measuring solution in the photometer, press the  key and read off the measured value.

12.3 Changing the submethod

To change the submethod (e.g. to determine the result in a different unit), press the  /  keys in the measuring menu of the selected test. The currently active submethod is indicated by the number marked with an underscore in the task bar. After selecting the desired submethod, its measuring range is briefly displayed. To display the measurement result after changing the submethod, perform the measurement again. The result is then displayed in the unit stored for the submethod. For some parameters only one submethod is available.

For detailed information on the individual tests, refer to the corresponding test instruction. Manual pages in pictogram form can be found on the Internet under www.mn-net.com and on the software DVD accompanying the photometer.

NOTICE

The first time a method is selected by entering the three-digit method number, submethod 1 is automatically called up. Afterwards, the last submethod used is called up with each subsequent call of a method.

12.4 Photometric basic functions

The basic photometric functions can be accessed by entering the appropriate method number 901 – 903 and 905 – 906 or by using the  /  keys from the basic functions list.

12.4.1 Methode 901: Measurement with factor

Select the desired wavelength using / and confirm with . Enter the factor via the keypad and confirm with . Perform the measurement by following the instructions on the display. Return to the method selection by pressing . The factor will be erased.

12.4.2 Methode 902: Measurement with standard

Select the desired wavelength using / and confirm with . Enter the standard concentration via the keypad and confirm with . Perform the measurement by following the instructions on the display.

Return to the method selection by pressing . The standard concentration will be erased

12.4.3 Methode 903: Measurement of absorbance

Select the desired wavelength using / and confirm with . Perform the measurement by following the instructions on the display. Return to the method selection by pressing .

12.4.4 Methode 905: Measurement of transmission

Select the desired wavelength using / and confirm with . Perform the measurement by following the instructions on the display. Return to the method selection by pressing .

12.4.5 Methode 906: Measurement of nephelometric turbidity

Perform the measurement by following the instructions on the display. Return to the method selection by pressing .

12.4.6 Turbidity calibration

The turbidity calibration can be performed via this menu and is started by pressing . One after the other the turbidity solutions from the test kit *NANOCONTROL NANOTURB* (REF 925702) are put into the photometer and measurement is started with . After the calibration is finish it can be confirmed with .

NOTICE

Loss of turbidity calibration: The turbidity calibration cannot be reset to factory settings after it has been performed. Cancel an accidentally triggered turbidity calibration by pressing the key.

12.5 Correction value determination

These instructions are only to be used in connection with the corresponding test instructions of the *NANOCOLOR®* cuvette tests.

The photometric analysis of water samples with inherent color or turbidity always requires determination of a correction value. Colors and turbidity cause an increased light absorption, thus leading to deviating results. Determination of correction values require individual procedures for every test. For example, it is not always sufficient to measure the intrinsic color of the sample without reagents and subtract it from the measured value. In many cases, the reagents change

the intrinsic color or turbidity of the sample. All changes to the sample during analysis, such as dilution or addition of chemicals that influence pH or redox potential, must be reproduced as in the original analysis.

With the PF-12^{Plus}, the measurement program for the correction value is started after the measurement of the (turbid or colored) sample (value A) by pressing . The instrument requests inserting and measuring the correction value (value B). The corrected measurement result is displayed and stored. For some tests it is necessary to measure an additional correction blank value.

Underlying procedure:

Determine measured value according to test instruction = A

Determine correction value according to special instruction = B

Analysis result = A – B

Exceptions: Methods in which decreasing absorbance is measured against a reagent blank value. Then the following applies: Analysis result = A + B. The corresponding special instructions point out this fact.

It is very important to subtract only values with equal dimensions (e.g. mg/L N; mg/L NH₄; mmol/m³; A).

If the correction factor is so low in the same matrix with several samples that it is negligible for the measurement problem at hand, the determination can be omitted. However, this only results from practical testing and is not evident beforehand!

12.6 Determination of correction values for NANOCOLOR® tube tests

When measuring the correction value, a clean, empty cuvette filled with distilled water is used as the zero value

Test	Tube for correction (value B)
0-01 Zirconium	Proceed as described in the test instructions for Test 0-01, but do not add NANOFIX R2, seal, mix.
0-02, 0-03, 0-04, 0-05, 0-06, 0-08 Ammonium 3 2000	Proceed as described in the test instructions for test 0-02 / 0-03 / 0-04 / 0-05 / 0-06 / 0-08, but do not add NANOFIX R2, close, mix.
0-07 AOX 3	Almost all colors and turbidities are destroyed under test conditions and do not interfere. Resistant colors and turbidities cause deviating results which cannot be circumvented.
0-09 Lead 5	The original test already contains a correction.
0-13 Boron 3	No correction possible.
0-14 Cadmium 2	Fill empty test tube with 4.0 mL sample, add 0.2 mL R2, close and mix.
0-15 Carbonate hardness 15	Open carbonate hardness test tube, add 4.0 mL sample solution, close, mix and adjust to zero (value B). Open test tube again, add NANOFIX R2, close and shake well. Measure after 2 min (analytical result = A – B).
0-16 Permanganate Index 10	Almost all colors and turbidities are destroyed under test conditions and do not interfere. COD resistant colors and turbidities cause deviating results which cannot be circumvented.
0-17 Chlorine / Ozone 2	Fill empty test tube with 4.0 mL sample for each test.
0-18 Chlorine dioxide 5	Fill empty test tube with 4.0 mL sample for each test.

Test	Tube for correction (value B)
0-19 Chloride 200	Open chloride test tube, add 1.0 mL sample solution and 1.0 mL distilled water, close, mix.
0-21 Chloride 50	Open chloride test tube, add 4.0 mL sample solution and 1.0 mL distilled water, close, mix.
0-24 Chromate 5	Fill empty test tube with 4.0 mL sample, add 0.2 mL R2, close, mix.
0-25 Easily liberated cyanide 04	No correction possible.
0-11, 0-12, 0-20, 0-22, 0-23, 0-26, 0-27, 0-28, 0-29, 0-30, 0-33, 0-36, 0-38 COD 40–60000	Almost all colors and turbidities are destroyed under test conditions and do not interfere. COD resistant colors and turbidities cause deviating results which cannot be circumvented.
0-31 Cyanide 08	Proceed as described in the test instructions for test 0-31, but add 0.5 mL distilled water instead of 0.5 mL R3.
0-32 Anionic surfactants 4	Proceed as described in the test instructions for test 0-32, but instead of 0.5 mL R2 add 0.5 mL distilled water.
0-34 Cationic surfactants 4	Proceed as described in the test instructions for test 0-34, but instead of 0.5 mL R2 add 0.5 mL distilled water.
0-35 DEHA 1	Open DEHA test tube, add 4.0 mL sample, close and mix.
0-37 Iron 3	No correction possible.
0-39, 0-77 TOC 30–300 Compact	No correction possible.
0-40 Fluoride 2	No correction possible.
0-41 Formaldehyde 8	Open formaldehyde test tube, add 2.0 mL sample solution, close and mix.
0-42 Zinc 6	Proceed as described in the test instruction for test 0-42, but instead of 0.2 mL R3 add 0.2 mL distilled water.
0-43 Hardness 20	Open hardness test tube, add 0.2 mL sample solution, close and mix.
0-44 Hardness Ca / Mg	No correction possible.
0-45 Potassium 50	Open potassium test tube, add 2.0 mL sample solution, close and mix.
0-46 Formaldehyde 10	Fill empty test tube with 2.0 mL distilled water, 2.0 mL sample solution and 1.0 mL R2, close and mix.
0-47 Nonionic surfactants 15	No correction possible.
0-48, 0-51 TOC 30–300 Neo	No correction possible.
0-49 Silver 3	Proceed as described in the test instruction for test 0-49, but instead of 0.5 mL R3 add 0.5 mL DMSO.
0-50 organic acids 3000	Proceed as described in the test instruction for test 0-50, but instead of 2.0 mL R5 add 2.0 mL distilled water.
0-52 Complexing agents 10	Fill empty test tube with 4.0 mL sample and 1.0 mL distilled water, close and mix (analytical result = A + B).
0-53 Copper 5	Proceed as described in the test instruction for test 0-53, but add 0.1 mL DMSO instead of 0.1 mL R2.
0-54 Copper 7	Fill empty test tube with 4.0 mL sample and 0.4 mL distilled water, add 0.2 mL R2, close and mix.
0-56 Molybdenum 40	No correction possible.

Test	Tube for correction (value B)
0-57 Hydrocarbons 300	The original test already contains a correction.
0-58 Manganese 10	Fill empty test tube with 4.0 mL sample solution, 0.5 mL distilled water and 0.5 mL R2, close and mix. Add 1 measuring spoon R3, close and shake vigorously.
0-59 total Chromium 2	Almost all colors and turbidities are destroyed under test conditions and do not interfere. For samples which are still colored or turbid after decomposition: Fill empty test tube with 4.0 mL digestion solution.
0-60 Sulfate MR400	Fill empty test tube with 1.0 mL distilled water. Add 4.0 mL sample.
0-62 Sulfate LR 200	Fill empty test tube with 1.0 mL distilled water. Add 4.0 mL sample.
0-63 Sulfate HR 1000	Fill empty tube with 3.7 mL distilled water and 1 mL sample.
0-64, 0-65 Nitrate 8-50	Open nitrate test tube, add 0.5 mL sample and 0.5 mL 2-propanol, close, mix. Blank value for correction: Open nitrate test tube, add 0.5 mL distilled water and 0.5 mL 2-propanol, close, mix.
0-66 Nitrate 250	Open nitrate test tube, add 0.2 mL sample and 0.5 mL 2-propanol, close, mix. Blank value for correction: Open nitrate test tube, add 0.2 mL distilled water and 0.5 mL 2-propanol, close, mix.
0-67 total Kjeldahl Nitrogen TKN 16	No correction possible.
0-68 Nitrite 2	Fill empty test tube with 4.0 mL sample, add 0.2 mL R2, close, mix.
0-69 Nitrite 4	Open nitrite test tube, add 4.0 mL sample solution, close, mix.
0-70 POC 200	Proceed as described in the test instruction for test 0-70, but instead of 1.0 mL R2 add 1.0 mL distilled water.
0-71 Nickel 4	Proceed as described in the test instruction for test 0-71, but instead of 1.0 mL R2 add 1.0 mL NaOH 14 %.
0-73 Sulfide 3	Fill empty test tube with 0.5 mL sulfuric acid 50 %, add 1 measuring spoon R2 and 4.0 mL sample solution, close and shake gently. Add 200 µL R3, close, mix.
0-73 Sulfide 3	Fill empty test tube with 0.5 mL sulfuric acid 50 %, add 1 measuring spoon R2 and 4.0 mL sample solution, close and shake gently. Add 200 µL R3, close, mix.
0-74 Phenol-Index 5	Fill empty test tube with 1.0 mL R2 and add 4.0 mL sample solution, close and mix. Add 1 NANOFIX R2, close and mix.
0-55, 0-76, 0-80, 0-81, 0-95 ortho- and total- Phosphate LR 1-45	Proceed as described in the test instruction for test 0-55/0-76/0-80/0-81/0-95, but instead of R4 add 0.2 mL distilled water, close, mix.
0-79 ortho- and total Phosphate 50	Proceed as described in the test instruction for test 0-79, but instead of R3 add 1.0 mL sulfuric acid 20 %, close, mix.
0-82 Oxygen 12 8-22, 8-25 BSB _s	The original tests already contain a correction.
0-83, 0-88, 0-92 total Nitrogen TNb 22/60/220	Almost all colors and turbidities are destroyed under test conditions and do not interfere. For samples which are still colored or turbid after decomposition, correction values are determined as described above for test 0-64.

Test	Tube for correction (value B)
0-84 Residual hardness 1	Open Residual hardness test tube, add 5.0 mL sample solution, close and mix.
0-85 Starch 100	No correction possible.
0-89 Sulfite 10	Open test tube, add 4.0 mL sample solution and 0.2 mL distilled water, close and mix.
0-90 Sulfite 100	Fill empty test tube with 0.2 mL R2, 4.0 mL sample solution and 1.0 mL distilled water, close, mix. Analytical result = A + B
0-91 Thiocyanat 50	Fill empty test tube with 4.0 mL sample.
0-96 Zinc 4	Fill empty test tube with 4.0 mL sample, add 0.2 mL R2, close and mix.
0-97 Tin 3	Proceed as described in the test instruction for test 0-97, add 1.0 mL ethanol instead of R4.
0-98 Aluminum 07	Proceed as described in the test instruction for test 0-98, add 0.5 mL distilled water instead of R3.
8-38 Ethanol 1000	Open ethanol test tube, add 4.0 mL R1 and 0.5 mL sample solution (consider dilution), mix, add 2 drops R3, close and mix.
8-59 Methanol 15	Open methanol test tube, add 3.0 mL R1 and 1.5 mL sample solution (consider dilution), mix, add 2 drops R3, close and mix.
8-71 Peroxide 2	Fill empty test tube with 4.0 mL sample.
8-90 Sludge activity TTC 150	No correction possible.

Table 1: Correction value determination for *NANOCOLOR*[®] tubes tests

13. NTU-Check

The PF-12^{Plus} enables a warning of interfering turbidity when measuring a test kit. To activate/deactivate the NTU-Check, enter 958 as method in the start screen. Use the keys  /  to deselect the check and confirm your selection with .

Enter the desired warning limit using the keypad and confirm with .

NOTICE

When this function is activated, a nephelometric turbidity measurement is performed parallel to each measurement. MACHEREY-NAGEL recommends the use of 10 NTU as a warning limit.

14. Quality control

Each institution must implement its own QC policy.

To check the performance of the test kit and instrument used, we recommend measuring several standard solutions with known concentrations within the measuring range.

The photometric accuracy of the measurement of the PF-12^{Plus} can be checked using the *NANOCONTROL* *NANOCHECK* 2.0 color solutions (REF 925703) provided by MACHEREY-NAGEL.

15. Memory

All PF-12^{Plus} measurements are stored in the instrument memory. The measured value memory can be read out using the free MACHEREY-NAGEL data export software and the data can be transferred to the PC in this way. The PF-12^{Plus} itself allows a comprehensive memory selection and search for measurement results in the instrument memory.

The memory can be accessed by pressing the key  when the instrument is in the method call menu. The icon-based menu for the instrument memory appears. The icons have the following meanings and functions:

15.1 Display memory

Confirming the “Display memory” icon with the  key opens the entire memory content. Use the  /  keys to scroll through the individual entries.

15.2 Search in memory

By confirming the icon “Search memory” with  the menu for memory selection opens. The user can select data with a specific date, location or method number. The search criteria are reset every time the menu is selected. Select the search criteria with  /  and confirm with . Enter the appropriate criteria with the keypad. By confirming with  the search criteria are saved. Select the option “Find” with  /  and confirm with  to start the search. If no criteria are defined, the complete memory is shown and you can scroll through the data sets with  / .

15.3 Delete memory

Confirming the “Clear memory” icon by pressing the  key will clear the entire memory of the device. Before deleting the memory, a security query appears which must be confirmed again with  before it can be executed.

15.4 Output from memory


The measured value memory can be output to a terminal program via the USB interface of the PF-12^{Plus}. By selecting the “Output memory” icon and confirming with , the entire memory content is sent via the interface of the instrument.

16. Special methods

16.1 User defined methods

The PF-12^{Plus} allows the user to program up to 50 user-defined methods. These methods can be defined as linear (with one factor) or non-linear (functions up to the fourth degree + Ln-functions). For these methods all options as for the pre-programmed methods are available, e.g. location, date, dilution, data storage.


Select

Press **PROG** to enter the special method menu from the method selection menu.

Create

With **PROG** and for example 16 method P16 is selected.


When selecting an unassigned method number, only the “Create method” selection is available. After confirming with **OK** the menu-guided entry of the method characteristics is carried out. The letters and numbers for method name, filter selection, factor, decimal places, measuring range, unit and dimension can be entered via the keyboard.

In the “Factor” query, the factors F0 to F4 and FLn can be entered in sequence, which perform the calculation in the form of the following formula:

$$y = F^0 + F^1 \times A^1 + F^2 \times A^2 + F^3 \times A^3 + F^4 \times A^4 + FLn \times \ln A \quad \text{with } A = \text{Absorbance.}$$

By pressing the **ENT** key before entering the respective factor, a change of sign takes place. It is always possible to correct and re-enter an incorrect entry by using the **↑**/**↓** keys. It is also possible to abort the programming at any time with **ESC**.

When selecting an already created method, after pressing the **PROG** key and entering 2 digits, several options are available for selection which can be scrolled through with the **↑**/**↓** keys.

Run method

By confirming with **OK** the selected measuring program is executed.


Edit method

Select the option “Edit method?” by pressing **↑**/**↓** and confirm with **OK** to edit entries menu-guided.

NOTICE

Starting with firmware version >3.0.1, it is no longer possible to edit methods. To edit a method, create it as a new method and then delete the existing one.

03.06.2014 09:30
User method: P01
Delete method?
↑↓ OK

Delete

Select the option “Delete method?” by pressing / and confirm with  to delete this method after confirming an additional safety query.

To create and manage the user defined methods one can also use the free-of-charge *NANOCOLOR*[®] data export software from MACHEREY-NAGEL. This can be downloaded via www.mn-net.com. For all other information, refer to the software manual, which is available from the same link.

16.2 Sample numbering

To differentiate between samples, the PF-12^{Plus} assigns a sample number to every measurement (starting with 0001). Subsequent measurements are numbered accordingly. You have the option to change the sample number after every measurement and to assign it to the current result. To do so, press  before the next measurement. An entry option appears to enter the 4 digit sample number via the keypad. Once you have entered the sample number, it is displayed together with all other sample information at the bottom of the display. The following measurements are then automatically numbered. The sample numbering is carried out on the basis of a daily counter. With each new day the PF-12^{Plus} starts to count with the sample number 0001 and counts upwards for each further measurement.

16.3 Identification of sample location

To differentiate between sampling locations, you have the option to enter a sample location after performing a measurement (up to 12 characters). Press  key to open the entry option for the sample location. Once you have entered the sample location, it is displayed together with all other sample information at the bottom of the display. The sample location is added to the result in the memory automatically. If you have entered additional sample locations before, you can press to choose or edit the sample location from a list of up to 20 previously defined sample locations with the keys /. For editing the sample locations choose the appropriate entry with the / keys and replace the existing letters with the keypad. Deleting single digits can be accomplished by typing the key  two times.

16.4 Sample dilution

During chemical analyses samples often have to be diluted to fit within a certain measuring range. To avoid manual calculations, you can enter the dilution into the PF-12^{Plus}, which will then re-calculate the result, taking the dilution into account.

Example:

Expected result: between 80 and 200 mg/L

Measuring range photometer: 0.1 – 10.0 mg/L

Necessary dilution: (200 : 10 mg/L): at least 1:20 (1 + 19)

Recommended dilution: 1:25 to 1:50 to get within the middle of the measuring range

Dilutions are entered as 1 part sample plus X parts distilled water as to create distinct readings.

You can enter dilutions between 1 + 1 and 1 + 999.

Press **F3** before the next measurement to open the entry option for the dilution. The screen for entering the dilution appears. Once you have entered the dilution, it is displayed together with all other sample information at the bottom of the display. The result is immediately re-calculated using the entered dilution.

NOTICE

The entered dilution only impacts the current result. In the case of large dilutions, e.g. 1 + 999, it is possible that the result is displayed in a different dimension as to enable proper display in the screen.

16.5 Preprogrammed special methods

```
Method:
- Special methods
  Settings
  Basic functions
  Tube tests
↑↓
```

Besides the MACHEREY-NAGEL NANOCOLOR® and VISOCOLOR® ECO test kits, the PF-12^{Plus} is also equipped with preprogrammed special methods.

These can be chosen via the list entry "Special methods" with the keys **↑**/**↓**. Confirming with **OK** starts the chosen method.

17. Maintenance and cleaning

NOTICE

The unit must be switched off and disconnected from the power supply for all cleaning work.

WARNING

Fire hazard: Do not use flammable detergents or solvents to clean the unit.

WARNING

Incorrect measurements: A dirty cuvette slot can lead to contamination of the inserted cuvettes and incorrect results. Check the cleanliness of the cuvette slot at regular intervals and during maintenance. Clean the cuvette slot in case of contamination (see chapter 17.1 Cleaning of the device, page 65).

WARNING

Dispose of chemicals, waste and used cuvettes according to national, regional and local regulations.

NOTICE

Do not clean the device with acetone or similar products.

NOTICE

Make sure that there is no cuvette and no cuvette adapter in the cuvette slot.

NOTICE

Only use accessories and spare parts approved by the manufacturer. If parts not approved by the manufacturer are used, the operator is responsible for the conformity of the device and the warranty will be voided.

17.1 Cleaning of the device

The housing of the photometer can be cleaned with a damp cloth. Wipe off splashes immediately. The cuvette slot can be cleaned and dried with a soft, lint-free cotton cloth. In general, photometers and cuvettes must always be kept clean. The cuvettes are cleaned with a lint-free cloth.

18. Errors, sources and solutions

NOTICE

Malfunction or damage of the device due to incorrect maintenance. Maintenance of the device may only be carried out by qualified personnel.

Depending on the operating status, different messages can be displayed. The source of the error can be either an operating error or a malfunction of the device.

In case of recurring errors contact MACHEREY-NAGEL.

Error	Source	Solution
AD-overflow	Too much light at detector	Contact MACHEREY-NAGEL for calibration.
 Remove cuvette!	If, when determining the air value, the temperature deviates by $\pm 10\text{ }^{\circ}\text{C}$ from the last measurement or $dE > \pm 10\text{ mE}$, it cannot be determined with certainty whether a cuvette is in the PF-12 ^{Plus} .	Remove cuvette and confirm with 
Display stays dark	No power supply	Check if batteries / rechargeable batteries are inserted or USB cable is plugged in.
Instrument is turned on, but doesn't react	Instrument crashed	Turn instrument off, wait 10 s, turn instrument on.

Error	Source	Solution
When powered from the PC, the device always restarts	The power supply via the USB interface of the computer is not sufficient.	Use another USB interface or a USB hub with external power supply. The power supply must meet the specifications in the technical data of the PF-12 ^{Plus} .
RTC error	The battery of the real time clock is defective or empty.	Set the date and time again and then restart the device. If the date and time do not show the set values, the battery of the real-time clock must be replaced. Contact MACHEREY-NAGEL. Use only the battery recommended by MACHEREY-NAGEL for replacement.
Battery/Accu empty!	The batteries or the inserted battery pack are empty. The power supply is not sufficient to perform a measurement.	Insert new batteries or charge the battery pack.

19. Service

19.1 Spare parts, accessories and consumables

NOTICE

Only use accessories and spare parts approved by the manufacturer. If parts not approved by the manufacturer are used, the operator is responsible for the conformity of the device and the warranty will be voided.

WARNING

Risk of injury: Injury to persons and damage or malfunction of the photometer or equipment when using parts not approved by the manufacturer. The replacement parts in this section are approved by the manufacturer.

Description	REF
Rugged case with foam inlay for PF-12 ^{Plus} , accessories and VISOCOLOR® ECO test kits	914351.1
Calibration cuvette	916908
USB-cable A/B for data export and power supply	919687
USB mains adapter for PF-12 ^{Plus} with Euro plug (type C)	919220
Optional battery pack	919201
Test tubes (20 pieces; 16 mm OD)	91680
Operating instruction	919252
VISOCOLOR® ECO Test instructions for photometric determinations	931504
NANOCONTROL NANOCHECK 2.0	925703
NANOCONTROL NANOTURB	925702

19.2 Disposal

NOTICE

Disposal via public disposal systems is not permitted. Contact your local MACHEREY-NAGEL representative.

Disposal according to EU Directive 2012/19/EU. In accordance with the EU Directive 2012/19/EU MACHEREY-NAGEL takes back the old device and disposes it free of charge.

19.3 Warranty, liability and complaints

The warranty for this device is 24 months from date of purchase. The original invoice serves as proof and must be presented when making a claim. In case of improper handling and/or maintenance of the device, the warranty expires. It does not cover defects that are due to a power supply other than the external power supply included in the delivery.

The warranty is limited to the repair of defective parts or - at MACHEREY-NAGEL's discretion - the delivery of a faultless replacement unit. The warranty period of 24 months shall not be affected by any warranty claim. There is no right of withdrawal. Further claims are excluded. These include in particular all claims for damages resulting from consequential or indirect damages. In addition, our general terms and conditions of sale and delivery apply in the currently valid version as printed on all price lists.

19.4 Contact

MACHEREY-NAGEL GmbH & Co. KG

If you still have questions or need technical assistance after reading the operating instruction, contact:

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MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany