MACHEREY-NAGEL

URYXXON® 500

User Manual

[en](#) [de](#) [es](#) [it](#) [ro](#)

Rapid Tests

MACHEREY-NAGEL

www.mn-net.com

Contents

Instruction for use	3
Gebrauchsanweisung	35
Instrucciones de uso	69
Istruzioni per l'uso	101
Instrucțiuni de utilizare	133

MACHEREY-NAGEL

Reflection photometer URYXXON® 500

Instruction for use

Rapid Tests

MACHEREY-NAGEL

www.mn-net.com

1. Introduction	6
1.1 Intended Purpose	6
1.2 Intended user	6
1.3 Indications for use	6
1.4 Patient population	6
1.5 System description	6
1.5.1 Measuring principle	6
1.5.2 Functional principle URYXXON® 500	6
1.6 Product labeling	6
1.7 Technical data	7
1.8 Symbols	8
2. Safety precautions	9
2.1 Power supply	9
2.2 Biohazard	9
2.3 Protective clothing	9
2.4 Inappropriate application	9
2.5 Housing damage	9
2.6 Cable break	9
2.7 Transport	9
3. Instrument set-up	10
3.1 Package content	10
3.2 Set-up location	10
4. Initial operation	11
4.1 Insert the transport rack	11
4.2 Power supply	11
4.3 Connections/interface	12
4.4 Insert printer paper	12
4.5 Switching on the instrument	12
4.6 Calibration	13
4.7 Switching off the device	13
5. Quick start	13
6. Operation and user guidance	14
6.1 Touch screen	14
6.2 Operating the user interface	14
6.2.1 Activity-buttons	14
6.2.2 Option-fields	14
6.2.3 Edit-buttons	14
6.2.4 Input-menu	14
6.3 Stand-by mode	14
7. Main menu	15
8. Settings	15
8.1 Strip	15
8.1.1 Units	15
8.1.2 Customized order	15
8.2 Interface	16
8.2.1 HW settings	16
8.2.1.1 Baud rate	16
8.2.1.2 Bit length	16
8.2.1.3 Parity	16
8.3 Printer	16
8.4 Sediment options	17
8.4.1 Enter sediment	17
8.4.1.1 Edit sediment parameters	17
8.4.2 Filter	17
8.4.2.1 Edit filter	17
8.4.3 Print after	18
8.4.4 Print sediment settings	18
8.5 Language	18
8.6 Date/Time	18
8.7 User mode	18
8.7.1 Single user	19
8.7.2 User / Admin	19
8.7.3 Multi user	19
8.7.3.1 Edit users	19
8.8 Customization	20
8.8.1 Clarity	20
8.8.2 Color	20
8.8.3 Welcome line	20
8.8.4 Printer header	20

8.9 Sound	20
8.10 Print settings	20
8.11 Autostart	20
9. Enter sediment	21
9.1 Microscope	21
9.2 Print	22
9.3 Send to a PC	22
10. Memory	22
10.1 Show all	22
10.1.1 Print	22
10.1.2 Send to a PC	22
10.1.3 Microscope	22
10.2 Search	23
10.2.1 Date	23
10.2.2 Filter	23
10.2.3 Action	24
11. Check mode	24
11.1 Internal check	24
11.2 Control measurement	24
11.3 Control strip	24
11.4 Stored data	24
12. General information for measurements	25
12.1 Urine sample	25
12.2 Test strip preparation	25
13. Measurement	25
13.1 Sample and patient information	25
13.2 Single measurement	26
13.3 Successive measurements	26
13.4 Measurements with worklist	26
13.4.1 Create a worklist	26
13.4.2 Edit worklist	27
13.4.3 Measuring with the worklist	27
13.4.4 Options	27
13.4.4.1 Load list	27
13.4.4.2 Print list	27
13.4.4.3 Delete list	27
13.4.5 Emergency measurement	27
14. Quality control	28
14.1 Check mode	28
14.1.1 Internal check	28
14.1.2 Control measurement	28
14.1.3 Control strip measurement	28
14.1.4 Stored data	28
15. Connection of external devices	29
15.1 Computer	29
15.2 Barcode reader / Keyboard	29
15.3 Connection to a LIS	29
16. Cleaning and maintenance	29
16.1 Daily cleaning	29
16.2 Weekly cleaning	29
17. Service	30
17.1 Spare parts and accessories	30
17.2 Consumables	30
17.3 Error messages	30
17.4 Disposal	31
17.5 Analytical Performance Characteristics	31
17.5.1 Table of results	31
17.5.2 Precision	31
17.5.3 Assay reportable range	32
17.5.4 Cut-off concentration	32
17.6 Changing the cut-off points	32
17.7 Comparative study	32
17.8 Interference study	33
17.9 Warranty	33
17.10 General usage comments	33
17.11 Literature	33
17.12 CE Certification	34
17.13 Contact	34
17.14 Version History	34

1. Introduction

This instruction for use is provided for the URYXXON® 500 models REF 930080 and 930080.XXX.

The URYXXON® 500 is a reflection photometer for the semi-quantitative evaluation of URYXXON® Stick 10 urine test strips. Measured values may be displayed, printed and transferred to a computer. The URYXXON® 500 is designed for *in vitro* diagnostic use (IVD) and should be used by healthcare professionals, only.

Equipped with a touch screen and intuitive user guidance, the URYXXON® 500 assures a quick, efficient and convenient work flow. Therefore, the URYXXON® 500 is the ideal instrument for modern urine analysis.

1.1 Intended Purpose

The URYXXON® 500 is a reflection photometer for the evaluation of urine test strips. The URYXXON® 500 is intended for use with compatible test strips from MACHEREY-NAGEL. The *in vitro* analysis system is an aid to diagnoses in the areas of kidney function, urinary tract infections, metabolic diseases (e.g. diabetes mellitus) and liver function. The URYXXON® 500 is intended for use by professional users. The URYXXON® 500 enables integration into laboratory information and management systems (LIMS) and supports bidirectional data transfer (from / to the device).

1.2 Intended user

The URYXXON® 500 is intended for use by trained professionals with experience in human urine sampling and is not suitable for self-testing.

1.3 Indications for use

For indication regarding the indications for use refer to the instructions for use of the used urine test strips.


1.4 Patient population

There are no restrictions regarding the patient population with MACHEREY-NAGEL urine test strips.

1.5 System description

1.5.1 Measuring principle

The test strip is placed on a transport tray rack and is moved below a movable measuring head by the transport rack.

The strip is illuminated with four LED and a detector registers the intensity of light reflected by the test strip. Using an internal calibration, the intensity of the reflected light is correlated to the concentration of the analyte. By comparing the remission values with a range table (so-called cut-off limits) the measurement result is determined. Whenever samples are strongly alkaline, a density correction is automatically conducted.

1.5.2 Functional principle URYXXON® 500

A measurement is started by placing a strip on the transport rack. If the Autostart-feature is turned off, the measurement is started by pressing the start panel in the display. The result is shown, stored in the memory, printed out and released via the interfaces after the measurement has been completed. After ten minutes the instrument will go into stand-by. Touching the screen will reactivate the instrument. All user inputs are performed via the touch-screen (see chapter 6. Operation and user guidance, page 14).

1.6 Product labeling

Exemplary name plate on the device



Exemplary sticker on the box

1.7 Technical data

General info	Size	24 cm x 28 cm x 15 cm (9.5" x 11" x 5.9")
	Weight	3.9 kg (10.5 lb)
	Measurements/h	350
	Power supply	External power pack 100–250 V (Type SB-076 A0F-11), Input 47–63 Hz, Output 7.5 V 6 A or GST40A07 / Globtek Inc GTM96600-4507.5T3.
	Transport condition	-10 °C–45 °C, max. 80 % humidity (non-condensing)
	Operation condition	10 °C–40 °C, max. 80 % humidity (non-condensing), max. 2000 m
	Storage condition	10 °C–40 °C, max. 80 % humidity (non-condensing)
Measurement	Principal	Reflection photometer
	Light source	4 LEDs, wave length: 470, 550, 610, 910 nm
	Detector	Photodiode
User Interface	Touch screen	W x H = 10 x 7.5 cm Color display Resolution 320 x 240 Pixel
	User menu	Settings and user guidance in 8 languages
	Printer	Thermo transfer
Interface	1 x USB	
	2 x RS232	9600 Baud, 8, N, 1 uni-/bidirectional mode (base setting)
	1 x PS/2	Keyboard / barcode scanner
Memory	Capacity	500 measurements
		20 control measurements
Accessories (optional)	Keyboard	Data input
	Barcode scanner	Data input
	Serial cable	PC connection
Protection class	Equipment:	Protection Class III
	Power Supply:	Protection Class I
Overvoltage category	Equipment:	Category II
	Power supply:	Category I
Degree of pollution		2

1.8 Symbols

The following table explains the symbols that are used on the device, the power supply, the package as well as the instruction for use. For the explanation of the symbols used on the strips and control solutions, refer to the respective instructions for use.

	Caution
	Non-sterile
	Keep away from sunlight
	Biological risk
	Keep dry
	Temperature limits
	Do not use if package is damaged and consult instructions for use
	Consult instructions for use or consult electronic instructions for use
	Do not dispose of the device in common household waste.
	Manufacturer
	Catalogue number
	Batch code
	Serial Number
	Date of Manufacture
	Use-by date
	<i>In vitro</i> diagnostic medical device
	CE Mark
	Direct current
	Do not reuse
	Device for near patient testing
	Device not for self testing
	Contains <n> tests

2. Safety precautions

BE SURE TO READ THE FOLLOWING SAFETY INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE.

Failure to follow these instructions may result in serious injury to the operator, malfunction or damage to the equipment.

Keep this operating instruction in a safe place for future reference. Follow the safety notes and instructions in the operating instruction and observe the stickers and notices on the device.

Do not work on internal parts of the unit. Non-compliance will invalidate any warranty claim.

Use of the hazard warnings:

DANGER

Indicates an imminent or potentially hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Indicates an imminent or potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may result in malfunction or damage of the instrument. Information that clarifies facts in the text and that requires special attention.

2.1 Power supply

To ensure safe functioning of the instrument, only use the included power pack!

WARNING

Risk of electric shock. Only use the included power supply Type SB-076 A0F-11. Use of other power supplies can damage the device and harm the user.

WARNING

Only use the device with protective grounding.

2.2 Biohazard

Using the instrument makes contact to urine inevitable. Therefore, use protective gloves whenever measuring with the instrument and be cautious when handling urine samples.

CAUTION

Biological Hazard: Urine residues

Urine and used test strips pose a risk of infection. Always wear protective gloves when handling and disposal. Dispose of used test strips according to the regulations for handling potentially infectious material.

2.3 Protective clothing

We recommend to always wear gloves and a lab coat as protective clothing.

2.4 Inappropriate application

The URYXXON® 500 is designed for the professional evaluation of URYXXON® Stick 10 urine test strips. Inappropriate use, or use by untrained staff terminates any warranty claim for the user with respect to MACHEREY-NAGEL.

2.5 Housing damage

If the housing is damaged, the instrument has to be sent in for repair. Correct functioning and measuring results are not guaranteed any more.

2.6 Cable break

In the case of cable break, immediately stop using the instrument and replace the power pack.

WARNING

Risk of electric shock. Cable break can lead to exposure to conduction cables. Stop using the device immediately and remove the power supply from the outlet if possible.

2.7 Transport

To transport the instrument, the cardboard box from the initial shipment including the protective foamed plastic is ideal. Therefore, we recommend to keep the box. In case you do not have the original box anymore, use a large cardboard box and plenty of soft cushioning material, so that the instrument cannot shift during transport.

Make sure to remove the transport rack before you ship the instrument.

3. Instrument set-up

3.1 Package content

Open the box carefully with a sharp instrument. Make sure not to damage the content. Take out the instrument and all other parts carefully. Check the package, instrument and all accessories for visible damage. If a part is damaged, contact your distributor or the technical support (see chapter 17.13 Contact, page 34).

The following picture shows the complete package content. Check the shipment for completeness. If a part is missing, contact your distributor or the technical support as well.

NOTICE

Keep the original box and the packaging material from the initial shipment to optimally protect the instrument in case of a return shipment.

3.2 Set-up location

Only use the instrument in a suitable location. It should be positioned on a dry, clean, level and horizontal surface. For detailed operation conditions please refer to the table in chapter 1.7 Technical data, page 7.

Make sure to protect the instrument from drastic changes in temperature. Avoid positioning the instrument close to windows, hot plates, ovens etc. To ensure correct functioning and reliable results, the instrument should not be exposed to direct sunlight or other bright, focused light sources.

NOTICE

When operating outside these limits, the correct functioning of the device is not guaranteed.

- ① URYXXON® 500
- ② Transport rack
- ③ Thermal printer paper rolls
- ④ Instruction for use
- ⑤ Serial cable RS232
- ⑥ Power cable
- ⑦ Power adapter

4. Initial operation


This chapter describes the first steps to a successful start-up. The correct execution of the steps described here must always be ensured by the user.

4.1 Insert the transport rack


NOTICE

Only start the instrument once the transport rack has been inserted and all power connections are complete.

Before turning on the instrument for the first time, the transport rack has to be inserted. Remove the black transport tray (8) from the instrument. To do so, pull on the black handle on the transport tray's right side.



Before you insert the transport rack, make sure there is no packaging material blocking the mechanics. Take the transport rack (2) and align its feed holes with the corresponding rubber pins on the transport panel.

The mounting direction of the transport rack is determined by the asymmetrical feed holes / pins.

Push the transport tray gently back into the instrument.

4.2 Power supply

CAUTION

Multiple hazards: Only qualified personnel should perform the instructions described in this step.

WARNING

Risk of electric shock. Only use the included power supply Type SB-076 A0F-11. Use of other power supplies can damage the device and harm the user.

CAUTION

Risk of injury due to the danger of tipping over due to slipping. Place the device on a flat surface. Do not stack the device.

WARNING

Take care not to overload the mains socket. There is a risk of overload and fire. Make sure that the power cord is not damaged. Check the suitability of the power source used for the equipment. Only connect the power supply to grounded outlets.

WARNING

Risk of injury: Familiarize yourself with the device before working with it and read this document carefully. Do not use the device unless you have received instruction in its use.

NOTICE

Do not switch on the device until all data and power cables have been fully assembled and connected.


The instrument requires a power supply of 7.5 V DC, which is supplied by the power pack ((Type SB-076 A0F-11) 100–250 V, 47 · 63 Hz) or GST40A07 / Globtek Inc GTM96600-4507.5T3.

Connect the power cable (6) to the transformer (7), and then connect the transformer cable to the DC IN port (a) at the back of the instrument. Connect the cable to a grounded outlet.

4.3 Connections / interface

Apart from the power connection (1), the instrument has 4 additional interfaces (1 x USB (2), 1 x PS/2 (3), 2 x RS232 (4)) to connect different devices or a laboratory information system (LIS) (see above).

4.4 Insert printer paper

NOTICE

Only use special thermal printer paper!

To open the printer paper cover, press the little button above the touch screen. Move the printer paper cover backwards.

Unroll appr. 10 cm (4") of printer paper and place the printer paper into the chamber, so that the paper unrolls downwards and towards the touch screen.

Close the printer paper cover with an audible click, so that the printer paper comes out from the chamber.

4.5 Switching on the instrument

You can find the On / Off switch (5) on the back of the instrument. Turn on the instrument by pressing the switch (see pictures below). Right after switching on, the instrument performs a self-test. The touch screen lights up and the start screen appears. The little light above the transport tray turns green, thus signalling operational readiness. If a serious error is detected during the self-test, no measurements can be performed. In this case, contact the MACHEREY-NAGEL technical service.

On

Off

4.6 Calibration

The instrument performs an automatic, internal calibration before each measurement.

4.7 Switching off the device

Always switch off the device after use by pressing the on / off switch **5**. If you do not use the device for a longer period of time, disconnect the device from the power outlet.

5. Quick start

NOTICE

The instrument features multiple different settings, options and functions. We recommend to read the instruction for use carefully and to make individual adjustments to ensure the instrument functions at its best.


NOTICE

Only use consumables within their expiration date.

Directly after turning on the instrument, it is ready and can perform measurements. In order to start a measurement, dip a test strip with all test pads into a urine sample. Draw the long edge of the strip along the rim of the sample container and press the same edge briefly to an absorbent towel.





Dip test strip

Remove excess urine

Place the urine test strip straight on the transport table under the strip detector (a) as shown in the figure. The measurement starts automatically as soon as the detector has recognized the strip. To do this, push the strip up to the front edge of the transport table.

To enter additional information for the sample, choose the activity-button before placing the strip on the tray. A menu appears in which you can enter the patient ID, sequence no., as well as pre-defined values for clarity and color.

The result is automatically printed and saved after the measurement.

6. Operation and user guidance

The instrument is designed for easy and intuitive use. Nevertheless, read the following paragraphs carefully to familiarize yourself with the instrument's functions and user guidance.

6.1 Touch screen

NOTICE

Do not touch the touch screen with pointed or hard items.

All operations are carried out with a touch-sensitive screen. An external keyboard can be connected via the PS/2 interface (see chapter 15. Connection of external devices, page 29).

6.2 Operating the user interface

Note some basic rules for operating the instrument's user interface:

- All menu items that you can access are framed in blue. That means, only buttons with a blue frame can be opened, edited or activated by touch. Depending on the button, a new menu opens up, or the respective activity is carried out.
- Entries in large, blue font can be accessed directly. They can be edited via the input-menu .
- Entries in small, black font cannot be changed.

The options to access functions can be sub-divided into four categories:

6.2.1 Activity-buttons

Blue-framed activity-buttons include a symbol/inscription and trigger a certain activity. The following symbols are used:

- escape/cancel
- confirm
-  – print information
-  – enter additional sediment information
-  – send data to a connected PC/LIS
-  – scroll in a list up/down
-  – browse from one screen/entry to another

6.2.2 Option-fields

Option-fields consist of blue-framed choice-buttons. Choice-buttons are circles. If a choice-button is active, it is filled with blue. Within an option field, only one choice-button can be active.

Touch a choice-button to activate it. Press to confirm and save your changes. If you leave menus that include option-fields by pressing , none of your changes will be adopted.

6.2.3 Edit-buttons


Edit-buttons have blue frames and enable you to enter information, e.g. patient ID or sequence no. Sometimes, edit-buttons already include information. In this case, the font is large and blue, thus indicating the possibility to edit the information. To edit/enter information, press the button. An input-menu with alphanumeric input-buttons appears. By touching the numbers/letters in the input-menu, you can enter information into the edit-button.

6.2.4 Input-menu

The input-menu appears whenever you want to make an entry. Above the input-buttons you will find a window which shows your new entry.

You can enter upper case letters, lower case letters or numbers. To switch between the options, press the button that shows the desired form of entry (abc/ABC/123). Initially, the instrument is opting for capital letters. Once you make a change, the instrument remembers your last choice. If you open the input-menu again, the former choice is still active.

For new entries touch a letter/number. To enter a different letter/number than the first one on the button, press the button repeatedly.

You can delete an entry by pressing . Press to confirm and save your changes. If you leave the input-menu by pressing , none of your changes will be adopted.

6.3 Stand-by mode


NOTICE

If the device is activated from the stand-by mode, it will always go to the start screen. Any changes in menus or changes in settings that have not been saved before the instrument went into stand-by mode will be lost.

After 10 minutes without any operation the device goes into stand-by mode and the screen saver appears. By pressing the touch screen you can activate the system again and get directly to the start screen.


7. Main menu

Once the instrument has been switched on and performed the self-test, the start screen appears.

The start screen is the starting point for measurements and all other settings. You can access the main *Menu* from the start screen, which allows you to change the instrument settings, manage saved data, add sediment information to your results and perform quality control measurements.


To enter the main *Menu*, press  in the upper left corner of the start screen. The main *Menu* with multiple options appears.

Press the desired option to display additional choices.


8. Settings


In the *Settings* menu, you can adapt the instrument to your requirements. Touch the activity-buttons to make changes.

8.1 Strip

In the *Strip* menu, you can change the units and the order in which the parameters are displayed.

Press  to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.1.1 Units


The option-field *Units* allows you to choose the units for the different parameters. Three different units (Conv, SI, ARB) and two mixed entries (Conv + ARB, SI + ARB) are available. Choose your desired unit(s) by pressing the respective choice-button.

You can find additional information concerning the different units in the table of results (see chapter 17.5 *Analytical Performance Characteristics*, page 31).

8.1.2 Customized order



The option-field *Customized order* enables you to change the order of the parameters when results are printed and displayed. By default, the order is identical to the order of the test pads on the URYXXON® Stick 10 and the choice-button *off* is enabled.

If you want to change the order, press the choice-button *on*. The adjacent activity-field  is activated. Once you can press it, a new menu appears, displaying the existing parameters. Alongside are scroll-buttons   to reach all parameters in the list. The current order of the parameters is displayed on the right half of the screen.

To change the orders of parameters, press the activity-button in the lower part of the screen. The current order disappears and only the number 1: remains.

By pressing the parameters on the right side of the screen you can define a new order for all parameters. You can pick the parameters either by pressing the parameter directly, or by scrolling through the list using the scroll-buttons and pressing the activity-button . Parameters that have already been selected are not available anymore. They are listed on the right side of the screen with a preceding position number.

If you are not satisfied with your selection; press . The selection on the right hand side is discarded and you can enter your personal order again.


Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.2 Interface

NOTICE

Changes in the interface settings directly influence data transfer. This is especially important when the instrument is connected to a LIS.

In the *Interface* menu you can change the interface settings.


You can decide whether the interface communication should be uni- or bidirectional. To do so, press the respective choice-button in the option-field. Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

NOTICE

If bidirectional communication is enabled, the receiver has to be configured in the same way. If this is not the case the device will report a communication error.

8.2.1 HW settings

In the *HW settings* menu, you can carry out additional hardware adjustments. Once you press the activity-button , a new menu with multiple option-fields appears. Within the menu, the bold values show the default settings.

Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.2.1.1 Baud rate

You have five different choices for the baud rate (1200, 2400, 4800, 9600, 19200). Press the respective choice-button to change the baud rate.

8.2.1.2 Bit length


You have two different choices for the bit length (7 Bit / 8 Bit). Press the respective choice-button to change the bit length.

8.2.1.3 Parity

You have three choices for the parity (none, even, odd). Press the respective choice-button to change the parity.


8.3 Printer

The *Printer* menu allows you to turn the printer on and off. Press the respective choice-button in the option-field. Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.4 Sediment options

In the *Sediment options* menu, you have the option to enter additional sediment parameters for every result. This function enables you to manage both urine chemistry results and microscopic findings for individual samples.

NOTICE

To add sediment information to a sample, go to the main menu and choose "Enter sediment" (see chapter 9. *Enter sediment*, page 21).


8.4.1 Enter sediment

The option-field **Enter sediment** allows you to turn the enter sediment function on or off. To do so, press the respective choice-button. To edit the available sediment parameter, press the activity-button **Edit sediment parameters**. You can only edit the sediment parameters, if the choice-button **On** is activated. Press **✓** to confirm and save your changes. If you leave the menu by pressing **✕**, none of your changes will be adopted.

8.4.1.1 Edit sediment parameters


The instrument can store up to 30 different sediment parameters. To label the sediment parameters, you can use up to 10 characters each. For every parameter, you can define 5 categories or increments. To label these categories you can use up to 8 characters each.

The *Edit sediment parameters* menu consists of 6 edit-buttons, in which you can enter individual information. In the upper left corner of the screen you can see the number of the parameter you are currently editing.

The first edit-button in the center of the screen is meant to assign a name to the sediment parameter. To enter the parameter's name, press the empty edit-button. The input-menu **abc** appears and you can enter the parameter. Confirm your entry by pressing **✓**. The new parameter has been adopted.

You can enter up to 5 increments (e.g. 0–5, 6–15, 16–30, ...) or categories (e.g. small, medium, large, ...) for every parameter into the edit-buttons on the right side of the screen. Press the edit-buttons and the input-menu **abc** appears in which you can enter your increments or categories. Confirm your entry by pressing **✓**. The new category has been adopted.

To define more parameters, use the browse-arrows **◀▶** in the bottom of the screen. The sediment parameters you have created in the menu can now be used to enter sediment information in the *Enter sediment* menu (see chapter 9. *Enter sediment*, page 21). Press **✓** to confirm and save your changes. If you leave the menu by pressing **✕**, none of your changes will be adopted.


8.4.2 Filter

If the filter is activated, you can enter specific filter criteria to select results, for which sediment information have to be entered. The results are selected to enter sediment findings in the *Enter sediment* menu (see chapter 9. *Enter sediment*, page 21).


By default, the choice-button **off** is activated in the option field and filter criteria are not installed. To activate the filter and to define filter criteria, press the choice-button **on**. The adjacent activity-button **Edit filter** is activated.

8.4.2.1 Edit filter


To get to the *Edit filter* menu, press **Edit filter**. In the menu, you can define the results that require sediment findings. On the left side of the screen, you find a list with all parameters found on the test strips. You can determine a filter criterion for every parameter. Such a criterion is always one possible measurement value. Thus, you can define for which results sediment findings have to be added. (e.g. Blood – 10 / Ket – 100 / ...)

You can choose the parameters either by using the scroll arrows **↑↓** and the activity-button **Select**, or by directly touching the parameter. Once you have chosen a parameter, a new list appears on the right side of the screen.

The list on the right shows the different gradations for the respective parameter. Choose your desired value by pressing the respective button. The value is displayed on the right side of the screen.

If you want to delete a filter criterion you have already entered, select the criterion in the list on the left side of the screen and press the activity-button . The parameter is removed from the filter criteria on the right side. If you want to enter a new criterion for the parameter, select the parameter again on the left side of the screen and make a new selection.

After you have defined your filter criteria, all results that are equal to, or exceed one of the criteria will be selected for sediment findings. The activity-button in the main menu gets you to all the results that match your filter criteria.

Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

NOTICE

The "Enter sediment" function and the sediment filter only list results that were measured after the function has been set / enabled. If you activate or change the sediment function after measurements have been performed, those measurements will not be captured by the new filter criteria. Therefore, the previous measurements cannot be accessed via the "Enter sediment" menu. If the filter is not activated, all results will be selected for sediment findings. However, the memory function enables you to enter sediment findings for all results as well (see chapter 10. Memory, page 22).

8.4.3 Print after

The option-field **Print after** enables you to determine at which point results should be printed. By default, the choice-button **Measurement** is activated. This means all results will be printed immediately after the measurement has been performed, irrespective of the need for additional sediment information.


If you want results to be printed after possible sediment findings have been entered, press the choice-button **Sediment**. A result that is listed for additional sediment findings is now only printed once the sediment findings have been entered. Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.4.4 Print sediment settings

The printer icon  in the middle of the screen allows you to print all sediment settings. The basic settings for sediment entry, sediment filter, as well as all defined sediment parameters will be printed.

8.5 Language

You can choose from 8 different languages in the language menu. Press the choice-button for the respective language.


Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.6 Date/Time

You can adjust the date and time format and change the date and time in the *Date/Time* menu.

Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

Time

You can set the hours and minutes separately. Press the edit-buttons to change the hours and minutes. The input-menu appears so that you can make changes.

Time format

The time format can be set to 12 h or 24 h. Press the respective choice-button. Depending on the format, you can adjust the activity-button as well.

Date

The date is arranged by year, month and day. Press the edit-buttons to change the date. The input-menu appears so that you can make changes.

Date format

You have three different date format options: (yyyy-mm-dd, mm-dd-yyyy, dd-mm-yyyy).


Press the respective choice-button and the date will be displayed accordingly.

NOTICE

Verify the date and time after starting the instrument for the first time.

8.7 User mode

The *User mode* menu allows you to adjust the settings for access rights and user administration. You have three different modes to choose from, enabling you to define passwords and create a group for multiple users.


8.7.1 Single user

The choice-button **single user** is set by default. In the single user mode, every user can use the instrument without restrictions to access rights. If the choice-button **single user** is activated, both the activity-buttons and are inactive.

8.7.2 User/Admin


The choice-button **User/Admin** allows you to set up one administrator for the instrument. If the Admin mode is active, only the administrator can access and change the instrument settings. All other users are excluded from the settings. The Admin mode is protected by a password. To set-up the administrator, press the choice-button **User/Admin**.

Then press the activity-button . The input-menu appears and you can enter a password. Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.


As soon as you have chosen the "User / Admin" Mode and defined a password, the password is requested whenever you turn on the instrument. To sign in as an administrator, choose the activity-button after the device has been turned on. You have 10 seconds to press the button, which is indicated by a timer. Once you press the activity-button the input-menu appears and you can enter the password. The instrument starts and you can access and change the settings normally.

If a different user than the administrator wants to use the instrument, entering the password is not necessary. Press , or wait for 10 seconds until the timer has run out. The instrument then starts automatically. You can use the device for normal measurements, but cannot access the instrument settings.


8.7.3 Multi user

The *multi user* mode allows you to set-up up to 10 users with individual passwords. Only the first user in the user group has administrator rights. All other users cannot access the instrument settings. If the *multi user* mode is set, all print-outs show the user who performed the measurement. To set-up multiple users, press the choice-button **Multi-user** and press the activity-button .

8.7.3.1 Edit users


Within the *Edit user* menu, the middle of the screen shows an empty list with places for new users. The first place is underlined and designated for the administrator. The asterisk on the edge of the place indicates the current user.

To initially set-up a new user, first create the user with administration rights. Go to the first entry with the scroll-arrows and press the activity-button . The input-menu appears and you can enter a name for the administrator. After confirming the name with , you have to enter a password. Choose a password and confirm your entry with . The *Edit users* menu appears again and shows the name of the administrator in the first place. To add more users, choose an empty place and repeat the steps above.

To delete a user from the list, use the scroll-arrows to move to the respective place. Then press the activity-button . After having set-up the users, leave the menu via . Confirm the activation of the user mode *Multi user* with .


As soon as you have created different users and the instrument is set in the "Multi user" mode, you have to choose a user every time the instrument is turned on. Once you turn on the instrument, a login screen appears. The screen contains a list with all users that have been created. Use the scroll-arrows to choose the desired user and confirm your choice with . The input-menu appears to enter the password. Confirm your entry with and the instrument starts as usual. If you leave the input-menu via , the login screen appears again.

8.8 Customization

The *Customization* menu allows you to assign 4 different values for **Clarity** and **Color**, as visual sample information. Additionally, you can define a welcome line and a printer header.

Once you open the menu, the first screen allows you to enter values for clarity and color. The browse-arrows allow you to switch to a second screen in which you can enter welcome line and printer header.

Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.8.1 Clarity

Clarity is a visual urine parameter, which you can assign to any sample as optional patient information (see *chapter 13.1 Sample and patient information, page 25*). There are four different options. By default, they are defined as “Clear”, “Turbid”, “Cloudy” and “Dark”. To create new names, press the edit-buttons for clarity. The input-menu appears and you can enter a new name (10 characters max.). Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.


8.8.2 Color

Color is a visual urine parameter, which you can assign to any sample as optional patient information (see *chapter 13.1 Sample and patient information, page 25*). There are four different options. By default, they are defined as “Yellow”, “Red”, “Brown” and “Bright”. To create new names, press the edit-buttons for color. The input-menu appears and you can enter a new name (10 characters max.). Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.8.3 Welcome line

As a welcome line, you can enter an individual message that always appears in the initial print-out when the instrument is switched on.

By default, there is no set welcome line. To create your own welcome line, press the edit-button. The input-menu appears and you can enter two lines for the new welcome line (24 characters in 2 lines max.). Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.


8.8.4 Printer header

The printer header is the heading for every print out. By default, there is no set printer header. Press the edit-button to create your own header. The input-menu appears and you can enter two lines for the new printer header (24 characters in 2 lines max.). Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.9 Sound

The *Sound* menu allows you to enable or disable the operating sounds. Press the respective choice-button.


Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

8.10 Print settings


You can print an overview over all settings. Only the sediment settings can be printed separately (see *chapter 8.4 Sediment options, page 17*). To start printing, press the activity-button .

8.11 Autostart

The *Autostart* menu allows you to enable/disable the autostart function. Press the respective choice-button.

Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted. If the autostart function is active, the device automatically recognizes a strip that is placed on the transport tray and starts the measurement. If the autostart function is inactive, you will find the new activity-button in the upper right corner of the start screen. In this case, the instrument only starts a measurement if you press the activity-button .


9. Enter sediment


The *Enter sediment* menu only shows urine chemistry results for which an additional entry of sediment findings is necessary. You can determine freely, which results need to receive additional sediment findings as well as the actual sediment parameters (see *chapter 8.4 Sediment options, page 17*).

NOTICE

To add sediment parameters, you have to activate the Enter sediment function and define sediment parameters first.

Press the activity-field . A new screen appears which shows the results of samples that still need additional sediment findings. You can set the filter for these results under sediment settings (see *chapter 8.4 Sediment options, page 17*). Multiple options for the samples are available.


NOTICE

The sediment filter only filters results that were measured after the filter has been set and activated. Already existing results will not be filtered afterwards



9.1 Microscope

To enter sediment findings for the displayed result, press the activity button . A new screen appears, in which you can choose from a list of sediment parameters on the left side of the screen (see *chapter 8.4 Sediment options, page 17*).

Choose your desired sediment parameter with the scroll-arrows and confirm your selection with the activity-button . Alternatively, you can directly press the sediment parameter.

As soon as you have chosen one sediment parameter, a list with the sediment parameter's gradation appears. Press the respective gradation, which is then inserted in the upper right of the screen. You can enter more sediment information for the result in the same way.

To delete sediment parameters you have already chosen, press the respective parameter in the list on the left side of the screen. Press the activity-button **Clear**, which will remove the parameter from the result's sediment findings.

Press **✓** to confirm and save your changes. The result is processed depending on the settings (printed, sent) and removed from the enter sediment list. If you leave the menu by pressing **X**, none of your changes will be adopted.

To get to the next result, use the browse-arrows **◀▶** at the bottom right of the screen.

9.2 Print

Press the printer icon to print the result in the display. Thus, you can print results that already have been printed for a second time.


9.3 Send to a PC

The current result is sent to a connected computer (see chapter 15. Connection of external devices, page 29).

10. Memory


The instrument saves results of up to 500 measurements. By default after each measurement the result is saved automatically. After 500 measurements, new data will overwrite the eldest saved dataset.

The *Memory* menu allows you to manage the saved data. To get to the memory, press the activity-button **Memory**. Three additional activity-buttons, **Show all**, **Search** and **Delete** appear.

10.1 Show all

Press the activity-button **Show all** to show the whole memory. Once you press the button, you get to the result screen, which automatically displays the latest result. You see the sequence number and the patient ID, date and time of the measurement, as well as the measuring values. Positive values are always in red font and marked with an asterisk. If specified, you also see the additional parameters color and clarity, as well as any added sediment parameters.

To get to other results, you can browse through the memory using the browse-arrows **◀▶** in the bottom left. The following options are available for the results:

10.1.1 Print

The current result is printed.

10.1.2 Send to a PC

The current result is sent to a connected computer (see chapter 15. Connection of external devices, page 29).

10.1.3 Microscope

You can add supplementary sediment findings for any saved value (see chapter 9. Enter sediment, page 21).

NOTICE

If you change any values, you might have to re-send them to a connected LIS.

10.2 Search

The activity-button allows you to search systematically for specific results.

Search memory	
Date	All
Filter	All
Action	Display
<input type="button" value="X"/> <input type="button" value="✓"/>	

You can edit date, filter and action to create targeted search queries. Choose the desired options via the buttons. Press to confirm and save your changes. If you leave the menu by pressing , none of your changes/selections will be adopted.

10.2.1 Date

In order to choose results from a specific time frame, press the activity-button . Three new activity-buttons appear , and .

Search memory	
Date	All
Filter	Today
Action	Day
<input type="button" value="X"/> <input type="button" value="✓"/>	

The activity-button does not filter a certain time frame and is set by default. only chooses results from the current day. If you have chosen or the memory selection screen appears again and your selection for the date is displayed right next to the button.

The activity-button enables you to choose a concrete day. After pressing the activity-button , the current date is displayed between two browse-arrows . Use the arrows to choose a specific date.

Search memory	
Date	Day <input type="button" value="◀"/> 2009.10.21 <input type="button" value="▶"/>
Filter	All
Action	Display
<input type="button" value="X"/> <input type="button" value="✓"/>	

10.2.2 Filter

The activity-button allows you to select results with specific properties from the memory. You can choose from multiple filter criteria.

Search memory		
Date	All	Not printed
Filter	ID	Not sent
Action	Seq.No.	Positive
	Seq.No. Range	Sediment
		Error
<input type="button" value="X"/> <input type="button" value="✓"/>		

All

Selects all data for the respective time frame.

ID

Selects a specific patient ID. After you press the ID, an empty edit-button appears. To enter the patient ID, touch the button and enter the ID via the input-menu . Press , to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted. The ID can optionally also be entered with a barcode scanner.

Search memory	
Date	All
Filter	ID <input type="text" value=""/>
Action	Display
<input type="button" value="X"/> <input type="button" value="✓"/>	

Seq.No.

Selects a specific sequence number. After you press the Seq. No., an edit-button with four "0" appears. To enter a new number, touch the button and enter the number via the input-menu. Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

Search memory	
Date	All
Filter	Seq.No. <input type="text" value="0000"/>
Action	Display
<input type="button" value="X"/> <input type="button" value="✓"/>	

Seq.No. range

Selects a specific sequence number range. After you press the Seq.No., two empty edit-buttons with four "0" appear. To enter the sequence numbers, touch the buttons and enter the numbers via the input-menu . Press , to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

Search memory	
Date	All
Filter	Seq.No. range <input type="text" value="0000"/> - <input type="text" value="0000"/>
Action	Display
<input type="text" value="✕"/> <input type="text" value="✓"/>	

Not printed

Selects all results that have not yet been printed.

Not sent

Selects all results that have not yet been sent.

Positive

Selects all results with at least one positive value.

Sediment

Selects all results that already contain sediment information.

Error

Selects all results with error messages.

10.2.3 Action

The activity-button determines which action will be performed for the results that have been specified before using the date and filter functions. Press the button and all action options appear. Press , to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

Search memory	
Date	Display
Filter	Print
Action	Send
	Print list
	Delete
<input type="text" value="✕"/> <input type="text" value="✓"/>	

Display

Displays the selected results.

Print

Prints the selected results.

Send

Sends the selected results to a PC/LIS.

Print list

Prints the selected patient IDs and sequence numbers only.

Delete (Filter selection)

Deletes the results selected with the search function.

Delete (complete memory)

The activity-button within the memory-menu deletes the entire memory. To do so, press the button and confirm the following security question.

NOTICE

Saved data of control measurements will not be deleted with this function.

NOTICE

Make sure the memory does not contain anymore relevant data, before you delete results. Once data has been deleted from the instrument memory, it cannot be recovered. We therefore recommend sending the data to a PC first.

11. Check mode

The *Check mode* menu allows you to access the quality control functions (see chapter 14. Quality control, page 28).

Check mode	
Settings	Internal check
Enter sediment	Control measurement
Memory	Control strip
Check mode	Stored data
<input type="text" value="✕"/>	

11.1 Internal check

The self-test is repeated and the measuring values are displayed.

11.2 Control measurement

This button initiates a control measurement with artificial urine.

11.3 Control strip

This button initiates the measurement of a colored control strip as an inspection equipment monitoring tool. For ordering of the colored control strips contact MACHEREY-NAGEL or your local distributor.

11.4 Stored data

This button allows you to access the saved results of the last 20 control measurements. Press the button to display the results. You can either print or send the data.

12. General information for measurements

Before the actual testing procedure with the instrument can be started, some preparations concerning the test strip are necessary. Adhere to this information for every testing procedure and for every test as to ensure correct and safe results.

12.1 Urine sample

- Use fresh urine (within 4 h after sample collection)
- Do not centrifuge sample
- Mix sample well
- Cooled urine has to be warmed to room temperature (15–30 °C)

12.2 Test strip preparation

CAUTION

Biological Hazard: Urine residues

Urine and used test strips pose a risk of infection. Always wear protective gloves when handling and disposal. Dispose of used test strips according to the regulations for handling potentially infectious material.

1. Take one test strip out of the box.
2. Do not touch test pads with fingers.
3. Dip the test strip for 1 second into the sample. Make sure all test pads are immersed completely.
4. Draw the edge of the strip along the rim of the sample container to remove excess urine.
5. Touch the long edge of the strip briefly (1 second) to an absorbent paper towel to remove excess urine between the test pads.
6. The test strip is ready to be measured.

13. Measurement

13.1 Sample and patient information

Right after the instrument has been switched on, it is ready to perform measurements. As soon as the start screen appears, the green light above the transport tray indicates readiness.

The right half of the start screen with the heading **Worklist** is initially not activated and is only used as soon as patient IDs are imported (see *chapter 13.4 Measurements with worklist, page 26*).

You can add supplemental information concerning the patient and sample. To do so, press the activity-button **Edit** in the lower left of the screen. A new menu appears, in which you can add the supplemental information.

Press to confirm and save your changes. If you leave the menu by pressing , none of your changes will be adopted.

ID

Press the edit-button for ID. The input-menu appears and you can edit the patient ID.

Seq.No

This edit-button contains the current sequence number for the next measurement. The sequence number always starts at 0001 and is set back after every change in dates. To change the sequence number, press the edit-button containing the number. The input-menu appears and you can edit the sequence number.

Clarity


This option-field allows you to add visual information about the sample's clarity. Press the choice-buttons to select one of the four different options. You can set the description for the four parameters in the settings (see *chapter 8.8 Customization, page 20*). By default, the clarity parameters are set to "Clear", "Turbid", "Cloudy" and "Dark".

Color

This option-field allows you to add visual information about the sample's color. Press the choice-buttons to select one of the four different options. You can set the description for the four parameters in the settings (see *chapter 8.8 Customization, page 20*). By default, the color parameters are set to "Yellow", "Red", "Brown" and "Bright".

13.2 Single measurement

The green light above the transport tray indicates that the instrument is ready to measure. If the autostart function is disabled, you have to initiate a measurement by pressing in the upper right of the start screen (see chapter 8.11 Autostart, page 20).

Place the test strip on the left side of the instrument onto the black transport tray. The light briefly flashes green, thus indicating automatic strip detection. Right after, the transport rack moves to take the strip into the instrument. During that time, the light is red, because no other strip should be placed on the tray while a strip is drawn into the housing.

The measurement of the strip is automatic. The measuring time for one single strip is about 1 minute. When the measurement is completed, the results are printed. The result is saved automatically and can be found in the instrument's memory. You can also send the result to a PC (see chapter 15. Connection of external devices, page 29).

NOTICE

If you choose "print after Sediment entry" in the sediment settings, then results that require sediment findings will be printed after you have entered the sediment findings, not right after the measurement has been performed.

13.3 Successive measurements

NOTICE

Always enter patient and sample information before you place the strip on the transport tray. Once the strip is drawn into the housing, information cannot be edited anymore.

Generally, successive measurements follow the same principles as single measurements.

Place the test strip on the left side of the instrument onto the black transport tray. The light briefly flashes green, thus indicating automatic strip detection. Right after, the transport rack moves to take the strip into the instrument. During that time, the light is red, because no other strip should be placed while the strip is drawn into the housing. The measurement of the strip is automatic.

As soon as the strip has been drawn into the housing, the light is green again and the display shows a new sequence number. Enter necessary supplemental information as usual, then place the next strip onto the transport tray. The light flashes briefly and the next strip is drawn in.

Thus, you can measure any number of strips successively.

13.4 Measurements with worklist


With the instrument, you can create a list of samples and their corresponding patient IDs and then measure those consecutively. The worklist is displayed in the right half of the start screen. As soon as you have imported a patient ID with a barcode scanner, keyboard or your LIS, the IDs are added to the worklist automatically.

In the lower part of the screen, the activity-button is always activated. Once you press the activity-button the activity-button appears, which allows you to load a worklist from your LIS at any time.


13.4.1 Create a worklist

In order to manually create a worklist, you have to connect a barcode scanner or a keyboard first.

Import the samples' barcoded patient IDs with a barcode scanner into the worklist, or use a keyboard to enter patient IDs. Make sure to import the data in the same order as you want to measure the samples. The first imported ID appears on the left side of the start screen in the ID field. All other imported IDs are listed in the worklist area on the right side of the screen. You can always add new IDs to an existing worklist, even if you already started with your measurements. Those additional IDs will be added to the end of the existing worklist.

You have the option to load a worklist from your LIS if it is connected to the instrument. To do so, press the activity-button underneath the worklist field. Then press the activity-button . The instrument can now receive data from your LIS. The first ID will be placed in the left side of the start screen as next sample ID. The other IDs will be listed in the worklist. You can always add additional IDs to the existing list with a barcode scanner or a keyboard.



NOTICE

You can only transfer worklist data from your LIS to an instrument, if the LIS is programmed and set accordingly. Such an adjustment cannot be made by the device alone.

13.4.2 Edit worklist

NOTICE
 You can only edit information for samples, as long as they are in the worklist. Therefore, always edit samples before you place the respective strip on the transport tray. As soon as the strip has been drawn into the housing, the data cannot be edited anymore.

As for every sample, you can edit the worklist entries and add additional information as well. To do so, press on an ID in the worklist, or choose an ID via the scroll-arrows   and press . The edit screen for measurements appears and you can add or edit ID, Seq.No., clarity and color, respectively (see chapter 8.8 Customization, page 20).

The activity-button  allows you to delete an entry from the worklist. If you press , the entry will be placed at the end of the list.

13.4.3 Measuring with the worklist

NOTICE
 It is very important to place the strips onto the transport rack according to the order defined by the worklist, so that the correct strips will be measured for the displayed ID.

When you measure with the worklist, patient IDs are already set for every sample. The sample's ID which will be measured next is always in the left side of the start screen under ID.


Dip the test strip into the sample. To start measuring, place the test strip on the left side of the instrument onto the black transport tray. The light briefly flashes green, thus indicating automatic strip detection. Right after, the transport rack moves to take the strip into the instrument. During that time, the light is red, because no other strip should be placed while the strip is drawn into the housing. The measurement of the strip is automatic.

Make sure to always place the strips onto the transport rack in the order corresponding to the worklist entries. Once the strip has been drawn in, a new sequence number appears and the ID on the left side is replaced with the first ID from the worklist on the right. The worklist on the right moves up by one entry.


If necessary, add supplemental information and place the next strip onto the transport rack. Thus, you can work your way to the end of the list.

Within a worklist, you cannot change the sequence numbers, since the individual measurements are pre-defined by the entries in the worklist.

13.4.4 Options

On the right side of the screen underneath the worklist, you can find the activity-button . These options include three choices.

13.4.4.1 Load list

The activity-button  allows you to load an existing worklist from your LIS (see chapter 13.4.1 Create a worklist, page 26).


13.4.4.2 Print list

The activity-button  allows you to print all IDs in the worklist, including the respective sequence number and supplemental information. Press the button to start printing.

13.4.4.3 Delete list

The activity-button  removes all entries in the worklist irrevocably. In order to delete the list, press the activity-button  and confirm the following security question with . The results of samples from the worklist that have already been measured remain in the instrument's memory.

13.4.5 Emergency measurement


You can slip in an emergency measurement while you are working with a worklist. Press the activity-button  in the lower left part of the screen. The menu for patient and sample information appears, and you can enter patient ID, etc. The new ID is placed as the next sample in the left part of the screen, while the other samples all move back one position. Then, perform the measurement as usual.

14. Quality control


Every facility has to implement its own QC policy. To ensure the instrument performs properly, it has to be checked on a regular basis. In addition to the internal self-test, we therefore recommend control measurements with known positive and known negative control solutions (e.g. Medi-Test Control REF 93038).

14.1 Check mode


To get to the check mode, press the activity-button . You find the check mode as the last option in the main menu. Within the menu *Check mode* there are four different options.

14.1.1 Internal check

Once you press the activity-button , the instrument starts a self-test. During the process, both the LEDs and the instrument's mechanics are checked for their operational capabilities. After the internal check, the instrument prints values for the intensity of the LEDs. If you have trouble with your instrument, you can send the LED values to our technical service to get initial feedback on your problem.

14.1.2 Control measurement

The activity-button  initiates quality control measurements with known positive and known negative quality control solutions. Use your normal test strips for the measurement. To start the control measurement, press the activity-button. First, you have to measure a strip with negative control solution. A new screen appears asking you to "Insert NEG strip". As soon as you have placed the respective strip, the instrument asks to "Insert POS strip" and you have to insert a strip with positive control solution. Then, the message "Please wait" is displayed until both strips have been measured. Once the measurements are completed, the results are printed and the start-screen appears.

Perform a normal measurement with negative and positive control solution (e.g. Medi-Test Control REF 93038). Preparation and the measurement itself are carried out exactly like a normal measurement (see *chapter 13. Measurement, page 25*). The results are printed and saved.

To evaluate the results, read the instructions of the control solution manufacturer. In case you are using Medi-Test Control solutions from MACHEREY-NAGEL, refer to the reference table provided in the instructions for evaluation.



We recommend performing tests with control solutions:

- In the beginning of every day
- When a new batch of strips is used
- When results are questionable
- Whenever the operators change


14.1.3 Control strip measurement

NOTICE


Do not touch the test fields of the control strip with your fingers.

The activity-button  initiates an inspection equipment monitoring process with the help of colored control strips. Only use the colored control strips provided by MACHEREY-NAGEL for this purpose. Empty the the waste container and clean the transport tray before starting the control strip measurement (see *chapter 16. Cleaning and maintenance, page 29*). To start the measurement with the colored control strip, press the activity-button. Follow the instruction and measure the colored control strip. Then, the message "Please wait" is displayed until the control strip has been measured. The result of the measurement is displayed and printed and can now be compared to the values given in the instruction leaflet of the colored control strips. Leave the menu via . In case the values do not match the ranges given in the instruction leaflet, contact MACHEREY-NAGEL or your local distributor.

14.1.4 Stored data

The activity-button  allows you to review the saved results of your control measurements. The results are displayed in exactly the same way as for normal measurements. However, here you cannot enter sediment information.

You can print and send the data for control measurements with the icons in the lower part of the screen.

15. Connection of external devices

15.1 Computer

For detailed information about how to connect your instrument to a PC, refer directly to the technical service.

15.2 Barcode reader / Keyboard

NOTICE

Switch off the device before connecting further devices. Do not switch on the device until all data and power cables have been fully set up and connected. Install a core ferrite before using further devices.

WARNING

Be sure to read the safety instructions carefully before using the external device. (barcode reader /keyboard).

To connect a barcode scanner or keyboard, turn off the instrument. Connect the scanner /keyboard via the PS/2 port with the instrument. When using a barcode reader or keyboard, a split core Ferrite has to be installed on the respective cable (10.5 mm, 1 MHz, 1 GHz, 270 ohm, STAR-TEC Series). Turn on the instrument. The scanner /keyboard will be recognized automatically and can be used immediately.

15.3 Connection to a LIS

To connect the instrument to an existing LIS, you will need technical support. To do so, contact your LIS provider.

16. Cleaning and maintenance

CAUTION

Biological Hazard: Urine residues

Urine and used test strips pose a risk of infection. Always wear protective gloves when handling and disposal. Dispose of used test strips according to the regulations for handling potentially infectious material.

Generally, the instrument has to be turned off and disconnected from the power supply, before you start cleaning it.

16.1 Daily cleaning

NOTICE

Do not clean the transport tray in a dish-washer!

The transport tray and the waste container have to be cleaned and emptied every day as to ensure the correct functioning of the instrument.

Carefully remove the black transport tray from the housing. To do so, pull on the black handle on the transport tray's right side. Empty the waste container. Rinse the table with clear water, afterwards clean the transport tray with 70–85 % ethanol or isopropanol. Make sure that the transport tray is dry before it is reinserted into the device.

16.2 Weekly cleaning

NOTICE

Never use aggressive, organic solvents such as acetone for cleaning the instrument. We recommend 70 % isopropanol for cleaning.

Housing, transport rack and touch screen should be cleaned at least once a week.

Wipe the housing with a damp cloth. Use 70–85 % ethanol or isopropanol. Pay attention to a sufficient exposure time. Transport plate, the surface underneath the transport table and the display can be cleaned in the same way. Make sure that no water penetrates the housing or gets into the printer.

We also recommend cleaning the transport rack on a regular basis in the same way as the housing. Make sure not to use strongly alkaline cleaning agents on the transport rack.

17. Service

17.1 Spare parts and accessories

Item	REF
Transport tray	930751
Transport rack	930752
Power pack	930753
Power cable EU	930754.EU
Power cable Great Britain	930754.GB
Printer paper	93071
PC connection cable (RS232)	919772
USB cable	919687
Instruction for use German / English	930081
Barcode scanner	93074

17.2 Consumables

The following test strips can be evaluated with the URYXXON® 500:

Item	REF
Medi-Test URYXXON® Stick 10	93068

Measurements should be regularly monitored with positive and negative control solutions (Medi-Test Control, REF 93038).

17.3 Error messages

NOTICE

If any of the problems persist, contact your local MACHEREY-NAGEL representative.

The instrument has several error codes which are displayed in case specific problems occur. All error codes are displayed as a number in the touch screen. Measurements with an error will be stored in the instruments memory including the respective error code.

The following table lists the existing error codes, probable causes and possible corrective activities:

Error (code)	Cause	Solutions
E-02	Wrong strip / strip is bent	Strip was not identified as URYXXON® Stick 10 → check strip → use only URYXXON® Stick 10 → Check that the used test strips are not bent.
E-04	Dry strip	Some parameters were detected to be off limit, which may indicate (partially) dry test pads → dip all test pads into the sample completely
E-06	Communication error	A connected device cannot work with the instrument's data protocol (only bi-directional communication) → only connect external devices that can communicate with the instrument → adjust the hardware settings
E-07	Strip not aligned	Strip is not parallel to the side of the transport tray → insert strip carefully → do not insert strip while transport rack is moving and the light is red → empty the waste container, which may cause jams if it is full → do not move transport tray during measurements → make sure the transport rack is placed correctly on the transport plate
E-08	Ref-Pad Error	Wrong reference pad value → turn the instrument off and on again → perform internal check → make a control measurement and check the result's accuracy
E-09	External light source	Stray light is causing problems during measurement → re-position the instrument so it is not exposed to direct sunlight or very bright light sources
E-11	Mechanical error	Mechanical resistance or a foreign object are disturbing the measurement → turn off the instrument, remove transport tray and carefully remove any objects obstructing the mechanics → do not move the transport tray during measurements
E-13	Waste container is full	Too many strips are in the waste container → remove transport tray and empty waste container
Transport tray cannot be removed	The instrument has been switched of during transport rack movement.	Switch on the instrument and wait until the self test is finished. The transport tray can now be removed again.

17.4 Disposal

NOTICE

Disposal using public waste disposal facilities is not permitted. In case of disposal, contact your MACHEREY-NAGEL representative.

Waste disposal according to EU Directive 2012/19/EU. In compliance with national legal regulations (EU Directive 2012/19/EU), MACHEREY-NAGEL disposes the old instruments.

17.5 Analytical Performance Characteristics

17.5.1 Table of results

Param.	CONV	SI	ARB
BLD	NEG	NEG	NEG
	10 Ery/ μ L	10 Ery/ μ L	+
	50 Ery/ μ L	50 Ery/ μ L	++
	250 Ery/ μ L	250 Ery/ μ L	+++
UBG	NORM	NORM	NORM
	2 mg/dL	35 μ mol/L	+
	4 mg/dL	70 μ mol/L	++
	8 mg/dL	140 μ mol/L	+++
BIL	NEG	NEG	NEG
	1 mg/dL	17 μ mol/L	+
	2 mg/dL	35 μ mol/L	++
	4 mg/dL	70 μ mol/L	+++
PRO	NEG	NEG	NEG
	30 mg/dL	0.3 g/L	+
	100 mg/dL	1 g/L	++
NIT	NEG	NEG	NEG
	POS	POS	+
	KET	NEG	NEG
GLU	NEG	NEG	NEG
	NORM	NORM	NORM
	50 mg/dL	2.8 mmol/L	+
	150 mg/dL	8.3 mmol/L	++
pH	≥ 500 mg/dL	≥ 27.8 mmol/L	+++
	5	5	5
	6	6	6
	6.5	6.5	6.5
	7	7	7
SG	8	8	8
	9	9	9
	1.000	1.000	1.000
	1.005	1.005	1.005

Param.	CONV	SI	ARB
	1.010	1.010	1.010
	1.015	1.015	1.015
	1.020	1.020	1.020
	1.025	1.025	1.025
	1.030	1.030	1.030
LEU	NEG	NEG	NEG
	25 Leu/ μ L	25 Leu/ μ L	+
	75 Leu/ μ L	75 Leu/ μ L	++
	500 Leu/ μ L	500 Leu/ μ L	+++

Meaning of the used abbreviations:

BLD - Blood, UBG - Urobilinogen, BIL - Bilirubin, PRO - Protein, NIT - Nitrite, KET - Ketone, GLU - Glucose, SG - specific gravity of the urine (density), LEU - Leukocytes

17.5.2 Precision

The precision was evaluated by using the commercially available control solution Medi-Test Control, manufactured by MACHEREY-NAGEL, Dueren, Germany. This Control is used for daily routine system check of the URYXXON® 500 urinalysis system. It is available as a 2 level set (Medi-Test Control N = level 1 (N): negative-normal, Medi-Test Control P = level 2 (P): positive values). Repeatability and reproducibility amount to 100 % in all cases.

The results are summarized in the following table.

Precision / Repeatability within series

	BIL	BLD	GLU	KET	LEU	NIT	PRO	UBG	pH	SG
Total measurements	40	40	40	40	40	40	40	40	40	40
True positive	20	20	20	20	20	20	20	20	20*	20*
False positive	0	0	0	0	0	0	0	0	n/a	n/a
False negative	0	0	0	0	0	0	0	0	n/a	n/a
True negative	20	20	20	20	20	20	20	20	20*	20*

Precision / Repeatability batch-to-batch

	BIL	BLD	GLU	KET	LEU	NIT	PRO	UBG	pH	SG
Total measurements	60	60	60	60	60	60	60	60	40	40
True positive	30	30	30	30	30	30	30	30	20*	20*
False positive	0	0	0	0	0	0	0	0	n/a	n/a
False negative	0	0	0	0	0	0	0	0	n/a	n/a
True negative	30	30	30	30	30	30	30	30	20*	20*

Precision / Reproducibility day-to-day

	BIL	BLD	GLU	KET	LEU	NIT	PRO	UBG	pH	SG
Total measurements	100	100	100	100	100	100	100	100	100	100
True positive	50	50	50	50	50	50	50	50	50*	50*
False positive	0	0	0	0	0	0	0	0	n/a	n/a
False negative	0	0	0	0	0	0	0	0	n/a	n/a
True negative	50	50	50	50	50	50	50	50	50*	50*

* pH and density do not provide negative / positive interpretation. Obtained values were in accordance to target values of Control N and P.

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17.5.4 Cut-off concentration

The Cut-Off-Concentration (Medical decision point) has been defined as the concentration, where minimum 50% of the results are positive.

Parameter	Cut-off concentration URYXXON® 500
Blood	approx. 4 Ery/μL
Urobilinogen	approx. 2.0 mg/dL
Bilirubin	approx. 0.8 mg/dL
Protein	approx. 15 mg/dL
Nitrite	approx. 0.05 mg/dL
Ketones	approx. 10 mg/dL
Glucose	approx. 35 mg/dL
Leucocytes	approx. 10 Leu/μL

17.6 Changing the cut-off points

NOTICE

Changes in the cut-off points are a major intervention for the instrument's measuring procedure. Changes in the cut-off points will alter the measuring results. Such results have to be verified. If you do not have the capability of verifying such new measurements, we strongly recommend not changing the cut-off points. MACHEREY-NAGEL cannot guarantee any measurements made with changed cut-off points (see *chapter 17.10 General usage comments, page 33*).

It is possible to change the instrument's pre-defined cut-off points. If you want to make changes to the cut-off points, contact MACHEREY-NAGEL directly.

17.7 Comparative study

In a comparative study 240 urine samples were evaluated with urine test strips on the URYXXON® 500 and similar competitor instruments.

The results of the study are summarized in the following table.

The exact agreement for parameters blood, urobilinogen, bilirubin, protein, nitrite, ketones, glucose, pH and leukocytes amounts to ≥ 80 % in every case. For parameter density agreement ± 1 bin amounts to 88 % and 43 % when compared to two different competitor products. Agreement ± 2 bins amounts to 100 % and 85 %, respectively. Overall average exact agreement is 91 % and 86 % when compared to competitor products. Overall agreement ± 1 bin amounts to 100 % and 98 %.

Parameter	Comparison		Comparison	
	URYXXON® 500 – Competitor 1 ***		URYXXON® 500 – Competitor 2 ***	
	exact agreement	agreement ± 1 bin	exact agreement	agreement ± 1 bin
Blood	84 %	100 %	86 %	100 %
Urobilinogen	94 %	100 %	98 %	100 %
Bilirubin	91 %	100 %	84 %	100 %
Protein	92 %	100 %	90 %	100 %
Nitrite	92 %	N/A	97 %	N/A
Ketones	99 %	100 %	97 %	100 %
Glucose	95 %	100 %	94 %	100 %
pH	81 %	100 %	87 %	99 %
Density	88 %*	100%**	43 %*	85 %**
Leukocytes	90 %	100 %	85 %	99 %
overall average	91 %	100 %	86 %	98 %

* this value represents agreement ± 1 bin

** this value represents agreement ± 2 bins

*** similar competitor devices with similar intended purpose (reflectometric device for automatic read-out of urine test strips)

Diagnostic specificity, diagnostic sensitivity, negative predictive value and positive predictive value for parameters blood, urobilinogen, bilirubin, protein, nitrite, ketones, glucose and leukocytes are summarized in the following table. Overall, diagnostic specificity is 95 % and diagnostic sensitivity is 98 %. The average negative predictive value is 98 %; the average positive predictive value is 86 %. Values for the parameters pH and density cannot be calculated because these parameters do not allow negative/positive interpretation.

Parameter	Diagnostic specificity	Diagnostic sensitivity	NPV	PPV
Blood	100 %	91 %	87 %	100 %
Urobilinogen	96 %	100 %	100 %	75 %
Bilirubin	89 %	100 %	100 %	73 %
Protein	88 %	100 %	100 %	83 %
Nitrite	90 %	100 %	100 %	73 %
Ketones	99 %	100 %	100 %	97 %
Glucose	100 %	99 %	99 %	100 %
pH	N/A	N/A	N/A	N/A
Density	N/A	N/A	N/A	N/A
Leukocytes	94 %	97 %	99 %	86 %
overall average	95 %	98 %	98 %	86 %

Likelihood ratios for parameters BLD, UBG, KET, GLU and LEU were calculated as LR+ > 10 and LR- < 0.1, which demonstrates convincing diagnostic evidence. For parameters BIL, PRO and NIT, likelihood ratios were calculated as LR+ 5-10, which demonstrates high diagnostic evidence. LR- was <0.1 for these parameters. Calculation of LR+ and LR- is not possible for pH and SG as diagnostic sensitivity and diagnostic specificity are not applicable for these parameters.

17.8 Interference study

The influence of possible interfering substances is described for each test parameter:

Blood: Normal concentrations of ascorbic acid (≤ 40 mg/dL) do not influence the test result. False-positive reactions can be caused by residues of cleaning agents which contain peroxide or other cleaning agents, as well as menstrual blood.**

Urobilinogen: The detection is inhibited by higher concentrations of formaldehyde (> 30 mg/dL). Nitrite concentrations > 2.5 mg/dL and prolonged exposure of the urine to light can lead to low or false-negative values.* Overly high or false-positive results can be caused by dyes (e.g. betanin) or medications excreted in the urine.**

Bilirubin: The detection is inhibited by higher concentrations of ascorbic acid (> 40 mg/dL) and nitrite (> 2.5 mg/dL). Prolonged exposure of the urine to light can lead to low or false-negative values.* Excreted dyes (e.g. betanin) and medications (e.g. phenazopyridine > 0.1 mg/dL) can simulate a positive result as well as urine indican at a concentration of > 40 mg/dL.**

Protein: False-positive findings can occur in the case of extremely alkaline urine (pH > 9), disinfectant residues (e.g. benzalkonium chloride > 12.5 mg/dL) in the urine container or in the presence of quinine (> 50 mg/dL).**

Nitrite: False-negative results can occur in the case of antibiotic therapy, and in the case of an overly low nitrate level in the urine as a result of low-nitrate food or severe dilution (diuresis). Microbes without the ability to form nitrite can also be present. A false-positive reaction color can be caused by phenazopyridine (> 0.1 mg/dL) or dyes (e.g. betanin) excreted in the urine.**

Ketone: Phthalein compounds in concentrations of > 75 mg/dL may result in false positive results.**

Glucose: Normal concentrations of ascorbic acid (≤ 40 mg/dL) do not influence the test result.* False-positive reactions can be caused by oxidising cleaning agents in the sample container.**

pH: The dye betanin has no influence on the evaluation of the pH test field up to the highest tested concentration of 500 mg/dL.

Density: In the case of elevated protein excretion (> 500 mg/dL), the density values determined are too low.

Leukocytes: An attenuated reaction can be expected if preparations with nitrofurantoin (> 2 mg/dL) or phenazopyridine (> 0.2 mg/dL) are taken.* Formaldehyde (as a preservative, > 30 mg/dL) and dyes (e.g. betanin) can lead to a false-positive reaction. In the case of specimens from female patients, a false-positive reaction can be simulated by vaginal discharge.**

* Interference study with urine with pathological finding (first positive scale value).

** Interference study with urine without pathological findings.

17.9 Warranty

The warranty period is 24 months. The original commercial invoice serves as the guarantee card and must be presented in case of a claim during the given period. The warranty is not valid in cases of improper handling and the warranty does not cover defects caused by the external power supply. The warranty is limited to repairs of defective parts or, in our sole discretion, replacement with a brand new and defect-free instrument.

Any warranty claims made or alleged do not influence the total warranty period of 24 months. The right of withdrawal does not exist. Further claims are excluded and these especially include claims for compensation caused by direct or indirect damages. Additionally our current Terms and Conditions for the Sale of Goods as printed on price lists shall apply.

17.10 General usage comments

Only use the instrument as instructed to analyse URYXXON® Stick 10 test strips. Do not open the instrument or make any other unauthorized modifications. The instrument is a highly sensitive measuring device with high accuracy for the reflectometric evaluation of URYXXON® Stick 10 urine test strips. During the production process, first all optical components are adjusted with specially designed tools. Subsequently, technical experts calibrate the gauge head. Unauthorized modifications or opening the instrument in an inappropriate way can result in maladjustments of the gauge head, or similar damage, which in turn can lead to wrong results.

NOTICE

Due to the reasons stated above, MACHEREY-NAGEL declines any responsibility for the correct functioning of the instrument if it has been opened or used improperly and objects to any warranty claims made in such a case. MN's warranty and liability obligations as to product(s) or to any party on account thereof shall be void and of no effect if the product(s) or any part thereof have been subject to unauthorized change, abuse, misuse, accident or alteration, unauthorized application or installation not in compliance with the product(s) specifications as set forth in the product leaflets/instructions hereof ("unauthorized action") and such unauthorized action resulted in the non-conformance or defect of the product(s).

NOTICE

Note that any serious incident that has occurred in relation to the product shall be reported immediately to the manufacturer and the competent authority of the member state in which the incident occurred. European vigilance contact points: http://ec.europa.eu/health/md_sector/contact_en

17.11 Literature

1. L. Rodrigues, N. S. Leite-de-Lima, C. Landes, J. G. C. Luz, Changes in admission laboratory tests in patients with maxillofacial fractures and the influence of dento-alveolar trauma, *Dent Traumatol*, 2020, 36, 291 – 297.
2. M. J. C. Schot, S. van Delft, A. M. J. Kooijman-Buiting et al, Analytical performance, agreement and user-friendliness of six point-of-care testing urine analysers for urinary tract infection in general practice, *BMJ Open*, 2015, (5):e006857.

17.12 CE Certification

The CE mark declares that the product complies with the harmonization legislation of the European Community listed below: European Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS 2)

This device complies with the following regulation:
- 2017/746- IVD-Regulation (EU)

17.13 Contact

MACHEREY-NAGEL GmbH & Co. KG

Valenciennes Str. 11

52355 Düren · Germany

Tel.: +49 24 21 969-0

info@mn-net.com · www.mn-net.com

Technical support: If you have any questions about the function or use of the device, or if you suspect a defect, contact the technical support of MACHEREY-NAGEL:

Tel.: +49 24 21 969-333

support@mn-net.com

17.14 Version History

Instruction for use URYXXON® 500, V1 | March 2010

Instruction for use URYXXON® 500, V2 | March 2022

Instruction for use URYXXON® 500, V2.1 | January 2023

Instruction for use URYXXON® 500, V2.2 | January 2025

Reason for revision: Implementation of revised device labeling and contact information.

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17.11 Literaturhinweise

1. L. Rodrigues, N. S. Leite-de-Lima, C. Landes, J. G. C. Luz, Changes in admission laboratory tests in patients with maxillofacial fractures and the influence of dento-alveolar trauma, *Dent Traumatol*, 2020, 36, 291–297.
2. M. J. C. Schot, S. van Delft, A. M. J. Kooijman-Buiting et al, Analytical performance, agreement and user-friendliness of six point-of-care testing urine analysers for urinary tract infection in general practice, *BMJ Open*, 2015, (5):e006857.

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The screenshot shows the 'Date/Time' configuration screen. At the top, it says 'Date/Time'. Below that, there are four rows of controls:

- Time:** Two input boxes containing '15' and '17', followed by a dropdown menu showing 'pm'.
- Time format:** Two radio buttons, '24h' (selected) and '12h'.
- Date:** Three input boxes containing '2009', '10', and '21'.
- Date format:** Three radio buttons, 'yyyy-mm-dd' (selected), 'mm-dd-yyyy', and 'dd-mm-yyyy'.

 At the bottom, there are two buttons: a back button with an 'X' and a confirm button with a checkmark.

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The screenshot shows the 'User mode' configuration screen. At the top, it says 'User mode'. Below that, there are three radio buttons: 'Single user', 'User/Admin' (selected), and 'Multi user'. Underneath, there are two input fields: 'Password' and 'Edit users'. At the bottom, there are two buttons: a back button with an 'X' and a confirm button with a checkmark.

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This screenshot is identical to the previous one, but the 'Password' input field is highlighted with a blue border, indicating it is the active field.

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The screenshot shows the 'Login' screen. At the top, it says 'Login'. Below that, there is a numeric keypad with digits 0-9. At the bottom, there is an 'Admin' button. At the very bottom, there are two buttons: a back button with an 'X' and a confirm button with a checkmark.

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