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
US +1 888 321 62 24 sales-us@mn-net.com

A011653/xxx

Solid reagents for the oxidative digestion of heavy metal samples and total phosphorus in a heating block or microwave

Contents:

- 30 g NanOx Metal digestion reagent (AR)
- 1 measuring spoon 85 mm orange
- 30 g NanOx Metal neutralization reagent (NR)
- 1 measuring spoon 70 mm white
- 1 measuring spoon 85 mm black

Hazard warnings:

Information on hazards can be found on the outer label and in the safety data sheet. The safety data sheet can be downloaded at www.mn-net.com/SDS.

NANOCOLOR® tests for use with NanOx Metal:

Test	Name	Note
0-98	Aluminum 07	Only possible with microwave digestion
1-02	Aluminum	Only possible with microwave digestion
0-14	Cadmium 2	
0-24	Chromate 5	Digestion without neutralization (steps 1 – 6 only)
1-25	Chromate	Digestion without neutralization (steps 1 – 6 only)
1-51	Cobalt	
0-53	Copper 5	
0-54	Copper 7	Digestion without neutralization (steps 1 – 6 only) Use double the amount of R2 when performing the test.
1-53	Copper	Digestion without neutralization (steps 1 – 6 only) Use double the amount of R2 when performing the test.
0-37	Iron 3	Digestion without neutralization (steps 1 – 6 only)
1-28	Iron LR	
1-36	Iron	Digestion without neutralization (steps 1 – 6 only)
0-71	Nickel 4	
1-62	Nickel	
0-55	Phosphate 45	
0-76	Phosphate 1	
0-79	Phosphate 50	Digestion without neutralization (steps 1 – 6 only)
0-80	Phosphate 15	
0-81	Phosphate 5	
0-42	Zinc 6	
0-96	Zinc 4	
1-95	Zinc	
0-01	Zirconium 100	

Principle:

Oxidative digestion for the determination of total phosphorus and complex-bound metals or metal ions that are present in an oxidation state in which they would not be detectable without digestion, e.g., chromium(III).

Interferences:

Samples that consume large amounts of oxidizing agent (e.g., COD values above 1000 mg/L O₂) pose a risk of incomplete digestion.

Procedure:

You will need the following accessories for digestion in the heating block:

- NANOCOLOR® heating block
- NANOCOLOR® reaction tubes 16 mm OD (REF 91680)
- NANOCOLOR® piston pipette (REF 916909)
- NANOCOLOR® pipette tips (REF 916916)

1. Open the reaction tube. Pipette 6 mL of sample into the tube.
2. Add 1 level orange measuring spoon of NanOx Metal digestion reagent AR.
3. Close the tube and shake vigorously.
4. Heat in a heating block for 30 minutes at 120 °C or for 1 hour at 100 °C.
5. Remove the tube from the heating block. Allow to cool to room temperature.
If the liquid is cloudy or not colorless, repeat the process (steps 2 – 5).
6. Swirl the test tube once and then open it.
7. Check for peroxide (e.g., with QUANTOFIX® Peroxide 25, REF 91319).
If peroxide is still detectable, repeat steps 3 – 7.
8. Add 3 level white measuring spoons of NanOx Metal neutralization reagent NR.
The pH value must be between 3 and 7. Add more neutralizing reagent if necessary.
9. Close the tube and shake vigorously.

The digested solution can now be used as a sample for subsequent analysis.

Procedure for XL digestion (e.g., for multiple determinations and standard tests):

You will need the following accessories for XL digestion in the heating block:

- Heating block NANOCOLOR® C2 M (REF 919350.1)
- NANOCOLOR® reaction tubes 22 mm OD (REF 91622)
- NANOCOLOR® piston pipette (REF 916909)
- NANOCOLOR® pipette tips (REF 916916)

1. Open the reaction tube. Pipette 17 mL of sample into the tube.
2. Add 3 level orange measuring spoons of NanOx Metal digestion reagent AR.
3. Close the reaction tube and shake vigorously.
4. Heat in a heating block for 30 minutes at 120 °C or for 1 hour at 100 °C.
5. Remove the tube from the heating block. Allow to cool to room temperature.
If the liquid is cloudy or not colorless, repeat the process (steps 2 – 5).
6. Swirl the test tube once and then open it.
7. Check for peroxide (e.g., with QUANTOFIX® Peroxide 25, REF 91319).
If peroxide is still detectable, repeat steps 3 – 7.
8. Add 3 level white measuring spoons of NanOx Metal neutralization reagent NR.
The pH value must be between 3 and 7. Add more neutralizing reagent if necessary.
9. Close the tube and shake vigorously.

The digested solution can now be used for subsequent analysis. For standard tests, use 17 mL instead of 20 mL of sample and multiply the measurement result by the corresponding dilution factor (x 1.18).

Notes:

If heavily soiled and after repeated use, clean the digestion vessel with a bottle brush and rinse with distilled water.

The instructions for use for pressure digestion in the microwave can be found on our website.

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