

## User manuals

NucleoSpin® DNA Forensic

740840 / .10 /  
.50 / .250

Dear valued customer,

Thank you for choosing MACHEREY-NAGEL for your application. We have attached a short protocol for your review.

To obtain the best results we recommend to follow the detailed protocol available online, especially when you are a first time user of this kit.

All important links to the above mentioned product are listed in this leaflet.

Important information regarding product components, specifications, safety instructions, and processing protocols, can all be found on the product website and accessed easily via the QR code.

## QR-Code product website



[qr.mn-net.com/qr/\(241\)740840](https://qr.mn-net.com/qr/(241)740840)

Use the following QR code or the link below for direct access to the user manual.

## QR-Code user manual



[qr.mn-net.com/qr/\(IFU\)740840](https://qr.mn-net.com/qr/(IFU)740840)
















It is strongly recommended to read the detailed protocol section of the user manual if using the kit for the first time. However, experienced users may refer to the protocol at a glance. The protocol at a glance is designed to be used only as a supplemental tool for quick reference while performing the purification procedure.

We are constantly improving our products and we reserve the right to make changes or additions to protocols. Please check for updated revisions for previously downloaded manuals.

This leaflet does not replace the full manual!












## Protocol at a glance

## Protocol for the isolation of genomic DNA from forensic samples

<b>1 Lyse sample</b>		Add 20 $\mu$ L Liquid Proteinase K, 5 $\mu$ L TCEP solution (14 mg/mL) and 450 $\mu$ L Buffer FOL  Mix 56 °C, 1 h
<b>2 Adjust DNA binding conditions</b>		Add 580 $\mu$ L FOB Vortex
<b>3 Bind DNA</b>		Load 600 $\mu$ L lysate
		1 min at 11,000 x g.
		Load remaining lysate 1 min at 11,000 x g.
<b>4 Wash with FOW1</b>		Add 400 $\mu$ L FOW1 1 min at 11,000 x g.
		
<b>5 Wash with FOW2</b>		Add 400 $\mu$ L FOW2
1 <sup>st</sup> wash		1 min at 11,000 x g.
<b>6 Wash with FOW2</b>		Add 400 $\mu$ L FOW2
2 <sup>nd</sup> wash		1 min at 11,000 x g.
<b>7 Dry silica membrane</b>		2 min at 11,000 x g.
		
<b>8 Elute DNA</b>		Add 50 – 100 $\mu$ L Buffer FOE. 1 min RT
		1 min at 11,000 x g.

## Protocol at a glance

## Protocol for the isolation of genomic DNA from bone or teeth samples

<b>Sample preparation</b>		<p>Mill up to 150 mg of bone or teeth to a fine powder.</p> <p>Transfer up to 150 mg of bone powder into a 2 mL microcentrifuge tube.</p>
<b>1 Pre-lyse samples</b>		<p>Add 200 µL phosphate buffer (0.5 M EDTA/0.25 M PO<sub>4</sub><sup>3-</sup>, pH 8), 700 µL Buffer T1 and 20 µL Proteinase K.</p> <p>Mix 56 °C shaking overnight.</p> <p>Afterwards 48 h at 2–8 °C on a shaking incubator.</p>
<b>2 Lyse sample</b>		<p>Mix, Add 800 µL Buffer B3</p> <p>Mix</p> <p>10 min, 70 °C</p> <p>Mix</p> <p>10 min 5,000 x g.</p> <p>Transfer 800 µL of the supernatant to a new 2 mL microcentrifuge tube.</p>
<b>3 Adjust binding conditions</b>		<p>Add 1160 µL Binding Buffer FOB</p> <p>Mix</p>
<b>4 Bind DNA to NucleoSpin® DNA Forensic Colum</b>		<p>Load 600 µL lysate</p> <p>1 min 11,000 x g</p> <p>Discard the flowthrough and place the column back in the Collection Tube.</p>
<b>5 Repeat binding step</b>		<p>Load 600 µL</p> <p>1 min 11,000 x g</p> <p>Discard the flowthrough and place the column back in the Collection Tube.</p>
<b>6 Load the remaining lysate to the column</b>		<p>1 min at 11,000 x g</p> <p>Discard the flowthrough and place the column in a new Collection Tube (provided)</p>
<b>7 Wash with FOW1</b>		<p>Add 400 µL Buffer FOW1 to the column and centrifuge for 1 min at 11,000 x g. Discard the flowthrough and place the column back into the Collection Tube.</p>
<b>8 Wash with FOW2 (1<sup>st</sup>)</b>		<p>Add 400 µL Buffer FOW2</p> <p>1 min 11,000 x g</p>
<b>9 Wash with FOW2 (2<sup>nd</sup>)</b>		<p>Add 400 µL Buffer FOW2</p> <p>1 min 11,000 x g</p>
<b>10 Dry silica membrane</b>		<p>2 min 11,000 x g</p>
<b>11 Elute DN</b>		<p>Add 100 µL Buffer FOE</p> <p>1 min at room temperature</p> <p>1 min 11,000 x g</p>

## QR-Code Safety Data Sheet

Please find a digital version of the safety data sheets by following the link below:



[www.mn-net.com/sds](http://www.mn-net.com/sds)



We strongly recommend to carefully read the detailed protocol section of the product's user manual. If you have any questions about the protocol or product, please contact our Technical Support.

## Contact MN

### Germany and international

MACHEREY-NAGEL GmbH & Co. KG  
Valenciennner Str. 11 · 52355 Düren · Germany  
Tel.: +49 24 21 969-0  
Toll-free: 0800 26 16 000 (Germany only)  
E-mail: [info@mn-net.com](mailto:info@mn-net.com)

### Technical support Bioanalysis

Tel.: +49 24 21 969-333  
E-mail: [support@mn-net.com](mailto:support@mn-net.com)

### USA

MACHEREY-NAGEL Inc.  
924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA  
Toll-free: 888 321 6224 (MACH)  
E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

### France

MACHEREY-NAGEL SAS  
1, rue Gutenberg – BP135 · 67720 Hoerdt Cedex · France  
Tel.: +33 388 68 22 68  
E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €

Siret 379 859 531 00020 · RCS Strasbourg B379859531  
N° intracommunautaire FR04 379 859 531

### Switzerland

MACHEREY-NAGEL AG  
Hirsackerstr. 7 · 4702 Oensingen · Switzerland  
Tel.: +41 62 388 55 00  
E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)