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BioFix® Nitrification Inhibition Test *N-Tox*

Method:

Amperometric method of determining the specific retardation effect of the substances contained in a sample on nitrite oxidation in connection with microbial nitrification with the aid of defined bacterial strains, preferably of the genus *Nitrobacter*. The results are specified as % retardation of oxygen consumption in the sample solution in comparison to a non-inhibited control solution.

Range:	0–100 % inhibition
Reaction time:	10 min
Reaction temperature:	room temperature

Contents of reagent set:

5 reaction vessels	1 bottle, white, containing 5 mL <i>N-Tox</i> R3
1 bottle, natural colour, containing 100 mL <i>N-Tox</i> R1	1 bottle, black, containing 100 mL <i>N-Tox</i> R4
1 test tube containing 2 mL <i>N-Tox</i> R2	

Safety precautions:

This rapid test does not involve any hazardous substances subject to compulsory marking requirements. According to information sheet B006 1/92 ZH 1/346 issued by the employers' liability insurance association for the chemical industry (Germany), the employed nitrifying microorganisms belong to risk group I, i. e. no risk applies to humans or vertebrates.

Storage:

Reagent *N-Tox* R2 must be stored ice-cold at -20 ± 2 °C. The remainder of BioFix® Nitrification Inhibition Test *N-Tox* can also be stored ice-cold together with reagent *N-Tox* R2 or in a cool state at +2 °C to +8 °C. Note the printed use-by date. Reagent R2 should only be removed from the ice compartment, defrosted and homogenized directly prior to use. When carrying out several tests, always return to the ice compartment for interim storage.

Please read the information below carefully before carrying out a test and be sure to observe the instructions!

- In order to carry out testing in the correct manner, the following materials and equipment are required: Starter Set for BioFix® Nitrification Inhibition Tests (REF 970101), CHROMAFIL® membrane filter, single-use, unsterile, pore size 45 µm (REF 91652), oxygen measuring device, magnetic stirrer, stand with clamps, piston pipettes with tips.
- For optimum testing, the **oxygen content of the sample** or the sample dilution should be **> 8.0 mg/L O₂**. If necessary, enrich with oxygen before testing by briefly aerating or vigorously shaking the sample (dilution).
- The **mixing speed of the magnetic stirrer** significantly influences the amount of oxygen consumption! Therefore, carry out all measurements of a test series with a constant mixing speed. Optimal speed: 200 to 400 rpm
- Adjust temperature of reagents **R1, R3 and R4 to room temperature**, prior to carrying out the tests! This is best achieved by removing these reagents from the refrigerator on the day before analysis.
- Shake reagents R1 and R4 briefly before use!
- First remove reagent R2 from the freezer immediately before use, fully defrost and finally immediately store below freezing again.** Also store below freezing in between multiple tests.
- Reagents **R2 and R3** are to be drawn up using the microliter syringe **without air bubbles!**
- It is essential that the microliter syringe be rinsed after each pipetting operation by drawing up distilled water several times!**
- Rinse the electrode adapter, the reaction vessel and the oxygen electrode after **each** measuring operation, in order to avoid contamination with foreign matter when measuring the next sample.
- The reaction vessels contained in the test kit can be reused after rinsing with distilled water.

Procedure:

It is recommended to carry out a new control measurement after at least every second sample measurement (determination of sample's oxygen self-consumption or determination of the sample's nitrification inhibition) during a measuring sequence!

Owing to the principle of oxygen measurement that this quick test is based on as the measurement size, **samples with oxygen consuming contents** can falsify the measurement results and conceal potential nitrification inhibition effects. To compensate for such effects, an additional calculation of the oxygen self-consumption of the sample or sample dilution (step 4) must be carried out in addition to a control measurement (step 2) and the actual measurement of nitrification inhibition of the sample or sample dilution (step 5), which is taken into consideration during the calculation of the end result.

Recommended Measuring Sequence with Calculation of the Sample's Oxygen Self-Consumption:

C → SC-S → NI-S → C → SC-S → NI-S → C → SC-S → NI-S → C →
 C = Control measurement (step 2)

Recommended Measuring Sequence without Calculation of the Sample's Oxygen Self-Consumption:

C → NI-S → NI-S → C → NI-S → NI-S → C → NI-S → NI-S → C →
 C = Control measurement (step 2)

SC-S = Self-consumption of the sample or sample dilution (step 4)

NI-S = Nitrification inhibition of the sample or sample dilution (step 5)

Step 1: Preparation of nitrite-oxidizing bacteria (Reagent *N-Tox* R2)

First take the reagent *N-Tox* R2 from the freezer immediately before the test, **fully (!)** defrost and homogenize by shaking briefly. Immediately after adding to the control and sample preparations, replace the reagent *N-Tox* R2 in the freezer!

*Note: The bacteria in reagent *N-Tox* R2 also naturally lose activity with increased age even if stored correctly. This does not influence the measurement result of the test in % inhibition, however, since the relationship between possible inhibition of the sample preparation in comparison to the oxygen consumption in the control always remains constant even with increased bacteria activity.*

Step 2: Control

- Fit electrode adapter to O₂ electrode*. Provide reaction vessel with stirring bar and place on magnetic stirrer.
- Shake reagent *N-Tox* R1 vigorously for 30 s** to enrich the oxygen content, then fill the reaction vessel with **reagent *N-Tox* R1** until it overflows, **avoiding air bubbles** (approx. 7–8 mL).
- Carefully insert O₂ electrode with fitted electrode adapter into the reaction vessel until it is sealed **air-tight and free of air bubbles** (surplus medium overflows). **Stir** this preparation on the magnetic stirrer for approx. **1–2 min** until the meter indicates a constant oxygen level. **Make sure there are no air-bubbles in the entire test mixture!**
- In the meantime, wet the entire surface of the electrode adapter and in particular, the two small injection holes with distilled water (e. g. using a squeeze bottle), to ensure an air-tight seal.
- Remove prepared reagent *N-Tox* R2 (see Step 1) from the ice compartment, **fully** defrost and homogenize by shaking briefly. Then **add 100 µL of reagent *N-Tox* R2** via the small injection hole of the electrode adapter, using a microliter syringe. Return reagent *N-Tox* R2 to the ice compartment immediately.
- After 2 minutes** (and before adding reagent *N-Tox* R3 !), **read** and record the **oxygen concentration** in the control **O_{CO}** (O₂ concentration of control at time t = 0 min).
*Important: Rinse microliter syringe thoroughly by drawing up distilled water several times in the course of the 2 minutes and before reagent *N-Tox* R3 is pipetted!*
- Now add **100 µL of reagent *N-Tox* R3 immediately** via the small injection hole of the electrode adapter, using a microliter syringe.
- After 10 minutes, read** and record the **oxygen concentration** in the control **O_{C10}** (O₂ concentration of control at time t = 10 min).

* For precise information on fitting the electrode adapter to the oxygen electrode, please refer to the instruction leaflet enclosed with the "Starter Set for BioFix® Nitrification Inhibition Tests" (REF 970101).

Step 3: Sample Preparation

- In case of substantially turbid samples, first carry out non-sterilizing preliminary filtration using standard (folded) filters, by centrifuging or comparable methods.
- Set the pH value of the sample to **pH 7.8 ± 0.2**, using 0.1 N NaOH or 0.1 N HCl.
- Fine filter 20 mL sample** using CHROMAFIL® membrane filters, single-use, unsterile, pore size 0.45 µm (REF 91652) to remove the largest parts of the naturally present microflora, which can falsify the results.
- Place **16 mL finely filtered sample solution** and **4 mL reagent *N-Tox* R4** in a suitable container (e. g. suitably sized beaker).
- Mix this mixture well and shake hard for 30 s to aerate. Then use as **sample solution** for calculating the oxygen self-consumption of the sample (step 4) or for calculating the nitrification inhibition (step 5).

Step 4: Oxygen self-consumption of the sample

*Note: For **not** oxygen self-consuming samples, step 4 can be skipped. While evaluating the test results and calculating the nitrification inhibition, the variable $\Delta O_{SC-S} = 0!$*

- Fit electrode adapter to O₂ electrode*. Provide reaction vessel with stirring bar and place on magnetic stirrer.
- Fill the reaction vessel with the **sample solution, pH 7.8** from step 3 to the overflow point **without air bubbles** (approx. 7–8 mL).
- Carefully insert O₂ electrode with fitted electrode adapter into the reaction vessel until it is sealed **air-tight and free of air bubbles** (surplus medium overflows). **Stir** this preparation on the magnetic stirrer for approx. **1–2 min** until the meter indicates a constant oxygen level. **Make sure there are no air-bubbles in the entire test mixture!**
- In the meantime, wet the entire surface of the electrode adapter and in particular, the two small injection holes with distilled water (e. g. using a squeeze bottle), to ensure an air-tight seal.
Important: No reagent *N-Tox* R2 is added!
- After 2 minutes** (and before adding reagent *N-Tox* R3!), **read** and record the **oxygen concentration** in the sample preparation **O_{SC-S0}** (O₂ concentration of sample preparation for determination of oxygen self-consumption at time t = 0 min).
- Now **add 100 µL of reagent *N-Tox* R3 immediately** via the small injection hole of the electrode adapter, using a microliter syringe.
- After 10 minutes, read** and record the **oxygen concentration** in the sample preparation **O_{SC-S10}** (O₂ concentration of sample preparation for determination of oxygen self-consumption at time t = 10 min).

Step 5: Nitrification inhibition of the sample

- Fit electrode adapter to O₂ electrode*. Provide reaction vessel with stirring bar and place on magnetic stirrer.
- Fill the reaction vessel with the **sample solution, pH 7.8** from step 3 to the overflow point **without air bubbles** (approx. 7–8 mL).
- Carefully insert O₂ electrode with fitted electrode adapter into the reaction vessel until it is sealed **air-tight and free of air bubbles** (surplus medium overflows). **Stir** this preparation on the magnetic stirrer for approx. **1–2 min** until the meter indicates a constant oxygen level. **Make sure there are no air-bubbles in the entire test mixture!**
- In the meantime, wet the entire surface of the electrode adapter and in particular, the two small injection holes with distilled water (e. g. using a squeeze bottle), to ensure an air-tight seal.
- Remove reagent *N-Tox* R2 from the ice compartment, **fully** defrost and homogenize by shaking briefly. Then **add 100 µL of reagent *N-Tox* R2** via the small injection hole of the electrode adapter, using a microliter syringe. Return reagent *N-Tox* R2 to the ice compartment immediately.
- After 2 minutes** (and before adding reagent *N-Tox* R3!), **read** and record the **oxygen concentration** in the sample preparation **O_{NI-S0}** (O₂ concentration of sample preparation for determination of nitrification inhibition at time t = 0 min).
- Now **add 100 µL of reagent *N-Tox* R3 immediately** via the small injection hole of the electrode adapter, using a microliter syringe.
- After 10 minutes, read** and record the **oxygen concentration** in the sample preparation **O_{NI-S10}** (O₂ concentration of sample preparation for determination of nitrification inhibition at time t = 10 min).

Step 6: Evaluation

- Oxygen consumption in the control: $\Delta O_C = O_{CO} - O_{C10}$
 - Oxygen self-consumption of the sample: $\Delta O_{SC-S} = O_{SC-S0} - O_{SC-S10}$
 - Oxygen consumption in sample preparation for determination of nitrification inhibition: $\Delta O_{NI-S} = O_{NI-S0} - O_{NI-S10}$
- #: $\Delta O_C \cdot C = 0$, if there is no oxygen self-consumption of the sample

Corrected oxygen consumption of the sample preparation taking oxygen self-consumption in consideration: $\Delta O_S = \Delta ONI-S - \Delta O_{SC-S}$

Final result:

% inhibition of nitrite oxidation = $[(\Delta O_C - \Delta O_S) : \Delta O_C] \times 100$

For evaluation purposes, we recommend using the evaluation sheet provided on the back of these instructions for use. Please feel free to copy this sheet for your own personal use.

* For precise information on fitting the electrode adapter to the oxygen electrode, please refer to the instruction leaflet enclosed with the "Starter Set for BioFix® Nitrification Inhibition Tests" (REF 970101).

Interpretation of Results:

The inhibition values obtained can be interpreted in the following way:

0–10 % inhibition	Sample does not inhibit nitrification
10–20 % inhibition	Sample may inhibit nitrification
20–80 % inhibition	Sample does inhibit nitrification
> 80 % inhibition	Dilute sample and test again

As a further method for comparing and estimating the obtained inhibition data at different sample times or from different locations, which is independent of the standards, the result can also be given as **G_{NI}-Value** according to results of other biotoxicity tests (e. g. luminescent bacteria test according to DIN EN ISO 11348). This is defined as the reciprocal of the first dilution of the sample, that causes a nitrification inhibition of less than 20 %.

Example: If 35 % inhibition is obtained at a 1:2 dilution of a sample and a 15 % inhibition at a 1:4 dilution, in this case the G_{NI}-Value would be G_{NI} = 4.

Disposal:

The polystyrene box with the empty square bottles and test tubes can be disposed of with normal domestic refuse without any problems. The test preparations can be drained away by pouring into a sink. Restrictions may arise if the sample contains harmful or toxic substances which are subject to special disposal requirements. The user is responsible for disposing of such test preparations in accordance with the applicable rules and regulations.

Date of analysis: _____ Name of analyst: _____

Control (Step 2)			Sample Measurement								Results
			Sample Name	Self-Consumption (Step 4)			Nitrification Inhibition (Step 5)			Corrected O ₂ -Consumption	
O _{C0} [mg/L]	O _{C10} [mg/L]	ΔO _C [mg/L] = O _{C0} - O _{C10}		O _{SC-S0} [mg/L]	O _{SC-S10} [mg/L]	ΔO _{SC-S} [mg/L] = O _{SC-S0} - O _{SC-S10}	O _{NI-S0} [mg/L]	O _{NI-S10} [mg/L]	ΔO _{NI-S} [mg/L] = O _{NI-S0} - O _{NI-S10}	ΔO _S [mg/L] = ΔO _{NI-S} - ΔO _{SC-S}	Inhibition [%] = [(ΔO _C - ΔO _S) : ΔO _C] x 100

O_{C0}: Oxygen concentration [mg/L] in control at the beginning of the test at time t₀
 O_{C10}: Oxygen concentration [mg/L] in control at the end of the test after t = 10 min incubation
 ΔO_C: Oxygen consumption [mg/L] in control after t = 10 min incubation
 O_{SC-S0}: Oxygen concentration [mg/L] in sample preparation for determination of oxygen self-consumption at the beginning of the test at time t₀
 O_{SC-S10}: Oxygen concentration [mg/L] in sample preparation for determination of oxygen self-consumption at the end of the test after t = 10 min incubation
 ΔO_{SC-S}: Oxygen consumption [mg/L] in sample preparation for determination of oxygen self-consumption after t = 10 min incubation
 O_{NI-S0}: Oxygen concentration [mg/L] in sample preparation for determination of nitrification inhibition at the beginning of the test at time t₀
 O_{NI-S10}: Oxygen concentration [mg/L] in sample preparation for determination of nitrification inhibition at the end of the test after time t = 10 min incubation
 ΔO_{NI-S}: Oxygen consumption [mg/L] in sample preparation for determination of nitrification inhibition after time t = 10 min incubation

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