

MACHEREY-NAGEL

User manual



Purification of His-tag proteins

- Protino® Ni-IDA 150 Packed Columns
- Protino® Ni-IDA 1000 Packed Columns
- Protino® Ni-IDA 2000 Packed Columns
- Protino® Ni-IDA Resin
- Protino® Ni-IDA 96 Ni-IDA

February 2023 / Rev. 10

Contact MN

Germany and international

MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0
Toll-free: 0800 26 16 000 (Germany only)
E-mail: info@mn-net.com

Technical Support Bioanalysis

Tel.: +49 24 21 969-333
E-mail: support@mn-net.com

USA

MACHEREY-NAGEL Inc.
924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA
Toll-free: 888 321 6224 (MACH)
E-mail: sales-us@mn-net.com

France

MACHEREY-NAGEL SAS
1, rue Gutenberg – BP135 · 67720 Hoerdt Cedex · France
Tel.: +33 388 68 22 68
E-mail: sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €
Siret 379 859 531 00020 · RCS Strasbourg B379859531 · N° intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG
Hirsackerstr. 7 · 4702 Oensingen · Switzerland
Tel.: +41 62 388 55 00
E-mail: sales-ch@mn-net.com

Table of contents

1	Components	4
1.1	Contents and storage	4
1.2	Additional materials to be supplied by user	5
2	Product description	7
2.1	The basic principle	7
2.2	About this user manual	7
3	Product description	9
3.1	Specifications	9
3.2	Purification under native and denaturing conditions	10
3.3	Binding capacity of Protino® Ni-IDA	11
3.4	Culture size	14
3.5	Binding, washing, and elution	17
3.6	Compatibility of reagents	18
4	Safety instructions	20
4.1	Disposal	20
5	Purification of polyhistidine-tagged proteins from <i>E.coli</i> under native conditions	21
5.1	Preparation of buffers for purification under native conditions	21
5.2	Preparation of cleared lysates under native conditions	22
5.3	Protino® Ni-IDA Packed Columns – purification under native conditions	23
5.4	Protino® Ni-IDA Resin – gravity-flow column chromatography under native conditions	26
5.5	Protino® Ni-IDA Resin – batch gravity-flow purification protocol under native conditions	27
5.6	Protino® Ni-IDA Resin – batch purification protocol under native conditions	28
5.7	Protino® Ni-IDA Resin – medium pressure column chromatography under native conditions	29
5.8	Protino® 96 Ni-IDA – purification under native conditions	30
6	Purification of polyhistidine-tagged proteins from <i>E.coli</i> under denaturing conditions	32
6.1	Preparation of buffers for purification under denaturing conditions	32
6.2	Cell extract preparation under denaturing conditions	33
6.3	Protino® Ni-IDA Packed Columns – purification under denaturing conditions	34
6.4	Protino® Ni-IDA Resin – gravity-flow column chromatography under denaturing conditions	35
7	Cleaning, recharging, and storage	36
8	Appendix	37
8.1	Troubleshooting	37
8.2	Ordering information	39
8.3	Product use restriction / warranty	40

1 Components

1.1 Contents and storage

Protino® Ni-IDA 150 Packed Columns		
REF	10 preps 745150.10	50 preps 745150.50
Protino® Ni-IDA 150 Packed Columns	10	50
8 x LEW Buffer	5 mL	30 mL
4 x Elution Buffer	8 mL	20 mL
User manual	1	1

Protino® Ni-IDA 1000 Packed Columns		
REF	5 preps 745160.5	50 preps 745160.50
Protino® Ni-IDA 1000 Packed Columns	5	50
8 x LEW Buffer	30 mL	140 mL
4 x Elution Buffer	8 mL	100 mL
Plastic Washer	5	8
User manual	1	1

Protino® Ni-IDA 2000 Packed Columns		
REF	5 preps 745170.5	25 preps 745170.25
Protino® Ni-IDA 2000 Packed Columns	5	25
8 x LEW Buffer	30 mL	140 mL
4 x Elution Buffer	20 mL	100 mL
Plastic Washer	5	8
User manual	1	1

Kit contents (continued)

Protino® Ni-IDA Resin				
REF	745210.5	745210.30	745210.120	745210.600
Protino® Ni-IDA Resin	5 g	30 g	120 g	600 g
User manual	1	1	1	1

Protino® 96 Ni-IDA		
REF	1 x 96 preps 745300.1	4 x 96 preps 745300.4
Protino® 96 Ni-IDA Plate	1	4
8 x LEW Buffer	2 x 30 mL	140 mL
4 x Elution Buffer	2 x 20 mL	100 mL
User manual	1	1

Storage conditions

All kit components can be stored at 15–25 °C and are stable up to one year.

1.2 Additional materials to be supplied by user

Reagents

- Lysozyme
- **Protino® Ni-IDA 150/1000/2000 Packed Columns**
Purification under native conditions: Kits already contain buffer stock solutions that have to be prepared according to the instructions, section 5.3.
Purification under denaturing conditions: Denaturing Solubilization Buffer, Denaturing Elution Buffer, additional LEW Buffer (sodium phosphate, sodium chloride, urea, and imidazole). For buffer compositions refer to section 6.1.
- **Protino® Ni-IDA Resin**
Purification under native conditions: LEW Buffer, Elution Buffer (sodium phosphate, sodium chloride, imidazole). For buffer compositions refer to section 5.1.
Purification under denaturing conditions: LEW buffer, Denaturing Solubilization Buffer, Denaturing Elution Buffer, (sodium phosphate, sodium chloride, urea, and imidazole). For buffer compositions refer to section 6.1.
- **Protino® 96 Ni-IDA**
Purification under native conditions: Kits already contain buffer stock solutions that have to be prepared according to the instructions, section 5.8.

Consumables

- Appropriate centrifugation / collection tubes
- Protino® Columns for gravity-flow column chromatography using Protino® Ni-IDA Resin
For column IMAC using **Protino® Ni-IDA Resin** we generally recommend gravity-flow procedure. For this MACHEREY-NAGEL offers **Protino® Columns 14 mL** and **35 mL**. **Protino® Columns** are empty polypropylene columns with an inserted filter frit. Separate frits for covering the column bed are also included. Protino® Columns are available with volume capacities of 14 mL and 35 mL (see ordering information). They can be used to retain up to 1.4 g and 3.5 g of Protino® Ni-IDA Resin, respectively. These maximum amounts of resin correspond to a protein binding capacity of 14 mg and 35 mg respectively (for 6 x His-GFPuv, concentration 2 mg/mL). For detailed information on binding capacity please also refer to sections 3.3 and 3.4).

Table 1: Protino® Columns to be used with Protino® Ni-IDA Resin

	Volume capacity [mL]	Max. amount of Protino® Ni-IDA Resin per column [g]	Protein binding capacity ¹ [mg]
Protino® Columns 14 mL	14	1.4	28
Protino® Columns 35 mL	35	3.5	70

Equipment

- Appropriate centrifuge, sonicator

¹ Protein binding capacity refers to 6 x His-GFPuv.

2 Product description

2.1 The basic principle

Protino® Ni-IDA products enable fast and convenient purification of recombinant polyhistidine-tagged proteins by immobilized metal ion affinity chromatography (IMAC). Protino® Ni-IDA is a dry silica-based resin precharged with Ni^{2+} ions. Binding of protein is based on the interaction between the polyhistidine tag of the recombinant protein and immobilized Ni^{2+} ions. The chelating group of Protino® Ni-IDA is based on IDA (iminodiacetic acid), which enables strong and efficient binding of target protein onto the IMAC matrix (Figure 1). In contrast to traditional IDA matrices, Protino® Ni-IDA shows an optimized, low density of IDA ligands, which is created by a special manufacturing process. This non-saturating surface concentration of IDA eliminates non-specific interactions of contaminating proteins with the adsorbent. As a result, Protino® Ni-IDA ensures higher target protein purity.

IDA is a threedentate chelator which occupies three of the six binding sites in the coordination sphere of the Ni^{2+} ion. The remaining three coordination sites are usually occupied by water molecules and can be exchanged with histidine residues of the recombinant protein (Figure 1).




Figure 1 Protino® Ni-IDA – Structure of IDA in complex with Ni^{2+}

2.2 About this user manual

For quick orientation in this user manual please follow the corresponding cross-reference given below.

Table 2: Protocol guide

Product	Application	Page(s)
Protino® Ni-IDA Packed Columns	Gravity flow column chromatography	23, 32
Protino® Ni-IDA Resin	Gravity flow column chromatography	26, 35
	Batch binding (in combination with gravity flow column chromatography)	27
	Batch purification	28
	Medium pressure column chromatography (FPLC™)	29
Protino® 96 Ni-IDA	Gravity flow column chromatography	30

Protino® Ni-IDA Packed Columns:

Experienced users who are performing the purification of His tagged proteins using **Protino® Ni-IDA Packed Columns** may refer to the protocol-at-a-glance instead of this user manual (see section 5.3.1). The protocol-at-a-glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure. First-time users are strongly advised to read this user manual.

The **Protino® Ni-IDA Packed Columns** protocols in this manual are organized as follows: The culture volumes and volumes of the respective buffers used for a particular column size are highlighted. Each procedural step is arranged like the following example (taken from section 5.3):

150	1000	2000
-----	------	------

3 Column equilibration

Equilibrate Protino® Ni-IDA Packed Columns with **1 x LEW Buffer**. Allow the column to drain by gravity

0.4 mL	4 mL	12 mL
--------	------	-------

Protino® Ni-IDA 150 / 1000 Packed Columns are designed to fit into most 15 mL conical centrifuge tubes (e.g., BD Falcon REF 352097) for convenient fraction collection.

For example, if you are using **Protino® Ni-IDA 150 Packed Column** you are requested to refer to the white boxes. These boxes indicate the volume of bacterial culture or buffer to be used. The respective buffer is highlighted in **bold type** within the instruction.

Referring to the a.m. example there has to be used 400 µL of LEW Buffer for column equilibration when using **Protino® Ni-IDA 150 Packed Columns**.

All technical literature is available on the internet at www.mn-net.com.

Please contact Technical Service regarding information about changes of the current user manual compared to previous revisions.

3 Product description

3.1 Specifications

Table 3: Specifications
Specifications Protino® Ni-IDA Packed Columns

	150 ¹	1000 ¹	2000 ¹
Protein capacity (6 x His-GFPuv, ~ 32 kDa)	800 µg (40 mg resin)	5 mg (250 mg resin)	10 mg (500 mg resin)
Application	Gravity flow columns		
Physical form	Ready-to-use columns, filled with dry matrix precharged with Ni ²⁺		

Specifications Protino® Ni-IDA Resin

Protein capacity (6 x His-GFPuv, ~ 32 kDa)	20 mg/g resin 10 mg/mL bed volume
Application	Batch Gravity flow column FPLC™
Physical form	Dry matrix, precharged with Ni ²⁺
Max. pressure	145 psi (10 bar)

Specifications Protino® 96 Ni-IDA

Protein capacity (6 x His-GFPuv, ~ 32 kDa)	1 mg/well (50 mg resin/well)
Application	Gravity flow chromatography
Physical form	Ready-to-use 96-well plates, filled with dry matrix precharged with Ni ²⁺

¹ The nomenclature of the Protino® Ni-IDA Packed Columns is independent from the binding capacity but is to reflect the dimension of the column such as small, medium, and large.

The maximum column capacities of the respective sizes are as follows:

150: 1 mL; 1000: 3 mL; 2000: 6 mL

Table 3: Specifications**Specifications Protino® Ni-IDA Packed Columns, Resin, and 96 kits**

Matrix	Macroporous silica
Density	0.5 g/mL (1 g resin corresponds to 2 mL bed volume)
Chelating group	IDA (iminodiacetic acid)
Mean particle size	90 µm
pH stability	Since silica is susceptible to hydrolysis at high pH, buffers with pH > 8.4 should not be used)
Storage	≤ 25 °C
Recommended imidazole concentration for load/wash	0 mM
Recommended imidazole concentration for elution	≤ 250 mM

- Protino® Ni-IDA products enable routine purification of recombinant polyhistidine-tagged proteins under native or denaturing conditions.
- Although designed for the purification of polyhistidine-tagged proteins from *E.coli*, Protino® Ni-IDA products can also be used for the purification of polyhistidine-tagged proteins from other expression systems including insect cells, mammalian cells, and yeast.
- The capacity of Protino® Ni-IDA (see Table 3) was determined by using polyhistidine-tagged green fluorescent protein (6 x His-GFPuv, ~ 32 kDa) expressed in *E.coli*. Capacities will vary for each His-tagged protein.
- For research use only.

3.2 Purification under native and denaturing conditions

This manual describes methods for the preparation of cell extracts from *E.coli* and procedures for the purification of polyhistidine-tagged recombinant proteins using Protino® Ni-IDA.

If recombinant proteins are expressed in *E.coli* ideally the target proteins remain soluble in the cytoplasm. However, especially proteins that are highly expressed accumulate in insoluble aggregates, which are called inclusion bodies. For solubilization of inclusion bodies buffers containing large amounts of denaturants are used. This manual includes instructions for isolation of soluble proteins (purification under native conditions, see section 5) as well as insoluble proteins from inclusion bodies (purification under denaturing conditions, see section 6).

In general for purification of polyhistidine-tagged proteins, the bacterial cells are disrupted using lysozyme in combination with sonication. After centrifugation, soluble target protein is found in the supernatant while inclusion bodies remain in the pellet. The clear supernatant can directly be subjected to further purification using Protino® Ni-IDA Packed Columns or Protino® Ni-IDA Resin under native conditions (see section 5). In case of massive formation of inclusion bodies

the target protein is extracted from the pellet using a denaturant (8 M urea) and further purified using protocols for the purification under denaturing conditions (see section 6). If the distribution of the recombinant protein is unknown it is recommended to perform SDS-PAGE analysis using the crude cell extract prior to centrifugation and the clear supernatant after centrifugation. While the crude cell extracts will contain both soluble and insoluble target protein, only soluble target protein is found in the supernatant.

3.3 Binding capacity of Protino® Ni-IDA

3.3.1 General information

The binding capacity of Protino® Ni-IDA strongly depends on the **characteristics** of the polyhistidine-tagged protein, for example amino acid composition, molecular weight, 3-D structure, oligomerization properties, etc. Furthermore, the absolute yield also depends on the total **amount** and **concentration** of the target protein in the sample which in turn directly correlate with the expression level and the cell density of the expression culture. **Therefore binding capacity will vary for each polyhistidine-tagged protein and has to be determined for each expression experiment.**

The binding behaviour of any polyhistidine-tagged protein to Protino® Ni-IDA can be examined by calculating the amount of protein that is eluted as a function of the amount of protein that has been loaded (see Figure 2). Please note that the resulting graph will vary in dependence on characteristics and concentration of the individual His-tag protein. The binding curve can be divided in three stages:

- 1. Stage of maximum recovery.** At this stage the loaded protein is bound to the resin nearly quantitatively and can be eluted nearly quantitatively, too (the binding curve is almost linear, see Figure 2, • eluted His-GFPuv).
- 2. Stage of increasing yield / decreasing recovery.** At this stage the binding curve becomes non-linear and finally binding approaches saturation. The protein yield increases with further increasing amount of loaded protein.
- 3. Stage of maximum yield / minimum recovery.** When loading excess protein, the available binding sites of the resin are saturated. The amount of eluted protein reaches a maximum.

The **binding capacity** for each individual protein can be defined as the yield, at which the binding curve changes from the stage of maximum recovery to the stage of increasing yield / decreasing recovery. This point is an optimal compromise between protein load and recovery and will vary for each individual protein.

3.3.2 Binding capacity

The binding capacity of Protino® Ni-IDA is exemplified using the green fluorescent protein (6 x His-GFPuv, ~ 32 kDa) at a concentration of 2 mg/mL.

Please note that different recombinant proteins may show a different binding behaviour.

Figure 2 shows a plot of the amount of eluted 6 x His-GFPuv against the amount of loaded 6 x His-GFPuv. The binding curve can be divided in three stages:

1. Stage of maximum recovery:
< ~ 20 mg 6 x His-GFPuv load/g resin
2. Stage of increasing yield / decreasing recovery:
> ~ 20 mg 6 x His-GFPuv load/g resin
3. Stage of maximum yield / minimum recovery:
> ~ 70 mg 6 x His-GFPuv load/g resin

Under the above mentioned conditions the binding capacity of Protino® Ni-IDA for 6 x His-GFPuv is approximately 20 mg protein per g of resin (see arrow, Figure 2). At this point the protein recovery is > 80 %. Consequently the following amounts of 6 x His-GFPuv have to be loaded: For optimal recovery: load ~ 20 mg protein per 1 g of Protino® Ni-IDA Resin, for maximum yield: load ~ 75 mg protein per 1 g of Protino® Ni-IDA Resin.




Figure 2 Binding behaviour of 6 x His-GFPuv to Protino® Ni-IDA Resin

Gravity flow columns packed with 40 mg of Protino® Ni-IDA Resin were loaded with increasing volumes of an E.coli lysate containing 6 x His-GFPuv (protein concentration 2 mg/mL). After washing with 640 µL LEW Buffer the target protein was eluted with 960 µL Elution Buffer. Yield (left axis) and recovery (right axis) of 6 x His-GFPuv are plotted versus the amount of loaded protein. For convenient analysis the values are converted to mg 6 x His-GFPuv per 1 g resin.

Recovery rates and yield can be increased by using samples containing higher concentrated polyhistidine-tagged protein (6 x His-GFPuv). Figure 3 shows that the yield of purified polyhistidine-tagged protein is not only depending on the total amount of target protein loaded on the column (also see Figure 2) but also on its concentration in the lysate. Consequently the concentration of target protein in the sample should be as high as possible.




Figure 3 Binding behaviour of 6 x His-GFPuv to Protino® Ni-IDA at different concentrations of the polyhistidine-tagged protein in the sample

Recombinant 6 x His-GFPuv was expressed in *E.coli*. The concentration of the target protein in the culture reached 100 mg/L. 1 g cells were lysed in 2 mL LEW Buffer according to section 3.4 in order to obtain a highly concentrated lysate. The concentration of 6 x His-GFPuv in the lysate was 8 mg/mL. Gravity flow columns packed with 40 mg of Protino® Ni-IDA Resin were loaded with increasing volumes (amounts) of the lysate. After washing, the target protein was eluted with Elution buffer. The yield of 6 x His-GFPuv is plotted against the amount of loaded target protein. The same test was performed using diluted lysates with concentrations of 0.5 and 2 mg/mL. For convenient analysis the values are converted to mg 6 x His-GFPuv per 1 g resin.

Please note: The higher the protein concentration in the sample and the higher the total amount of protein loaded on Protino® Ni-IDA Packed Columns or Resin, the higher will be the absolute yields. For example if loading ~ 75 mg 6 x His-GFPuv (concentration: 8 mg/mL) per g Protino® Ni-IDA Resin, a maximum yield of ~ 50 mg/g can be obtained.

3.4 Culture size

As outlined above, the protein yield depends on various parameters. However, some recommendations on protein load and culture size can be given as a starting point.

Note that yield and protein load are exemplified for the 6 x His-GFPuv (~ 32 kDa) and may vary from protein to protein.

- Use rather high concentrations of the target protein in the sample.
- For maximum yield use an excess amount of polyhistidine-tagged protein in the loaded sample. For example apply up to 75 mg of anticipated 6 x His-GFPuv per 1 g of Protino® Ni-IDA Resin.
- For maximum recovery use up to 20 mg of 6 x His-GFPuv per 1 g of Protino® Ni-IDA Resin.

The concentration of the polyhistidine-tagged protein in the culture may vary from < 1 mg/L up to 200 mg/L depending on cell density and expression level. It is recommended to determine the protein concentration for each expression experiment, for example via SDS-PAGE. On average, 250 mL of culture will produce approximately 1 g of pelleted, wet cells.

- Transfer the cell lysate from a 200 – 750 mL (high expression at 100 mg/L) or 2000 – 7500 mL (low expression at 10 mg/L) *E.coli* culture to 1 g of Protino® Ni-IDA Resin.

In order to obtain highly concentrated lysates, lyse wet cells in 2–5 mL LEW Buffer per 1 g wet mass. The volume of LEW Buffer should be adjusted according to the amount of polyhistidine-tagged protein in the culture. For example, 1 g cells may be resuspended in 2–5 mL LEW Buffer if a protein is expressed at 50–200 mg/L. For cultures with lower target protein content 1 g cells should be resuspended in 2 mL of LEW Buffer.

For recovering polyhistidine-tagged protein from *E.coli* cultures we recommend treatment with lysozyme in combination with sonication. If you are purifying recombinant protein from eukaryotic cells, treat the cells with an appropriate buffer containing a mild detergent (Sambrook *et al.*, 1989).

Table 4: Determination of culture and buffer volume requirements

Concentration of HisTag protein in the culture	Results in	protein load [mg]	Recommended <i>E.coli</i>		Native LEW Buffer ² (resuspension) [mL]	Denaturing		
			culture volume [mL]	pelet wet mass ¹ [g]		LEW Buffer ² (resuspension) [mL]	LEW Buffer ² (washing of IB ⁴) [mL]	DS Buffer ³ (lysis of IB ⁴) [mL]
Protino® Ni-IDA Resin, 1 g								
high, ~ 100 mg/L	Recovery _{max}	20	200	0.8	4	4	8	
	Yield _{max}	75	750	3	15	15	30	
low, ~ 10 mg/L	Recovery _{max}	20	2000	8	16	16	80	
	Yield _{max}	75	7500	30	60	60	300	
Protino® Ni-IDA 150 Packed Columns (40 mg Resin)								
high, ~ 100 mg/L	Recovery _{max}	0.8	8	0.04	0.2	0.2	0.4	
	Yield _{max}	3	30	0.12	0.6	0.6	1.2	
low, ~ 10 mg/L	Recovery _{max}	0.8	80	0.32	0.64	0.64	3.2	
	Yield _{max}	3	300	1.2	2.4	2.4	12	
Protino® Ni-IDA 1000 Packed Columns (250 mg Resin)								
high, ~ 100 mg/L	Recovery _{max}	5	50	0.2	1	1	2	
	Yield _{max}	18.8	188	0.75	15	15	7.5	
low, ~ 10 mg/L	Recovery _{max}	5	500	2	4	4	20	
	Yield _{max}	18.8	1880	7.5	24	24	75	
Protino® Ni-IDA 2000 Packed Columns (500 mg Resin)								
high, ~ 100 mg/L	Recovery _{max}	10	100	0.4	5	5	4	0.8
	Yield _{max}	37.5	375	1.5	30	30	15	4.8
low, ~ 10 mg/L	Recovery _{max}	10	1000	4	8	8	40	8
	Yield _{max}	37.5	3750	15	48	48	150	48

¹ On average, 250 mL of culture will produce approximately 1 g of pelleted, wet cells.² Lysis-Equilibration-Wash Buffer³ Denaturing Solubilization Buffer⁴ Inclusion Bodies




Figure 4 Required culture volumes and LEW Buffer volumes for maximum protein recovery or maximum yield in dependence on protein expression level. **The volumes noted below are only exemplary and are shown for 1 g of Protino® Ni-IDA Resin (purification under native conditions).** Please use these recommendations as a starting point to evaluate optional purification results. Note that purification conditions have to be optimized for each individual polyhistidine-tagged protein.

3.5 Binding, washing, and elution

In contrast to traditional IDA matrices, Protino® Ni-IDA shows an optimized, low density of IDA ligands, which is created by a special manufacturing process. This non-saturating surface concentration of IDA eliminates non-specific interactions of contaminating proteins with the adsorbent. Since virtually no contaminating host proteins bind to Protino® Ni-IDA, stringent washing procedures are generally not necessary. Therefore LEW Buffer, which is used for lysis, equilibration, and washing, does not contain any imidazole. Bound polyhistidine-tagged protein can competitively be eluted by adding imidazole. The recommended Elution Buffer contains 250 mM imidazole in order to recover even strong binding, multimeric proteins with more than one polyhistidine tag (also see buffer compositions section 5.1 and 6.1). However, as shown in Figure 5, depending on the protein, elution may be equally effective in the presence of much lower imidazole concentrations. If, for example, the stability or integrity of the target protein in 250 mM imidazole is a concern the concentration of imidazole in the eluent may readily be reduced.




Figure 5 Purification of polyhistidine-tagged GFPuv using Protino® Ni-IDA and Ni-IDA Agarose

Recombinant GFPuv was expressed in E.coli, lysed, loaded onto each gravity flow column, and eluted by a stepwise imidazole gradient. Eluted fractions were analyzed by SDS-PAGE. Pure polyhistidine-tagged protein can be eluted from Protino® Ni-IDA (left panel) at much lower imidazole concentrations than from Ni-IDA Agarose (right panel). In addition, Ni-IDA Agarose releases contaminating proteins from 10 mM to 100 mM imidazole.

M = Marker proteins, CL = Cleared lysate.

3.6 Compatibility of reagents

Buffer components that chelate metal ions, such as EDTA and EGTA, should not be used since they strip Ni^{2+} ions from the matrix.

Do not use buffers with $\text{pH} > 8.4$, since silica dissolves in solutions of high pH .

Table 5: Reagent compatibility chart

Reagent	Effect	Comments
Sodium phosphate	Used in LEW and Elution Buffer in order to buffer the solutions at $\text{pH} 8$	50 mM is recommended. The pH of any buffer should be adjusted to 8, although in some cases a pH between 7 and 8 can be used
Tris	Coordinates with Ni^{2+} ions, causing a decrease in capacity	10 mM may be used, sodium phosphate buffer is recommended
Sodium Chloride	Prevents ionic interactions and therefore unspecific binding	Up to 2 M can be used, at least 0.3 M should be used
Imidazole	Binds to immobilized Ni^{2+} ions and competes with the polyhistidine-tagged proteins	Should not be included in LEW Buffer
Urea	Solubilizes protein	Use 8 M for purification under denaturing conditions
GuHCl	Solubilizes protein	Up to 6 M can be used
β -mercaptoethanol	Prevents formation of disulfide bonds; Can reduce Ni^{2+} ions at higher concentrations	Up to 50 mM in samples has been used successfully in some cases
DTT, DTE	Can reduce Ni^{2+} ions at higher concentrations	Up to 10 mM in samples has been used successfully in some cases
Glutathione reduced	Can reduce Ni^{2+} ions at higher concentrations	Up to 30 mM in samples has been used successfully in some cases
Glycerol	Prevents hydrophobic interactions between proteins	Up to 50 % can be used
EDTA	Coordinates with Ni^{2+} ions, causing a decrease in capacity at higher concentrations	Not recommended, but up to 1 mM in samples has been used successfully in some cases

Table 5: Reagent compatibility chart

Ethanol	Prevents hydrophobic interactions between proteins	Up to 20 % can be used; Ethanol may precipitate proteins, causing low flow rates and column clogging
SDS	Interacts with Ni ²⁺ ions, causing a decrease in capacity	Not recommended, but up to 0.5 % in samples has been used successfully in some cases
Triton, Tween	Removes background proteins	Up to 2 % can be used

4 Safety instructions

When working with the **Protino® Ni-IDA** kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles).

For more information consult the appropriate Material Safety Data Sheets (MSDS available online at www.mn-net.com/msds).



4.1 Disposal

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

5.2 Preparation of cleared lysates under native conditions

- 1 Refer to Table 3, section 3.4 for detailed information on culture and buffer volume requirements**

Thaw the cell pellet from an *E.coli* expression culture on ice (if frozen). Resuspend **1 g of pelleted, wet cells** in **2–5 mL LEW Buffer** (for details see section 3.4). Pipette up and down, or use stirring until complete resuspension without visible cell aggregates. Perform this step on ice.

- 2 Add lysozyme** to a final concentration of 1 mg/mL. Stir the solution **on ice for 30 min**.
-

- 3** Sonicate the suspension on ice according to the instructions provided by the manufacturer (e.g., use 10 × 15 s bursts with a 15 s cooling period between each burst).

Carefully check samples' appearance after sonication. If the lysate is still viscous from incomplete fragmentation of DNA, add 5 µg/mL DNase and stir on ice for 15 min.

- 4** *Centrifuge the crude lysate at 10,000 × g for 30 min at 4 °C to remove cellular debris. Carefully transfer the supernatant to a clean tube without disturbing the pellet.*

If the supernatant is not clear, centrifuge a second time or filter through a 0.45 µm membrane (e.g., cellulose acetate) to avoid clogging of the IMAC column with insoluble material.

Store supernatant on ice.

Proceed to section 5.3, 5.4, 5.5, 5.6, 5.7, or 5.8.

5.3 Protino® Ni-IDA Packed Columns – purification under native conditions

5.3.1 Protocol-at-a-glance

Note: This section only refers to Protino® Ni-IDA 150, 1000, and 2000 Packed Columns.

		Protino® Ni-IDA Packed Columns		
		150	1000	2000
1	Cultivate and harvest cells	4,500– 6,000 x g 15 min at 4 °C	4,500– 6,000 x g 15 min at 4 °C	4,500– 6,000 x g 15 min at 4 °C
2	Preparation of working solutions (per column)			
	8 x LEW Buffer	0.3 mL	1.5 mL	2.5 mL
	+water	+2.1 mL	+10.5 mL	+17.5 mL
	= 1 x LEW Buffer	= 2.4 mL	= 12 mL	= 20 mL
	4 x Elution Buffer	0.25 mL	1.5 mL	3 mL
	+water	+0.75 mL	+4.5 mL	+9 mL
	= 1 x Elution Buffer	= 1 mL	= 6 mL	= 12 mL
3	Cell extract preparation			
	Refer to section 5.2.			
4	Column Equilibration			
	1 x LEW Buffer	320 µL	2 mL	4 mL
5	Binding	Load clarified lysate onto the column	Load clarified lysate onto the column	Load clarified lysate onto the column
6	Washing			
	1 x LEW Buffer	2 x 320 µL	2 mL	4 mL
7	Elution			
	1 x Elution Buffer	3 x 240 µL	3 x 1.5 mL	3 x 3 mL

5.3.2 Procedure

Note: Experienced users may refer to the protocol at a glance, page 23.

Protino® Ni-IDA Packed Columns

	150	1000	2000
1 Cultivate and harvest cells			
Harvest cells from an <i>E.coli</i> expression culture by centrifugation at 4,500–6,000 x g for 15 min at 4 °C . Remove supernatant. Store cell pellet at -20 °C if not processed immediately.			
2 Preparation of working solutions			
Prepare 1 x LEW (Lysis / Equilibration / Wash) Buffer and 1 x Elution Buffer by diluting the supplied stock solutions.			
<i>Note: If precipitate is observed in the stock solutions, warm and shake them to dissolve precipitate prior to diluting the buffers.</i>			
Mix 8 x LEW Buffer			
	0.3 mL	1.5 mL	2.5 mL
with deionized water			
	2.1 mL	10.5 mL	17.5 mL
to get a final volume of 1 x LEW Buffer sufficient for one column run.			
	2.4 mL	12 mL	10 mL
Mix 4 x Elution Buffer			
	0.25 mL	1.5 mL	3 mL
with deionized water			
	0.75 mL	4.5 mL	9 mL
to get a final volume of 1 x LEW Buffer sufficient for one column run.			
	2.4 mL	8 mL	12 mL
3 Cell Extract Preparation			
Refer to section 5.2. For detailed information on culture and buffer volumes for cell extract preparation also see Table 3, section 3.4.			

Protino® Ni-IDA Packed Columns

150

1000

2000

4 Column equilibration

Equilibrate Protino® Ni-IDA Packed Columns with **1 x LEW Buffer**. Allow the column to drain by gravity.

320 µL

2 mL

4 mL

Protino® Ni-IDA 150/ 1000 Packed Columns are designed to fit into most 15 mL conical centrifuge tubes (e.g., BD Falcon REF 352097) for convenient fraction collection.

5 Binding

Add the cleared lysate (see section 5.2) to the pre-equilibrated column and allow the column to drain by gravity.

6 Washing

Wash the column with **1 x LEW Buffer**. Allow the column to drain by gravity.

2 x 320 µL

2 x 2 mL

2 x 4 mL

7 Elution

Elute the polyhistidine-tagged protein in a new collecting tube by adding **1 x Elution Buffer**. Allow the column to drain by gravity

3 x 240 µL

3 x 1.5 mL

3 x 3 mL

*Note: Depending on protein characteristics 90% of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/ or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

5.4 Protino® Ni-IDA Resin – gravity-flow column chromatography under native conditions

For column IMAC using Protino® Ni-IDA Resin we generally recommend gravity-flow procedure. This protocol describes gravity-flow column chromatography using **Protino® Ni-IDA Resin** and **Protino® Columns** (see section 2.1).

Prepacked columns filled with 40 mg, 250 mg or 500 mg Protino® Ni-IDA Resin may readily be used (Protino® Ni-IDA 150, 1000, or 2000 Packed Columns, see section 5.3 and ordering information).

Note: When using other types of chromatography columns please note that the pore size of the filter frit should be around 50 µm to ensure appropriate flow rates.

1 Column preparation

Transfer the **appropriate amount** of **Protino® Ni-IDA Resin** to an empty **Protino® Column**. To achieve tight packing, gently tap the column on a hard surface until the bed height remains constant. Place a separate filter frit on top of the column bed by using a lab pen. Gently tap on the frit to ensure that there is no gap between column bed and filter frit.

1 g of Protino® Ni-IDA Resin will result in 2 mL bed volume.

The amount of resin required depends on the amount of polyhistidine-tagged protein to be purified. The binding capacity of Protino® Ni-IDA Resin varies from protein to protein. See section 3.3 for general guidelines.

2 Column equilibration

Equilibrate the column with **4 bed volumes** of **LEW Buffer**. Allow the column to drain by gravity.

3 Binding

Add the supernatant (cleared lysate, see section 5.2) to the pre-equilibrated column and allow the column to drain by gravity.

Apply at least 1.5 bed volumes of sample.

4 Washing

Wash the column twice with **4 bed volumes** of **LEW Buffer**. Allow the column to drain by gravity.

5 Elution

Elute the polyhistidine-tagged protein in three fractions. Add **3 × 3 bed volumes** of **Elution Buffer** and collect separately. Allow the column to drain by gravity.

*Note: Depending on protein characteristics 90% of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

5.5 Protino® Ni-IDA Resin – batch gravity-flow purification protocol under native conditions

Although we recommend gravity flow procedure, polyhistidine-tagged proteins may be purified by the following batch/gravity-flow protocol.

Note: Usually the yield is not significantly increased using time-consuming shaking.

1 Batch binding

Add the **appropriate amount** of **Protino® Ni-IDA Resin** directly to the cleared lysate (see section 5.2). Gently mix the material on an orbital shaker for 5 – 15 min. Do not use a magnetic stirrer to avoid generating fine particles through excessive physical force.

The amount of resin required depends on the amount of polyhistidine-tagged protein to be purified. The binding capacity of Protino® Ni-IDA Resin varies from protein to protein. See section 3.3 for general guidelines.

1 g of Protino® Ni-IDA Resin will result in 2 mL bed volume.

The length of time required for optimal binding will vary from protein to protein.

2 Transfer

Transfer the lysate-resin mixture to an empty chromatography column, for example Protino® Columns (see section 1.2 and ordering information). Let the resin settle by gravity flow.

3 Washing

Wash the column with **8 bed volumes** of **LEW Buffer**. Allow the column to drain by gravity.

4 Elution

Elute the polyhistidine-tagged protein in three fractions. Add **3 x 3 bed volumes** of **Elution Buffer** and collect separately. Allow the column to drain by gravity.

*Note: Depending on protein characteristics 90% of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

5.6 Protino® Ni-IDA Resin – batch purification protocol under native conditions

Although we recommend gravity flow procedure polyhistidine-tagged proteins may be purified by the following batch protocol.

1 Batch binding

Add the **appropriate amount** of **Protino® Ni-IDA Resin** directly to the cleared lysate (see section 5.2) filled in a centrifugation tube. Close the tube and mix the suspension gently, for example on an orbital shaker for 5 – 15 min.

The amount of resin required depends on the amount of polyhistidine-tagged protein to be purified. The binding capacity of Protino® Ni-IDA Resin varies from protein to protein. See section 3.3 for general guidelines.

1 g of Protino® Ni-IDA Resin will result in 2 mL bed volume.

The length of time required for optimal binding will vary from protein to protein.

Sediment the resin by gravity or centrifugation at **500 x g** for **1 min**. Carefully decant the supernatant and discard it.

2 Washing

Add **4 bed volumes** of **LEW Buffer** and mix for 5 min.

Sediment the resin by gravity or centrifugation at **500 x g** for **1 min**. Carefully decant the supernatant and dispose of it.

Repeat the washing step one or two more times (total wash **2–3 x 4 bed volumes** of **LEW Buffer**).

3 Elution

Add **3 bed volumes** of **Elution Buffer** and mix for 5 min.

Sediment the resin by gravity or centrifugation at **500 x g** for **1 min**.

Carefully decant or pipette the eluate in a new tube.

Repeat the elution step two more times (total elution **3 x 3 bed volumes** of **LEW Buffer**).

*Note: Depending on protein characteristics 90 % of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

5.7 Protino® Ni-IDA Resin – medium pressure column chromatography under native conditions

According to the physical stability of the Protino® Ni-IDA Resin it is ideally suited for medium pressure column chromatography. The rigid matrix can be run under high flow rates and high back pressure. Furthermore Protino® Ni-IDA Resin does not shrink or swell upon hydration.

As a starting point run columns at flow rates of 0.5–1.0 mL/min/cm². If the polyhistidine-tagged protein does not bind, further reduce the flow rate. Optimal flow rates have to be determined empirically, because dissociation rates vary widely from protein to protein.

1 Column preparation -slurry packing

Make a ~ 10 % (w/v) slurry of **Protino® Ni-IDA Resin** in degassed deionized water. Do not use a magnetic stirrer to avoid generating fine particles through excessive physical force. Slowly pour the suspension into the column. Avoid introducing air bubbles.

1 g of Protino® Ni-IDA Resin will result in 2 mL bed volume.

The amount of resin required depends on the amount of polyhistidine-tagged protein to be purified. The binding capacity of Protino® Ni-IDA Resin varies from protein to protein. See section 3.3 for general guidelines.

Allow the resin to settle. Insert and adjust top adapter and connect the column to the chromatography system according to the manufacturer's instructions. Avoid introducing air bubbles.

2 Column equilibration

Equilibrate the column with **4 bed volumes** of **LEW Buffer** or until the baseline at 280 nm is stable.

3 Binding

Apply the cleared lysate (see section 5.2) to the column.

4 Washing

Wash the column with **8 bed volumes** of **LEW Buffer** or until the baseline at 280 nm is stable.

Do not add imidazole to the LEW Buffer.

5 Elution

Elute the polyhistidine-tagged protein with **5–10 bed volumes** of **Elution Buffer** using a step gradient.

When monitoring protein elution note that imidazole absorbs at 280 nm.

*Note: Depending on protein characteristics 90 % of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

5.8 Protino® 96 Ni-IDA – purification under native conditions

1 Cultivate and harvest cells

E.coli cells may be cultivated in 96-well plates (1 mL culture volume), 24-well plates (5 mL culture volume), or any other appropriate cultivation vessel, if larger culture volumes are required.

Harvest cells by centrifugation at **2000 x g** for **15 min at 4 °C**. Store cell pellets at -20 °C or -70 °C for at least 1 h.

2 Preparation of working solutions

Prepare 1 x LEW (Lysis / Equilibration / Wash) Buffer and 1 x Elution Buffer by diluting the supplied stock solutions.

Note: If precipitate is observed in the stock solutions, warm and shake them to dissolve precipitate prior to diluting the buffers.

Mix **30 mL** of **8 x LEW Buffer**
with **210 mL** of **deionized water**
to get a **final volume** of **240 mL 1x LEW Buffer** sufficient for 96 preps.

Mix **20 mL** of **4 x Elution Buffer**
with **60 mL** of **deionized water**
to get a **final volume** of **80 mL** of **1 x Elution Buffer** sufficient for 96 preps.

3 Cell extract preparation

Use standard procedures for the preparation of cell extracts, such as lysozyme treatment, sonication, or detergent treatment. Note that optimal sample preparation steps have to be determined empirically depending on the characteristics of the of the polyhistidine-tagged protein and host organism.

Cell extracts from > 5 mL *E.coli* expression culture

For preparation of cell extracts from large volumes of *E.coli* expression culture (> 5 mL culture) refer to section 5.2.

Cell extracts from ≤ 5 mL *E.coli* expression culture

For preparation of cell extracts from up to 5 mL *E.coli* expression culture we recommend the following protocol as a starting point for further optimization.

Thaw cell pellets at room temperature. Resuspend each pellet in **1 mL** of **1 x LEW Buffer** containing **0.2 mg/mL lysozyme**. Incubate at room temperature for 30 min in a shaker. If the lysate is still viscous, add 15 U of Benzonase® per well, mix, and incubate at room temperature for 30 min. Benzonase® reduces lysate viscosity by rapidly hydrolysing DNA and RNA.

Centrifuge the crude lysate at **5,000 x g** for **30 min at 4 °C** to remove cellular debris.

If the supernatant is not clear, centrifuge a second time to avoid clogging of the Protino® 96 Ni-IDA Plate with insoluble material.

Store supernatant on ice.

4 Equilibration

Equilibrate the Protino® 96 Ni-IDA Plate by adding **500 µL** of **1x LEW Buffer** to each well. Allow the wells to drain by gravity.

5 Binding

Add the cleared lysate (supernatant, see step 3 or section 5.2) to the pre-equilibrated wells of the Protino® 96 Ni-IDA Plate. Allow the wells to drain by gravity.

Note: The lysate should not contain any imidazole, since most proteins do not bind to the resin in the presence of even low imidazole concentrations.

6 Washing

Wash each well twice by adding **800 µL** of **1x LEW Buffer**. Allow the wells to drain by gravity.

Note: Stringent washing procedures by using buffers with low imidazole concentrations are generally not necessary. Therefore LEW Buffer does not contain any imidazole.

7 Elution

Elute the polyhistidine-tagged protein in a suitable collecting plate (e.g., Rack of Tube Strips, see ordering information) by adding **3 × 250 µL** of **1x Elution Buffer**. Allow the wells to drain by gravity.

Note: To increase the concentration of the polyhistidine-tagged protein in the final eluate add 80 µL of Elution Buffer and discard the flow-through (this flow-through mainly contains LEW Buffer from the washing steps). Afterwards elute the protein by adding 3 × 250 µL of Elution Buffer.

*Note: Depending on protein characteristics 90% of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

6.2 Cell extract preparation under denaturing conditions

We recommend this protocol if expression leads to the formation of inclusion bodies. Cells are disrupted under native conditions using lysozyme together with sonication. After centrifugation the polyhistidine-tagged protein is extracted and solubilized from the pellet by using a denaturant (8 M urea). The extract obtained is clarified by centrifugation and applied to Protino® Ni-IDA Packed Columns or Protino® Ni-IDA Resin under denaturing conditions. Purification of polyhistidine-tagged proteins under denaturing conditions is similar to purification under native conditions except that the cell extract and buffers loaded on the column contain 8 M urea. For buffer compositions see section 6.1.

1 Isolation of inclusion bodies

Refer to **Table 3, section 3.4** for detailed information on culture and buffer volume requirements.

Thaw the cell pellet from an *E.coli* expression culture on ice (if frozen). Resuspend **1 g** of **pelleted, wet cells** in **5 mL LEW Buffer** (without denaturant) on ice (also see section 3.4). Pipette up and down, or use stirring until complete resuspension without visible cell aggregates.

Add **lysozyme** to a final concentration of 1 mg/mL. Stir the solution **on ice** for **30 min**.

Sonicate the suspension on ice according to the instructions provided by the manufacturer (e.g., use 10 × 15 s bursts with a 15 s cooling period between each burst).

Carefully check samples' appearance after sonication. If the lysate is still viscous from incomplete fragmentation of DNA, add 5 µg/mL DNase and stir on ice for 15 min.

Centrifuge the crude lysate at **10,000 x g** for **30 min at 4 °C** to collect the inclusion bodies. Discard supernatant. Keep pellet on ice.

2 Solubilization of inclusion bodies

Resuspend the pellet in **10 mL LEW Buffer** per g wet cells to wash the inclusion bodies.

Centrifuge the suspension at **10,000 x g** for **30 min at 4 °C**. Discard supernatant.

Resuspend the pellet in **2.0 mL Denaturing Solubilization Buffer per g wet cells** to solubilize the inclusion bodies. Homogenization or sonication may be necessary to resuspend the pellet. Dissolve the inclusion bodies by stirring on ice for 60 min.

Centrifuge at **10,000 x g** for **30 min at 20 °C** to remove any remaining insoluble material. Carefully transfer the supernatant to a clean tube without disturbing the pellet.

If the supernatant is not clear centrifuge a second time or filter through a 0.45 µm membrane (e.g., celluloseacetate) to avoid clogging of the IMAC column with insoluble material.

Save supernatant.

Proceed to section 6.3 or 6.4.

6.3 Protino® Ni-IDA Packed Columns – purification under denaturing conditions

Protino® Ni-IDA Packed Columns

150

1000

2000

1 Cell Extract Preparation

Refer to section 6.2. For detailed information on culture and buffer volumes for cell extract preparation also see Table 3, section 3.4.

2 Solubilization of inclusion bodies

Refer to section 6.2. For detailed information on culture and buffer volumes for cell extract preparation also see Table 1, section 3.4.

3 Column equilibration

Equilibrate Protino® Ni-IDA Packed Columns with **Denaturing Solubilization Buffer**. Allow the column to drain by gravity.

320 µL

2 mL

4 mL

Protino® Ni-IDA 150 / 1000 Packed Columns are designed to fit into most 15 mL conical centrifuge tubes (e.g., BD Falcon 352097) for convenient fraction collection.

4 Binding

Add the supernatant (solubilized protein, see section 6.2) to the pre-equilibrated column and allow the column to drain by gravity.

5 Washing

Wash the column with **Denaturing Solubilization Buffer**. Allow the column to drain by gravity.

6 Elution

Elute the polyhistidine-tagged protein in a new collecting tube by adding **Denaturing Elution Buffer**. Allow the column to drain by gravity.

3 × 240 µL

3 × 1.5 mL

3 × 3 mL

*Note: Depending on protein characteristics 90% of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

6.4 Protino® Ni-IDA Resin – gravity-flow column chromatography under denaturing conditions

For column IMAC using Protino® Ni-IDA Resin we generally recommend the gravity-flow procedure. This protocol describes gravity-flow column chromatography using Protino® Ni-IDA Resin and **Protino® Columns** (see section 1.2).

Prepacked columns filled with 40 mg, 250 mg or 500 mg Protino® Ni-IDA Resin readily be used (Protino® Ni-IDA 150, 1000, or 2000 Packed Columns, see section 6.3 and ordering information).

Note: When using other types of chromatography columns please note that the pore size of the filter frit should be around 50 µm to ensure appropriate flow rates.

1 Column preparation

Transfer the appropriate amount of **Protino® Ni-IDA Resin** to an empty **Protino® Column**. To achieve tight packing gently tap the column on a hard surface until the bed height remains constant. Place a separate filter frit on top of the column bed by using a lab pen. Gently tap on the frit to ensure that there is no gap between column bed and filter frit.

1 g of Protino® Ni-IDA Resin will result in 2 mL bed volume.

The amount of resin required depends on the amount of polyhistidine-tagged protein to be purified. The binding capacity of Protino® Ni-IDA Resin varies from protein to protein. See section 3.3 for general guidelines.

2 Column equilibration

Equilibrate the column with **4 bed volumes** of **Denaturing Solubilization Buffer**. Allow the column to drain by gravity.

3 Binding

Add the supernatant (solubilized proteome, see section 6.2) to the pre-equilibrated column and allow the column to drain by gravity.

4 Washing

Wash the column with **8 bed volumes** of **Denaturing Solubilization Buffer**. Allow the column to drain by gravity.

5 Elution

Elute the polyhistidine-tagged protein in three fractions. Add **3 x 3 bed volumes of Elution Buffer** and collect separately. Allow the column to drain by gravity.

*Note: Depending on protein characteristics 90% of the eluted protein can be found in the **first** elution fraction.*

Use protein assay and/or SDS-PAGE analysis to determine which fraction(s) contain(s) the majority of the polyhistidine-tagged protein.

7 Cleaning, recharging, and storage

Cleaning

After use, wash resin with 10 bed volumes of LEW Buffer and 10 bed volumes of deionized water. Wash with 2 bed volumes of 20 % ethanol and store resin in 20 % ethanol at 4 °C.

Recharging

Depending on the nature of the sample Protino® Ni-IDA Resin can be reused 3–5 times. Reuse should only be performed with identical polyhistidine-tagged proteins to avoid possible cross-contamination.

After the final elution step wash Protino® Ni-IDA Resin with the following solutions:

5 bed volumes of	100 mM EDTA, pH 8
5 bed volumes of	deionized water
5 bed volumes of	100 mM NiSO ₄
10 bed volumes of	deionized water

After equilibrating with LEW buffer the resin is ready for reuse.

Complete regeneration

If a complete regeneration is mandatory, wash resin with the following solutions:

2 bed volumes of	6 M GuHCl, 0.2 M acetic acid
5 bed volumes of	deionized water
3 bed volumes of	2 % SDS
5 bed volumes of	deionized water
5 bed volumes of	100 % EtOH
5 bed volumes of	deionized water
5 bed volumes of	100 mM EDTA pH 8
5 bed volumes of	deionized water
5 bed volumes of	100 mM NiSO ₄
10 bed volumes of	deionized water

8 Appendix

8.1 Troubleshooting

Problem	Possible cause and suggestions
Sample does not enter column bed	<i>Sample / lysate contains insoluble material</i>
	<ul style="list-style-type: none">If the sample is not clear use centrifugation or filtration (0.45 µm membrane) to avoid clogging of the IMAC column.
Protein does not bind to the resin	<i>Sample / lysate contains genomic DNA</i>
	<ul style="list-style-type: none">Lysate may remain viscous from incomplete shearing of genomic DNA after sonication. Add 5 µg/mL DNase and incubate on ice for 10 min.
Protein elutes with wash buffer	<i>Problems with vector construction</i>
	<ul style="list-style-type: none">Ensure that protein and tag are in frame.
Protein does not elute from column	<i>Incorrect binding conditions</i>
	<ul style="list-style-type: none">Check composition of buffers and verify pH 7–8. Ensure that there is no chelating or strong reducing reagent or imidazole present.
Protein does not elute from column	<i>Incorrect buffer composition</i>
	<ul style="list-style-type: none">Check composition of buffers and verify pH 7–8. Ensure that there are no chelating or strong reducing agents or imidazole present.
Protein does not elute from column	<i>Elution conditions are too mild.</i>
	<ul style="list-style-type: none">Increase concentration of imidazole.

Problem	Possible cause and suggestions
Unwanted proteins elute with polyhistidine-tagged protein	<i>Insufficient wash</i>
	<ul style="list-style-type: none">• Use larger volumes for washing step.
	<i>Binding and wash conditions are too mild</i>
	<ul style="list-style-type: none">• Add small amounts of imidazole (1 – 10 mM). Verify that the imidazole concentration is low enough to bind the polyhistidine-tagged protein.
	<i>Contaminating proteins and target protein are linked together via disulfide bonds</i>
	<ul style="list-style-type: none">• Add up to 30 mM β-mercaptoethanol to reduce disulfide bonds.
	<i>Contaminating proteins are proteolytic products of target protein</i>
	<ul style="list-style-type: none">• Perform cell lysis at 4 °C.• Include protease inhibitors.
	<i>Expression is too low</i>
	<i>Contaminating host proteins have a better chance to bind to the resin when only small amounts of target protein are present in the lysate. Very low amounts of polyhistidine-tagged protein are not able to replace the majority of contaminating proteins effectively.</i>
<ul style="list-style-type: none">• Increase expression level.• Increase amount of starting cell material.• Do not exceed recommended lysis volumes.	

8.2 Ordering information

Product	REF	Pack of
Protino® Ni-IDA Resin	745210.5	5 g
	745210.30	30 g
	745210.120	120 g
	745210.600	600 g
Protino® Ni-IDA 150 Packed Columns (contains 40 mg of resin each)	745150.10	10 preps
	745150 .50	50 preps
Protino® Ni-IDA 1000 Packed Columns (contains 250 mg of resin each)	745160.10	10 preps
	745160.50	50 preps
Protino® Ni-IDA 2000 Packed Columns (contains 500 mg of resin each)	745170.5	5 preps
	745170.25	25 preps
Protino® Ni-TED Resin	745200.5	5 g
	745200.30	30 g
	745200.120	120 g
	745200.600	600 g
Protino® Ni-TED 150 Packed Columns (contains 40 mg of resin each)	745100.10	10 preps
	745100.50	50 preps
Protino® Ni-TED 1000 Packed Columns (contains 250 mg of resin each)	745110.10	10 preps
	745110.50	50 preps
Protino® Ni-TED 2000 Packed Columns (contains 500 mg of resin each)	745120.5	5 preps
	745120.25	25 preps
Protino® Columns 14 mL (empty gravity flow columns for use with e.g. Protino® Ni-IDA Resin)	745250.10	10 columns
Protino® Columns 35 mL (empty gravity flow columns for use with e.g. Protino® Ni-IDA Resin)	745255.10	10 columns
NucleoBond® Rack Small (for Protino® Ni-TED /IDA 150 Packed Columns)	740562	1
NucleoBond® Rack Large (for Protino® Ni-TED /IDA 1000 /2000 Packed Columns)	740563	1
Rack of MN Tube Strips (1 rack, 12 strips with 8 tubes each, and 12 Cap Strips)	740637	1 set

8.3 Product use restriction / warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety instructions from the respective Material Safety Data Sheet of the product.

This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY-NAGEL literature. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY-NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.

Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:
MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-333
support@mn-net.com


Trademarks:

Benzonase® is a registered trademark of Merck KGaA, Darmstadt, Germany.


FPLC™ is a trademark of Amersham Biosciences Corp.

Protino® is a registered trademark of MACHEREY-NAGEL GmbH & Co KG

All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment). Regarding these products or services we can not grant any guarantees regarding selection, efficiency, or operation.



Plasmid DNA
Clean up
RNA
DNA
Viral RNA and DNA
Protein
High throughput
Accessories
Auxiliary tools



MACHEREY-NAGEL

www.mn-net.com

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

