


MACHEREY-NAGEL

User manual





















- NucleoSpin® totalRNA FFPE
- NucleoSpin® totalRNA FFPE XS

November 2023 / Rev. 9

Total RNA from FFPE samples

Protocol at a glance (Rev.09) – page 1

	NucleoSpin® totalRNA FFPE	NucleoSpin® totalRNA FFPE XS
1 Sample preparation	Insert FFPE section(s) in a microcentrifuge tube	Insert FFPE section(s) in a microcentrifuge tube
2 Deparaffinize sample	 1 mL Paraffin Dissolver 56 °C, 5 min Vortex hot sample 16,000 x g, 2 min	 1 mL Paraffin Dissolver 56 °C, 5 min Vortex hot sample 16,000 x g, 2 min
	 170 µL MLF 16,000 x g, 2 min	 140 µL MLF 16,000 x g, 2 min
	Remove Paraffin Dissolver	Remove Paraffin Dissolver
3 Lyse sample A Quick protocol (perform method 3A or 3B)	 15 µL Proteinase K Mix gently 56 °C, 15 min	 12 µL Proteinase K Mix gently 56 °C, 15 min
	 15 µL MKA Vortex 0 °C, 5 min 16,000 x g, 5 min	 12 µL MKA Vortex 0 °C, 5 min 16,000 x g, 5 min
	 Transfer sample 80 °C, 15 min	 Transfer sample 80 °C, 15 min
Lyse sample B Protocol for difficult to lyse cells (perform method 3A or 3B)	 15 µL Proteinase K Mix gently 56 °C, 90 min	 12 µL Proteinase K Mix gently 56 °C, 90 min
	 15 µL MKA Vortex 0 °C, 5 min 16,000 x g, 5 min	 12 µL MKA Vortex 0 °C, 5 min 16,000 x g, 5 min
	 Transfer sample	 Transfer sample
4 Adjust binding conditions	 500 µL MX Vortex RT, 1 min	 400 µL MX Vortex RT, 1 min

Total RNA from FFPE samples

Protocol at a glance (Rev.09) – page 2

















	NucleoSpin® totalRNA FFPE		NucleoSpin® totalRNA FFPE XS	
5 Bind RNA		Load sample 16,000 x g, 15 s		Load sample 16,000 x g, 15 s
6 Wash and dry silica membrane	1st 	700 µL MW2 16,000 x g, 15 s	1st 	400 µL MW2 16,000 x g, 15 s
	2nd 	250 µL MW2 16,000 x g, 1 min	2nd 	200 µL MW2 16,000 x g, 1 min
7 Optional: Digest DNA		50 µL rDNase RT, 15 min		25 µL rDNase RT, 15 min
		100 µL MX RT, 1 min 16,000 x g, 15 s		50 µL MX RT, 1 min 16,000 x g, 15 s
	1st 	700 µL MW2 16,000 x g, 15 s	1st 	400 µL MW2 16,000 x g, 15 s
	2nd 	250 µL MW2 16,000 x g, 1 min	2nd 	200 µL MW2 16,000 x g, 5 min
8 Elute highly pure RNA		30–50 µL RNase-free H ₂ O RT, 1 min 16,000 x g, 1 min		5–30 µL RNase-free H ₂ O RT, 1 min 16,000 x g, 1 min

Table of contents

1 Components	4
1.1 Kit contents	4
1.2 Reagents, consumables, and equipment to be supplied by user	6
1.3 About this user manual	6
2 Product description	7
2.1 The basic principle	7
2.2 Kit specifications	8
3 Storage conditions and preparation of working solutions	10
4 Safety instructions	12
4.1 Disposal	12
5 Protocols	13
5.1 NucleoSpin® totalRNA FFPE	13
5.2 NucleoSpin® totalRNA FFPE XS	17
5.3 DNA digestion in the RNA eluates	21
6 Appendix	22
6.1 Troubleshooting	22
6.2 Ordering information	25
6.3 References	26
6.4 Product use restriction/warranty	27

1 Components

1.1 Kit contents

NucleoSpin® totalRNA FFPE			
REF	10 preps 740982.10	50 preps 740982.50	250 preps 740982.250
Paraffin Dissolver	15 mL	60 mL	3 × 100 mL
Lysis Buffer MLF	10 mL	10 mL	50 mL
Precipitation Buffer MKA	1 mL	1 mL	10 mL
Binding Buffer MX	13 mL	60 mL	250 mL
Reaction Buffer for rDNase	7 mL	7 mL	30 mL
Wash Buffer MW2 (Concentrate)*	6 mL	25 mL	100 mL
RNase-free H ₂ O	13 mL	13 mL	13 mL
rDNase, RNase-free (lyophilized)*	1 vial (size A)	1 vial (size C)	5 vials (size C)
Liquid Proteinase K	0.6 mL	0.8 mL	5 mL
NucleoSpin® RNA Columns (light blue rings)	10	50	250
Collection Tubes (2 mL)	10	50	250
Collection Tubes (1.5 mL)	10	50	250
User manual	1	1	1

* For preparation of working solutions and storage conditions see section 3.

Kit contents *continued*

NucleoSpin® totalRNA FFPE XS			
REF	10 preps 740969.10	50 preps 740969.50	250 preps 740969.250
Paraffin Dissolver	15 mL	60 mL	3 × 100 mL
Lysis Buffer MLF	10 mL	10 mL	50 mL
Precipitation Buffer MKA	1 mL	1 mL	10 mL
Binding Buffer MX	13 mL	60 mL	250 mL
Reaction Buffer for rDNase	7 mL	7 mL	30 mL
Wash Buffer MW2 (Concentrate)*	6 mL	25 mL	100 mL
RNase-free H ₂ O	13 mL	13 mL	13 mL
rDNase, RNase-free (lyophilized)*	1 vial (size A)	1 vial (size C)	3 vials (size C)
Liquid Proteinase K	0.2 mL	0.8 mL	3 × 1.25 mL
NucleoSpin® RNA FFPE XS Columns (light blue rings)	10	50	250
Collection Tubes (2 mL)	10	50	250
Collection Tubes (1.5 mL)	10	50	250
User manual	1	1	1

* For preparation of working solutions and storage conditions see section 3.

1.2 Reagents, consumables, and equipment to be supplied by user

Reagents

- 96–100 % ethanol (preferably undenaturated ethanol)

Consumables

- 1.5 mL microcentrifuge tubes (for sample lysis)
- Disposable pipette tips

Equipment

- Manual pipettors
- Centrifuge for microcentrifuge tubes
- Vortex mixer
- Thermal heating-block (adjustable to 56 °C and 80 °C)
- Personal protection equipment (lab coat, gloves, goggles)

1.3 About this user manual

It is strongly recommended that first time users of **NucleoSpin® totalRNA FFPE (XS)** kit read the detailed protocol sections of this user manual. Experienced users, however, may refer to the Protocol at a glance instead. The Protocol at a glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure.

All technical literature is available on the internet at www.mn-net.com.

Please contact Technical Service regarding information about changes of the current user manual compared to previous revisions.

2 Product description

Formalin-fixed and paraffin-embedded (FFPE) tissue is commonly used in histopathological analysis. Recently, there is more and more interest in also investigating DNA modifications, RNA expression or miRNA profiles of old, archived FFPE samples. However, fixation, embedding and storage lead to crosslinking and fragmentation of RNA. Especially crosslinks cannot be detected by standard quality control assays, such as gel electrophoresis, spectrophotometry or microfluidic analysis, but the efficiency of enzymatic reactions is significantly reduced, for example in RT-PCR. Standard RNA purification procedures do not remove these chemical modifications and therefore result in low RNA yield or poor downstream application performance. The **NucleoSpin® totalRNA FFPE (XS)** procedure implements buffers and procedural steps to efficiently decrosslink nucleic acids and yield high quality RNA for the most demanding applications.

2.1 The basic principle

The **NucleoSpin® totalRNA FFPE (XS)** kits provide a convenient, reliable, and fast method to isolate RNA from formalin-fixed, paraffin-embedded (FFPE) tissue. Paraffin Dissolver (patent pending) enables an efficient deparaffinization of the sample.

The tissue sample is then heat incubated with Proteinase K to digest the fixed tissue, release nucleic acids, and gently remove crosslinks. Optimal binding conditions for even small RNA (e.g., miRNA) are adjusted and the lysate is applied to the **NucleoSpin® RNA Column/ NucleoSpin® RNA FFPE XS Column**. RNA is bound to the silica membrane. Residual DNA remaining on the membrane is removed by convenient on-column rDNase digestion. Washing steps remove salts, metabolites, and macromolecular cellular components. Pure RNA is finally eluted in a small volume of RNase-free water, yielding highly concentrated RNA.

Nucleic acid preparation using **NucleoSpin® totalRNA FFPE (XS)** can be performed at room temperature. The eluate, however, should be treated with care. RNA is very sensitive to trace contaminations of RNases, often found on general lab ware, fingerprints, and dust. To ensure RNA stability keep RNA frozen at -20 °C for short-term or -70 °C for long-term storage.

2.2 Kit specifications

NucleoSpin® totalRNA FFPE (XS) is recommended for the isolation of total RNA from formalin-fixed, paraffin-embedded (FFPE) tissue samples, typically as thin sections (approx. 3–20 µm thickness). Formalin-fixed samples which are not embedded in paraffin can also be used as sample material by omitting the deparaffinization steps.

The sample size can be up to ~10 sections (1–20 µm) of FFPE. The amount of embedded tissue can be up to 50 mg for **NucleoSpin® totalRNA FFPE** or up to 5 mg for **NucleoSpin® totalRNA FFPE XS** (1 × 10 µm section with 1 cm² tissue is approximately 1 mg).

RNA yield strongly depends on sample type, quality, and amount. Furthermore, the procedures of tissue sampling, post sampling delay before fixation, fixation time, embedding, and storage conditions have a high impact on RNA quality and yield. For more details see, for example, Chung J.Y. et al. (2008); van Maldegem F. et al. (2008); von Ahlfen S. et al. (2007); Castiglione F. et al. (2007); Leyland-Jones B.R. et al. (2008).

RNA concentration: RNA can be eluted highly concentrated and ready to use in a small volume of 30–50 µL (**NucleoSpin® totalRNA FFPE**) or even 5–30 µL (**NucleoSpin® totalRNA FFPE XS**).

RNA size distribution: RNA isolated from formalin-fixed, paraffin-embedded tissue shows size distribution from 15 to 5,000 bases. Often short sized RNA from ca. 100–300 bases predominate, especially when the sample material is old. However, samples which were subjected to good tissue fixation, embedding, and storage conditions can yield RNA even larger than 5,000 bases.

RNA integrity: RNA Integrity Numbers (RIN) according to Agilent 2100 Bioanalyzer assays depend on sample type and quality. In general the quality of RNA extracted from FFPE samples is poor. Typical RIN of RNA isolated with NucleoSpin® totalRNA FFPE (XS) kits are in range of 2–6.

rDNase is supplied with the kit for a convenient removal of DNA by on-column digestion. For more demanding downstream applications, DNA can also be digested in solution as described in section 5.2.

Table 1: Kit specifications at a glance

Parameter	NucleoSpin® totalRNA FFPE	NucleoSpin® totalRNA FFPE XS
Technology	Silica membrane technology	Silica membrane technology
Format	Mini spin columns	Mini spin columns – XS design
Sample material	Up to ~10 sections with up to 50 mg of tissue	Up to ~10 sections with up to 5 mg of tissue
Typical yield	Strongly depends on sample quality and amount	Strongly depends on sample quality and amount
Elution volume	30 – 50 µL	5 – 30 µL
Preparation time	70 min/6 preps (90 min including optional rDNase digest)	70 min/6 preps (90 min including optional rDNase digest)
Use	For research use only	

3 Storage conditions and preparation of working solutions

Storage conditions:

- Store lyophilized rDNase at 4 °C upon arrival (stable until: see package label).
- All other kit components should be stored at room temperature (15 – 25 °C) and are stable until: see package label. Storage at lower temperatures may cause precipitation of salts. If precipitation occurs, incubate the bottle for several minutes at about 30 – 40 °C and mix well until the precipitate is redissolved.

Before starting any **NucleoSpin® totalRNA FFPE (XS)** protocol prepare the following:

- **RNase-free rDNase:** Add the indicated volume of Reaction Buffer for rDNase to the rDNase vial and incubate for 1 min at room temperature. Gently swirl the vial to completely dissolve the rDNase. Be careful not to mix the enzyme too vigorously as rDNase is sensitive to mechanical agitation. Dispense into aliquots and store at -20 °C. The frozen working solution is stable for at least 6 months. Do not freeze/thaw the aliquots more than three times.
- **Wash Buffer MW2:** Add the indicated volumes of 96 – 100% ethanol to the MW2 concentrate. Stored at room temperature, the buffer is stable for at least one year.

NucleoSpin® totalRNA FFPE			
REF	10 preps 740982.10	50 preps 740982.50	250 preps 740982.250
Wash Buffer MW2 (Concentrate)	6 mL Add 24 mL 96 – 100 % ethanol	25 mL Add 100 mL 96 – 100 % ethanol	100 mL Add 400 mL 96 – 100 % ethanol
RNase-free rDNase (lyophilized)	1 vial (size A) Add 0.75 mL Reaction Buffer for rDNase	1 vial (size C) Add 3 mL Reaction Buffer for rDNase	5 vials (size C) Add 3 mL Reaction Buffer for rDNase to each vial

NucleoSpin® totalRNA FFPE XS			
REF	10 preps 740982.10	50 preps 740982.50	250 preps 740982.250
Wash Buffer MW2 (Concentrate)	6 mL Add 24 mL 96 – 100 % ethanol	25 mL Add 100 mL 96 – 100 % ethanol	100 mL Add 400 mL 96 – 100 % ethanol
RNase- free rDNase (lyophilized)	1 vial (size A) Add 0.75 mL Reaction Buffer for rDNase	1 vial (size C) Add 3 mL Reaction Buffer for rDNase	3 vials (size C) Add 3 mL Reaction Buffer for rDNase to each vial

4 Safety instructions

When working with the **NucleoSpin® totalRNA FFPE (XS)** kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles). For more information consult the appropriate Material Safety Data Sheets (MSDS available online at www.mn-net.com/msds).



The waste generated with the **NucleoSpin® totalRNA FFPE (XS)** kit has not been tested for residual infectious material. A contamination of the liquid waste with residual infectious material is highly unlikely due to strong denaturing lysis buffer and Proteinase K treatment but it cannot be excluded completely. Therefore, liquid waste must be considered infectious and should be handled and discarded according local safety regulations.

4.1 Disposal

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

5 Protocols

5.1 NucleoSpin® totalRNA FFPE

Before starting the preparation:

Check that rDNase and Buffer MW2 were prepared according to section 3.

Set incubator(s) to 56 °C (for paraffin melting and lysis step) and 80 °C (for decrosslink step).

Please note that lysis step 3A is the standard method for most common sample materials while lysis step 3B is utilized for difficult to lyse cells

1 Sample preparation

Insert FFPE section(s) in a microcentrifuge tube (not provided with the kit).

2 Deparaffinize sample

Add **1 mL Paraffin Dissolver** to the sample.

Incubate **5 min** at **56 °C** (to melt the paraffin).

Vortex the hot sample.

Make sure that paraffin completely melts during the heat incubation step and mix well after melting to completely dissolve the paraffin.

Centrifuge sample for **2 min** at **16,000 x g**.



+ 1 mL Paraffin Dissolver

56 °C, 5 min

Vortex hot sample



16,000 x g, 2 min



***Attention:** Check Buffer MLF prior to use! If a white precipitate is visible, heat the buffer for several minutes at 30–40 °C until the precipitate is completely dissolved and mix thoroughly.*

Add **170 µL Buffer MLF**. Do not mix!

Centrifuge sample for **2 min** at **16,000 x g**.

Remove and discard Paraffin Dissolver by pipetting it off.

***Note:** Slight residues of Paraffin Dissolver do not affect the following steps.*



+ 170 µL MLF



16,000 x g, 2 min

Remove Paraffin Dissolver

3 A) Lyse sample – method A

(perform method 3 A or 3 B) "Quick protocol"

Add **15 µL Proteinase K**.

Mix by gently shaking or pipetting up and down. Do not vortex!

Incubate for **15 min** at **56 °C** to lyse sample tissue.

If tissue is still visible, continue incubation until sample is digested.



+ 15 µL Proteinase K

Mix gently

56 °C, 15 min

Add **15 µL Buffer MKA** and vortex briefly.

Incubate for **5 min** on **ice**.

Centrifuge for **5 min** at **16,000 x g**.



+ 15 µL MKA
Vortex
0 °C, 5 min



16,000 x g,
5 min

Transfer clear supernatant to a new 1.5 mL microcentrifuge tube (not provided).

Incubate at **80 °C** for exactly **15 min**.

Shorter incubation or lower temperature leads to incomplete decrosslinking. Longer incubation might lead to degradation of RNA.



Transfer sample
80 °C, 15 min

3 B) Lyse sample – method B

(perform method 3A or 3B) "Protocol for difficult to lyse samples"

Add **15 µL Proteinase K**.

Mix by gently shaking or pipetting up and down. Do not vortex!

Incubate for **90 min** at **56 °C** to lyse sample tissue.

If tissue is still visible, continue/increase incubation time up to overnight until sample is digested.



+ 15 µL Proteinase K
Mix gently
56 °C, 90 min

Add **15 µL Buffer MKA** and vortex briefly.

Incubate for **5 min** on **ice**.

Centrifuge for **5 min** at **16,000 x g**.



+ 15 µL MKA
Vortex
0 °C, 5 min



16,000 x g,
5 min

Transfer clear supernatant to a new 1.5 mL microcentrifuge tube (not provided).



Transfer sample

4 Adjust binding conditions

Add **500 µL Buffer MX** and mix by vortexing (2 × 5 s).

Incubate for **1 min** at **room temperature**.



+ 500 µL MX
Vortex
RT, 1 min

5 Bind RNA

Place a **NucleoSpin® RNA Column** in a new **Collection Tube (2 mL)**.

Load sample onto the column and centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



Load sample

**16,000 x g,
15 s**

6 Wash and dry silica membrane**1st wash**

Add **700 µL Buffer MW2** to the NucleoSpin® RNA Column.

Centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



+ 700 µL MW2

**16,000 x g,
15 s**

2nd wash

Add **250 µL Buffer MW2** to the NucleoSpin® RNA Column.

Centrifuge for **1 min** at **16,000 x g** to dry the membrane completely.

If the flowthrough in the collection tube has touched the NucleoSpin® RNA Column after 2nd wash, discard flowthrough and centrifuge again.



+ 250 µL MW2

**16,000 x g,
1 min**

7 Optional: Digest DNA

Add **50 µL rDNase** directly onto the silica membrane of the NucleoSpin® RNA Column.

Incubate at **room temperature** for **15 min**.



+ 50 µL rDNase

RT, 15 min

Add **100 µL Buffer MX**.

Incubate for **1 min** at **room temperature**.

Centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



+ 100 µL MX

RT, 1 min



**16,000 x g,
15 s**

1st wash

Add **700 µL Buffer MW2** to the NucleoSpin® RNA Column.

Centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



+ 700 µL MW2



**16,000 x g,
15 s**

2nd wash

Add **250 µL Buffer MW2** to the NucleoSpin® RNA Column.

Centrifuge for **1 min** at **16,000 x g** to dry the membrane completely.

If the flowthrough in the collection tube has touched the NucleoSpin® RNA Column after 2nd wash, discard flowthrough and centrifuge again.



+ 250 µL MW2



**16,000 x g,
1 min**

8 Elute highly pure RNA

Place the NucleoSpin® RNA Column in a new **Collection Tube (1.5 mL)**.

Add **30 µL** (for high concentration) or **50 µL** (for high yield) RNase-free H₂O to the column.

Incubate for **1 min** at room temperature.

Centrifuge for **1 min** at **16,000 x g**.

Keep the eluted RNA on ice or freeze at **-20 °C** (short-term storage) or **-70 °C** (long-term storage).



**+ 30–50 µL
RNase-free
H₂O**

RT, 1 min



**16,000 x g,
1 min**

5.2 NucleoSpin® totalRNA FFPE XS

Before starting the preparation:

Check that rDNase and Buffer MW2 were prepared according to section 3.

Set incubator(s) to 56 °C (for paraffin melting and lysis step) and 80 °C (for decrosslink step).

Please note that lysis step 3A is the standard method for most common sample materials while lysis step 3B is utilized for difficult to lyse cells

1 Sample preparation

Insert FFPE section(s) in a microcentrifuge tube (not provided with the kit).

2 Deparaffinize sample

Add **1 mL Paraffin Dissolver** to the sample.

Incubate **5 min** at **56 °C** (to melt the paraffin).

Vortex the hot sample.

Make sure that paraffin completely melts during the heat incubation step and mix well after melting to completely dissolve the paraffin.

Centrifuge sample for **2 min** at **16,000 x g**.



+ 1 mL Paraffin Dissolver

56 °C, 5 min

Vortex hot sample



16,000 x g, 2 min



***Attention:** Check Buffer MLF prior to use! If a white precipitate is visible, heat the buffer for several minutes at 30–40 °C until the precipitate is completely dissolved and mix thoroughly.*

Add **140 µL Buffer MLF**. Do not mix!

Centrifuge sample for **2 min** at **16,000 x g**.



+ 140 µL MLF



16,000 x g, 2 min

Remove and discard Paraffin Dissolver by pipetting it off.

***Note:** Slight residues of Paraffin Dissolver do not affect the following steps.*

Remove Paraffin Dissolver

3 A) Lyse sample – method A

(perform method 3 A or 3 B) "Quick protocol"

Add **12 µL Proteinase K**.

Mix by gently shaking or pipetting up and down. Do not vortex!



+ 12 µL Proteinase K

Mix gently

Incubate for **15 min** at **56 °C** to lyse sample tissue.

If tissue is still visible, continue incubation until sample is digested.

56 °C, 15 min

Add **12 µL Buffer MKA** and vortex briefly.

Incubate for **5 min** on **ice**.

Centrifuge for **5 min** at **16,000 x g**.



+ 12 µL MKA

Vortex

0 °C, 5 min



**16,000 x g,
5 min**

Transfer clear supernatant to a new 1.5 mL microcentrifuge tube (not provided).

Incubate at **80 °C** for exactly **15 min**.

Shorter incubation or lower temperature leads to incomplete decrosslinking. Longer incubation might lead to degradation of RNA.



**Transfer
sample**

80 °C, 15 min

3 B) Lyse sample – method B

(perform method 3A or 3B) "Protocol for difficult to lyse samples"

Add **12 µL Proteinase K**.

Mix by gently shaking or pipetting up and down. Do not vortex!

Incubate for **90 min** at **56 °C** to lyse sample tissue.

If tissue is still visible, continue / increase incubation time up to overnight until sample is digested.



**+ 12 µL
Proteinase K**

Mix gently

56 °C, 90 min

Add **12 µL Buffer MKA** and vortex briefly.

Incubate for **5 min** on **ice**.

Centrifuge for **5 min** at **16,000 x g**.



+ 12 µL MKA

Vortex

0 °C, 5 min



**16,000 x g,
5 min**

Transfer clear supernatant to a new 1.5 mL microcentrifuge tube (not provided).



**Transfer
sample**

4 Adjust binding conditions

Add **400 µL Buffer MX** and **mix** by vortexing (2 × 5 s).

Incubate for **1 min** at **room temperature**.



+ 400 µL MX

Vortex

RT, 1 min

5 Bind RNA

Place a **NucleoSpin® RNA FFPE XS Column** in a new **Collection Tube (2 mL)**.

Load sample onto the column and centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



Load sample

**16,000 x g,
15 s**

6 Wash and dry silica membrane

1st wash

Add **400 µL Buffer MW2** to the **NucleoSpin® RNA FFPE XS Column**.

Centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



+ 400 µL MW2

**16,000 x g,
15 s**

2nd wash

Add **200 µL Buffer MW2** to the **NucleoSpin® RNA FFPE XS Column**.

Centrifuge for **1 min** at **16,000 x g** to dry the membrane completely.

If the flowthrough in the collection tube has touched the NucleoSpin® RNA FFPE XS Column after 2nd wash, discard flowthrough and centrifuge again.



+ 200 µL MW2

**16,000 x g,
1 min**

7 Optional: Digest DNA

Add **25 µL rDNase** directly onto the silica membrane of the **NucleoSpin® RNA FFPE XS Column**.

Incubate at **room temperature** for **15 min**.



**+ 25 µL rDNase
RT, 15 min**

Add **50 µL Buffer MX**.

Incubate for **1 min** at **room temperature**.

Centrifuge for **15 s** at **16,000 x g**.

Discard flowthrough and place the column back into the collection tube.



**+ 50 µL MX
RT, 1 min**

**16,000 x g,
15 s**

1st wash

Add 400 µL Buffer MW2 to the NucleoSpin® RNA FFPE XS Column.

Centrifuge for 15 s at 16,000 x *g*.

Discard flowthrough and place the column back into the collection tube.



+ 400 µL MW2

16,000 x *g*,
15 s

2nd wash

Add 200 µL Buffer MW2 to the NucleoSpin® RNA FFPE XS Column.

Centrifuge for 1 min at 16,000 x *g*.

Discard flowthrough and place the column back into the collection tube.

Centrifuge for 5 min at 16,000 x *g* to dry the membrane.



+ 200 µL MW2

16,000 x *g*,
1 min



16,000 x *g*,
5 min

8 Elute highly pure RNA

Place the NucleoSpin® RNA FFPE XS Column in a new **Collection Tube (1.5 mL)**.

Add **5 µL** (for high concentration) to **30 µL** (for high yield) **RNase-free H₂O** to the column.

Incubate for **1 min** at **room temperature**.

Centrifuge for **1 min** at **16,000 x *g***.

Keep the eluted RNA on ice or freeze at -20 °C (short-term storage) or -70 °C (long-term storage).



+ 5–30 µL
RNase-free
H₂O

RT, 1 min



16,000 x *g*,
1 min

5.3 DNA digestion in the RNA eluates

Comments on DNA removal:

Although the on-column rDNase digest in the standard protocol is very efficient, there are still certain applications which require even lower quantities of residual DNA.

For example, RT-PCR reactions with primers that do not differentiate between cDNA (derived from RNA) and contaminating genomic DNA.

DNA digestion in solution can efficiently degrade contaminating DNA. This requires stringent RNase control and optionally repurification of the RNA (in order to remove buffer, salts, DNase, and digested DNA).

The high quality, RNase-free, recombinant DNase (rDNase) provided with the kit facilitates such a digestion in solution.

A Digest DNA (Reaction setup)

Add **1/10 volume of rDNase** (dissolved in Reaction Buffer for rDNase) to the eluted RNA (e.g., add 3 µL enzyme to 30 µL RNA eluate).

B Incubate for 10 min at 37 °C.

C Inactivate rDNase

Incubate the sample for **5 min at 75 °C** to inactivate the rDNase. Put the sample on ice.

In most cases a further purification (in order to remove inactivated rDNase, buffer and salts) is not necessary. If nevertheless a repurification is required, NucleoSpin® RNA Clean-up XS is recommended (see section 6.2 ordering information).

6 Appendix

6.1 Troubleshooting

Problem	Possible cause and suggestions
Poor RNA quality or yield	<i>RNase contamination</i>
	<ul style="list-style-type: none"> • Create an RNase-free working environment. Wear gloves during all steps of the procedure. Change gloves frequently. Use sterile, disposable polypropylene tubes and filter tips. Keep tubes closed whenever possible during the preparation. Glassware should be oven-baked for at least 2 hours at 250 °C before use.
	<i>Poor sample quality</i>
	<ul style="list-style-type: none"> • Sample quality very much influences the obtainable RNA amount and quality. For aspects concerning sample harvest, fixation, embedding, and storage refer to: Castiglione F. et al. (2007), Chung J.Y. et al. (2008), Leyland-Jones B.R. et al. (2008), von Ahlfsen S. et al. (2007), von Maldegem F. et al. (2008).
	<i>Reagents not applied or restored properly</i>
	<ul style="list-style-type: none"> • Always dispense exactly the buffer volumes given in the protocols! • Always follow the given instructions closely, with specific attention paid to order and mode of mixing (shaking, vortexing, etc). • Add the indicated volume of 96 – 100 % ethanol to Buffer MW2 Concentrate and mix thoroughly. • Store kit components at room temperature. Storage at lower temperatures may cause salt precipitation. If precipitation occurs, incubate the bottle for several minutes at about 30 – 40 °C and mix well until the precipitate is redissolved. • Keep bottles tightly closed in order to prevent evaporation or contamination.
	<i>Ionic strength and pH influence A_{260} absorption as well as ratio A_{260}/A_{280}</i>
	<ul style="list-style-type: none"> • For absorption measurement, use 5 mM Tris pH 8.5 as diluent. Please also see: Manchester K.L. (1995) and Wilfinger W.W. et al. (1997).
	<i>Proteinase digestion time</i>
	<ul style="list-style-type: none"> • Depending of the nature of the sample, an optimal digestion time from 15 min to 3 hours has to be determined empirically. If tissue residues are still visible after 15 min continue the incubation for up to 3 hours. If a large portion of the sample still remains undigested, continue digestion overnight. An overnight incubation is not recommended if the tissue digested well within 3 hours.

Contami- nation of RNA with genomic DNA	<p><i>rDNase not active</i></p> <ul style="list-style-type: none">• Reconstitute and store lyophilized rDNase according to instructions given in section 3. <p><i>rDNase solution not properly applied</i></p> <ul style="list-style-type: none">• Pipette rDNase solution directly onto the center of the silica membrane and close the lid in order to press the solution into the membrane. <p><i>Too much cell material used</i></p> <ul style="list-style-type: none">• Reduce quantity of cells or tissue used.• Use larger PCR targets (e.g., > 500 bp) or intron spanning primers for RNA analysis.• Use support protocol for subsequent rDNase digestion in the eluate (section 5.3).
Suboptimal performance of RNA in downstream experiments	<p><i>Carry-over of ethanol or salt</i></p> <ul style="list-style-type: none">• Do not let the column flowthrough touch the column outlet after the second Buffer MW2 wash. Be sure to centrifuge at the recommended speed and time in order to remove ethanolic Buffer MW2 completely.• Check that Buffer MW2 has been equilibrated to room temperature before use. Washing at lower temperatures lowers efficiency of salt removal.• Depending on the robustness of the RT-PCR system used, RT-PCR might be inhibited if complete eluates are used as template for RT-PCR. Use less eluate as template. <p><i>Store isolated RNA properly</i></p> <ul style="list-style-type: none">• Eluted RNA should always be kept on ice for optimal stability since trace contaminations of omnipresent RNases (general lab ware, fingerprints, dust) will degrade the isolated RNA. For short term storage freeze at -20 °C, for long term storage freeze at -70 °C.

Silica abrasion from the membrane

Discrepancy
between A_{260}
quantification
values
and PCR
quantification
values

- Due to the typically low RNA content in very small samples and the resulting low total amount of isolated RNA, a RNA quantification via A_{260} absorption measurement is often hampered due to the low sensitivity of the absorption measurement. When performing absorption measurements close to the detection limit of the photometer, the measurement may be influenced by minor amounts of silica abrasion. In order to prevent incorrect A_{260} -quantification of small RNA amounts centrifuge the eluate for 30 s at $> 11,000 \times g$ and take an aliquot for measurement without disturbing any sediment. Alternatively, use a silica abrasion insensitive RNA quantification method (e.g., RiboGreen® fluorescent dye).

Measurement not in the range of photometer detection limit

Unexpected
 A_{260}/A_{280} ratio

- In order to obtain a reliable A_{260}/A_{280} ratio it is necessary that the initially measured A_{260} and A_{280} values are significantly above the detection limit of the photometer used. An A_{280} value close to the background noise of the photometer will cause non reliable A_{260}/A_{280} ratios.
-

6.2 Ordering information

Product	REF	Pack of
NucleoSpin® totalRNA FFPE	740982.10/.50/.250	10/50/250 preps
NucleoSpin® totalRNA FFPE XS	740969.10/.50/.250	10/50/250 preps
NucleoSpin® RNA XS	740902.10/.50/.250	10/50/250 preps
NucleoSpin® RNA Clean-up XS	740903.10/.50/.250	10/50/250 preps
NucleoSpin® RNA	740955.10/.50/.250	10/50/250 preps
NucleoSpin® RNA/Protein	740933.10/.50/.250	10/50/250 preps
NucleoSpin® TriPrep*	740966.10/.50/.250	10/50/250 preps
rDNase Set	740963	1
Paraffin Dissolver (blue)	740343.60	60 mL
Paraffin Dissolver	740968.25	25 mL
Collection Tubes (2 mL)	740600	1000

Visit www.mn-net.com for more detailed product information.

* DISTRIBUTION AND USE IN THE USA IS PROHIBITED FOR PATENT REASONS.

6.3 References

Castiglione F. et al. (2007): Real-time PCR analysis of RNA extracted from formalin-fixed and paraffin-embedded tissues: effects of the fixation on outcome reliability. *Appl Immunohistochem Mol Morphol* 15(3): 338–342.

Chung J. Y. et al. (2008): Factors in tissue handling and processing that impact RNA obtained from formalin-fixed, paraffin-embedded tissue. *Journal of Histochemistry & Cytochemistry*. 56(11): 1033–1042.

Leyland-Jones B. R. et al. (2008): Recommendations for collection and handling of specimens from group breast cancer clinical trials. *J. Clin. Oncol* 26(34): 5638–5644.

Manchester K.L. (1995): Value of A_{260}/A_{280} ratios for measurement of purity of nucleic acids. *Biotechniques* 19, 208–209.

von Ahlfen S. et al. (2007): Determinants for RNA quality from FFPE samples. *PLoS ONE*. Issue 12, e1261.

von Maldegem F. et al. (2008): Effects of processing delay, formalin fixation, and immunohistochemistry on RNA recovery from formalin-fixed paraffin-embedded tissue sections. *Diagn Mol Pathol*. 17(1): 51–58.

Wilfinger W.W. et al. (1997): Effect of pH and ionic strength on the spectrophotometric assessment of nucleic acid purity. *Biotechniques* 22, 474–481.

6.4 Product use restriction/warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety instructions from the respective Material Safety Data Sheet of the product.


This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY-NAGEL literature. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY-NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.


Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:
MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-333
support@mn-net.com



Plasmid DNA
Clean up
RNA
DNA
Viral RNA and DNA
Protein
High throughput
Accessories
Auxiliary tools



MACHERY-NAGEL

www.mn-net.com

MACHERY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

