DiaQuant®

Peracetic acid

de

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



en

Description:

DiaQuant® Peracetic Acid are test strips for the semi-quantitative detection of peracetic acid (PAA) (CH_3COOOH) in rinse solutions.

Measuring range:

5–50 mg/L peracetic acid (PAA)

Color gradation:

0 · 5 · 10 · 20 · 30 · 50 mg/L peracetic acid (PAA)

Pack content:

100 test strips

Reaction principle:

Peracetic acid oxidizes an aromatic amine forming a blue color.

General information:

Only remove the required number of test strips from the container, and then close the container securely again immediately afterwards. Do not touch the test field. Ideally, the color evaluation should be carried out by diffuse daylight. Artificial light can lead to difficulties or an incorrect color reading.

Hazard warnings:

This test does not contain hazardous substances that must be labelled.

Instructions for use:

1. Insert the test strip into the test solution for 1 s.
2. Shake off any excess liquid.
3. Wait 5 s.
4. Then immediately compare the test field with the color scale. Take the value which matches closest with the colored test field (reading accuracy: $\pm \frac{1}{2}$ colored field of the scale). The reaction color of the test field may change after the value has been taken. It is therefore crucial to evaluate the coloration within the prescribed timescale in order to achieve correct results. Ignore any color changes that occur after the reaction time (5 s).

Quality control:

To check the correct functioning of the test strips, use a peracetic acid solution with a concentration of 20 mg/L. For this purpose, first prepare a stock solution of 2000 mg/L by mixing 1.115 mL of 39% peracetic acid with distilled water to a volume of 250 mL in a volumetric flask. Then, dilute 5 mL of the stock solution in 500 mL of hydrogen peroxide solution (approx. 2000 mg/L) (= 20 mg/L PAA). Immediately perform the measurement with the test strip. If the control solution produces a negative result even after repeating the process, then the remaining unused test strips must be discarded. Even during a negative control (inserting a test strip into distilled water) no blue coloration may occur. Possible reasons for both faults may be that the use-by-date has been exceeded, the test strip container was left open for too long or the test strips were stored incorrectly.

Interferences:

If the solution to be tested has a pH value of between 2 and 9, the reaction will take place without interferences. Very strong acid solutions have to be buffered using sodium acetate; alkaline solutions have to be set to a pH value of 5–7 using citric acid. The presence of other strong oxidants such as halogens and halogen oxides will also lead to positive results.

Storage:

Avoid exposing the strips to sunlight and moisture. Keep container cool and dry (storage temperature 4–30 °C).

If correctly stored, the test strips may be used until the use-by-date printed on the packaging.

Additional information:

The test strip container stopper contains a non-toxic drying agent. If swallowed, drink plenty of water.

Disposal:

Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

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