MACHEREY-NAGEL

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1. α-acids (MEBAK-method)

1.1 Background

When beer is brewed, iso-alpha acids are formed from the initially insoluble alpha acids in the hops during boiling, which give the beer bitterness – a process called isomerization. The bitterness units (BU) in beer are derived from the alpha acid content and influence the taste, durability of the froth and have bacteriostatic properties.

1.2 Method

The bitter substances are extracted from an acidified sample using iso-octane. Interferences can be removed by rinsing the extract with acidic methanol. The concentration is then determined spectrophotometrically

1.3 Measuring range

0–60 mg/L α-acids

1.4 Cuvette and wavelenght

10 mm rectangular cuvette / 360 nm

1.5 Accessoires

- Centrifuge tubes with solvent-tight screw cap, 100–110 mL capacity
- Shaking machine
- Centrifuge, 3000 rpm
- Spectrophotometer, 360 nm, e.g. REF 919600.1
- 10 mm rectangular cuvettes QS, e.g. REF 91933

1.6 Reagents

- Hydrochloric acid, 6N
- Iso-octane (2,2,4-Trimethylpentan), spectroscopically pure
- Sodium sulphate, anhydrous
- Methanol
- Hydrochloric acid, 4N
- Sodium hydroxide, 6N, carbonate-free
- Acidic methanol solution: mix 64 mL methanol and 36 mL 4N hydrochloric acid (prepare daily)
- Alkaline methanol solution: make up 0.2 mL 6N sodium hydroxide with methanol to 100 mL (prepare daily)

1.7 Sample preparation

- Clarify wort or cloudy beer by centrifuging at 3000 rpm for 15 min (do not filter)
- Remove carbon dioxide from wort or beer without loss of foam.

1.8 Procedure

Blank preparation:

- 5 mL iso-octane into a 25 mL volumetric flask
- Fill up to the mark with alkaline methanol solution and mix

Main value:

- Pipette 50 mL of the sample tempered to 20 °C into a centrifuge tube
- Add 3.0 mL of 6N hydrochloric acid and 25 mL of iso-octane
- Close the centrifuge tube and shake mechanically for 30 min at optimum mixing intensity
- Centrifuge for 5 min at 3000 rpm to separate the phases and break the emulsion
- remove the lower aqueous phase by aspiration with a pipette and discard. Add enough sodium sulphate to the iso-octane phase so that it is clear after brief vigorous shaking
- Pipette 10 mL of the iso-octane phase into a 25 mL mixing cylinder
- Add 10 mL of acidic methanol solution, mix the contents by rotating the cylinder 100 times
- Pipette 5 mL of the supernatant clear iso-octane phase into a 25 mL volumetric flask
- Fill up to the mark with alkaline methanol solution and mix
- Measure the absorbance of the iso-octane solution against the blank at 360 nm

1.9 Measurement

Call up and measure method 3–40

1.10 Literature

Brautechnische Analysemethoden, volume II, 4. edition, 2002, ISBN 3-9805814-5-4, chapter 2.18.2, page 116ff

2. Anthocyanogens acc. to Harris and Ricketts

2.1 Background

Anthocyanogens are phenolic compounds that enter the beer from hops and malt depending on the processing. This can lead to changes in the color, taste (stability), foam or chemical-physical stability of the beer. Polymerization or condensation can lead to undesirable by-products (e.g. proteins).

2.2 Method

First, the anthocyanogens are adsorbed onto polyamide. The resulting adsorbate is dissolved in a butanol-hydrochloric acid solution and heated to produce a red solution, which is then measured spectrophotometrically.

2.3 Measuring range

1 – 100 mg/L Anthocyanogens

2.4 Cuvette and wavelength

10 mm rectangular cuvette / 550 nm

2.5 Accessories

- Shaking machine
- Centrifuge
- Mixing cylinder with ground-glass stopper, 50 mL
- Frit 1 G4
- Suction flask
- Test tubes with ground-glass stopper, 30 mL, graduation up to 25 mL
- Vacuum pump
- Spectrophotometer, 550 nm, e.g. REF 919600.1
- 10 mm rectangular cuvettes OS, e.g. REF 91933

2.6 Reagents

- MN-polyamide SC 6
- Solution 1: *n*-butanol / 37 % hydrochloric acid 5 + 1 (V/V)
- Solution 2: Dissolve 120 mg iron(II)-sulphate ($\text{FeSO}_4 \times 7 \text{H}_2\text{O}$) in 100 mL solution 1

2.7 Sample preparation

- Wort and centrifuge green beers for 10 min at 3000 rpm

2.8 Procedure

Blank preparation:

- Pipette 10 mL dist. water (blank value) into a 50 mL mixing cylinder
- Treat blank as main sample from step 2.

Main value

- 5 mL beer (or wort) + 5 mL distilled water into a 50 mL mixing cylinder
- Rinse each 0.5 g polyamide powder with 10 mL distilled water into the mixing cylinder
- Shake the mixing cylinder mechanically for 40 minutes
- Filter each suspension through a 1 G4 frit and rinse twice with approx. 20 mL H₂O
- Suck the frits and polyamide powder dry, transfer the residue quantitatively into a test tube using a spatula and rinse with 15 mL of solution 1
- Add 0.5 mL of solution 2 to each and heat the test tube in a boiling water bath for 30 min (stir well with a glass rod during the first 5 min)
- Remove the glass rod, rinse with a little solution 1, temper the test tube to 20 °C, and fill up to 25 mL with solution 1.

2.9 Measurement

Call up and measure method 3–41

2.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.17.2, page 109ff

https://www.mebak.org/methode/b-590_42_111/anthocyanogenic-method-harris-and-ricketts/843

3. Beer color, photometric (MEBAK method)

3.1 Background

The color of the beer is also a quality feature that is checked during the brewing process. Pale beers have a low EBC value, dark beers have a higher value.

3.2 Method

The absorbance of the sample is measured at 430 nm and converted using a specified factor.

3.3 Measuring range

1 – 60 EBC units

3.4 Cuvette and wavelength

10 mm rectangular cuvette / 430 nm

3.5 Accessories

- Spectrophotometer, 430 nm (± 0.5 nm), e.g. REF 919600.1
- 10 mm rectangular cuvettes OS, e.g. REF 91933

3.6 Interferences

A spectrometric absorption curve does not reflect the color impression of the human eye, as light of the same intensity in different parts of the spectrum affects the eye differently. In addition, the absorbance curves at 430 nm are very steep, so that slight measurement errors can occur. There are also differences when comparing light-colored beers with diluted dark beers.

3.7 Sample preparation and procedure

- Dilute the sample so that the absorbance is within the linearity of the spectrophotometer
- Filter the sample through a membrane filter. Filtration can be omitted if the turbidity of the diluted sample is below 1 EBC turbidity unit.
- If necessary, clarify the sample by adding 0.1 % diatomaceous earth and filtering upstream of the membrane filtration.
- Transfer the sample to a 10 mm rectangular cell and measure against dist. Water (blank value).

3.8 Measurement

Call up and measure method 3 – 42

3.9 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.13.2, page 88ff

4. Bitterness units in wort (MEBAK method)

4.1 Background

The predominant bitter compounds in wort and beer are the so-called iso- α -acids. In wort, β - and δ -acids can also be detected. There are also oxidation products of hop bitterness, which also contribute to the bitter taste.

4.2 Method

The sample is acidified and then extracted with iso-octane. The bitter substances extracted in the iso-octane extract are measured spectrophotometrically.

4.3 Measuring range

20–60 BU

4.4 Cuvette and wavelengths

10 mm rectangular cuvette / 275 nm

4.5 Accessories

- Centrifuge tubes with solvent-tight caps, 35 mL
- Glass beads
- Shaking machine
- Centrifuge, 3000 rpm
- Spectrophotometer, 275 nm, e.g. REF 919600.1
- 10 mm rectangular cuvette QS, e.g. REF 919120

4.6 Reagents

- Hydrochloric acid, 6N
- Iso-octane (2,2,4-trimethylpentane), spectroscopically pure (absorbance measured in 10 mm QS rectangular cuvette at 275 nm against H₂O < 0.010) (e.g. Uvasol)

4.7 Sample preparation

- Clarify wort by centrifuging at 3000 rpm for 15 min (do not filter the sample)
- Remove carbon dioxide from the sample without loss of foam

4.8 Procedure

- Pipette 10 mL (or 5 mL wort + 5 mL distilled H₂O) of the sample tempered to 20 °C into a centrifuge tube
- Add 0.5 mL 6N hydrochloric acid, 20 mL iso-octane and 3 glass beads
- Close the centrifuge tube and shake mechanically for 15 min at 20 °C
- Centrifuge for 3 min at 3000 rpm to separate the phases and break the emulsion
- Measure the absorbance of the iso-octane extract in a 10 mm rectangular cuvette at 275 nm against iso-octane of the same quality (blank value)

4.9 Measurement

Call up and measure method 3–43

4.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.18.1, page 114ff

5. Bitterness units in beer (REF 918172)

5.1 Background

(See 4.1)

5.2 Method

(See 4.2)

5.3 Measuring range

2.0–80.0 BE

5.4 Cuvette and wavelengths

10 mm rectangular cuvette / 275 nm

5.5 Accessories

- Spectrophotometer, 275 nm, e.g. REF 919600.1
- 10 mm rectangular cuvettes QS, e.g. REF 919120
- pipette

5.6 Reagents

- Reagents from the test box

5.7 Procedure

- Pipette a portion of the iso-octane (upper phase) into the rectangular cuvette. Zero the photometer. Pipette back the iso-octane.
- Add 1 mL sample via the inner wall of the cuvette. Shake vigorously for 20 seconds, wait 1 min for phase separation.
- Pipette 2.5 mL of the iso-octane phase into the rectangular cuvette and measure.

5.8 Measurement

Call up and measure method 1 – 72

5.9 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.18.1, page 114ff

<https://www.mebak.org/methode/b-400-17-110/bittereinheiten/663>

6. Total carotenoids fraction 1 – 4 (MEBAK method)

6.1 Background

Carotenoids are plant color pigments that occur in many fruits and are responsible for their yellow to reddish color. In the human body, carotenoids are converted into vitamin A, among other things. The concentration of carotenoids is detected in drinks such as vitamin juices.

6.2 Method

The sample is mixed with Carrez solution and then extracted with acetone and petroleum ether. Measurements are made with a spectrophotometer.

6.3 Measuring range

0.1 – 8.0 mg/L carotenoids

6.4 Cuvette and wavelength

10 mm rectangular cuvette / 450 nm

6.5 Accessories

- Centrifuge tubes, 60 to 100 mL capacity
- Centrifuge
- Shaking funnel, 200 mL
- Round bottom flask, 250 mL
- Chromatography column: H = 25 – 30 cm, I. D. = 2 cm
- Rotary evaporator
- Photometer, wavelength 450 nm, e.g. REF 919600.1
- Glass cuvettes, 10 mm, e.g. REF 91933

6.6 Reagents

- Carrez I solution: aqueous potassium ferricyanide(II) solution, $K_4[Fe(CN)_6] \cdot 3 H_2O$, 150 g/L
- Carrez-II solution: aqueous zinc sulphate solution, $ZnSO_4 \cdot 7 H_2O$, 300 g/L acetone pur.
- Petroleum ether pure, boiling range 40 – 60 °C
- Sodium sulphate anhydrous, pure.
- Toluene p. a.
- Elution mixture A: Petroleum ether + toluene = 4 parts + 1 part
- Elution mixture B: Petroleum ether + toluene = 2 parts + 1 part
- Aluminum oxide 90 active, neutral (Merck No.: 101077), partially deactivated:
- Add 12 g H_2O to 100 g Al_2O_3 and shake until all lumps have disintegrated (keep sealed).

6.7 Sample preparation and procedure

A: Determination of total carotenoids

- Homogenize the beverage by shaking vigorously - Dilute concentrates to juice strength
- Pipette 5 to 50 mL of the beverage into the centrifuge tube, depending on the test material
- if necessary, add water to 50 mL
- add 1 mL each of Carrez I and Carrez II solution
- Mix the solution well with a glass rod and centrifuge after standing for approx. 2 minutes (approx. 5 min)
- Select the amount of test material used so that the extraction of the petroleum ether extract made up to 100 mL does not exceed 0.5 units (for example: Orange juice: 25 mL; orange nectar: 50 mL; passion fruit juice: 5 to 20 mL; passion fruit nectar: 25 mL; tangerine juice: 20 mL).
- Discard the supernatant, only weakly colored to colorless liquid after carefully decanting off the precipitate
- Extract the adsorbed carotenoids from the precipitate with 40 mL acetone. To do this, first swirl the precipitate together with the acetone vigorously with a glass rod and then centrifuge again for about 3 min (about 5 min) and pour the now bright yellow colored acetone into a shaking funnel into which 50 mL (measuring cylinder) of petroleum ether has previously been added
- Shake vigorously and separate the separating water phase (approx. 10 mL) from the organic phase after allowing it to stand for a short time
- Extract the precipitate treated with the acetone again with approx. 20 mL acetone in the manner described above to remove any carotenoids still adhering to it.
- after centrifugation, add this second extraction by carefully decanting the organic phase in the shaking funnel
- wash the organic phase with 50 mL H₂O by shaking briefly
- discard the water phase after the small emulsion has formed again
- For drying, place the carotenoid-containing petroleum ether phase in a centrifuge tube into which 2 g anhydrous sodium sulphate has previously been weighed.
- stir briefly with a glass rod and centrifuge
- carefully rinse the petroleum ether phase into a 100 mL volumetric flask
- Add approx. 30 mL of petroleum ether to the sodium sulphate precipitate adhering to the centrifuge tube to remove any carotenoids still adhering and centrifuge with a glass rod after stirring.
- combine this weakly colored petroleum ether with the first fraction and fill the volumetric flask to the mark with petroleum ether
- Measure the petroleum ether extract at 450 nm in 10 mm cuvettes against petroleum ether

B: Determination of the individual carotenoid fractions

Preparation of the chromatography column

- Slurry 28–30 g partially deactivated aluminum oxide in a beaker with approx. 30 mL petroleum ether
- Fill the chromatography column charged with a little petroleum ether and stuffed with a little glass wool and drain the petroleum ether to just above the aluminum oxide
- Make sure that the petroleum ether and the aluminum oxide do not exhibit any major temperature differences in order to ensure bubble-free filling of the column. The filling height should be 8 cm.

Fractionation of the carotenoids

- The petroleum ether extract obtained from the determination of the total carotenoids
- quantitatively into a 250 mL round bottom flask and evaporate petroleum ether on the rotary evaporator
- at about 40 °C in a vacuum
- Take up the residue immediately in a few milliliters of petroleum ether and transfer to the chromatography column using a pipette (e.g. Eppendorf).
- rinse with small portions (approx. 10 mL in total) of petroleum ether
- Avoid swirling up the aluminum oxide when adding the sample solution,
- allow to soak in straight
- Elute at an elution rate of about 2 drops per second
- Elute 4 fractions and collect them separately in an appropriate volumetric flask:

Fraction I (so-called carotene fraction): Petroleum ether

Fraction II (so-called cryptoxanthin ester fraction): Elution mixture A

Fraction III (so-called xanthophyll ester fraction): Elution mixture B

Fraction IV (apocarotene and others): Acetone

- Measure fractions against petroleum ether at 450 nm

6.8 Measurement

Call up and measure method 3 – 45

6.9 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 3.7.2.1, page 276ff

<https://www.mebak.org/methode/b-590-46-110/gesamtcarotin-photometrisch/851>

7. Total carotenoids in carrot juice (MEBAK method)

7.1 Background

Carotenoids are plant color pigments that occur in many fruits and are responsible for their yellow to reddish color. In the human body, carotenoids are converted into vitamin A, among other things. The concentration of carotenoids is detected in drinks such as vitamin juices.

7.2 Method

The carotenoids are extracted from the sample and the extract solution is washed several times. The spectrophotometric measurement then takes place.

7.3 Measuring range

0.1 – 8.0 mg/L carotenoids

7.4 Cuvette and wavelength

10 mm rectangular cuvette / 450 nm

7.5 Accessories

- Separating funnel 250 mL, amber glass
- 50 mL volumetric flask, amber glass
- 100 mL volumetric flask, amber glass
- 20 mL volumetric flask
- 50 mL volumetric flask
- Spectrophotometer, wavelength 450 nm, e.g. REF 919600.1
- Glass cuvettes, 1 cm path length, e.g. REF 91933

7.6 Reagents

- Acetone, p.a.
- *n*-hexane, p.a.
- Methanol, p.a.
- Extraction solution: acetone/*n*-hexane (50:50, v/v)
- NaCl solution, 10 %
- Sodium sulphate, anhydrous, p.a.

7.7 Procedure

- Pipette 10 mL of homogenized sample into a 50 mL volumetric flask and fill with H₂O
- (for ACE beverages: 10 mL in 20 mL volumetric flask)
- Add 5 mL of the diluted sample with 50 mL extraction solution to the separating funnel and shake
- Add 100 mL NaCl solution, shake, wait for phase separation and discard the usually not completely clear water phase
- if the aqueous phase is still clearly yellow in color, add 50 mL methanol to the sample solution, shake, add another 50 mL extraction solution and proceed as described above
- Wash the carotenoid-containing organic phase twice with 50 mL NaCl solution each time
- Place the organic phase in a beaker with 2–3 g Na₂SO₄, stir briefly and decant the hexane phase into a 100 mL volumetric flask (for ACE beverages: 50 mL volumetric flask)
- Wash carotenoids adhering to the Na₂SO₄ several times with 5 mL hexane until almost colorless (usually three times is sufficient), decant hexane each time and add to the sample, make up to the mark with hexane
- measure in 10 mm cuvettes loaded with 1 g Na₂SO₄ at 450 nm against *n*-hexane as soon as the salt has sunk

7.8 Measurement

Call up and measure method 3–46

7.9 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 3.7.2.2, page 280ff

<https://www.mebak.org/methode/b-590-46-110/gesamtcarotin-photometrisch/851>

8. Total polyphenols (MEBAK method)

8.1 Background

Depending on the processing, phenolic compounds from hops and malt enter the beer, which can lead to changes in the color, taste (stability), foam or chemical-physical stability of the beer. Polymerization or condensation can lead to undesirable by-products (e.g. proteins).

8.2 Method

Polyphenols react in an alkaline environment with iron(III) ions to form colored complexes. The concentration is measured spectrophotometrically.

8.3 Measuring range

1 – 800 mg/L Polyphenols

8.4 Cuvette and wavelength

10 mm rectangular cuvette / 600 nm

8.5 Accessories

- Centrifuge
- Spectrophotometer, 600 nm, e.g. REF 919600.1
- 10 mm rectangular cuvette OS, e.g. REF 91933

8.6 Reagents

- Carboxymethylcellulose-ethylenediaminetetraacetic acid solution (CMC-EDTANa):
 - c. Weigh in 10 g CMC (low viscosity) and 2 g EDTA-Na₂
 - d. Dissolve the substances in approx. 500 mL H₂O while stirring. solution with H₂O to 1 L, if necessary, clarify by centrifugation.
- Ammonium iron(III) citrate, 3.5 %:
 - c. Dissolve 3.5 g ammonium iron(III) citrate, green (16 % Fe) in H₂O to 100 mL; the solution must be completely clear and can be kept for about 1 week
 - d. Ammonia, diluted:
Dilute 1 part concentrated ammonia (d = 0.91) with 2 parts H₂O

8.7 Sample preparation

- Remove carbon dioxide from beer by shaking
- Clarify cloudy wort or beers by centrifugation

8.8 Procedure

Blank value preparation:

- Dose 10 mL sample (remove carbon dioxide from beer by shaking cloudy wort or clarify beer by centrifuging) into a 25 mL volumetric flask
- Add 8 mL CMC-EDTA solution and mix well
- Add 0.5 mL diluted ammonia solution and mix well
- Top up with H₂O to 25 mL and mix again.

Main Value:

- Mix well 10 mL test sample + 8 mL CMC-EDTA solution in
- a 25 mL volumetric flask
- Add 0.5 mL iron(III) solution and mix well again
- Add 0.5 mL diluted ammonia solution, mix thoroughly
- Fill up to 25 mL with H₂O and mix
- after 10 min measure absorbance in 10 mm rectangular cuvette at 600 nm against blank sample

Note:

- Always mix the solutions well after adding the individual reagents

8.9 Measurement

Call up and measure method 3–47

8.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.17.1, page 107ff

<https://www.mebak.org/methode/b-590-42-111/anthocyanogene-methode-harris-und-ricketts/843>

9. Copper Cuprethol (MEBAK method)

9.1 Background

Copper is contained in spraying agents and can therefore enter the brewing process via the raw materials. However, it is also found in equipment or pipes used in the brewing process. Heavy metal ions such as copper have a negative effect on the taste and colloidal stability of beer.

9.2 Method

The copper forms a color complex with diethanolamine and carbon disulfide. This is measured spectrophotometrically.

9.3 Measuring range

0.10–5.00 mg/L Cu

9.4 Cuvette and wavelength

10 mm rectangular cuvette / 445 nm

9.5 Accessories

- Spectrophotometer, 445 nm, e.g. REF 919600.1
- Cuvettes, 10 mm pathlength, e.g. REF 91933

9.6 Reagents

- Solution 1: Dissolve 4 g diethanolamine in 200 mL methanol (stable)
- Solution 2: Dissolve 1 g carbon disulphide in 200 mL methanol (stable)
- Solution 3: Mix 1 part by volume of methanol with 1 part by volume of solution 1 (stable)
- Solution 4: (cuprethol reagent): mix 3 parts by volume of solution 1 with 1 part by volume of solution 2
- Solution 2 (prepare this solution fresh every day)
- Acetate buffer, pH 4.6: Dissolve 105 g sodium acetate ($\text{CH}_3\text{CO}_2\text{Na} - 3 \text{H}_2\text{O}$) in approx. 800 mL H_2O , add 65 mL concentrated acetic acid and make up to 1 L with H_2O

9.7 Sample preparation

- 50 mL of beer free of carbon dioxide (allow beer foam to collapse)
- Pipette into 100 mL conical flask
- Add 25 mL acetate buffer and mix

9.8 Procedure

Blank preparation:

- Pipette 20 mL of the buffered beer solution into a 50 mL Erlenmeyer flask
- Add 2 mL of solution 3 and mix

Main value:

- Pipette 20 mL of the buffered beer solution into a 50 mL conical flask
- Add 2 mL of solution 4 and mix
- Measurement must be carried out within 10 min after addition of solution 4 in a 10 mm rectangular cuvette (against the blank value)

9.9 Measurement

Call up and measure method 3–48

9.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.29.4, page 152ff

<https://www.mebak.org/methode/b-420-37-009/kupfer-cuprethol/2261>

10. Nickel (MEBAK-method)

10.1 Background

Like copper and other heavy metals, nickel impairs the colloidal stability of beer as well as the foam stability. Nickel gets into the beer via beer foam stabilizers or contact with stainless steel.

10.2 Method

The sample is extracted and then mixed with dimethylglyoxime to form a colored complex. The measurement is carried out spectrophotometrically.

10.3 Measuring range

0.03–5.00 mg/L Ni

10.4 Cuvette and wavelength

10 mm rectangular cuvette / 530 nm

10.5 Accessories

- Spectrophotometer, 530 nm, e.g. REF 919600.1
- Cells, 10 mm path length, e.g. REF 91933
- Kjeldahl flask
- Separating funnel, 125 mL

10.6 Reagents

- Nitric acid, 65 % (d = 1.40)
- Perchloric acid, 70 % (d = 1.67)
- Hydroxylamine hydrochloride, 10 %
- Sodium citrate, 20 %: dissolve 20 g trisodium citrate-2-hydrate with H₂O to 100 mL
- Ammonia, 25 % (d = 0.91)
- Dimethylglyoxime, 1 % in ethanol
- chloroform
- Ammonia, diluted: Dilute 1 part 25 % ammonia with 49 parts H₂O
- Hydrochloric acid, 0.33 N
- Sodium tartrate, 20 %: dissolve 20 g sodium tartrate-2-hydrate with H₂O to 100 mL
- Potassium peroxodisulphate, 4 %
- Sodium hydroxide solution, 5 N
- Dimethylglyoxime: dissolve 1 g dimethylglyoxime in 5 mL 5 N sodium hydroxide solution and mix with H₂O to 100 mL
- Phenolphthalein, 1 % in ethanol

10.7 Sample preparation

- Condense 100 mL beer in a Kjeldahl flask to a large extent
- evaporate three times with 65 % nitric acid
- then evaporate twice with 70 % perchloric acid to destroy all organic substances
- Transfer the residue to a 125 mL separating funnel

10.8 Procedure

Blank value:

- 15 mL 0.33 N hydrochloric acid, 2 mL sodium tartate, 10 mL potassium peroxy-disulphate, 0.6 mL dimethylglyoxime sodium hydroxide solution and 2.5 mL 5 N sodium hydroxide solution in a 50 mL volumetric flask with distilled water. Fill up to the mark with distilled water.

Main value:

- Residue from preparation Add 2 mL hydroxylamine hydrochloride solution, 5 mL sodium citrate solution and 3 drops of phenolphthalein
- add 25 % ammonia until a pink coloration appears
- add a further 4 drops of ammonia and 2 mL of alcoholic dimethylglyoxime solution
- Bring to approx. 60 mL with H₂O
- Shake out three times with 5 mL chloroform for 30 seconds each time (any interfering copper compounds can be removed from the chloroform extract by shaking for one minute with 5 mL diluted ammonia solution (1 + 49)).
- Extract the chloroform phase with 10 mL of 0.33 N hydrochloric acid for 1 min.
- Extract again with 5 mL 0.33 N hydrochloric acid (nickel passes into the hydrochloric acid solution, making sure that no chloroform residue remains in the combined hydrochloric acid extracts)
- Collect the hydrochloric acid extracts in a 50 mL volumetric flask, add 2 mL sodium tartrate solution, 10 mL potassium peroxodisulphate solution, 0.6 mL dimethylglyoxime sodium hydroxide solution and 2.5 mL 5 N sodium hydroxide solution
- Make up to 50 mL with H₂O and leave to stand for at least 30 min (but no longer than 2 h)
- Measure the absorbance of the color solution in a 10 mm cuvette at 530 nm against the blank value

10.9 Measurement

Call up and measure method 3–49

10.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.29.6, page 156ff

<https://www.mebak.org/methode/b-590-72-111/nickel-ni-spektralphotometrisch/3212>

11. Thiobarbituric acid index (MEBAK method)

11.1 Background

The thiobarbituric acid index provides information on the thermal stress of malt and wort. During malt and wort preparation, increased thermal stress occurs at least twice, which can jeopardize aroma stability.

11.2 Method

Thiobarbituric acid is added to the sample and the resulting yellow coloration is measured spectrophotometrically.

11.3 Measuring range

1 – 100 TBI

11.4 Cuvette and wavelength

10 mm rectangular cuvette OS / 448 nm

11.5 Accessories

- Water bath, 70 °C
- Brown reaction tubes with ground joint, 20 mL or 25 mL
- Spectrophotometer, 448 nm, e.g. REF 919600.1
- 10 mm rectangular cuvettes OS, e.g. REF 91933

11.6 Reagents

- Acetic acid, 90 %:
Dilute 225 g acetic acid 100 % (glacial acetic acid) with H₂O to 250 g.
- Thiobarbituric acid, 0.02 mol/L:
Dissolve 0.288 g 2-thiobarbituric acid (M = 144.15 g/mol) in 100 mL volumetric flask with 90 % acetic acid while heating in a water bath and, after cooling to 20 °C, make up to the mark with 90 % acetic acid. Prepare fresh daily.
- Diatomaceous earth

11.7 Sample preparation

- Clarify turbid test solutions by filtration over diatomaceous earth.

Dilution

- Dilute condiments and beers 10-fold with H₂O
- Dilute condiments 5-fold with H₂O

Blank value

- Add 5 mL of 90 % acetic acid to 10 mL of the diluted sample, shake and proceed as for the main value.

11.8 Procedure

Main value:

- Add 5 mL of thiobarbituric acid to 10 mL of the diluted sample and shake.
- Place in a water bath at 70 °C for 70 min (avoid direct sunlight, make sure that the temperature in the bath only drops by 1 – 2 °C for a short time when adjusting the test tubes).
- After the reaction time has elapsed, quickly cool the test tubes down to 20 °C (strong running cold water or cold bath).
- Measure the yellow coloration immediately in 10 mm rectangular cuvettes at 448 nm against H₂O.

Note

- Follow the analysis instructions exactly, as they are of an empirical nature.

11.9 Measurement and result

- Method 3 - Call up 50 and measure
- Calculate the result: $TBI = (EH - EL) * F$
- EH = Absorbance of the main value
- EL = Absorbance of the blank value
- F = Dilution factor

11.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5-4, Chapter 2.4, page 35ff

<https://www.mebak.org/methode/b-400-05-111/thiobarbitursaeurezahl-tbz/646>

12. Vicinal diketones (REF 985010)

12.1 Background

Diacetyl (= 2,3-butanedione) is formed by oxidative decarboxylation of intermediate products during the primary fermentation of yeast. A concentration above 0.10 mg/L leads to a bad taste in beer.

12.2 Method

Diacetyl (and derivatives) react to form a quinoxaline derivative. The analysis is carried out spectrophotometrically at 335 nm.

12.3 Measuring range

0.015–0.600 mg/kg diacetyl

12.4 Cuvette and wavelength

Round cuvette / 335 nm

12.5 Accessories

- Accessories for steam distillation (see special instructions)
- Piston pipette for 4 mL and suitable pipette tip
- Spectrophotometer, 335 nm, e.g. REF 919600.1

12.6 Reagents

- Anti-foam emulsion (free of diketones)
- Reagents from the test box

12.7 Sample preparation

Steam distillation necessary, see special instructions

12.8 Procedure

Blind value:

- Pipette 4 mL sample into the cuvette
- Shake vigorously and zero the photometer

Main Value:

- After measuring the blank value, add 1 *NANOFIX*[®] to the solution
- Shake vigorously
- Wait 5 min and then measure

12.9 Measurement

Call up and measure method 0–10

12.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.23, page 134ff

<https://www.mebak.org/methode/b-420-21-111/vicinale-diketone-photometrisch/731>

13. Cyclamate (MEBAK method)

13.1 Background

Cyclamate is a sweetener that is added, for example, to energy-reduced/sugar-free non-alcoholic fruit-based drinks. However, only a maximum quantity of 250 mg/L is permitted.

13.2 Method

Cyclamate is converted into a dichloro compound and extracted with heptane. The measurement is carried out spectrophotometrically.

13.3 Measuring range

0–400 mg/L cyclamate

13.4 Cuvette and wavelength

10 mm rectangular cuvette QS/314 nm

13.5 Accessories

- Spectrophotometer, 314 nm, e.g. REF 919600.1
- Rectangular cuvette QS, 10 mm, e.g. REF 919120
- Separating funnel, 100 mL
- Round filter (e.g. MN 640w)

13.6 Reagents

- Heptane, p. a.
- Sulphuric acid, 50 % (v/v)
- Sodium hypochlorite, free chlorine content 1.5 % (*VISOCOLOR*[®] *ECO* free chlorine 2, REF 931016 can be used to check the free chlorine content)
- Sodium hydroxide solution, 200 g/L
- Sodium sulphate, anhydrous p. a.
- Sodium cyclamate or *N*-cyclohexylsulfamic acid
- Potassium hexacyanoferrate, p. a.
- Carrez I solution:
 - Dissolve 150 g potassium hexacyanoferrate(II) with H₂O to 1000 mL
- Zinc sulphate, p. a.
- Carrez-II solution:
 - 300 g zinc sulphate Dissolve to 1000 mL with H₂O
- Cyclohexylsulfamic acid - stock solution:
 - Dissolve approx. 1 g cyclohexylsulfamic acid to the nearest 0.1 mg in H₂O and place in
 - 1 L volumetric flask with H₂O to the mark
- Cyclohexylsulfamic acid - standard solution:
 - Pipette 5 mL, 10 mL, 20 mL and 30 mL of the stock solution into 100 mL volumetric flasks and fill to the mark with H₂O. The concentrations of the standard solutions are 50, 100, 200 and 300 mg/l

13.7 Sample preparation and performance

- The cyclamate content in the sample solution should not exceed 300 mg/L, dilute liquid foods accordingly.
- Weigh solid foodstuffs into volumetric flasks and bring cyclamate into solution by ex-traction with H₂O
- Clarify cloudy samples with Carrez I and Carrez II solution
- To do this, pipette 50 mL of sample solution into a 100 mL volumetric flask, add 1 mL each of Carrez-I and Carrez-II solution, fill to the mark with H₂O and filter
- Acidify 10 mL sample or standard solutions in the separatory funnel with 2 mL H₂SO₄ (50 %)
- Add 20.0 mL *n*-heptane and 5 mL sodium hypochlorite solution and shake vigorously for 1 min
- Drain off the lower aqueous phase and discard. Any emulsion remains in the separating funnel
- Wash the heptane phase twice with 10 mL NaOH solution (20 %) and three times with 50 mL H₂O chlorine-free solution (the heptane phase must no longer smell of chlorine). From the addition of the sodium hypochlorite solution until the chlorine-free heptane phase is
- chlorine-free heptane phase work quickly
- Dry the organic phase over 1 – 2 g Na₂SO₄ and filter through a black belt filter
- Measure the absorbance of the sample solution and the standard solutions at 314 nm in a 10 mm cuvette against *n*-heptane as a blank value

13.8 Measurement

Call up and measure method 3 – 52

13.9 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 3.11.4, page 302ff

14. Iron (MEBAK method)

14.1 Background

Like copper and other heavy metals, iron impairs the colloidal stability of beer as well as its foam stability. Iron can enter the beer in many ways, e.g. via the raw materials used, pipes, beer foam stabilizers or filtration aids.

14.2 Method

Fe^{3+} is reduced to Fe^{2+} and combined with ferrozine in a complex. The concentration is determined spectrophotometrically.

14.3 Measuring range

0–1 mg/L Fe

14.4 Cell and wavelength

40 mm rectangular cuvette / 560 nm

14.5 Accessories

- Analytical balance, readability 0.1 mg
- Pipettes, 0.1 mL, 2 mL, 5 mL
- Spectrophotometer, 560 nm, e.g. REF 919600.1
- 40 mm rectangular cuvette

14.6 Reagents

- Prepare all solutions with iron-free H_2O .
- Buffer solution, pH 4.3:
 - Dissolve 75 g ammonium acetate and 150 g conc. acetic acid in approx. 800 mL H_2O , check pH and make up to 1 L.
- Ferrozine reagent:
 - Dissolve 0.257 g ferrozine or ferrospectral in 50 mL buffer (solution will keep for 2 weeks)
- Ascorbic acid, 2.5 % (prepare fresh daily)
- Hydrochloric acid, conc.
- Iron(III) standard solution to determine the calibration curve:
 - Dissolve 863.4 mg ammonium iron(III) sulphate [$\text{NH}_4\text{Fe}(\text{SO}_4)_2 - 12 \text{H}_2\text{O}$] in H_2O in a 1 L volumetric flask and, after adding 0.1 mL conc. hydrochloric acid, make up to the mark with H_2O . After diluting 50 mL of this solution with H_2O to 1 L, a standard solution containing 5 mg/mL Fe^{3+} is obtained.

14.7 Sample preparation

- Remove the carbon dioxide from the beer and allow the foam to collapse completely.

14.8 Procedure

- Pipette 40 mL beer into a 50 mL volumetric flask and add 2 mL ferrozine reagent and 1 mL ascorbic acid solution.
- Fill up to the mark with H_2O .
- Prepare the blank in exactly the same way, but without adding the ferrozine reagent; prepare a separate blank for each beer.
- Measure the absorbance of the solution in a 40 mm rectangular cuvette at 560 nm against the corresponding blank value.

14.9 Measurement

Call up and measure method 3–53

14.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.29.3, page 149ff

<https://www.mebak.org/methode/b-590-71-111/eisen-fe-spektralphotometrisch/3208>

15. Free amino nitrogen (FAN, MEBAK method)

15.1 Background

Amino acids and other low-molecular nitrogen compounds have an influence on the fermentation process. They can lead to undesirable by-products during fermentation. In addition, reactions with reducing sugars can occur, the products of which have an influence on the aroma, color and redox potential of the beer.

15.2 Method

The sample is mixed with a color reagent containing ninhydrin and heated. The resulting coloration is measured spectrophotometrically.

15.3 Measuring range

0–400 mg/L N

15.4 Cuvette and wavelength

10 mm rectangular cuvette / 570 nm

15.5 Accessories

- Test tubes with ground glass stopper, 16 × 150 mm
- Variable pipette 1.0–5.0 mL
- Pipette tips for pipette
- Boiling water bath
- Water bath at 20 °C
- Spectrophotometer, 570 nm, e.g. REF 919600.1
- 10 mm rectangular cuvette OS, e.g. REF 91933

15.6 Reagents

- Color reagent: Dissolve 10.0 g disodium hydrogen phosphate ($\text{Na}_2\text{HPO}_4 - 12 \text{H}_2\text{O}$), 6.0 g potassium dihydrogen phosphate (KH_2PO_4), 0.5 g ninhydrin, 0.3 g fructose in H_2O and make up to 100 mL (This solution can be kept for 2 weeks in a dark bottle. The pH value must be 6.6–6.8).
- Dilution solution: Dissolve 2 g potassium iodate in 600 mL H_2O and add 400 mL 96 % ethanol
- Stock solution: Dissolve 107.2 mg glycine in 100 mL H_2O . Store this stock solution at 0 °C.
- Standard solution: Make up 1 mL stock solution to 100 mL with H_2O .
- This standard solution contains 2 mg/L amino nitrogen.

15.7 Sample preparation

- Dilute wort 100-fold, beer 50-fold (1–3 mg/L amino nitrogen)

15.8 Procedure

- Analyze sample, standard solution or blank in triplicate.
- Pipette 2 mL of the diluted sample or standard solution or H_2O into a test tube.
- Add 1 mL of color reagent and mix.
- Close the test tube loosely with a glass stopper to avoid evaporation losses.
- Heat in a boiling water bath for exactly 16 min, then cool in a water bath at 20 °C for 20 min.
- Add 5 mL dilution solution.
- Measure absorbance within 30 min in a 10 mm rectangular cuvette at 570 nm against the equally treated blank (H_2O + color reagent).

Correction for dark seasonings and beers (make three corrections):

- Transfer 2 mL of the diluted sample to a test tube
- add 1 mL H₂O instead of the color reagent and proceed as described above
- Measure against H₂O after adding 5 mL of dilution solution

15.9 Measurement

Call up and measure method 3–54

15.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.8.4.1, page 62ff

<https://www.mebak.org/methode/b-400-11-111/freier-amino-stickstoff-fan-ninhydrin-methode-spektral-photometrisch/655>

16. Total carotenoids Wesergold (MEBAK method)

16.1 Background

Carotenoids are plant color pigments that occur in many fruits and are responsible for their yellow to reddish color. In the human body, carotenoids are converted into vitamin A, among other things. The concentration of carotenoids is detected in drinks such as vitamin juices.

16.2 Method

The sample is mixed with THF, then filtered using acetone and measured spectrophotometrically.

16.3 Measuring range

12–60 mg/L Carotenoids

16.4 Cuvette and wavelength

10 mm rectangular cuvette / 450 nm and 550 nm

16.5 Accessories

- Spectrophotometer, 450 nm and 550 nm, e.g. REF 919600.1
- 10 mm rectangular cuvettes QS, e.g. REF 919120

16.6 Reagents

- Tetrahydrofuran (THF), stabilized 2,6-Di-tert-butyl-4-methylphenol (Merck 109731)
- Acetone

16.7 Sample preparation and Procedure

- Dilute carrot juice 1 + 4 with water (ACE-beverages/Multivitamin with 50 RDA without dilution)
- Measure 0.50 mL diluted sample into a test tube (pipette), add 0.5 mL THF solution, mix, leave to settle
- Add 4.0 mL acetone, mix and filter through the membrane filter
- Measure sample using an equally treated water portion as blank value at 450 nm (E1) and at 550 nm (E2)

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17. Iso- α -acids (MEBAK method)

17.1 Background

During beer brewing, iso-alpha acids are formed from the initially insoluble alpha acids in the hops during boiling, which give the beer bitterness - a process called isomerization. The bitterness units (BE) in beer are derived from the alpha acid content and influence the taste, durability of the head and have bacteriostatic properties. Iso-alpha acids are sensitive to light, which is why the beer is protected by dark bottles.

17.2 Method

In an acidified sample, the bitter substances are extracted using iso-octane. Interferences can be removed by rinsing the extract with acidic methanol. The concentration is then determined spectrophotometrically.

17.3 Measuring range

0–60 mg/L

17.4 Cuvette and wavelength

10 mm quartz glass cuvette / 255 nm and 360 nm

17.5 Accessories

- Centrifuge tubes with solvent-tight screw cap, 100–110 mL capacity
- Shaking machine
- Centrifuge, 3000 rpm
- Spectrophotometer, 255 and 360 nm, e.g. REF 919600.1
- Quartz cuvettes, 10 mm layer thickness, e.g. REF 919120

17.6 Reagents

- Hydrochloric acid, 6 N
- Iso-octane (2,2,4-trimethylpentane), spectroscopically pure (absorbance measured in 1 cm cuvette at 275 nm against H₂O < 0.005).
- Sodium sulphate, anhydrous
- methanol
- Hydrochloric acid, 4 N
- Sodium hydroxide, 6 N, carbonate-free
- Acidic methanol solution: mix 64 mL methanol and 36 mL 4 N hydrochloric acid (prepare fresh daily)
- Alkaline methanol solution: make up 0.2 mL 6 N sodium hydroxide with methanol to 100 mL (prepare fresh daily)

17.7 Sample preparation

- Centrifuge the wort or cloudy beer at 3000 rpm for 15 min clarify (do not filter the sample)
- Remove carbon dioxide from beer without loss of foam

17.8 Procedure

- Pipette 50.0 mL of the sample tempered to 20 °C into a centrifuge tube
- Add 3.0 mL of 6 N hydrochloric acid and 25.0 mL of iso-octane
- Close the centrifuge tube and shake mechanically for 30 min at optimum mixing intensity
- Centrifuge for 5 min at 3000 rpm to separate the phases and break the emulsion
- Remove the lower aqueous phase by aspiration with a pipette Add enough sodium sulfate to the remaining, usually somewhat emulsified iso-octane phase so that it is clear after brief, vigorous manual shaking
- Pipette 10.0 mL of this into a 25 mL mixing cylinder
- Add 10.0 mL of acidic methanol solution and shake the contents for 3 minutes
- Pipette 5.0 mL of the supernatant clear iso-octane phase into a 25 mL volumetric flask
- Fill up to the mark with alkaline methanol solution and mix the contents
- Measure the absorbance of the iso-octane solution against a blank solution of 5.0 mL iso-octane and 20.0 mL alkaline methanol solution at 255 nm and 360 nm

17.9 Measurement

Call up and measure method 3–56

17.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.18.2, page 116ff

18. Copper ZDBT (MEBAK method)

18.1 Background

Copper is contained in spraying agents and can therefore enter the brewing process via the raw materials. However, it is also found in equipment or pipes used in the brewing process. Heavy metal ions such as copper have a negative effect on the taste and colloidal stability of beer.

18.2 Method

Copper is bound into a complex using ZDBT and extracted with trichloroethane. The measurement is carried out spectrophotometrically.

18.3 Measuring range

0–3 mg/L copper

18.4 Cuvette and wavelength

10 mm rectangular cuvette / 435 nm

18.5 Accessories

- Centrifuge tubes, 100 mL, with solvent-tight caps Spectrophotometer, 435 nm, e.g. REF 919600.1
- Quartz cuvettes, 10 mm layer thickness, e.g. REF 91933

Before use:

- Shake tubes with 15 mL H₂O, 3 mL diluted sulphuric acid and 5 mL ZDBT reagent to remove copper, discard rinsing solution and allow tubes to drain.
- Test tubes, treated in the same way
- Centrifuge
- Spectrophotometer, 435 nm
- Cuvettes, 10 mm path length

18.6 Reagents

- H₂O, copper-free (e.g. extract ion exchanger water with 0.5 % ZDBT solution)
- Standard copper solution, 0.01 mg/mL: 3.9295 g copper sulphate (CuSO₄ × 5 H₂O)
- Dissolve with H₂O to 1 L and dilute 100-fold immediately before use
- 1,1,1 - Trichloroethane
- Zinc dibenzyl dithiocarbamate (ZDBT), 0.5 % in trichloroethane
- Sulphuric acid, diluted: 1 part conc. sulphuric acid + 3 parts H₂O
- Hydrogen peroxide, 30 %
- Hexanol-1 (antifoam agent)

18.7 Sample preparation

- Decarbonize the beer and allow the foam to collapse completely

18.8 Procedure

Blank value:

- 25 mL dist. H₂O instead of a sample and treat in the same way as the main value
- Measure against trichloroethane (blank value)

Main value:

- Pipette 25 mL beer into centrifuge tube
- Add 1 mL H₂O and 3 mL diluted sulphuric acid
- Place in boiling water bath for 30 min
- After cooling, add 5 or 10 mL of color reagent
- Close the jar and shake for 10 min on a shaker - Centrifuge for 2 min at 2500 rpm
- Carefully remove part of the lower phase with a pipette and allow it to flow slowly into the test tube held at an angle in order to separate any entrained aqueous phase on the wall.
- fill 10 mm cuvette from the test tube and measure at 435 nm absorbance against pure trichloroethane (main value)

Correction value:

- To correct absorbent beer ingredients extracted with trichloroethane, carry out the procedure with pure trichloroethane (without ZDBT) and measure against trichloroethane as a reference (correction value). Do not use glassware treated with ZDBT solution for this purpose, otherwise the results will be too high due to carry-over of ZDBT.

To determine the calibration curve:

- Pipette 25 mL of decarbonated beer into each of 6 centrifuge tubes
- add 0; 1.0; 2.0; 3.0; 4.0 and 5.0 mL of the diluted copper standard solution (corresponding to 0; 0.4; 0.8; 1.2; 1.6 and 2.0 mg copper per liter of beer)
- Add 1 mL H₂O₂ and 3 mL diluted H₂SO₄ and proceed as described above
- Calculate the conversion factor F ($F = \text{mg Cu}^{2+}/\Delta E$) from the known copper content of the calibration solutions in mg/L divided by the respective determined corresponding absorbance (reduced by the absorbance value of the preparation without copper addition) or determine from calibration curves

18.9 Measurement an result

Call up an measure method 3–57

Calculation of the result: $\text{Copper (mg/L)} = (EH - EB - EK) * F$

EH = Absorbance of the main value

EB = Absorbance of the blank value

EK = Absorbance of the correction value

F = Conversion factor from absorbance to mg/L

18.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.29.5, page 153ff

<https://www.mebak.org/methode/b-420-37-009/kupfer-cuprethol/2261>

19. Photometric iodine assay (MEBAK)

19.1 Background

Iodine is used for the detection of starch. Amylose and amylopectin are components of starch and are broken down into sugars and smaller polysaccharides during the brewing process. If hexagonal structures are still present in the starch, iodine is deposited and colors the solution dark purple.

19.2 Method

Starch and higher molecular weight starch degradation products are precipitated with ethanol and centrifuged. The residue is dissolved in phosphate buffer and then iodine solution is added. The resulting coloration is measured spectrophotometrically.

19.3 Measuring range

0–1 Iodine value

19.4 Cuvette and wavelength

40 mm rectangular cuvette / 578 nm

19.5 Accessories

- Centrifuge tubes with ground glass stopper, 100–110 mL capacity
- Shaking machine
- Pipettes, 0.5 mL, 2 mL, 10 mL, 20 mL, 40 mL
- Spectrophotometer, 578 nm, e.g. REF 919600.1
- 40 mm rectangular cuvette OS
- “Plumper” or plastic spatula

19.6 Reagents

- Ethanol, 95 %
- Iodine solution, 1N (stock solution)

19.7 Procedure

- Iodine solution, 0.02 N (prepare daily from the stock solution)
- Phosphate buffer solution, 0.1 M, pH 3.5: 0.1 M KH_2PO_4 solution with
- Adjust 0.1 M phosphoric acid to pH 3.5

Blank value:

- Pipette 10 mL phosphate buffer solution and 0.5 mL 0.02 N iodine solution into 40 mm rectangular cuvette and mix
- Measure absorbance at 578 nm against phosphate buffer solution (EI)

Main Value:

- Pipette 10.0 mL centrifuged wort or beer free of carbon dioxide into centrifuge tube
- Add 40.0 mL ethanol and shake mechanically for 10 minutes
- Centrifuge for 5 min at 2500 rpm
- Carefully and as completely as possible decant off the clear phase
- Dissolve the residue with 20.0 mL phosphate buffer solution by mechanical shaking for 10 min
- Centrifuge the solution for 5 min at 2500 rpm
- Pipette 2 mL of the supernatant and 8 mL of phosphate buffer solution into a 40 mm rectangular cell and measure at 578 nm against phosphate buffer solution (EZ)
- Add 0.5 mL of 0.02 N iodine solution, mix the contents immediately with the “Plumper”, measure after 30 sec (EH)

19.8 Measurement and result

Call up and measure method 3–58

Calculate the result: $(EH - EI - 0.952 * EZ) * 5$

EH = absorbance of the test solution from the main experiment

EI = absorbance of the iodine solution (blank test)

EZ = absorbance of the centrifugation supernatant

0.952 = Volume change factor (10.0 mL → 10.5 mL)

5 = Dilution factor (2 + 8 parts)

19.9 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5–4, Chapter 2.3.2, page 34ff

<https://www.mebak.org/methode/b-400-04-110/photometrische-iodprobe/645>

20. Total carbohydrates in beer (MEBAK)

20.1 Background

In terms of quantity, carbohydrates make up the largest proportion of nutrients in beer.

20.2 Method

The sample is hydrolyzed and dehydrated and then reacts with anthrone. The resulting coloration is measured with a spectrophotometer.

20.3 Measuring range

0–5 g/L carbohydrate

20.4 Cuvette and wavelengths

10 mm rectangular cuvette / 625 nm

20.5 Accessories

- Test tubes with ground glass stopper (20 × 150 mm)
- Water bath, 95 °C ± 0.5 °C
- Cold bath
- Spectrophotometer, 625 nm, e.g. REF 919600.1
- Glass cuvettes, 10 mm path length, e.g. REF 91933

20.6 Reagents

- Sulphuric acid, 85 vol %: Carefully add 850 mL conc. H₂SO₄, d = 1.84 g/mL, to 150 mL H₂O while cooling and adjust to 1 L after complete cooling to room temperature and mixing with H₂O
- Anthrone reagent: Dissolve 1 g anthrone (Merck 101468) in 85 % H₂SO₄ and make up to 1 L with the acid. Prepare fresh.
- *D*-glucose standard: Dry pure anhydrous *D*-glucose at reduced pressure in a vacuum drying oven for 4 h at 100 °C and 100 mbar.
- *D*-glucose stock solution: Dissolve 0.4 g *D*-glucose anhydrous in H₂O and make up to make up to 1 L.
- *D*-glucose standard solution: Dilute 10 mL stock solution with H₂O to 100 mL. Prepare fresh.

20.7 Sample preparation

- Decarbonize beer, allow foam to collapse - Dilute 2 mL beer with H₂O
- dilute to 500 mL

20.8 Procedure

- Dilute 2 mL beer with H₂O to 500 mL
- Pipette 3 mL of diluted beer into test tube, cool to 2–4 °C
- Add 10 mL of cooled (2–4 °C) anthrone reagent and mix well while cooling
- Place tubes in a water bath at 95 ± 0.5 °C for exactly 20 min - cool quickly to 20 °C
- measure absorbance at 625 nm against an equally treated blank, using 3 mL H₂O instead of the diluted beer sample (EP)
- run 3 glucose standard solutions with each series of analyses, using 3 mL of gluco-se standard solution instead of the diluted beer sample and treating it in the same way (ES)

20.9 Measurement and calculation

Call up and measure method 3–59

Calculation of the result: Total carbon dioxide rate g/100 mL = EP : ES

EP = absorbance of the sample (measured against blank value)

ES = mean value of the absorbance values of the glucose standard solutions (measured against blank value)

20.10 Literature

Brautechnische Analysenmethoden, Volume II, 4th edition, 2002, ISBN 3 9805814 5 – 4, Chapter 2.11, page 85ff

21. Hop Storage Index (HSI) in Whole Hops and Hop Products (MEBAK method)

21.1 Background

Hops and hop products are subject to aging processes which depend upon the conditions during storage and in the surrounding environment. To assess the freshness of hop products, analytical methods capable of detecting aging compounds are necessary. With fresh hops, the HSI is usually below 0.40. However, it can also be above this value under normal conditions when freshly harvested. When assessing hops and hop pellets, the HSI alone does not allow any qualitative statement to be made about the type of change. However, an increase in HSI (delta HSI) can serve as a characteristic of natural change during storage.

21.2 Method

Measurement of the absorption ratio at 275 nm to 325 nm of the hop extract in an alkaline methanol solution.

21.3 Measuring range

0.25 – 0.51 HSI

21.4 Cuvette and wavelength

10-mm rectangular cuvette QS / 275 nm and 325 nm

21.5 Accessories

- Scale, Accuracy 0.05 g
- Glass bottle, 250 mL, solvent-tight screw cap
- Volumetric flasks, 10, 25, 50, 100, 500 mL
- Pipettes, 1, 2, 5, 100 mL (use glass pipettes only!)
- Shaking machine
- Centrifuge (optional)
- Spectrophotometer (UV), z.B. REF 919600.1
- 10 mm rectangular cuvette QS, z.B. REF 919120

21.6 Reagents

- Methanol (absorbance at 275 nm in 10 mm rectangular cuvette against H₂O less than 0.060)
- bidistilled water
- Sodium hydroxide solution, 6 N
- Alkaline methanol (1 mL 6 N sodium hydroxide solution in 500 mL methanol; prepare fresh daily)
- Toluene (absorbance at 275 nm in 10 mm rectangular cuvette against H₂O less than 0.800 after dilution with alkaline methanol in a ratio of 1:100)

21.7 Sample preparation

- Grind the appropriate amount of whole hops or hop product with the hammer mill using a sieve with round holes
- Carefully brush out the hop powder adhering to the inside of the mill and add it quantitatively to the grist
- Always divide the entire available grist
- Insert the beaker into the MEBAK sample divider in an anticlockwise direction
- Put on the dust protection ring and check that the beakers are seated correctly
- Switch on the appliance and slowly pour the hop powder into the funnel
- If 80 g of hops were used, the contents of the beaker after the sample division will be approx. 10 g, which must be used completely for the subsequent analysis.

21.8 Procedure

- Weigh out 5 g of hop powder or ground hop pellets to the nearest 0.05 g in a 250 mL glass bottle
- Add 100 mL toluene
- Close the glass bottle tightly and shake vigorously for 30 minutes (e.g. at 240 rpm)
- Centrifuge or sediment for 5 min at 1.000 g until the solids have settled sufficiently for a clear aliquot to be taken. This should not take longer than 1 hour
- Dilute 5 mL of clarified solution with 100 mL of methanol
- Then further dilute a suitable amount of this solution with alkaline methanol (absorbance at 325 nm should be in the linear measuring range)
- Measure the absorbances at 325 and 275 nm against toluene diluted in the same way as the sample (measure quickly, avoid the influence of UV light)

21.9 Measurement

Call up and measure method 3–61

21.10 Literature

<https://www.mebak.org/en/methode/r-300-12-110/hop-storage-index-hsi-in-whole-hops-and-hop-products/2619>

22. Contact

MACHEREY-NAGEL GmbH & Co. KG

If you have any questions or require technical assistance, please contact us:

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MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Deutschland