

Special information

Version 1

Determination of chloride in concrete

Instructions	
Parameter	Chloride
Product group	VISOCOLOR® and NANOCOLOR®
Reagents and equipment	<p>NANOCOLOR® Chloride 50 (REF 985021), NANOCOLOR® Chloride 200 (REF 985019) or NANOCOLOR® Chloride (REF 91820)</p> <p>As alternative: VISOCOLOR® Chloride CL 500 (REF 915004) or VISOCOLOR® ECO Chloride (REF 931018)</p> <p>Digital piston pipette 1.0 – 5.0 mL, adjustable (REF 916909)</p> <p>Volumetric flask 100 mL with NS 12/21 and PE-stopper (REF 91683)</p> <p>Nitric acid 1 mol/L</p>
Sample preparation	<ol style="list-style-type: none"> 1. Use bore dust (e.g. produced by a metal drill, Ø 30 mm) or milled concrete. 2. Put 4.0 g of the dried concrete sample into a volumetric flask 100 mL. 3. Add 6.0 mL of nitric acid c = 1 mol/L into the volumetric flask, close and mix it. Wait until the silicic acid settles down (i.e. until a cloudy precipitate formed). 4. Fill volumetric flask with distilled water up to mark. 5. Fill the content into a beaker after the precipitation has settled down.
Analysis	<p>Perform the test according instruction leaflet.</p> <p>The read off measurement result from the titration syringe, colour chart or photometer display is translated into % chloride (referred to the concrete mass) regarding to the following equation:</p> <p>Calculation:</p> $\frac{mg / L Cl^{-}}{400} = \% Cl^{-}$

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service:	Tel.:	+49 24 21 969-333
	E-Mail:	support@mn-net.com
For orders (international)	E-Mail:	sales@mn-net.com