

# Diclofenac Sodium and Related Substances – Ph. Eur. monograph 1002

MACHEREY-NAGEL application department · Dr. H. R. Wollseifen, T. Kretschmer, R. Nuessgen

## Application benefits

- HPLC method with faster separations within allowable adjustments.
- Shorter runtimes
- Optimized system suitability

## MN products

### REF 760101.46

EC HPLC column (analytical), NUCLEODUR® C18 Gravity, 5 µm, 250x4.6 mm

### REF 763157.46

EC HPLC column (analytical), NUCLEOSHELL® RP 18, 5 µm, 250x4.6 mm

### REF 763156.46

EC HPLC column (analytical), NUCLEOSHELL® RP 18, 5 µm, 150x4.6 mm

### REF 763136.46

EC HPLC column (analytical), NUCLEOSHELL® RP 18, 2.7 µm, 150x4.6 mm

### REF 763134.46

EC HPLC column (analytical), NUCLEOSHELL® RP 18, 2.7 µm, 100x4.6 mm

### REF 702107

Screw closure, N 9, PP, yellow, center hole, Silicone white/PTFE red, 1.0 mm

### REF 702079

Screw neck vial, N 9, 11.6x32.0 mm, 1.5 mL, label, flat bottom, amber, silanized

## MN application numbers

HPLC: 129520    HPLC: 129550  
HPLC: 129530    HPLC: 129560  
HPLC: 129540

## Keywords

Diclofenac Sodium, Ph. Eur. monograph 1002, NUCLEOSHELL® RP18, NUCLEODUR® Gravity C18, L1, European Pharmacopeia

## Introduction

The Ph. Eur. Monograph 1002 describes the separation of Diclofenac Sodium from impurities. This work starts using a fully porous HPLC phase and shows the benefits using superficially porous particles. The method optimization was performed to achieve shorter run time and system suitability results within allowable adjustments.




Figure 1: Diclofenac Sodium

## Ph. Eur. Monograph 1002 method parameters

| Method Parameter       | Description   |
|------------------------|---|
| Test solution          | Dissolved 50.0 mg of the substance to be examined in the mobile phase and diluted to 50.0 mL with the mobile phase.   |
| Reference solution (a) | Diluted 2.0 mL of the test solution to 100.0 mL with the mobile phase. Diluted 1.0 mL of this solution to 10.0 mL with the mobile phase. Dissolved the contents of a vial of diclofenac for system suitability CRS (containing impurities A and F) in 1 mL of the mobile phase. |
| Reference solution (b) |   |
| Column size            | 250 x 4.6 mm  |
| Stationary phase       | End-capped octadecylsilyl silica gel for chromatography R (5 µm)  |
| Mobile phase           | Mix 34 volumes of a solution containing 0.5 g/L of phosphoric acid R and 0.8 g/L of sodium dihydrogen phosphate R, previously adjusted to pH 2.5 with phosphoric acid R, and 66 volumes of methanol R.  |
| Flow rate              | 1.0 mL/min  |
| Detection              | 254 nm  |
| Injection              | 20 µL   |
| Run time               | 1.6 times the retention time of diclofenac (about 25 min)   |
| Elution order          | Impurity A<br>Impurity F<br>Diclofenac  |
| System suitability     | Minimum resolution of 4.0 between peaks due to Impurity F and Diclofenac  |
| Reference solution (b) |   |

\* Diclofenac for system suitability CRS\* (Y0001635), and Diclofenac sodium CRS (S0765000) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

Table 1: Ph. Eur. Monograph 1002 Diclofenac Sodium Method Details

## Chromatographic methodology improvements

Figure 2: a




Figure 2: b




Figure 2: c




Figure 2: d




Figure 2: e




Figure 2: a: EC HPLC column (analytical), NUCLEODUR® C18 Gravity, 5 µm, 250x4.6 mm, b: EC HPLC column (analytical), NUCLEOSHELL® RP 18, 5 µm, 250x4.6 mm, c: EC HPLC column (analytical), NUCLEOSHELL® RP 18, 5 µm, 150x4.6 mm, d: EC HPLC column (analytical), NUCLEOSHELL® RP 18, 2.7 µm, 150x4.6 mm, e: EC HPLC column (analytical), NUCLEOSHELL® RP 18, 2.7 µm, 100x4.6 mm

## Results

| Method Parameter                 | Allowed Adjustments (isocratic elution)   | Method 1 (figure 2: a) | Method 2 (figure 2: b) | Method 3 (figure 2: c) | Method 4 (figure 2: d) | Method 5 (figure 2: e) |
|----------------------------------|---|------------------------|------------------------|------------------------|------------------------|------------------------|
| Mobile phase pH                  | ± 0.2 units   | As specified           | As specified           | As specified           | As specified           | As specified           |
| Concentration of salts in buffer | ± 10%   | As specified           | As specified           | As specified           | As specified           | As specified           |
| Composition of the mobile phase  | ± 30% of the minor solvent component relative or 2% absolute, whichever is the larger. No other component is altered by more than 10% absolute. | As specified           | As specified           | As specified           | As specified           | As specified           |
| Stationary phase                 | No change of C18 allowed  | NUCLEODUR® C18 Gravity | NUCLEOSHELL® RP 18     | NUCLEOSHELL® RP 18     | NUCLEOSHELL® RP 18     | NUCLEOSHELL® RP 18     |
| Particle size                    | – 50%   | 5 µm as specified      | 5 µm as specified      | 5 µm as specified      | 2.7 µm (– 46%)         | 2.7 µm (– 46%)         |
| Column length                    | ± 70%   | 250 mm as specified    | 250 mm as specified    | 150 mm (– 40%)         | 150 mm (– 40%)         | 150 mm (– 40%)         |
| Column internal diameter         | ± 25%   | 4.6 mm as specified    | 4.6 mm as specified    | 4.6 mm as specified    | 4.6 mm as specified    | 4.6 mm as specified    |
| Flow rate                        | ± 50%   | 1.0 mL as specified    | 1.0 mL as specified    | 1.0 mL as specified    | 1.0 mL as specified    | 1.0 mL as specified    |
| Column temperature               | ± 10 °C   | ambient as specified   | ambient as specified   | ambient as specified   | ambient as specified   | ambient as specified   |

| Method Parameter   | Allowed Adjustments (isocratic elution)   | Method 1 (figure 2: a) | Method 2 (figure 2: b) | Method 3 (figure 2: c) | Method 4 (figure 2: d) | Method 5 (figure 2: e) |
|--------------------|---|------------------------|------------------------|------------------------|------------------------|------------------------|
| Injection volume   | May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory | 20 µL as specified     | 20 µL as specified     | 20 µL as specified     | 20 µL as specified     | 20 µL as specified     |
| Detection [nm]     | No change permitted   | 254 nm as specified    | 254 nm as specified    | 254 nm as specified    | 254 nm as specified    | 254 nm as specified    |
| Run time           | 1.6 times the retention time of Diclofenac  | 52.2 min               | 26.2 min (- 49,8%**)   | 15.4 min (- 70.6%**)   | 18.3 min (- 64.9%**)   | 12.2 min (- 76.6%**)   |
| System suitability | Minimum resolution of 4.0 between peaks due to Impurity F and Diclofenac                                | 6.19                   | 10.30                  | 7.79                   | 9.79                   | 8.06                   |

\* European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques, \*\* runtime reduction in comparison to method 1.

## Conclusion


The fully porous NUCLEODUR® C18 Gravity, 5 µm, 250x4.6 mm HPLC column from MACHEREY-NAGEL fulfills all requirements of the Ph. Eur. monograph 1002. By using superficially porous NUCLEOSHELL® analytical columns the runtime of the method can be reduced by up to 76.6% (with NUCLEOSHELL® RP 18, 2.7 µm, 100x4.6 mm) compared to fully porous NUCLEODUR® silica gel, while keeping all method parameters well within the allowed adjustment range of the European Pharmacopoeia. We

were also able to improve the resolution as well as the peak intensity with NUCLEOSHELL® columns compared to the original method with fully porous silica gel.

MACHEREY-NAGEL shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. Information, descriptions, and specifications in this publication are subject to change without notice.

## The chromatography application database

- Free access to more than 3,000 application examples from HPLC, GC, TLC and SPE
- Free access to the MN application database: <https://chromaappdb.mn-net.com>



[www.mn-net.com](http://www.mn-net.com)

# MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG  
Valenciennner Str. 11  
52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com