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BOD₅-Nutrient Mixture (without N-allylthiourea)

Method:

The special chemical set **NANOCOLOR[®] BOD₅-Nutrient Mixture** is required for the photometric determination of BOD₅ with tube test **NANOCOLOR[®] BOD₅ - Test 8-22** (REF 985822) and **contains no nitrification inhibiting compound N-allylthiourea**. The specified concentrations of nutrients are in accordance with the recommendations of the German and European Standard Method **DIN EN 1899-H51**. Additional accessory for the determination of BOD₅ is the BOD₅ Accessories Set (REF 916918).

MACHEREY-NAGEL additionally offers a nutrient mixture containing the nitrification inhibiting compound **N-allylthiourea**, called "**BOD₅-Nutrient Mixture PLUS**" (REF 918995).

Contents of reagent set:

10 test tubes, clear glass with reagent R1, lyophilized
10 test tubes, brown glass with reagent R2, lyophilized

Hazard warning:

This reagent set does not contain any harmful substances which must be specially labelled as hazardous.

Procedure:

Always use **freshly prepared BOD₅-Nutrient Solutions!**

Requisite accessories: piston pipette with tips

Preparation of BOD₅-Nutrient Solutions R1 and R2:

1.	Open a clear glass vial with reagent R1 , add 10 mL aerated diluting water (see BOD ₅ -Accessories Set, REF 916918), close and shake to dissolve. Now you have BOD₅-Nutrient Solution R1 .
2.	Open a brown glass vial with reagent R2 , add 10 mL aerated diluting water (see BOD ₅ -Accessories Set, REF 916918), close and shake to dissolve. Now you have BOD₅-Nutrient Solution R2 .
3.	In connection with the BOD ₅ determination, add 0.25 mL BOD₅-Nutrient Solution R1 and 0.25 mL BOD₅-Nutrient Solution R2 for every 100 mL sample dilution. <i>For further information see package inserts of tube test NANOCOLOR[®] BOD₅ - Test 8-22 (REF 985822).</i>
<i>At the end of the day discard all BOD₅-Nutrient Solutions that are no longer needed!</i>	

Preparation of inoculated BOD₅-Nutrient Solutions R1 and R2:

If the sample itself does not contain a sufficient amount of appropriate microorganisms, so-called **inoculating water** is added to the BOD₅ test preparation via the BOD₅-Nutrient Solutions. This inoculating water can be obtained in the following ways:

- precipitated urban waste water from presettling tank outflow or mechanical inflow
- urban waste water from sewage disposal with max. 300 mg/L COD or 100 mg/L TOC
- surface water containing urban waste water

1.	Open a clear glass vial with reagent R1 , add 10 mL inoculating water , close and shake to dissolve. Now you have inoculated BOD₅-Nutrient Solution R1 .
2.	Open a brown glass vial with reagent R2 , add 10 mL inoculating water , close and shake to dissolve. Now you have inoculated BOD₅-Nutrient Solution R2 .
3.	In connection with the BOD ₅ determination, add 0.25 mL inoculated BOD₅-Nutrient Solution R1 and 0.25 mL inoculated BOD₅-Nutrient Solution R2 for every 100 mL sample dilution. <i>For further information see package inserts of tube test NANOCOLOR[®] BOD₅ - Test 8-22 (REF 985822).</i>
<i>At the end of the day discard all inoculated BOD₅-Nutrient Solutions that are no longer needed!</i>	

Note: The oxygen consumption in **control solutions with inoculated BOD₅-Nutrient Solutions** shall not exceed a value of 1.5 mg/L O₂ after 5 days of incubation at (20 ± 1) °C.

Examples of samples for which inoculated BOD₅-Nutrient Solutions need to be used:

- industrial, biological inactive waste water
- surface water containing no urban waste water
- deep frozen samples!
- control tests using a glucose / glutamine acid standard solution according to DIN EN 1899-1-H51 (**NANOCONTROL BOD₅-Standard**, REF 92582)!

Storage:

The reagent set is to be stored in dry conditions at room temperature. The expiry date printed on the kit is to be observed.

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