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en

IVD *In vitro* diagnostic medical device, suitable for self-testing

CE 0197

Instructions for use Medi-Test Glucose

REF 93001.005: 50 test strips in an aluminum box with desiccant stopper.

Intended use

Medi-Test Glucose urine test strips are test strips which are suitable for measuring the glucose content of human urine. The test strips yield semiquantitative results in the range of 0 – 1000 mg/dL. The test strips can be read by comparing them visually with the color scale. The results are used as a diagnostic aid for screening for diabetes mellitus and renal glycosuria, for example, during pregnancy. The urine test strips are suitable for self-testing.

Notes

Discuss any positive and any unclear result with your doctor. You should also do this if this test does not show any results but symptoms persist. If a false negative finding is suspected, repeat the test with a test strip from a previously unopened tin or ask your doctor. Your doctor will perform additional testing with the usual methods in his/her laboratory. In general, individual test strip results allow a definitive diagnosis and targeted therapy only in connection with other medical findings.

Use only urine samples which have stood for a maximum of 4 hours.

Use only clean containers for urine collection which are free of residues. For removal of excess liquid an absorbant paper tissue is required.

Substances which cause abnormal urine coloration can impair the analysis of the test strips. Do not touch the reaction zone. Always remove only the number of test strips needed. Immediately close the package tightly after removing test strips. Do not remove desiccant in stopper; do not damage stopper.

Do not use damaged test strips or tins.

Users with impaired color vision must be assisted by a person with normal color vision for the color comparison.

Store test strips out of the reach of children.

For use only outside of the body.

Do not reuse test strips!

Never make changes to the administration of prescribed medications based on a result from these test strips.

Quality control when used by professional staff

The test strips should be verified only with positive and negative control solutions. Medi-Test Control (REF 93038) is recommended as a control solution. The positive and negative controls should be performed when using a new test strip batch, and after 30 days in each case to verify the storage conditions. Each laboratory should determine its own target

values for adequate performance standards and review the test methods and sequences if these standards are not met.

Parameter

GLU **Glucose**

The detection is based on the glucose oxidase-peroxidase chromogen reaction. Except for glucose, no urine constituent which returns a positive reaction is known. Pathological glucose concentrations are indicated by a color change from green to blue-green. The test detects values starting from 30 mg glucose/dL urine. Yellow to pale green test fields should be evaluated as negative (or normal). The color comparison fields correspond to the following glucose concentrations:

neg. (yellow) · normal (yellow-green) · 50 · 150 · 500 · \geq 1000 mg/dL or

neg. (yellow) · normal (yellow-green) · 2.8 · 8.3 · 27.8 · \geq 55.5 mmol/L

Reactive substances*: Glucose oxidase 7 U; peroxidase 1 U; tetramethylbenzidine 96 µg.

Normal concentrations of ascorbic acid (\leq 40 mg/dL) do not influence the test result**.

False-positive reactions can be caused by oxidising cleaning agents in the sample container. Gentic acid does not interfere with the detection reaction.

In a comparison study (N = 125) with a reference method (visual evaluation of competitor test strips), a direct correspondence of 98 %, a diagnostic sensitivity of 100 % and a diagnostic specificity of 100 % was determined. LOT to LOT precision (N = 120), repeatability (N = 120), reproducibility (N = 120) and the negative/positive predictive value (N = 125) was determined to be 100 % in all cases. The positive likelihood ratio (LR+) is not applicable and the negative likelihood ratio (LR-) is 0. In a lay person study with 97 participants, an exact agreement of 96 % and an agreement of 100 % (target value +/- one bin) was obtained for three glucose levels (neg.; 50 and 150 mg/dL).

* quantity/cm² after impregnation.

** Interference study with urine with pathological finding (first positive scale value).

Shelf life

Protect test strips from sunlight and moisture. Store tin in a cool and dry location (storage temperature 4 – 30 °C). If stored properly, the test strips can be stored until the printed expiry date.

Disposal

Dispose of the used test strips taking applicable safety regulations into account.

Information on reporting obligation if incidents occur

We wish to point out that all serious incidents which occur in connection with the product must be reported to the manufacturer and the competent authority of the European member state or of the state in which the incident occurred. European vigilance contact points:
https://ec.europa.eu/health/md_sector/contact_en.

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Technical service

If you still have questions after reading the instructions or need technical assistance, please contact

MACHEREY-NAGEL GmbH & Co. KG

Valencienner Str. 11 · 52355 Düren · Germany

Tel.: +49 24 21 969-0 · Email: info@mn-net.com

Website: www.mn-net.com

Explanation of symbols

-  Declaration of Conformity
-  Device for self-testing
-  Please read instructions for use
-  In vitro diagnostic medical device
-  Manufacturer
-  Item number
-  Keep away from sunlight
-  Batch identification
-  Store in a dry place
-  Use by
-  Contains sufficient for <n> tests
-  Do not reuse
-  Permitted storage temperature range

Revision: 05/2023

Reason of revision: Addition of interfering substances and performance data.

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Assay procedure

Follow this sequence exactly:

- For control of dipping and reaction time use a timer accurate to the second. Note date and time.
- Open container. Remove one test strip only. Do not touch reaction zone/test field!
- Close container tightly immediately after removing test strip.
- Immerse test strip in the urine for approximately 1 second (fig. 1). The test field must be submerged.
- After removing the test strip from the urine sample, briefly dab the lateral edge onto an absorbent paper (fig. 2). Do not put the test strip down.
- Wait for 30 seconds (fig 3).
- Start the evaluation: On the container label a color sequence for the test field can be found. Assign the test field to a color value of its sequence. Find a match or select the color value that comes as close as possible to the test field (fig. 4).
- Note the result and submit it to your physician at your next visit.
- Finish the evaluation within 60 seconds after immersion at the latest. Observe this time limit, since the test field still changes its color intensity after this period.

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