MACHEREY-NAGEL

NANOCOLOR<sup>®</sup> TIC-Ex

- Handbuch
- Manual
- Manuel

<b>Deutsch</b> .....	<b>3</b>
<b>English</b> .....	<b>9</b>
<b>Français</b> .....	<b>15</b>
<b>Español</b> .....	<b>21</b>

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# CONTENTS

<b>1. Introduction</b> .....	<b>10</b>
<b>2. CE mark</b> .....	<b>10</b>
<b>3. Explanation of symbols</b> .....	<b>10</b>
<b>4. Safety instructions</b> .....	<b>10</b>
<b>5. Technical data</b> .....	<b>11</b>
<b>6. Setup and commissioning</b> .....	<b>11</b>
6.1 Working principle .....	11
6.2 Package content .....	11
6.3 Setting up the instrument .....	11
6.4 Operating elements .....	12
6.5 How to plug the instrument in .....	12
<b>7. Handling of the device</b> .....	<b>13</b>
<b>8. Cleaning of the device</b> .....	<b>13</b>
<b>9. Service</b> .....	<b>13</b>
9.1 Disposal .....	13
9.2 Warranty .....	13
9.3 Contact .....	14

## 1. Introduction

The *NANOCOLOR*<sup>®</sup> TIC-Ex is required during the sample preparation of MACHEREY-NAGEL TOC tube tests. For evaluation of the total organic carbon (TOC) prior to determination the inorganic carbon needs to be removed from the sample. The removal is accomplished by purging of air with help of *NANOCOLOR*<sup>®</sup> TIC-Ex through the prepared sample solution.

## 2. CE mark

The CE mark declares that the product complies with the harmonization legislation of the European Community listed below:

- European Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS 2)
- European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)
- European Directive 2014/30/EU on the harmonization of the laws of the member states relating to electromagnetic compatibility (EMC)
- European Directive 2014/35/EU on the harmonization of the laws of the member states relating to the making available on the market of electrical equipment designed for use within certain voltage limits (NSR)

## 3. Explanation of symbols

Description of the used symbols:

Caution: Follow the instructions in this section carefully to ensure the device does not become damaged.

Warning: The instructions in this section must be followed carefully to ensure that no further dangers or hazards exist, otherwise the device will very likely become damaged.

## 4. Safety instructions

Caution: The *NANOCOLOR*<sup>®</sup> TIC-Ex may only be used as specified by the manufacturer and only be correspondingly qualified personnel. Note all markings and signs mounted on the device or stated in the operating manual. Ignoring them can lead to personal injuries and/or damage to the device.

## 5. Technical data

Required power supply:	external mains adapter: Input 100 V–240 V; Output 9 V, 2.0 A
Dimensions:	350 mm x 300 mm x 235 mm
Weight:	4.95 kg (without power supply)
Operating environment (Operation)	5 °C–40 °C, 20%–80 % relative humidity
Operating environment (Transport, storage)	5 °C–40 °C, 20%–80 % relative humidity

## 6. Setup and commissioning

### 6.1 Working principle

The acidified sample solutions are put into the cuvette holder and placed in the premarked positions of the *NANOCOLOR*<sup>®</sup> TIC-Ex. The preinstalled pipette tips are dipped into the sample solutions by pressing down the lowering mechanism of the *NANOCOLOR*<sup>®</sup> TIC-Ex. After switching on the device, air is pumped through the sample solutions in the cuvettes for five minutes. Hereby the inorganic carbon is removed as carbon dioxide quantitatively. After elapsed time the process stops automatically. The pipette tips are removed from the solutions and can be disposed using the removing mechanism on top. Afterwards the sample solutions are treated following the next steps of the respective TOC test kit.

### 6.2 Package content

*NANOCOLOR*<sup>®</sup> TIC-Ex

Cuvette holder

Mains adapter 100–240 V, 50/60 Hz, 9 V  $\overline{\text{---}}$  + 3 adapters

User manual (this leaflet)

### 6.3 Setting up the instrument

To secure a problem-free sample purging process, please note the following remarks:

Attention:

Place the instrument on a hard and even surface.

Do not move the device during purging process.

Place the device in a dust free environment.

## 6.4 Operating elements

Operating element	Function
① Start button	Starts the air flow and the five minute timer
② Stop button	Stops air flow and resets the timer
③ LED	Shows device status
④ Removal mechanism	Mechanism to remove all pipette tips at once
⑤ Lowering mechanism	Mechanism for lowering and lifting of the purging unit
⑥ Cuvette holder	Holder for cuvettes during the purging process
⑦ Power plug	Connection of mains adapter
⑧ On/Off switch (I/O)	Switches device on and off


Fig. 1: Front view

Fig. 2: Rear view

## 6.5 How to plug the instrument in

Three adapters are provided for connecting the power pack to the available mains connection. The adapter matching the mains connection is plugged on to the power pack. After plugging the power pack cable into the jack “DC 9 V” (chapter 6.4 ⑦) and connecting the power pack to the power socket, the *NANOCOLOR*<sup>®</sup> TIC-Ex is ready for operation.

Warning: Electrical danger

The use of *NANOCOLOR*<sup>®</sup> TIC-Ex is only permitted to qualified personnel. Please consider the applicable safety instructions. Prior to each use make sure that the mains adapter and the power connection do not show any obvious defects, like cable breakage. In case of an issue the defective part has to be replaced to prevent usage.

## 7. Handling of the device

The handling of the **NANOCOLOR® TIC-Ex** is easy and intuitive. Turn on the device with the switch at the back of the instrument (chapter 6.4 ⑧). Assemble the needed amount of pipette tips (part of TOC tube tests) on the purging unit of the device. Place the prepared TOC sample tubes in the cuvette holder and insert the holder into the indentations of the **NANOCOLOR® TIC-Ex** housing. Insert the pipette tips into the samples by pushing down the purging unit (chapter 6.4 ⑤) with both hands until it is locked. The green LED (chapter 6.4 ③) indicates that the device is ready for use. To start the purging process, press the start button  (chapter 6.4 ①). The device starts pumping air through the samples for five minutes. The time is counted internally. The blinking orange LED (chapter 6.4 ③) indicates that the device is working. After five minutes an acoustic signal indicates the completed process; the air flow is stopped automatically. Lift up the purging unit by shortly pressing downwards. Proceed with the samples as described in the manual of the TOC tube test. The pipette tips can be removed by simultaneously pressing down the black handles on top (chapter 6.4 ④).

Pressing the stop button  (chapter 6.4 ②) will abort the process at any time. The device will directly stop the air flow and reset the timer. The red LED (chapter 6.4 ③) indicates that the process was aborted by the user. Pressing the stop button  a second time, restores the initial status. The green LED (chapter 6.4 ③) indicates that the device is ready for use again.

## 8. Cleaning of the device

Attention:

The instrument housing may be wiped with a wet cloth. Splashes and spillings of liquid need to be removed from the instrument immediately. The pipette tips are consumables and must be equipped freshly for each process (part of TOC tube tests).

## 9. Service

### 9.1 Disposal

Dispose in compliance with EU-directive 2012/19/EU (WEEE). In compliance with EU-directive 2012/19/EU (WEEE) MACHEREY-NAGEL disposes old instruments free of charge.

*Note: Disposal using public waste disposal facilities is not permitted. In case of disposal, please contact your local MACHEREY-NAGEL representative.*

### 9.2 Warranty

The warranty for this equipment has a duration of 24 months from the date of purchase. The original copy of the bill serves as a certificate and must be submitted in case of assertion of a warranty claim. The warranty expires in case of improper handling and/or maintenance of the equipment; it does not comprise defects due to the external power supply. The warranty is limited to the repair of faulty parts or - at our sole discretion - to the delivery of a faultless sub-

stitute. The warranty period of 24 months is not affected by claiming on the warranty during this period. There is no right of withdrawal. Further claims are excluded. Hereunto we count in particular all claims for damages evolving from consequential ages or indirect damages. Additionally the relevant version of our general sales terms and delivery conditions apply as printed on all price lists.

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MACHEREY-NAGEL GmbH & Co. KG
MACHEREY-NAGEL GmbH & Co. KG · Valenciener Str. 11 · 52355 Düren · Germany