

MACHEREY-NAGEL

User manual



■ NucleoTrap® CR

■ NucleoTrap®

October 2022 / Rev. 08

Contact MN

Germany and international

MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0
Toll-free: 0800 26 16 000 (Germany only)
E-mail: info@mn-net.com

Technical Support Bioanalysis

Tel.: +49 24 21 969-333
E-mail: support@mn-net.com

USA

MACHEREY-NAGEL Inc.
924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA
Toll-free: 888 321 6224 (MACH)
E-mail: sales-us@mn-net.com

France

MACHEREY-NAGEL SAS
1, rue Gutenberg – BP135 · 67720 Hoerdt Cedex · France
Tel.: +33 388 68 22 68
E-mail: sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €
Siret 379 859 531 00020 · RCS Strasbourg B379859531 · N° intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG
Hirsackerstr. 7 · 4702 Oensingen · Switzerland
Tel.: +41 62 388 55 00
E-mail: sales-ch@mn-net.com

PCR clean up, gel extraction

Protocol at a glance (Rev. 08)









	NucleoTrap®CR PCR clean up		NucleoTrap® Gel extraction
1 NucleoTrap®: Excise DNA fragment/ Solubilize gel slice			
NucleoTrap®CR: Adjust binding conditions	4 vol NT2/1 vol sample		300 µL NT1/100 mg gel
2 Bind DNA	10 µL silica matrix / 100 µL sample RT, 10 min 10,000 x g, 30 s	 	4 µL silica matrix/µg DNA 50 °C, 5 – 10 min 10,000 x g, 30 s
3 Wash silica matrix	1st 400 µL NT2 2nd 400 µL NT3 3rd 400 µL NT3 10,000 x g, 30 s 10,000 x g, 30 s 10,000 x g, 30 s	 	1st 500 µL NT2 2nd 500 µL NT3 3rd 500 µL NT3 10,000 x g, 30 s 10,000 x g, 30 s 10,000 x g, 30 s
4 Dry silica matrix	RT or 37 °C, 10 – 15 min		RT or 37 °C, 10 – 15 min
5 Elute DNA	25 – 50 µL NE RT, 10 – 15 min 10,000 x g, 30 s	 	25 – 50 µL NE RT, 10 – 15 min 10,000 x g, 30 s

Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Consumables and equipment to be supplied by the user	5
1.3	About this user manual	5
2	Product description	5
2.1	The basic principle	5
2.2	Kit specifications	6
2.3	Elution procedures	7
3	Storage conditions and preparation of working solutions	8
4	Safety instructions	9
4.1	Disposal	9
5	NucleoTrap [®] CR protocol – direct purification of PCR products	10
6	NucleoTrap [®] protocol – DNA extraction from agarose gels	12
7	Protocol for concentration, desalting, removal of enzymes, etc.	14
8	Appendix	15
8.1	Troubleshooting	15
8.2	Ordering information	16
8.3	References	16
8.4	Product use restriction/warranty	17

1 Components

1.1 Kit contents

NucleoTrap®CR		
REF	10 preps 740587.10	100 preps 740587
NucleoTrap®CR Suspension	100 µL	1000 µL
Buffer NT2	10 mL	100 mL
Wash Buffer NT3 (Concentrate)*	6 mL	25 mL
Elution Buffer NE**	13 mL	13 mL
User manual	1	1

NucleoTrap®		
REF	10 preps 740584.10	100 preps 740584
NucleoTrap® Suspension	100 µL	1000 µL
Buffer NT1	6 mL	60 mL
Buffer NT2	10 mL	100 mL
Wash Buffer NT3 (Concentrate)*	6 mL	25 mL
Elution Buffer NE**	13 mL	13 mL
User manual	1	1

*For preparation of working solutions and storage conditions see section 3.

** Composition of Elution Buffer NE: 5 mM Tris/HCl, pH 8.5

1.2 Consumables and equipment to be supplied by the user

Consumables:

- 96–100 % ethanol
- 1.5 mL microcentrifuge tubes

Equipment:

- Centrifuge for microcentrifuge tubes
- Manual pipettors and disposable tips
- Vortex mixer
- Heating-block
- Personal protection equipment (lab coat, gloves, goggles)

1.3 About this user manual

It is strongly recommended to read the detailed protocol sections of this user manual if using the NucleoTrap[®]CR/NucleoTrap[®] kits for the first time. However, experienced users may refer to the Protocol at a glance. The Protocol at a glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure.

All technical literature is available on the internet at www.mn-net.com.

Please contact Technical Service regarding information about changes of the current user manual compared to previous revisions.

2 Product description

2.1 The basic principle

With the NucleoTrap[®]CR/Trap[®] method, DNA binds in the presence of chaotropic salts (Buffer NT1 and Buffer NT2) to specially activated silica particles (matrix). Buffer NT1 contains additional components in order to dissolve agarose gel slices. Afterwards, the NucleoTrap[®]CR/Trap[®] matrix is added to the binding mixtures. Contaminations like salts and soluble macromolecular components are removed by a simple washing step with ethanolic Wash Buffer NT3. Pure DNA is finally eluted under low ionic strength conditions with slightly alkaline Elution Buffer NE (5 mM Tris-Cl, pH 8.5).

2.2 Kit specifications

- The NucleoTrap[®]CR kit is designed for direct purification of PCR products.
- The NucleoTrap[®] kit is designed for the purification of DNA from TAE / TBE agarose gels.
- In contrast to the NucleoTrap[®] matrix, the NucleoTrap[®]CR matrix will not bind DNA fragments < 100 bp due to a larger pore size of the silica matrix.
- Standard as well as low melting agarose gels can be used.
- The prepared DNA fragments can be used directly in applications like automated fluorescent DNA sequencing, PCR, or any kind of enzymatic manipulation.

Table 1: Kit specifications at a glance

Parameter	NucleoTrap [®] CR	NucleoTrap [®]
Technology	Silica matrix	Silica matrix
Format	Silica bead suspension	Silica bead suspension
Sample material	< 400 µL PCR reaction mixture	< 200 mg agarose gel
Use	For research use only	For research use only
DNA fragments from agarose gels	-	++
Desalting, removal of enzymes, nucleotides and/or labeling reagents like biotin or radioactive ATP, etc.	++	+
Direct purification of amplified DNA	++	-
Fragment size	100 bp–approx. 50 kbp	20 bp–approx. 50 kbp
Typical recovery	70–80 %	50–90 %
A ₂₆₀ /A ₂₈₀	1.7–1.9	1.7–1.9
Elution volume		20–50 µL
Preparation time		60 min/6 preps
Binding capacity		6 µg/10 µL matrix

-not recommended +possible ++optimal

2.3 Elution procedures

For the elution of DNA one of the following solutions can be used: Buffer NE (supplied)/TE buffer, pH 8.5/distilled water, pH 8.5.

If water is used, the pH should be checked and adjusted to pH 8–8.5 since deionized water usually exhibits a pH below 7. Furthermore, absorption of CO₂ Leads to a decrease in pH of unbuffered solutions.

Note: EDTA in TE buffer may cause problems in subsequent reactions. See Table 2 for correlation between fragment size and typical recoveries for purification of 1–5 µg of PCR fragments (for gel extraction, recoveries are approximately 10 % lower).

Table 2: DNA recovery with NucleoTrap[®]CR / NucleoTrap[®]

Fragment length	NucleoTrap [®] CR	NucleoTrap [®]
20 bp	0 %	50 %
40 bp	0 %	68 %
120 bp	68 %	78 %
200 bp	76 %	85 %
520 bp	80 %	87 %
2.5 kbp	81 %	88 %
5.3 kbp	80 %	86 %
8.7 kbp	76 %	80 %
19.4 kbp	74 %	74 %

3 Storage conditions and preparation of working solutions

Attention: Buffers NT1 and NT2 contain chaotropic salts! Wear gloves and goggles!

The NucleoTrap[®]CR/NucleoTrap[®] kits should be stored at room temperature and are stable until: see package label.

Before starting any NucleoTrap[®]CR/NucleoTrap[®] protocol prepare the following:

Wash Buffer NT3: Add the indicated volume of 96–100% ethanol to Wash Buffer NT3 Concentrate.

NucleoTrap [®] CR		
REF	10 preps 740587.10	100 preps 740587
Wash Buffer NT3 (Concentrate)	6 mL Add 24 mL ethanol	25 mL Add 100 mL ethanol

NucleoTrap [®]		
REF	10 preps 740584.10	100 preps 740584
Wash Buffer NT3 (Concentrate)	6 mL Add 24 mL ethanol	25 mL Add 100 mL ethanol

4 Safety instructions

When working with the **NucleoTrap®CR/NucleoTrap®** kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles). For more information consult the appropriate Material Safety Data Sheets (MSDS available online at <http://www.mn-net.com/msds>).



4.1 Disposal

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

5 NucleoTrap[®]CR protocol – direct purification of PCR products

Before starting the preparation:

- Check if Wash Buffer NT3 was prepared according to section 3.

1 Adjust DNA binding conditions

Add **4 volumes** of **Buffer NT2** to **1 volume** of sample (e.g., 400 μ L Buffer NT2 and 100 μ L PCR reaction mixture).



**+4 vol NT2 per
1 vol sample**

For sample volumes < 100 μ L adjust the volume of the reaction mix to 100 μ L using TE buffer (pH 7.5).

Note: If the volume of the PCR reaction mixture is > 100 μ L, the volumes of Buffer NT2 and NucleoTrap[®]CR Suspension must be increased proportionally. Example: a volume of 150 μ L reaction mixture needs 600 μ L of Buffer NT2, and 15 μ L NucleoTrap[®]CR Suspension to adjust proper binding conditions.

2 Bind DNA

Vortex the NucleoTrap[®]CR Suspension thoroughly until a homogeneous mixture results. Add **10 μ L** of **NucleoTrap[®]CR Suspension** to each **100 μ L** of **reaction mixture**.



**10 μ L
silica matrix
RT, 10 min**

Incubate the mixture for **10 min** at **room temperature** and vortex briefly every 2–3 min.



**10,000 x g,
30 s**

Centrifuge the sample at **10,000 x g** for **30 s** and discard the supernatant.

3 Wash silica matrix

1st wash

Add **400 μ L Buffer NT2** to the pelleted silica matrix and vortex briefly for resuspension of the pellet. Centrifuge for **30 s** at **10,000 x g** and remove the supernatant completely.



+400 μ L NT2



**10,000 x g,
30 s**

2nd wash

Add **400 μ L Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 s** at **10,000 x g** and remove the supernatant completely.



+400 μ L NT3



**10,000 x g,
30 s**

3rd wash

Add **400 μ L Buffer NT3** to the pelleted silica matrix and vortex briefly. Centrifuge for **30 s** at **10,000 x g**. Remove the supernatant and centrifuge the pellet again briefly. Remove residual Buffer NT3 completely.



+400 μ L NT3



**10,000 x g,
30 s**

4 Dry silica matrix

Dry the pelleted silica matrix at **room temperature** or at **37 °C** for **10–15 min**.

**RT or 37 °C,
10–15 min**

It is not recommended to dry the sample by vacuum since over-dried pellets lead to lower recoveries.

Residual ethanol from Buffer NT3 would inhibit subsequent reactions and has to be removed in this step.

5 Elute DNA

Add **25–50 µL Buffer NE** to the silica matrix. Resuspend the pellet by vortexing. Incubate the mixture at **room temperature** for **10–15 min**. Vortexing the mixture 2–3 times during incubation is recommended. Centrifuge the sample at **10,000 x g** for **30 s** and transfer the DNA containing supernatant to a clean tube (not provided). Repeating this step will increase the yield by approximately 10%.



+25–50 µL NE

**10,000 x g,
1 min**

Yield of larger fragments (> 5–20 kbp) can be increased by performing the incubation at 55 °C.

6 NucleoTrap® protocol – DNA extraction from agarose gels

Before starting the preparation:

- Check if Wash Buffer NT3 was prepared according to section 3.
- Set heating block to 50 °C.

1 Excise DNA fragment / Solubilize gel slice

Take a clean scalpel to excise the DNA fragment from agarose gel. Excise gel slice containing the fragment carefully to minimize the gel volume. Determine the weight of the gel slice and transfer it to a clean tube (not provided).



For each **100 mg** agarose gel add **300 µL NT1**.

For gels containing > 2% agarose, double the volume of Buffer NT1.

Note: If the weight of the gel slice is > 100 mg, the volume of Buffer NT1 must be increased proportionally. Example: a 150 mg gel slice (< 2% agarose) needs 450 µL Buffer NT1.



+300 µL NT1

2 Bind DNA

Vortex the NucleoTrap® Suspension thoroughly until a homogeneous mixture results. **For each µg of DNA** add **4 µL** of the **NucleoTrap® Suspension**, but at least 10 µL.

Incubate sample at 50 °C until the gel slice is dissolved (5–10 min). Vortex the sample briefly every 2–3 min until the gel slice is dissolved completely.

Centrifuge for **30 s** at **10,000 x g** and discard supernatant.



4 µL silica matrix/µg DNA

**50 °C,
5–10 min**



**10,000 x g,
30 s**

3 Wash silica matrix**1st wash**

Add 500 µL Buffer NT2 to the pelleted silica matrix and vortex briefly for resuspension of the pellet. Centrifuge for 30 s at 10,000 x g and remove the supernatant completely.



+500 µL NT2

10,000 x g,
30 s**2nd wash**

Add 500 µL Buffer NT3 to the pelleted silica matrix and vortex briefly. Centrifuge for 30 s at 10,000 x g and remove the supernatant completely.



+500 µL NT3

10,000 x g,
30 s**3rd wash**

Add 500 µL Buffer NT3 to the pelleted silica matrix and vortex briefly. Centrifuge for 30 s at 10,000 x g. Remove the supernatant and centrifuge the pellet again briefly. Remove residual Buffer NT3 completely.



+500 µL NT3

10,000 x g,
30 s**4 Dry silica matrix**

Dry the pelleted silica matrix at **room temperature** or at **37 °C** for **10–15 min**.

RT or 37 °C,
10–15 min

It is not recommended to dry the sample by vacuum since over-dried pellets lead to lower recoveries.

Residual ethanol from Buffer NT3 would inhibit subsequent reactions and has to be removed in this step.

5 Elute DNA

Add **25–50 µL Buffer NE** to the silica matrix. Resuspend the pellet by vortexing. Incubate the mixture at **room temperature** for **10–15 min**. Vortexing the mixture 2–3 times during incubation is recommended. Centrifuge the sample at **10,000 x g** for **30 s** and transfer the DNA containing supernatant to a clean tube (not provided). Repeating this step will increase the yield by approximately 10%.



+20–50 µL NE

10,000 x g,
30 s

Yield of larger fragments (> 5–20 kbp) can be increased by performing the incubation at 55 °C.

7 Protocol for concentration, desalting, removal of enzymes, etc.

Before starting the preparation:

- Check if Wash Buffer NT3 was prepared according to section 3.

1 Adjust DNA binding conditions

Add **4 volumes Buffer NT2** to **1 volume** of DNA containing **sample** (e.g., 400 μL Buffer NT2 and 100 μL reaction mixture).



**+4 vol NT2 per
1 vol sample**

2 Bind DNA

Vortex the NucleoTrap[®]CR / NucleoTrap[®] Suspension thoroughly until a homogeneous mixture results. For **each μg of DNA** add **4 μL of silica matrix, but at least 10 μL** .

Incubate the mixture for **10 min** at **room temperature** and vortex briefly every 2–3 min.

Centrifuge for **30 s** at 10,000 $\times g$ and discard supernatant.

Important note: Be aware of the NucleoTrap[®] Suspension binding fragments down to 20 bp (see Table 2, section 2.3).



**4 μL silica
matrix/ μg DNA**

**50 °C,
5–10 min**



**10,000 $\times g$,
30 s**

Continue with **section 5, step 3**.

8 Appendix

8.1 Troubleshooting

Problem	Possible cause and suggestions
Incomplete lysis of agarose slices	<i>High concentration of agarose</i>
	<ul style="list-style-type: none"> Use doubled volumes of Buffer NT1 for highly concentrated agarose gels.
	<i>Wrong buffer</i>
	<ul style="list-style-type: none"> Buffer NT2 cannot be used for gel dissolution.
	<i>Time and temperature</i>
	<ul style="list-style-type: none"> Check incubation temperature. Depending on the weight of gel slice, incubation (section 6, step 2) can be prolonged up to 20 min. Vortex every 2 min and check integrity of the gel slice. Heavy weight gel slices may be quenched or crushed before addition of Buffer NT1.
	<i>Reagents not applied properly</i>
	<ul style="list-style-type: none"> Add indicated volume of 96 – 100 % ethanol to Wash Buffer NT3 Concentrate and mix well before use.
No DNA yield	<i>Insufficient drying of the NucleoTrap®CR/NucleoTrap® silica matrix</i>
	<ul style="list-style-type: none"> Ethanol Wash Buffer NT3 has to be removed quantitatively before elution. Prolong the drying time up to 30 min. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots).
	<i>Isolation of large DNA fragments</i>
	<ul style="list-style-type: none"> Add room-temperature Elution Buffer NE and incubate at 55 °C for 10 – 15 min.

Problem	Possible cause and suggestions
Suboptimal performance of DNA in sequencing reactions	<i>Carry-over of ethanol/ethanolic Buffer NT3</i>
	<ul style="list-style-type: none"> Make sure to dry the silica matrix in order to achieve complete removal of ethanolic Buffer NT3 after the washing step. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots).
	<ul style="list-style-type: none"> Buffers other than Buffer NE, for example TE buffer (Tris/EDTA), were used for elution of DNA. Note: EDTA may inhibit sequencing reactions. In this case it is recommended to re-purify DNA and elute in Buffer NE or water.
	<i>Not enough DNA used for sequencing reaction</i>
	<ul style="list-style-type: none"> Quantitate DNA by agarose gel electrophoresis before setting up sequencing reactions.
	<i>NucleoTrap®CR or NucleoTrap® particles were not removed quantitatively</i>
	<ul style="list-style-type: none"> Centrifuge the eluate again and transfer the supernatant to a new tube.

8.2 Ordering information

Product	REF	Pack of
NucleoTrap®CR	740587.10	10 preps
	740587	100 preps
NucleoTrap®	740584	10 preps
	740584	100 preps
NucleoTrap® Suspension	740568	100 preps
Buffer NT1	740596.100	2 × 50 mL
Buffer NT2	740597	2 × 50 mL
Buffer NT3 Concentrate (for 100 mL Buffer NT3)	740598	20 mL
Collection Tubes (2 mL)	740600	1000

Visit www.mn-net.com for more detailed product information.

8.3 References

Vogelstein B., and D. Gillespie. 1979. Preparative and analytical purification of DNA from agarose. Proc. Natl. Acad. Sci. USA 76: 615–619.

8.4 Product use restriction/warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety instructions from the respective Material Safety Data Sheet of the product.

This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY NAGEL literature. No other statements or representations, written or oral, by MACHEREY NAGELS employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.


Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:
MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-333
support@mn-net.com

Patents / Trademarks:

CR is patented by Roche Diagnostics
NucleoTrap[®] is a trademark of MACHEREY-NAGEL GmbH & Co KG



MACHERY-NAGEL

www.mn-net.com

MACHERY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

