REF 945006 / 945007

12.23

de

BioFix® Lumi Leuchtbakterien

nach DIN EN ISO 11348-3: 2009-05, DEV L53

Methode:

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in accordance with DIN EN ISO 11348-3: 2009-05, DEV L53

Method:

Determination of the acute toxicity on **freeze-dried** luminescent bacteria in accordance with DIN EN ISO 11348-3. The measuring unit is the natural luminescence of the used microorganism *Vibrio fischeri* NRRL B-11177. The test system measures the light output of the luminescent bacteria after they have been challenged by a sample and compares it to the light output of a control (reagent blank) that contains no sample. A difference in light output between the sample and the control is attributed to the effect of the sample on the organisms.

Field of application:

Communal and industrial sewage water, surface water, diluted eluates of soil, waster and abandoned polluted areas.

Measuring range:

0 – 100 % inhibition

Indication of the results:

- % **Inhibition** of the light output in the sample in comparison to an uninhibited control
- **G_i-value** (in accordance with DIN EN ISO 11348-3: **LID Value**): reciprocal value of the first dilution of a sample for which the inhibition of the light output amounts to less than 20 %.
- **EC_{xx}**: concentration of a sample which causes an inhibition of the light output of exactly xx % (e.g. EC₅₀ value: concentration of the sample which causes an inhibition of the light output of 50 %):
- **TU values** ("Toxicity units" in accordance with the definition of the United States Environmental Protection Agency U.S. EPA): 100 divided by EC50 value.

Contents:	REF 945006
sufficient for	400 determinations 20 vials with freeze-dried luminescent bacteria in accordance with DIN EN ISO 11348-3 1 flask with 250 mL „BioFix® Lumi Medium for freeze-dried luminescent bacteria“ (in accordance with DIN EN ISO 11348-3)“
	REF 945007
sufficient for	200 determinations 10 vials with freeze-dried luminescent bacteria in accordance with DIN EN ISO 11348-3 1 flask with 125 mL „BioFix® Lumi Medium for freeze-dried luminescent bacteria“ (in accordance with DIN EN ISO 11348-3)“

Safety precautions:

This test kit does not contain any hazardous substances for which labelling is obligatory. The luminescent microorganism *Vibrio fischeri* NRRL B-11177 has never appeared as a pathogen up to now. In accordance with the code of practice B006 1/92 ZH 1/346 of the employer's liability insurance association of the chemical industry *Vibrio fischeri* is assigned to the risk group 1 i.e. there is no risk for humans and vertebrates.

Required accessories:

„BioFix® Lumi Diluent“ (REF 945601), Glass cuvettes 50 x 12 mm (REF 916912)

Storage:

The vials with freeze-dried BioFix® Lumi luminescent bacteria can be stored until the expiration date indicated on the packaging at temperatures of -20 ± 2 °C. The „BioFix® Lumi Medium for freeze-dried luminescent bacteria“ can be stored in a thawed out state in the refrigerator at temperatures of between +2 °C and +8 °C until the expiration date which is printed on it.

The freeze-dried luminescent bacteria should only be removed from the freezer compartment immediately prior to their reactivation. Reactivated luminescent bacteria should be immediately used! An intermediate storage of the reactivated luminescent bacteria should only take place in the refrigerator at temperatures of between +2 °C and +8 °C.

With the advancing storage period of the reactivated luminescent bacteria there can be a reduction of the natural light intensity and a shift of the sensitivity spectrum. The refreezing of reactivated luminescent bacteria cannot be recommended and is not subject to the guarantee obligations of MACHEREY-NAGEL.

Reactivation and dosage of freeze-dried BioFix® Lumi luminescent bacteria:

The reactivation and dosage of freeze-dried BioFix® Lumi luminescent bacteria is always carried out in the same way irrespective of the type of test. Further details to the test procedures can be derived from the instructions for use of the corresponding tests and test systems.

Step 1: Reactivation

1. Remove a deep frozen vial with BioFix® Lumi luminescent bacteria from the deep-freeze compartment and the precooled flask with „BioFix® Lumi Medium for freeze-dried luminescent bacteria“ from the refrigerator.
2. As quickly as possible add **6 mL** „**BioFix® Lumi Medium for freeze-dried luminescent bacteria**“ (“shock thawing”)
3. Dissolve the luminescent bacteria by shaking the vial several times
4. Following this, **add further 5 mL** of the cooled „**BioFix® Lumi Medium for freeze-dried luminescent bacteria**“.
5. Prior to further use of the reactivated luminescent bacteria store the solution for an interim period for 5 minutes in the refrigerator at a temperature of between +2 °C and +8 °C to stabilise it.

Step 2: Dosage

1. Transfer **0.5 mL** of the **reactivated luminescent bacteria suspension** to each provided cuvette.
2. Further details to the test procedure can be derived from the respective instructions for use.

Analytical quality assurance:

The observance of certain validity criteria is required by DIN EN ISO 11348-3. The observance of the required validity criteria is guaranteed by MACHEREY-NAGEL by means of the enclosed test certificate.

For the purposes of self-regulation and the inspection of the correct functionality of the test system a control measurement can be carried out by the user on location by means of the standard solutions stated within the standard method. The necessary information which is required to do this relating to such items as standard solutions and test concentrations can be derived from the enclosed test certificate.

pH-value:

According to the DIN EN ISO 11348-3 no correction of the pH is necessary if the pH value of the sample is between 6.0 and 8.5. If the pH value of the sample is not within the stated range then it must set to a value of 7.0 ± 0.2 pH-values of less than 6.0 and greater than 8.5 can lead to pH-related inhibitions of luminescence.

Interferences:

Non or hardly soluble or volatile substances or compounds which are difficult to dissolve, or substances which react with the diluent or the test suspension or which change their state during the test can falsify the test result or make it difficult to reproduce the test results.

Strong coloured (especially red or brown colored) or turbid samples can cause losses of light due to light adsorption or scattering. In accordance with the procedure described in DIN EN ISO 11348-3, annex A, these effects can be minimized by an appropriate sample treatment (centrifugation or sedimentation).

Oxygen is required for the bioluminescence (> 3 mg/L) and can be limiting in case of samples with a strong oxygen consumption. This can finally result in an inhibition of light emission.

If the water sample contains easily biodegradable nutrients a loss of the light intensity which is not caused by contaminants can occur.

As the test organism *Vibrio fischeri* is a marine organism the use of seawater samples can lead to an increase of luminescence masking potential inhibitory effects of the sample (see DIN EN ISO 11348 annex D).

Salt concentrations of the sample higher than 30 g/L NaCl or concentrations of other compounds with similar osmolality can in conjunction with the salts in the test solutions lead to hyperosmotic effects. To prevent these effects the overall salt concentration should not exceed the osmolality corresponding to 35 g/L NaCl.

Chlorine containing samples effect the test procedure so that a dechlorination is mandatory, e.g., by addition of a 1 % sodium thiosulfate solution.

Disposal:

Luminescent bacteria and test preparations can be disposed of without difficulty via the sink. Restrictions can arise if the sample contains hazardous or toxic substances which are subject to waste disposal regulations. The user is him/herself responsible for the orderly disposal of such test preparations in accordance with the respectively valid guidelines and regulations.

Literature:

DIN EN ISO 11348-3:2009-05

Determination of the inhibitory effect of water samples upon the light emission of *Vibrio fischeri* (luminescent bacteria test); part 3; procedure with freeze-dried bacteria.

AQS Data Sheet for Quality Assurance in Water, Waste Water and Silt Tests by the Federal Work Community Water of Germany (LAWA) P-9/8

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