REF 945025

05.23

de

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BioFix® Lumi Luminescent Bacteria, liquid-dried according to DIN EN ISO 11348-2

Method:

Determination of the acute toxicity of **liquid-dried** luminescent bacteria according to DIN EN ISO 11348-2. The measuring unit is the natural light output (luminescence) of the microorganism *Vibrio fischeri* NRRL B-11177. Obtained will be the level of inhibition of the light intensity in the sample in comparison to an uninhibited control. The level of inhibition of the light intensity in the sample is a measure of its toxicity.

Area of application:

Communal and industrial waste water, surface waters, aqueous eluates from soils, waste and polluted sites

Measurement range:

0 – 100 % inhibition

Specifications of the results:

- % **Inhibition** of the light intensity in the sample in comparison to an uninhibited control
- **G_L-Value** (according to DIN EN ISO 11348-2: **LID-Value** = Lowest Ineffective Dilution): Reciprocal value of the first dilution of a sample that causes an inhibition of the light intensity of less than 20 %.
- **ECXX**: Concentration of a sample that causes an inhibition of the light intensity of exactly xx % (e.g. EC50-Value: Concentration of the sample that causes an inhibition of the light intensity of 50 %)
- **TU-Values** (Toxicity Units according to the definition of the American Environmental Protection Agency U.S. EPA): 100 divided by EC50-Value

Content of the test kit:

10 vials with liquid-dried luminescent bacteria according to DIN EN ISO 11348-2
 10 vials with „BioFix® Lumi Reactivation Solution“ according to DIN EN ISO 11348-2, portioned ready-to-use
 1 flask with ready-to-use sodium chloride standard solution

Hazard warnings:

The test kit does not contain any designated hazardous materials.
 The luminescent bacteria strain *Vibrio fischeri* NRRL B-11177 has never been identified as a pathogen. According to the data sheet B006 1/92 ZH 1/346 of the society of the chemical industry, *Vibrio fischeri* is classed in risk group 1, i.e. it does not pose any risk to humans or vertebrates.

Storage:

The test reagents can be kept at -20 ± 2 °C until the date indicated on the packaging. First remove the liquid-dried luminescent bacteria from the freezer and thaw out immediately before the reactivation.
 Reactivated luminescent bacteria should be used up within at least 4 hours! Intermediate storage of reactivated luminescent bacteria must only be in undiluted state (see step 1) in a fridge at +2 °C to +8 °C.
 Prolonged storage time of the reactivated luminescent bacteria can lead to a reduction in the natural light intensity and an alteration of the sensitivity spectrum. Re-freezing reactivated luminescent bacteria is not recommended and is not subject to the guarantee obligations of MACHEREY-NAGEL.

Reactivation and dosage of the liquid-dried BioFix® Lumi luminescent bacteria:

The **reactivation** of liquid-dried BioFix® Lumi Luminescent Bacteria is always done using the same method regardless of the test. The **dosage** of the resuspended luminescent bacteria depends on the optional test method of the luminescent bacteria test according to DIN EN ISO 11348-2 without measuring GL1 (step 2.1) or with measuring GL1 (step 2.2).

Important note:

Further details for performing luminescent bacteria toxicity tests can be found in the instructions and manuals of the luminometer used!

Step 1: Reactivation

1. Remove a deep frozen vial of „BioFix® Lumi Reactivation Solution“, thaw out, and maintain a temperature of **+15 °C** in a suitable thermoblock.
2. Remove a deep frozen vial with **BioFix® Lumi luminescent bacteria** from the freezer immediately before the reactivation and **add 1 mL „BioFix® Lumi Reactivation Solution“**
3. Defrost and dissolve the luminescent bacteria in a water bath (room temperature) for **2 min.**
4. Then place the vial with the resuspended luminescent bacteria in a suitable thermoblock and maintain a temperature of **+15 °C for 15 min.**

Step 2.1: Dosage for measurement without GL1

1. **Mix** the resuspended luminescent bacteria **with the remaining „BioFix® Lumi Reactivation Solution“**.
2. Transfer **0.5 mL** of the **reactivated luminescent bacteria suspension** to each provided cuvette.
3. Maintain at +15 °C for 15 min.
4. Start the test by **measuring the initial light intensity I₀**.

Details for further performance of the test procedure has to be found in the instructions and manuals of the luminometer used.

or

Step 2.2: Dosage for measurement with GL1

1. Add 2 mL of „BioFix® Lumi Reactivation Solution“ to the resuspended luminescent bacteria and mix well.
2. Place **0.2 mL** of this suspension in the **test cuvettes for GL1 calculation** (control and sample mixtures).
3. Mix the remaining resuspended luminescent bacteria with the remaining „BioFix® Lumi Reactivation Solution“.
4. In the remaining test cuvettes provided, transfer **0.5 mL** of the **reactivated luminescent bacteria suspension**.
5. Maintain at **+15 °C for 15 min.**
6. Start the test by **measuring the initial light intensity I₀**.

Details for further performance of the test procedure has to be found in the instructions and manuals of the luminometer used.

Analytical quality assurance:

In the DIN EN ISO 11348-2, it is required that specific criteria are observed for validity. The accompanying MACHEREY-NAGEL test certificate guarantees that the required criteria have been observed.

To control and check the correct functionality of the test system, the user may also carry out a control measurement on-site with the standard solutions listed in the standard. The information such as standard substances and test concentrations necessary to do this can be found in the accompanying test certificate.

A test independent of the standard according to DIN EN ISO 11348-2 can be carried out using the accompanying ready-to-use sodium chloride standard solution. Mixing 0.5 mL of this standard solution with 0.5 mL resuspended luminescent bacteria (in accordance with step 2.1) results in a 40 % to 60 % inhibition of the light intensity after 30 min incubation at + 15 °C.

pH-Value:

According to the DIN EN ISO 11348-2, no pH correction is necessary with a sample pH value of 6.0 to 8.5. If the sample's pH value is outside this range, it must be adjusted to a value of 7.0 ± 0.2. pH values below 6.0 or greater than 8.5 can lead to pH induced light inhibitions

Interferences:

Non-dissolved, difficult to dissolve or volatile compounds or substances, which react with the diluent or the test suspension, or change their state during the test, can falsify the test result or affect the ability to reproduce the test results.

With strong coloured or cloudy samples, light losses can occur due to absorption or light scattering effects. By using colour correction cuvettes (REF 940006), these interferences can be compensated following the method described in the standard DIN EN ISO 11348-2, Appendix A. In very oxygen deficient samples, the oxygen necessary for the bioluminescence (> 0.5 mg/L) may be lacking. This can cause light inhibitions, which are not caused by toxicity.

If the water sample is exposed to organic substances with easily bio-degradable nutrients, a loss in light intensity may occur which is not caused by toxic matter. NaCl concentrations in the sample of less than 15 g/L (1.5 %) or greater than 50 g/L (5 %) or their osmolarity equivalents can cause osmotically induced light inhibitions.

Disposal:

It is perfectly safe to dispose of luminescent bacteria and test mixtures down the drain. There may be restrictions if the sample itself contains substances harmful to health or are poisonous and subject to specific disposal requirements. The user is personally responsible for the proper disposal of such test mixtures according to the relevant guidelines and regulations.

Literature:

DIN EN ISO 11348-2: Water quality - Determination of the inhibitory effect of water samples on the light emission of *Vibrio fischeri* (Luminescent bacteria test); Part 2: Method using liquid-dried bacteria AQS Data Sheet for Quality Assurance in Water, Waste Water and Silt Tests by the Federal Work Community Water of Germany (LAWA) P-9/5

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