


MACHEREY-NAGEL

User manual



■ Porablot PVDF

December 2023/ Rev. 04

Contact MN

Germany and international

MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0
Toll-free: 0800 26 16 000 (Germany only)
E-mail: info@mn-net.com

Technical Support Bioanalysis

Tel.: +49 24 21 969-333
E-mail: support@mn-net.com

USA

MACHEREY-NAGEL Inc.
924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA
Toll-free: 888 321 6224 (MACH)
E-mail: sales-us@mn-net.com

France

MACHEREY-NAGEL SAS
1, rue Gutenberg – BP135 · 67720 Hoerdt Cedex · France
Tel.: +33 388 68 22 68
E-mail: sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €
Siret 379 859 531 00020 · RCS Strasbourg B379859531 · N° intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG
Hirsackerstr. 7 · 4702 Oensingen · Switzerland
Tel.: +41 62 388 55 00
E-mail: sales-ch@mn-net.com

Table of contents

1	Introduction	4
1.1	About this user manual	4
2	Membrane handling and preparation	5
3	Buffers and solutions	6
4	Safety instructions	6
4.1	Disposal	6
5	Protein transfer	7
5.1	Electrotransfer of proteins	7
5.2	Protein dot blots	8
6	Detection	9
6.1	Blocking	9
6.2	Binding of antibodies	9
6.3	Detection by ¹²⁵ I-Labeled Antibodies	10
6.4	Immuno detection	10
7	Appendices	11
7.1	References	11
7.2	Appendix 1	11
7.3	Appendix 2	12
7.4	Appendix 3	13
8	Ordering information	14
9	Warranty	15

1 Introduction

Porablot PVDF is a naturally hydrophobic fluoropolymer membrane especially designed for Western transfer, protein binding assay, and protein sequencing applications. This pure white microporous solid phase support exhibits strong hydrophobic interactions with a wide range of proteins. The immobilized proteins can be used directly for protein sequencing or amino acid analysis and can be visually detected with all common staining reagents including Coomassie Blue, Amido Black, Ponceau S, and Indian Ink. Due to the hydrophobic nature of **Porablot PVDF** membrane it should be prewetted with methanol before equilibrating with the appropriate transfer buffer.

Electrophoretic transfer of proteins from a polyacrylamide gel to a membrane was first described by Towbin*. The protein (Western) transfer technique is a simple, rapid and sensitive method which is used to identify individual proteins in complex mixtures. The procedure can be divided into three parts:

- (I) Polyacrylamide gel electrophoresis of protein molecules
- (II) Electrotransfer of proteins from the polyacrylamide gel to membranes
- (III) Detection of electrophoretically transferred proteins

The versatility of the protein transfer procedure also allows a variety of gel systems to be used together with a wide range of detection systems making this method useful for a number of applications. In this guide we provide general procedures for the electrotransfer and detection of proteins to **Porablot PVDF** membranes and also include a section on protein blots.

The following procedures have been developed to take advantage of the unique properties of these membranes, and have given optimum results in our hands.

For research use only.

1.1 About this user manual

It is recommended to read the instructions of this user manual carefully before use. All technical literature is also available on the internet at www.mn-net.com.

Please contact Technical Service regarding information about changes of the current user manual compared to previous or updated revisions.

* Towbin, H. Staehelin, T., and Gordon, J. (1979), Proc. Nat. Acad.Sci. USA. 76:4350

2 Membrane handling and preparation

This membrane is mechanically very strong and resistant towards tearing or cracking, removal from the gel is particularly easy. The membrane should only be handled by the edges using gloves or forceps to prevent membrane contamination. Either scissors or a sharp scalpel must be used to cut the membrane.

The hydrophobic nature of **Porablot PVDF** membrane requires that it is prewet with alcohol before equilibrating with transfer buffer; the procedure is as follows: (a modified procedure is used for dot blots).

-
- 1 Cut the membrane to the desired size then lay it on the surface of 80–100 % (v/v) methanol or ethanol and leave for 3–5 seconds then fully submerge. The membrane will become translucent as it wets.
 - 2 Rinse the membrane with distilled water, totally immerse it in the transfer buffer and equilibrate for 5 minutes.
 - 3 The membrane should not be allowed to dry until blocked with a protein containing solution which renders the membrane hydrophilic.
-

3 Buffers and solutions

Electrophoresis Sample Buffer for proteins:

0.125 M Tris-HCl, pH 6.8 containing: 2 % (w/v) (SDS); 10 % (v/v) glycerol; 2 % (v/v) 2-mercaptoethanol; 0.01 % (w/v) Bromophenol Blue

Transfer Buffer:

25 mM Tris base, 0.192 M glycine, 20 % (v/v) methanol, pH 8.3

PBS (pH 7.4):

40 mM disodium hydrogen phosphate; 8 mM sodium dihydrogen phosphate, 150 mM sodium chloride

Casein Blocking Solution:

0.5 % (w/v) casein and 0.05 % Tween 20 in PBS

Heat the solution to 60 °C while stirring to dissolve the casein. Cool and filter through a 0.45 µm membrane filter to remove aggregates. Use a freshly prepared solution.

4 Safety instructions

When working with the **Porablot PVDF** kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles). For more information consult the appropriate Material Safety Data Sheets (MSDS available online at www.mn-net.com/msds).



4.1 Disposal

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

5 Protein transfer

We recommend to immobilize proteins on the membrane via electrotransfer. Slot blotting or spotting onto PVDF membranes can be performed on pre wetted membranes, but is not recommended. The optimum amount of protein for detection usually varies between 1 and 10 µg per band.

5.1 Electrotransfer of proteins

Note: This procedure should be carried out in conjunction with manufactures instructions for the use of a given electrotransfer apparatus. This is especially important, if a semi-dry blotting device is used.

Tank blotting

Use the whole gel or cut the gel into sections to be transferred.

- 1 Cut membrane to the size of gel or gel sections.

- 2 Gel treatment: Equilibrate gel in transfer buffer. Gels may be stained before transfer with Coomassie Blue, or after transfer with Fast Green, Amido Black, or any other appropriate stain. In this case soak the gel for one hour in the transfer buffer.

- 3 Assemble each gel section into a "sandwich" as follows:
 1. Saturate blotting papers (e.g., MN 218 B) with transfer buffer.
 2. Place gel on one sheet of filter paper.
 3. Ensure that the membrane is saturated with transfer buffer.
 4. Lay the wet membrane on the gel.
NOTE: Roll a clean pipet over the paper and membrane to remove trapped air bubbles.
 5. Place the other sheet of wet filter paper over the membrane. Complete the „sandwich“ by placing between two blotting papers (e.g., MN 218 B). Secure „sandwich“ between two plastic grids.
NOTE: Operation of the cooling system is necessary in order to optimise transfer results.

- 4 Fill electrotransfer apparatus with transfer buffer. Insert assembly into transfer apparatus with membrane positioned between the gel and anode. Connect heat exchanger to tap water.

- 5 Connect the tank to the power supply. Start the transfer following the manufacture's recommendations for current and transfer time.
NOTE: Large proteins (MW > 100 kD) may require longer transfer times.

- 6 The membrane can be rinsed in transfer buffer to remove excess gel fragments after transfer.

5.2 Protein dot blots

- 1 Cut the required membrane to the desired size taking appropriate precautions to prevent membrane contamination. If required use a pencil to draw a grid of 1 cm squares to act as a guide when loading samples.
 - 2 Prewet the membrane by laying it on the surface of 80 % (w/v) methanol in distilled water then totally submerge; leave the membranes in the 80 % methanol until prepared to proceed with step 3.
 - 3 Remove the membrane from the 80 % methanol and drain off the excess so that the membrane remains damp. Place the damp membrane on absorbent paper which has been moistened with 80 % methanol. Using a micropipette immediately apply the protein spots using a volume of 1 μ L with a protein concentration of between 1 – 10 mg/mL. After the protein spots have been applied the membrane can either be placed directly in the blocking solution or allowed to dry on a non-absorbent surface for 5 minutes at room temperature then briefly rewetted in 80 % methanol before placing in the blocking solution. Proceed with the blocking and staining procedure (section 5).
-

6 Detection

After electrotransfer, unlabeled proteins bound to the membrane may be detected by a variety of methods including immunological detection procedures which utilise ^{125}I labeled or enzyme conjugated antibody. If "direct" staining of all the separated protein is required then **Porablot PVDF** membrane can be stained with either Coomassie Blue or Indian Ink following the procedure given in Appendix 1.

6.1 Blocking

This procedure recommends the use of 0.5 % casein and 0.05 % Tween 20 in PBS as the blocking agent. This has provided consistently good results in both monoclonal and polyclonal antibody detection systems. It should be emphasised that the dilution of detecting antibodies (primary antibody and antibody enzyme conjugates) will affect the level of background. If background is still a problem this can be further reduced by diluting the detecting antibody in 0.1 to 0.5 % casein in PBS and filtering the dilution through a 0.45 μm membrane filter. Alternative blocking agents such as 5 % (w/v) bovine serum albumin in PBS 5 % (w/v), gelatine in PBS or 5 % nonfat dry milk may be tried.

-
- 1 Remove membrane from gel surface and place into a heat sealable plastic bag or an appropriate container.
-
- 2 Add at least 10 mL of 0.5 % casein or other appropriate blocking solution per 100 cm^2 membrane surface area to the bag or container. Incubate at room temperature (or up to 40 °C) for 30–60 minutes with constant agitation.
-

6.2 Binding of antibodies

-
- 1 Remove the blocking solution from the plastic bag or container.
-
- 2 Dilute first antibody in PBS which may contain 0.1–0.5 % casein. Refer to dilution recommended by the manufacturer or optimise dilution by experiment. If ^{125}I labelled antibodies are used, the counts should not exceed 5×10^4 cpm/mL.
-
- 3 Place 2 mL of antibody solution per 100 cm^2 of membrane in the bag or container.
-
- 4 Agitate the bag or container on a rotary shaker (250 rpm) for 30 min to 1 hour at room temperature, then remove unbound antibody as follows:
-
- 5 Remove membrane from the bag or container.
-
- 6 Briefly dip the membrane in PBS containing 0.05 % (v/v) Tween 20.
-
- 7 Place membrane in a larger plastic bag or container. Add 100 mL PBS, 0.05 % Tween 20 per 100 cm^2 membrane surface area.
-
- 8 Agitate bag or container on a rotary shaker (100 rpm) for 5 minutes at room temperature, then discard the buffer.
-
- 9 Repeat above steps (7 and 8) three times with the final wash being PBS alone.
-

6.3 Detection by ¹²⁵I-Labeled Antibodies

This detection procedure has previously been described by Gershoni*.

- 1 Blot membrane dry.
 - 2 Autoradiograph to a two-day exposure using an intensifying screen, at -70°C.
-

6.4 Immuno detection

Binding of secondary antibodies

- 1 Dilute the appropriate antibody-enzyme conjugate to the manufacturer's recommended dilution in PBS containing 0.05 % Tween 20 which may contain 0.1 – 0.5 % casein.
 - 2 Place the membrane in a fresh bag or container.
 - 3 Place 2 mL of conjugate solution per 100 cm² membrane into the bag.
 - 4 Agitate the bag or container on a rotary shaker (250 rpm) for 1 hour at room temperature.
 - 5 Wash 2 times 5 min with PBS containing 0.05 % Tween 20.
-

Detection

- 1 Prepare the substrate solution (see Appendix 3).
- 2 Equilibrate the membrane with detection buffer.
- 3 Immediately add the substrate solution to the membrane in the dish. (Use 2 mL of substrate solution per 100 cm² of membrane area)
- 4 Gently agitate for 1 – 2 minutes or until bands become visible on the membrane.
- 5 Wash the membrane in distilled water, blot gently and dry at 80°C for 1 – 2 minutes. Store protected from the light.

NOTE: If a permanent record is desired, densitometry reading or a photograph should be taken as soon as possible due to instability of the chromophore.

* Gershoni, J.M. and Palade, G.E. (1982), Anal. Biochem. 124:396

7 Appendices

7.1 References

- Gershoni, J.M. and Palade, G.E. (1982), *Anal. Biochem.* 124:396
- Towbin, H. Staehelin, T., and Gordon, J. (1979), *Proc. Nat. Acad.Sci. USA.* 76:4350

7.2 Appendix 1

Direct staining of proteins on Porablot PVDF membrane

Reagents:

<i>Coomassie Blue R</i>	Prepare a 0.1 % (w/v) methanol, 10 % (v/v) acetic acid in deionised water. Destain: 10 % (v/v) acetic acid, 40 % (v/v) methanol in deionised water.
<i>Amido Black</i>	Prepare a 0.5 % (w/v) solution in 50 % (v/v) methanol, 5 % (v/v) acetic acid in deionised water. Destain: 5 % (v/v) acetic acid, 50 % (v/v) methanol in deionised water.
<i>Ponceau S</i>	Dilute 20 mL of concentrate to 200 mL with deionised water. Destain: 5 % (v/v) acetic acid in deionised water
<i>Black Indian Ink</i>	Prepare a 0.05 % (v/v) solution in the following buffer: 0.15 M NaCl, 10 mM Na ₂ HPO ₄ , pH 7.2 containing 0.5 % (v/v) Triton X-100 (PBS-TRITON).

Staining Procedures:

<i>Coomassie Blue</i>	After transfer incubate the membrane for 30 minutes using 40 mL stain/100 cm ² of membrane; destain for 15 minutes or until background is eliminated.
<i>Amido Black</i>	After transfer incubate the membrane for 3 – 15 minutes using 40 mL stain / 100 cm ² of membrane; destain for 1 – 10 minutes.
<i>Ponceau S</i>	After transfer incubate the membrane for 10 minutes using 40 mL stain/100 cm ² of membrane; destain for 1 – 10 minutes.
<i>Indian Ink</i>	After transfer wash the membrane twice for 5 minutes in PBS-TRITON; incubate the membrane in the diluted Indian Ink for 20 minutes; destain using several changes of deionised water.

7.3 Appendix 2

Example antibody detection systems

Monoclonal antibody systems

- 1 Sample protein. Load 0.5 μ L of a 2.5 % (w/v) total rat brain homogenate in sample buffer per 5 mm well, on a 10 % SDS-PAGE gel.

 - 2 Primary antibody. Mouse monoclonal antibody – anti- β -tubulin diluted 1 / 100 or 1 / 1000 in blocking solution. After transfer, incubate the membrane with primary antibody for 1 hour at RT.

 - 3 Detection: directly linked polyclonal enzyme conjugate
 - (a) Goat anti-mouse IgG-alkaline phosphatase diluted 1 / 1000 incubated with membrane for 30 minutes at room temperature.or
 - (b) Goat-anti-mouse IgG peroxidase diluted 1 / 1000 incubated with membrane for 30 minutes at room temperature.
-

Detection: two stage system:

Rabbit anti-mouse immunoglobulin (Ig) diluted 1 / 1000 incubated with membrane for 30 minutes at room temperature.

- (a) Goat anti-rabbit Ig alkaline phosphatase linked diluted 1 / 1000 for 30 minutes at room temperature. or
- (b) Goat anti-rabbit Ig peroxidase linked diluted 1 / 1000 for 30 minutes at room temperature.

Polyclonal antibody system

- 1 Sample protein. Load 5 μ l of human serum albumin (HSA) in sample buffer to give 0.5 or 1 μ g per 5 mm well, on a 10 % SDS-PAGE gel.

 - 2 Primary antibody. Non-affinity purified Goat anti-HSA polyclonal antibody diluted 1 / 10,000 in blocking solution. After transfer, incubate the membrane with primary antibody for 1 hour at room temperature.
-

Detection:

Directly linked polyclonal enzyme conjugate

- (a) Rabbit anti-goat peroxidase linked diluted 1 / 1000 incubated with membrane for 1 hour* at room temperature. or
- (b) Rabbit anti-goat alkaline phosphatase linked diluted 1 / 1000 incubated with membrane for 1 hour* at room temperature.

* 30 minutes can be used.

7.4 Appendix 3

Example substrate systems

Alkaline phosphatase

Nitro Blue Tetrazolium (NBT)-5-bromo-4-chloro-3-indolyl phosphate (BCIP) substrate

-
- | | |
|---|---|
| 1 | NBT stock solution: 75 mg/mL NBT in 70 % (v/v) dimethylformamide (DMF) |
| 2 | BCIP stock solution: 50 mg/mL of BCIP in 70 % (v/v) DMF |
| 3 | Tris substrate buffer: 0.1 M Tris, 0.1 M NaCl, 50 mM MgCl ₂ , pH 8.5 |
-

Immediately before use, add:

- 33 µL NBT stock solution and
- 25 µL BCIP stock solution to
- 7.5 mL Tris substrate buffer

Peroxidase Diminobenzidene (DAB)

DAB stock solution: 0.278 g DAB (1.3 mM) per 100 mL in 50 mM phosphate buffer pH 7.2. Store frozen in aliquots. Dilute 1 in 10 before use.

Immediately before use add 6 µL H₂O₂ (30 % v/v) per 10 mL diluted DAB solution giving a final concentration of 0.02 % (v/v).

Alternative substrates:

- (a) 3-Amino-9-ethylcarbazole (AEC)
- (b) 4-chloronaphthol (4CN)

8 Ordering information

Description	Size	Pack of	Pore size [μm]	REF
Nitrocellulose membranes				
Porablot NCP	0.3 × 3 m	1 roll	0.45	741280
Nitrocellulose membranes with supporting tissue				
Porablot NCL	0.3 × 3 m	1 roll	0.45	741290
Porablot NCL	200 × 200 mm	10 sheets		741291
PVDF membranes				
Porablot PVDF	0.25 × 3 m	1 roll	0.20	741260

9 Warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety instructions from the respective Material Safety Data Sheet of the product.


This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY-NAGEL literature. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY-NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.


Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:
MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-333
support@mn-net.com



Plasmid DNA
Clean up
RNA
DNA
Viral RNA and DNA
Protein
High throughput
Accessories
Auxiliary tools



MACHERY-NAGEL

www.mn-net.com

MACHERY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

