REF 945002 / 945003

04.23

de

BioFix® Lumi Leuchtbakterien

nach DIN EN ISO 11348-3: 2009-05, DEV L53

Methode:

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in accordance with DIN EN ISO 11348-3: 2009-05, DEV L53

Method:

Determination of the acute toxicity on **freeze-dried** luminescent bacteria in accordance with DIN EN ISO 11348-3. The measuring unit is the natural luminescence of the used microorganism *Vibrio fischeri* NRRL B-11177. The test system measures the light output of the luminescent bacteria after they have been challenged by a sample and compares it to the light output of a control (reagent blank) that contains no sample. A difference in light output between the sample and the control is attributed to the effect of the sample on the organisms.

Field of application:

Waste water, aqueous extracts and leachates, fresh water (surface water and ground water), seawater, brackish water, eluates of sediments (fresh water brackish water and seawater), interstitial pore water, and single substances dissolved in water.

Measuring range:

0 – 100% inhibition

Indication of the results:

- % **Inhibition** of the light output in the sample in comparison to an uninhibited control
- **G_{ib}-value** in accordance with DIN EN ISO 11348-3, annex B) : reciprocal value of the first dilution of a sample for which the inhibition of the light output amounts to less than 20 %.
- **EC₅₀**: concentration of a sample which causes an inhibition of the light output of exactly xx % (e.g. EC50 value: concentration of the sample which causes an inhibition of the light output of 50 %):
- **TU-values** "Toxicity units" in accordance with the definition of the United States Environmental Protection Agency U.S. EPA): 100 divided by the EC50 value.

Contents:	REF 945002
sufficient for	2000 determinations 20 vials with freeze-dried luminescent bacteria in accordance with DIN EN ISO 11348-3 1 flask with 25 mL „BioFix® Lumi Reconstitution solution“
	REF 945003
sufficient for	1000 determinations 10 vials with freeze-dried luminescent bacteria in accordance with DIN EN ISO 11348-3 1 flask with 15 mL „BioFix® Lumi Reconstitution solution“

Safety precautions:

This test kit does not contain any hazardous substances for which labelling is obligatory. The luminescent microorganism *Vibrio fischeri* NRRL B-11177 has never appeared as a pathogen up to now. In accordance with the code of practice B006 1/92 ZH 1/346 of the employer's liability insurance association of the chemical industry *Vibrio fischeri* is assigned to the risk group 1 i.e. there is no risk for humans and vertebrates.

Required accessories:

BioFix® Lumi Medium for freeze-dried luminescent bacteria" (REF 945608), glass cuvettes 50 × 12 mm (REF 916912)

Storage:

The vials with freeze-dried BioFix® Lumi luminescent bacteria can be stored until the expiration date indicated on the packaging at temperatures of -15 to -21 °C. The "BioFix® Lumi Reconstitution solution" can be stored in an unfrozen state in the refrigerator at temperatures of between +2 °C and +8 °C until the expiration date (see label).

The freeze-dried luminescent bacteria should only be removed from the freezer compartment immediately prior to their reactivation. Reactivated luminescent bacteria should be used within 4 hours! An intermediate storage of the reactivated luminescent bacteria in the refrigerator at temperatures between +2 °C and +8 °C is only acceptable **after dilution (see Step 2 of Procedure B)**. Upon storage of the reactivated luminescent bacteria a reduction of the natural light emission and a shift of the sensitivity spectrum can be possible. It is not recommended to refreeze the reactivated luminescent bacteria and is not subject to the guarantee obligations of MACHEREY-NAGEL.

Reactivation and dosage of freeze-dried BioFix® Lumi luminescent bacteria:

The reactivation and dosage (procedure A or procedure B) of freeze-dried BioFix® Lumi luminescent bacteria is always carried out in the same way irrespective of the type of test

Important Note:

Further details to the test procedures can be derived from the instructions for use of the corresponding tests, test systems and handbooks of the luminometers used.

Step 1: Reactivation

1. Remove a deep frozen vial with BioFix® Lumi luminescent bacteria from the deep-freeze compartment and the precooled flask with „BioFix® Lumi Reconstitution solution“ from the refrigerator.
2. As quickly as possible **pour 1 mL** precooled (+2 °C to +8 °C) „BioFix® Lumi Reconstitution solution“ to the BioFix® Lumi luminescent bacteria („shock thawing“).
3. Dissolve the luminescent bacteria by shaking the vial several times.
4. Fill the solution of reactivated luminescent bacteria in an appropriate glass cuvette (50 × 12 mm, REF 916912) and prior to further use store the solution for an interim period for 5 minutes at a temperature of +2 °C to +8 °C to stabilise it (e. g., in the Reagent well of the luminometer Microtox® Modell 500).

Step 2: Dosage

- Procedure A:**
1. Transfer 0.5 mL of precooled „BioFix® Lumi Medium for freeze-dried luminescent bacteria“ (REF 945608) to each provided cuvette.
 2. After that add **10 µL reactivated, undiluted luminescent bacteria suspension** to each cuvette.
 3. Incubate for **15 min** at +15 °C.
 4. Start the test by **measuring the initial luminescence I₀**.

Important! After the measurement of the initial luminescence and prior to the addition of the control or sample solution it is important to keep the cuvettes at +15 °C.

Further details to the test procedure can be derived from the respective instructions for use of the luminometers used.

- Procedure B:**
1. Mix in a beaker of appropriate size **reactivated, undiluted luminescent bacteria suspension** and precooled (+2 °C to +8 °C) „BioFix® Lumi Medium for freeze-dried luminescent bacteria“ (REF 945 608) with a **ratio of 1 + 50** (e. g., 0.5 mL reactivated luminescent bacteria plus 25 mL medium)
 2. Transfer **0,5 mL** of this **diluted, reactivated luminescent bacteria suspension** to each provided cuvette.
 3. Incubate for **15 min** at +15 °C temperieren.
 4. Start the test by **measuring the initial luminescence I₀**.

Important! After the measurement of the initial luminescence and prior to the addition of the control or sample solution it is important to keep the cuvettes at +15 °C.

Further details to the test procedure can be derived from the respective instructions for use of the luminometers used.

Analytical quality assurance:

The compliance of certain validity criteria is required by DIN EN ISO 11348-3. The compliance is guaranteed by MACHEREY NAGEL by means of the enclosed test certificate.

According to the above mentioned norm the sensitivity of the luminescent bacteria has to be controlled by the end user. The required information such as reference substances of concentrations can be derived from the enclosed test certificate

pH-value:

According to the DIN EN ISO 11348-3 no correction of the pH is required if the pH value of the sample is between 6.0 and 8.5. It may be noted that pH-values of less than 6.0 and greater than 8.5 can lead to pH-related inhibitions of luminescence.

Interferences:

Non or hardly soluble or volatile substances or compounds which are difficult to dissolve, or substances which react with the diluent or the test suspension or which change their state during the test can falsify the test result or make it difficult to reproduce the test results.

Strong coloured or turbid samples can cause losses of light due to light adsorption or scattering. In accordance with the procedure described in DIN EN ISO 11348-3, annex A, these effects can be minimized by an appropriate sample treatment (centrifugation or sedimentation) or by using colour correction cuvettes (REF 940006).

Oxygen is required for the bioluminescence (> 3 mg/L) and can be limiting in case of samples with a strong oxygen consumption. This can finally result in an inhibition of light emission.

If the water sample contains easily biodegradable nutrients a loss of the light intensity which is not caused by contaminants can occur.

As the test organism *Vibrio fischeri* is a marine organism the use of seawater samples can lead to an increase of luminescence masking potential inhibitory effects of the sample (see DIN EN ISO 11348 annex D).

Salt concentrations of the sample higher than 30 g/L NaCl or concentrations of other compounds with similar osmolality can in conjunction with the salts in the test solutions lead to hyperosmotic effects. To prevent these effects the overall salt concentration should not exceed the osmolality corresponding to 35 g/L NaCl.

Chlorine containing samples effect the test procedure so that a dechlorination is mandatory, e.g., by addition of a 1 % sodium thiosulfate solution.

Disposal:

Luminescent bacteria and test preparations can be disposed of without difficulty via the sink. Restrictions can arise if the sample contains hazardous or toxic substances which are subject to waste disposal regulations. The user him/herself is responsible for the orderly disposal of such test preparations in accordance with the respectively valid guidelines and regulations.

Literature:

DIN EN ISO 11348-3:2009-05

Determination of the inhibitory effect of water samples upon the light emission of *Vibrio fischeri* (luminescent bacteria test); part 3; procedure with freeze-dried bacteria.

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