

Customer Information

14.10.2022, Version 1, FHS

Importing the test NANOCOLOR® ECO Iron to the NANOCOLOR® VIS II and UV/VIS II


Overview

Photometer	NANOCOLOR® VIS II and UV/VIS II
Availability	14.10.2022
Test numbers	4-02 (16°mm), 4-52 (50 mm)
Required accessories	USB drive

Important note

Although importing of the method files does not affect the data and the costumers special methods on the instrument, MACHEREY-NAGEL recommends precautions to avoid data loss. Perform a backup of the instrument database before the import of the method file, via *Main Menu* → *Settings* → *Backup* → *USB drive* → *Backup*. Furthermore, it is recommended to save user defined methods via *Methods* → *Special Methods* → *User Methods* → *List* → *Export* to a USB drive prior to the import process. Detailed data backup instructions can be found in the manual of the respective photometer.

Test import

- Download the method ZIP-file via the link above and extract the method files on your computer.
- Copy the .mtd files to the upper level of an USB drive.
- Connect the USB drive to the photometer.
- Enter the list of the user methods via *Methods* → *Special Methods* → *User Methods* → *List*.
- Select the -icon and import the method file for the respective photometer.
- The successful import is indicated by a message on the screen and must be confirmed with *OK*. ATTENTION: The test is not displayed in the open list view after importing!
- The NANOCOLOR® ECO Manganese test kit can be selected by entering the number 402 (round tube) or 452 (rectangular cuvette) in the *Test number* menu or via inserting of the barcoded test tubes into the photometer.

Modifications

- Import the method files VIS_II_4_02.mtd and VIS_II_4_52.mtd for the implementation of the NANOCOLOR® ECO Iron (REF 976002) on the NANOCOLOR® VIS II.
- Import the method file UVVIS_II_4_02.mtd and UVVIS_II_4_52.mtd for the implementation of the NANOCOLOR® ECO Iron (REF 976002) on the NANOCOLOR® UV/VIS II.

Contact

For further questions, please do not hesitate to contact us:

Technical Support and Customer Service	Tel.:	+49 24 21 969-333
	E-Mail:	support@mn-net.com
For orders (international)	E-Mail:	sales@mn-net.com