

Customer information

23rd of September 2021, Version 3, RaKr

Firmware update

Spectrophotometers *NANOCOLOR*[®] VIS and *UV*/*VIS*

Overview	
Instrument	Spectrophotometer <i>NANOCOLOR</i> [®] VIS and <i>UV</i> / <i>VIS</i>
Firmware version	6.0.8 (<i>NANOCOLOR</i> [®] VIS) 8.0.8 (<i>NANOCOLOR</i> [®] <i>UV</i> / <i>VIS</i>)
Price	free of charge
Availability	from 23 rd of September 2021
Available on	https://www.mn-net.com/de/nanocolor-device-updates
Required accessories	USB stick (REF 919119)



Important notes

For *NANOCOLOR*[®] *UV*/*VIS* we do recommend an update of your instrument, if you are using the firmware version 8.0.0 or 8.0.1! A bug in the data memory can lead to a wrong documentation of the chemical form (e.g. PO₄-P) of measurements.

The firmware update for the photometer *NANOCOLOR*[®] *UV*/*VIS* and VIS will perform an automatic system reset of the device when updating from version 7.X.X or 5.X.X, respectively. That means all measurement results and user defined methods will be lost during the update. Make sure to save these data with the help of our *NANOCOLOR*[®] PC software prior to the update.

If you use the *NANOCOLOR*[®] PC Software for spectrophotometers together with your photometer, a device update from versions < 6.0.0 (*NANOCOLOR*[®] *UV*/*VIS*) and < 4.0.1 (*NANOCOLOR*[®] VIS) also requires an update of the PC software to version 7.0.0.0. Otherwise, the data export from the photometer to the PC is not possible.

After the successful update of the photometer, a zero calibration is required. Use the calibration cuvette delivered with the photometer and follow the display instructions for zero calibration.

Update procedure

- Download the update ZIP file via the link above and unpack the ZIP file on your computer.
- Copy the .ahx update file to the upper level of a USB stick and connect the USB stick to the instrument.
- Start the instrument. The update file will be detected automatically and the update can be started by pressing the "Update" icon (The update process may take several minutes. Do not switch off the instrument during this process!).
- Follow the instructions on the screen.

Changelog

- Implementation of new test *NANOCOLOR*[®] easily liberated cyanide 04 (REF 985025)
- Bugfixes

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service

Tel.: +49 24 21 969-332

E-Mail: csc@mn-net.com

For orders (international)

E-Mail: sales@mn-net.com