

Customer information

05.09.2019, version 2, RaKr

Firmware update for filter photometers *NANOCOLOR*[®] 500 D, 400 D and PF-12

Overview

Photometer	Filter photometers <i>NANOCOLOR</i> [®] 500 D, 400 D and PF-12
Firmware version	2.70 (<i>NANOCOLOR</i> [®] 500 D) 3.69 (<i>NANOCOLOR</i> [®] 400 D) 2.68 (Compact photometer PF-12)
Price	free of charge
Availability	from 09 th September 2019
Available on	http://www.mn-net.com/tabid/11162/default.aspx
Required accessories	<ul style="list-style-type: none"> Update software for filter photometers Connection cable for the respective device



Important notes

The firmware update for the filter photometers *NANOCOLOR*[®] 500 D, 400 D and PF-12 deletes the memory of the device. Make sure to save data prior to the update with our PC-software.

Update process

Download the update software for filter photometers from our website via the above mentioned link and follow the instructions on the respective page for updating the device.

Modifications

- Including new *NANOCOLOR*[®] test kits Sulfate HR 1000 (REF 985063) and Iron LR (REF 918128)
- New sub-method 0-155 for the test kit *NANOCOLOR*[®] carbonate hardness 15 (REF 985015): 20 - 270 mg/L CaCO₃
- Bug fixes

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service

Tel.: +49 24 21 969-332

E-Mail: csc@mn-net.com

For orders (international)

E-Mail: sales@mn-net.com
www.mn-net.com


MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com