

Customer information

Rev. 0

QUANTOFIX® Relax Cleaning instructions

Introduction

Thanks to the QUANTOFIX® Relax's design, there is no direct contact between the test strip and the detector during measurement. However, the device can still be damaged by contamination over time. When measuring dairy products or other foods, residues may remain on the carrier, which can be unhygienic.

This guide will teach you how to quickly and easily clean the QUANTOFIX® Relax during everyday lab use.



Instruction

1. Turn off the device

Turn off the device without manually retracting the test strip holder.



2. Remove the test strip holder

Remove the test strip holder from the carriage by pulling it forward. This may require a little force, especially at first.



3. Clean the test strip holder

Clean the plastic holder with warm soapy water or ethanol. It can also be cleaned in the dishwasher.



4. Install the test strip holder

Slide the test strip holder onto the carriage until it is completely retracted into the housing. This may require a little force, especially at first.



5. Turn on the device

Turn on the device, and the carriage will move into its operating position as usual. The recesses in the test strip holder and the carriage should be aligned so that the sensor underneath is exposed.



Video E-Training



Watch the step-by-step instruction video here:



Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service

Tel.: +49 24 21 969 333
E-Mail: support@mn-net.com

For orders (international)

E-Mail: sales@mn-net.com