de *NANOCOLOR*® TOC 300 Compact

REF 985077

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en *NANOCOLOR*® TOC 300 Compact REF 985077

Overview

The test is suitable for the photometric determination of TOC.

The test is suitable for wastewater and production water.

Measuring range: 20 – 300 mg/L TOC (method 0771 / 0772)

- Number of tests: 19
- Wavelength for photometric determination: 436 nm
- Shelf life: 12 months
- Reaction time: 60/120 minutes
- Storage temperature: 2 – 8 °C
- Storage conditions: protected from sunlight, upright.

Method

TOC determination is carried out in two steps:

1. Removal of all inorganic carbon (TIC)

2. Digestion of total organic carbon (TOC) and detection of the resulting carbon dioxide by means of an indicator. Analogus to DIN EN 1484.

Interferences

The following contaminants do not interfere with the test up to the indicated concentrations. The cumulative effect of different interfering ions has not been tested.

Data in mg/L:

- TIC: 10000
- Cl^-: 10000
- Ca^{2+}, Mg^{2+}: 20000
- $S_2O_5^{2-}$, S^{2-}, SO_3^{2-}: 10000
- NO_2^-: 5000
- NH_4-N: 3000

The method is suitable for the analysis of seawater after 1 + 9 dilution.

Reagents and accessories

Contents of reagents set:

- 19 test tubes R1
- 19 test tubes R2
- 1 tube with blanc solution
- 19 plastic tips for TIC-Ex
- 19 threaded couplings
- 20 adhesive labels 10 × 38 mm

Required devices:

- MACHEREY-NAGEL photometer
- MACHEREY-NAGEL heating block
- Protective covering with bores (REF 919309)
- *NANOCOLOR*® TIC-Ex (REF 916993)
- Digital piston pipette 100 – 1000 µL (REF 91677) with pipette tips (REF 91667)

Standards

- *NANOCONTROL* TOC 300 (REF 92578)

Sampling and preparation

See DIN EN ISO 5667-3-A21.

Adjust to pH 1 – 12 prior to analysis.

Quality control

The measurement of a blank value and a standard is recommended before every measuring series as quality control measure.

LOT-specific certificates are available at *www.mn-net.com*.

Procedure

1. Open test tube 1. Pipette 0.3 mL of sample into test tube
2. Insert open cuvette in the TIC-Ex and remove TIC for 5 min (after removal blow out the fluid residues in the tips by briefly actuating the TIC-Ex)
3. Open cuvette 2 and immediately close it with the threaded coupling
4. Immediately screw cuvette 1 tightly with cuvette 2 (hold cuvette combination vertically. Ensure that there is no film of liquid in the neck of cuvette 1)
5. Insert cuvette combination into the heating block (blue indicator solution on top). Heat in the thermoblock for 1 hour at 100 °C [Method 0771]. For extended digestion, heat for 2 hours at 100 °C [Method 0772].
6. Remove cuvette combination from the heating block. Cool to room temperature. The cuvettes must not be placed in a water bath or held under running water to speed up cooling.
7. Clean outside of blank tube
8. Measure blank value
9. Clean outside of test tube 2
10. Invert the cuvette combination (with cuvette 2 facing down)
11. Measure

Disposal

Information regarding disposal can be found in the safety data sheet. You can download the SDS from *www.mn-net.com/SDS*.

MACHEREY-NAGEL

www.mn-net.com

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com
FR +33 388 68 22 68 sales-fr@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com
US +1 888 321 62 24 sales-us@mn-net.com

Notes

The required digestion time depends on the sample matrix and is typically 1 hour. For samples that are difficult to oxidize or digest, a digestion time of 2 hours is recommended.

To increase the accuracy, it is recommended to perform the measurement with reagent blank value. For the blank value, use COD-free water (REF 918993) as sample.

When using other photometers, make sure measurements are possible in test tubes (16 mm OD) and calibrate the method.

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from *www.mn-net.com/SDS*.

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