



# NucleoSpin® 8 / 96 Plasmid product metamorphosis

## Stream lined product group

4 instead of 7 kits with optimised contents



### Product changes

Previous REF	Previous name		New REF	New name
740620/.5	NucleoSpin® Multi-8 Plasmid	replaced by	740621/.5	NucleoSpin® 8 Plasmid
740621/.5	NucleoSpin® Multi-8 Plus Plasmid	changes to	740621/.5	NucleoSpin® 8 Plasmid
740625.1/.4/.24	NucleoSpin® Multi-96 Plus Plasmid	changes to	740625.1/.4/.24	NucleoSpin® 96 Plasmid
740730/.1	NucleoSpin® Robot-8 Plasmid	replaced by	740621/.5	NucleoSpin® 8 Plasmid
		alternative	740461.4	NucleoSpin® 8 Plasmid Core Kit*
740461.4	NucleoSpin® 8 Plasmid Core Kit	stays	740461.4	NucleoSpin® 8 Plasmid Core Kit
740708.2/.4/.24	NucleoSpin® Robot-96 Plasmid	replaced by	740625.1/.4/.24	NucleoSpin® 96 Plasmid
		alternative	740616.4	NucleoSpin® 96 Plasmid Core Kit*
740616.24	NucleoSpin® Robot-96 Plasmid Core Kit	changes to	740616.4	NucleoSpin® 96 Plasmid Core Kit

- „Multi“ means intended for manual use, „Robot“ means intended for automated use with increased buffer volumes as required by robotic platforms.
- All new kits can be used manually as well as automated on liquid handling platforms (increased buffer volumes).
- „Plus“ means Filter Strips/Plates are included.
- All new kits include Filter Strips/Plates.
- „Core kits“ do not contain accessory plates (Culture Plate, MN Wash Plate, Elution Plate) and no Buffer AW.
- Additional accessory plates and buffers can be ordered separately.
- \* please check kit contents, NucleoSpin® 8 Plasmid Core Kit and NucleoSpin® 96 Plasmid Core Kit do not contain Buffer AW.

**For detailed information about product changes:  
Please see the reverse side.**

### Technical Service

If you need any further information, please contact us  
by e-mail: [pm-htp@mn-net.com](mailto:pm-htp@mn-net.com)  
or phone: +49 24 21 969-270  
or via your local distributor.

## Detailed information about product changes

Prev. REF	Previous name	New REF	New name	Prev. contents	New contents	
740621	<b>NucleoSpin® Multi-8 Plus Plasmid</b> (12x8)	740621	<b>NucleoSpin® 8 Plasmid</b> (12x8)	Buffer A1	30 mL	<b>50 mL</b>
				Buffer A2	30 mL	<b>50 mL</b>
				Buffer A3	40 mL	<b>50 mL</b>
				Buffer AW	60 mL	<b>100 mL</b>
				Buffer A4 (Concentrate)	50 mL	<b>80 mL</b>
				Buffer AE	20 mL	<b>25 mL</b>
				RNase A	12 mg	<b>20 mg</b>
				NucleoSpin® Plasmid Bindings Strips	12	12
				NucleoSpin® Plasmid Filter Strips	12	12
				MN Wash Plate (incl. 6 paper sheets)	1	1
				Rack of Tube Strips including Cap Strips	1	1
				Culture Plate including Gas-permeable Foil	not included	<b>1</b>
				Protocol	1	1
740621.5	<b>NucleoSpin® Multi-8 Plus Plasmid</b> (60x8)	740621.5	<b>NucleoSpin® 8 Plasmid</b> (60x8)	Buffer A1	2 x 75 mL	<b>250 mL</b>
				Buffer A2	2 x 75 mL	<b>250 mL</b>
				Buffer A3	2 x 100 mL	<b>250 mL</b>
				Buffer AW	2 x 125 mL	<b>500 mL</b>
				Buffer A4 (Concentrate)	3 x 100 mL	<b>2 x 200 mL</b>
				Buffer AE	100 mL	<b>125 mL</b>
				RNase A	2 x 30 mg	<b>100 mg</b>
				NucleoSpin® Plasmid Bindings Strips	60	60
				NucleoSpin® Plasmid Filter Strips	60	60
				MN Wash Plate (incl. 6 paper sheets)	5	5
				Rack of Tube Strips including Cap Strips	5	5
				Culture Plate including Gas-permeable Foil	not included	<b>5</b>
				Protocol	1	1
740625.1	<b>NucleoSpin® Multi-96 Plus Plasmid</b>	740625.1	<b>NucleoSpin® 96 Plasmid</b>	Buffer A1	35 mL	<b>50 mL</b>
				Buffer A2	35 mL	<b>50 mL</b>
				Buffer A3	2 x 35 mL	<b>50 mL</b>
				Buffer AW	75 mL	<b>100 mL</b>
				Buffer A4 (Concentrate)	3 x 20 mL	<b>80 mL</b>
				Buffer AE	25 mL	<b>25 mL</b>
				RNase A	15 mg	<b>20 mg</b>
				NucleoSpin® Plasmid Binding Plate	1	1
				NucleoSpin® Plasmid Filter Plate	1	1
				Culture Plate including Gas-permeable Foil	1	1
				Self-adhering PE Foil	5	<b>not included</b>
				MN Wash Plate (incl. 6 paper sheets)	1	1
				Elution Plate U-bottom (incl. 1 PE Foil)	1	1
Protocol	1	1				
740625.4	<b>NucleoSpin® Multi-96 Plus Plasmid</b>	740625.4	<b>NucleoSpin® 96 Plasmid</b>	Buffer A1	150 mL	<b>200 mL</b>
				Buffer A2	150 mL	<b>200 mL</b>
				Buffer A3	200 mL	<b>200 mL</b>
				Buffer AW	240 mL	<b>400 mL</b>
				Buffer A4 (Concentrate)	200 mL	<b>2 x 160 mL</b>
				Buffer AE	75 mL	<b>2 x 100 mL</b>
				RNase A	60 mg	<b>80 mg</b>
				NucleoSpin® Plasmid Binding Plate	4	4
				NucleoSpin® Plasmid Filter Plate	4	4
				Culture Plate including Gas-permeable Foil	4	4
				Self-adhering PE foil	20	<b>not included</b>
				MN Wash Plate (incl. 6 paper sheets)	4	4
				Elution Plate U-bottom (incl. 1 PE Foil)	4	4
Protocol	1	1				
740625.24	<b>NucleoSpin® Multi-96 Plus Plasmid</b>	740625.24	<b>NucleoSpin® 96 Plasmid</b>	contents are 6x the contents of 740625.4		
740461.4	<b>NucleoSpin® 8 Plasmid Core Kit</b> (48 x 8)			Buffer A1	200 mL	no change
				Buffer A2	200 mL	
				Buffer A3	200 mL	
				Buffer A4 (Concentrate)	2 x 100 mL	
				Buffer AE	2 x 50 mL	
				RNase A	80 mg	
				NucleoSpin® Plasmid Bindings Strips	48	
				NucleoSpin® Plasmid Filter Strips	48	
				Protocol	1	
740616.24	<b>NucleoSpin® Robot-96 Plasmid Core Kit</b> (6 x 4 x 96)	740616.4	<b>NucleoSpin® 96 Plasmid Core Kit</b> (4 x 96)	Buffer A1	6 x 200 mL	200 mL
				Buffer A2	6 x 200 mL	200 mL
				Buffer A3	6 x 200 mL	200 mL
				Buffer A4 (Concentrate)	6 x (4 x 100 mL)	<b>2 x 100 mL</b>
				Buffer AE	6 x (2 x 50 mL)	2 x 50 mL
				RNase A	6 x 80 mg	80 mg
				NucleoSpin® Plasmid Binding Plate	24	4
				NucleoSpin® Plasmid Filter Plate	24	4
				MN Wash Plate (incl. 6 paper sheets)	24	<b>not included</b>
				Protocol	6	1



[www.mn-net.com](http://www.mn-net.com)

**MACHEREY-NAGEL**



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

**Germany and international:**  
Tel.: +49 24 21 96 90  
Fax: +49 24 21 96 91 99  
E-mail: [info@mn-net.com](mailto:info@mn-net.com)

**Switzerland:**  
MACHEREY-NAGEL AG  
Tel.: +41 62 388 55 00  
Fax: +41 62 388 55 05  
E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

**France:**  
MACHEREY-NAGEL EUROL  
Tel.: +33 388 68 22 68  
Fax: +33 388 51 76 88  
E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

**USA:**  
MACHEREY-NAGEL Inc.  
Tel.: +1 484 821 0984  
Fax: +1 484 821 1272  
E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

