

MACHEREY-NAGEL

Direct Lysis and PCR from human and animal samples

Bioanalysis



NucleoType Mouse PCR and NucleoType Blood PCR

- Designed for genotyping from tail, ear, tissue, blood, and even hair samples
- Easy handling with the Blood Transfer Tool
- Inhibitor Removal Pearls simplify challenging blood sample processing

30 %
Discount!!!

MACHEREY-NAGEL

www.mn-net.com



Since 1911

Fast genotyping from multiple sample types

Introducing the new NucleoType product line

NucleoType Mouse PCR

- Simple 5 minute sample preparation
- For use with tail clips, ear punches, and even humanely collected animal hair or fur
- Direct PCR protocol available



Product at a glance

Technology	Simple sample preparation suitable for hot start PCR (optional direct PCR)
Preparation time	Lysis 2 min (18–25 °C); Enzyme inactivation 3 min (98 °C); PCR cycling 30–90 min (cycler and target size dependent)


NucleoType Blood PCR

- Novel Blood Transfer Tool enables easy and fast blood transfer
- Inhibitor Removal Pearls eliminate challenges from difficult blood samples
- For use with whole blood and blood cards

Product at a glance

Technology	Direct HotStart PCR for blood sample material
Preparation time	Preparation per sample: < 1 min; PCR cycling: 30–90 min (cycler and target size dependent)

Challenging blood samples



Kit	Cat. No.	Preps	List price \$	Promo price \$*
NEW! NucleoType Blood PCR	743201.100	100	103.00	72.10
NEW! NucleoType Blood PCR	743201.500	500	395.00	276.50
NEW! NucleoType Mouse PCR	743200.100	100	73.50	51.45
NEW! NucleoType Mouse PCR	743200.500	500	292.00	204.40



* Valid until February 28, 2019, only in the USA.

Please visit our website for additional product information:

www.mn-net.com/NucleoType

Promo code: **NUCLEOTYPE**

Quote "NUCLEOTYPE" when ordering by e-mail, phone, fax or through our WebShop.

www.mn-net.com

MACHERY-NAGEL



MACHERY-NAGEL Inc. · 2850 Emrick Blvd. · Bethlehem, PA 18020 · USA

US:
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com

DE / International:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

CH:
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

FR:
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

