

NucleoMag[®] Blood 200 µL Prefilled Plates IsoPure[™] Mini

Protocol details

Application	DNA purification from blood
Kit	NucleoMag [®] Blood 200 µL Prefilled Plate
REF	744510
Protocol name	NMBlood200PFP
Instructions	Please refer to the user manual NucleoMag [®] Blood 200 µL Chapter 6 Protocol for prefilled plates



Instant protocol transfer via QR-code

Procedure	
1	Connect the scanning device (included) to the instrument
2	Activate the instrument scanning software
3	Scan the QR-code for instant protocol transfer
4	Confirm protocol transfer on the instrument



Instrumentation

Product	Specification	REF
IsoPure [™] Mini	Automated nucleic extraction system for MACHEREY-NAGEL's NucleoMag [®] kits enabling parallel processing of up to 16 samples	747000
Android [™] tablet	Android [™] tablet with IsoPure [™] Mini App for simple protocol design and transfer	747001

Disclaimer

Information

MACHEREY-NAGEL GmbH & Co. KG makes every effort to include accurate and up-to-date information within this publication; however, it is possible that omissions or errors might have occurred. MACHEREY-NAGEL GmbH & Co. KG cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. For technical details and detailed procedures of the specifications provided in this document please contact your MACHEREY-NAGEL representative. This publication may contain reference to applications and products which are not available in all markets. Please check with your local sales representative. All mentioned trademarks are protected by law. All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment). Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or support by MACHEREY-NAGEL GmbH & Co. KG. Any views or opinions expressed herein by the authors' do not necessarily state or reflect those of MACHEREY-NAGEL. NucleoMag[®] is a registered trademark of MACHEREY-NAGEL GmbH & Co. KG, Düren, Germany. IsoPure[™] is a brand of Accuris Instruments.