

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 640 mf
Basis weight	85 g/m ²
Thickness	0.19 mm / 7.5 mil
Filtration speed	27 s
Surface	Smooth
Color	White
Material	Cellulose
Ash content	< 0.01 %
Average retention capacity	4–12 µm
Characteristics / Application	Ashless, fat-free, medium fast filtration, e.g. for sulfides, lead chromate, calcium oxalate

*These data are typical values and must not be regarded as specification.