

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Ascorbic acid

Parameter*	Value*
REF	91314
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Ascorbic acid
Gradation	0 • 50 • 100 • 200 • 300 • 500 • 700 • 1000 • 2000 mg/L Vitamin C
Method	Phosphomolybdenum blue
Color change	yellow → green-blue
Instrumental measuring range	25–1000 mg/L Vitamin C

*These data are typical values and must not be regarded as specification.

