

## Data sheet

## Filter paper

| Parameter*                    | Value*                                                            |
|-------------------------------|-------------------------------------------------------------------|
| Grade                         | MN 960                                                            |
| Basis weight                  | 185 g/m <sup>2</sup>                                              |
| Thickness                     | 0.47 mm / 18.5 mil                                                |
| Filtration speed              | 14 s                                                              |
| Surface                       | Smooth                                                            |
| Color                         | White                                                             |
| Material                      | Cellulose                                                         |
| Characteristics / Application | Fast filtration, e.g. for filtration of wines, fruit juices, etc. |

\*These data are typical values and must not be regarded as specification.