

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 439
Basis weight	66 g/m ²
Thickness	0.42 mm / 16.5 mil
Filtration speed	13 s
Maximum temperature	550 °C / 1022 °F
Surface	Smooth
Color	White
Material	Borosilicate glass micro fibers
Average retention capacity	1.5 µm
Characteristics / Application	Glass fiber filter without binder, suitable for determination of suspended solids according to DIN EN 872:2005

*These data are typical values and must not be regarded as specification.